

PSM+



PRODUCT CATALOGUE 2025

PATIENT MONITORING & DIAGNOSTIC CABLES AND ACCESSORIES

Shenzhen Parts-Med Technology Co.,Ltd

Add: C, 2/F, Lilan Technology Building, Jiazitang Community, 518107 Fenghuang Street, Guangming District, Shenzhen City, Guangdong Province, China.

Web: www.parts-med.com

Email: marketing@Parts-Med.com

Copyright Statement

Copyright Statement: The copyright of this work belongs to Parts-Med technology. And has been registered for copyright protection. Without prior written consent from Parts-Med, no entity or individual is allowed to copy, modify, adapt, translate, transmit, or plagiarize all or part of this work in any way or do any act that may infringe upon the copyright of the Parts-Med. All rights reserved.



Declarations & Terms

Product Updates : The listed products might be modified or discontinued.

Quotation : Quotations are on Ex-Factory based (Shenzhen in China) in major currencies of USD or Euro. Parts-Med may adjust the prices yearly. The proper discount might be available depending on the purchasing volume and items of the individual orders.


Payment : Bank Transfer and Letter of Credit are recommended. PayPal, West Union and Credit card are also acceptable. Customers will have to pay the bank fees at both sides.

Lead time : Approximately 2-3 weeks upon the approved orders. The exact delivery date should be finally confirmed by the sales representatives depending on the quantities and business seasons.

Return & Exchange : Returning or exchanging are possible within 10 days after the goods are received, on the condition that quality problems are found and approved by Parts-Med. Problematic products could be repaired or exchanged when they are returned to Parts-Med factory, provided that failures are found within the warranty period. If no faults are found by Parts-Med or items are out of warranty, the customers should pay for the shipping cost.

Legal Affairs : The original brands or trademarks indicated in the catalogues are for the purpose of referring to the compatibility with Parts-Med products. They are proprietarily owned by the respective manufacturers. Customers of Parts-Med are advised to be aware of the patent issues Parts-Med products might be involved with in the marketed areas, and are responsible for the relevant legal affairs.

Parts-Med Tech The Provider of Medical Cables & Accessories



About Parts-Med Tech

Parts-Med Tech is a leading medical supplies manufacturer specializing in patient monitoring supplies and accessories. Parts-Med Tech have adhered to the customer-oriented service concept and dedicated ourselves to provide a one-stop shop solution at a full range, quality reliable and multi-compatible products to the customers with various patient monitoring accessories including EEG electrodes, ECG cables, SpO2 sensor, Temp probe, NIBP cuff, fetal probe, etc.. Meanwhile, Parts-Med Tech possess internationally approved certifications allowing Parts-Med Tech to enter main markets worldwide to serve customers with professional, fast and reliable services.

Accessories for multi-parameter monitors

1. ECG Monitoring Cables & Lead Wires

- 1. 1 ECG Cables & Lead Wires Cross Reference
- 1. 2 Major Types or Brands ECG Lead Wires

2. Telemetry ECG Cables And Holter ECG Cables

3. Disposable ECG Monitoring Lead Wires

- 3. 1 Disposable ECG Lead Wires
- 3. 2 Disposable ECG Radiolucent Lead Wires

4. Diagnosis EKG Cables & Accessories

- 4. 1 EKG Trunk Cables
- 4. 2 EKG Lead Wires
- 4. 3 EKG Accessories

5. EEG & EMG Electrodes

6. SPO2 Sensor & Ext-Cable

- 6. 1 Reusable Spo2 Sensor
- 6. 2 Single patient use Spo2 Sensor
- 6. 3 Spo2 Adapter cable

7. Temperature Sensors

8. Disposable ECG Monitoring Lead Wires

- 8. 1 Disposable ECG Lead Wires
- 8. 2 Disposable ECG Radiolucent Lead Wires

9. Diagnosis EKG Cables & Accessories

- 9. 1 EKG Trunk Cables
- 9. 2 EKG Lead Wires
- 9. 3 EKG Accessories

10. EEG & EMG Electrodes

Compliance with

ANSI/AAMI EC53
IEC 60601-1
ISO 10993-1
ISO 10993-5
ISO 10993-10

Certification



NMPA

(National Medical Products Administration)

1. ECG Monitoring Cables & Lead Wires

1.1. ECG Cables & Lead wires Cross reference



Specification

Cable material	TPU Jacket, Light gray
Main Cable length	2.5m
Connector & Yoke	Reliable, secure connection & Golden plated pins
Lead wire length	0.9m














Compatibility

Major brands ECG or multi-parameter monitoring system




ECG Cable & lead wire	Pictures	REF# examples	Reference to REF#
One piece Fixed Cable (Shielded)		A05S018-034A	A: Product category(A= One piece Fixed Cable) 05: Leads (03=3-lead; 05=5-lead; 10=10-lead) S: Snap end C: Clip end 018: GE -11pin Plug code 034: Representative length3.4m) A: AHA code (I= IEC code)
Trunk Cable (Shielded)		B05018-MQ-025AI	B: Product category(B= cable yoke type) 05: Leads (03=3-lead; 05=5-lead; 10=10-lead) 018: GE -11pin Plug code MQ: Yoke type (MQ= GE Multi-Link type; NK=NK BR type...) 025: Representative length(2.5m) A: AHA code (I=IEC code; A=AHA+IEC)
Lead wire (Shielded)		C05SMQ-90A	C: Product category (C= Reusable Lead wire) 05: Leads (03=3-lead; 05=5-lead; 10=10-lead) S: Snap end C: Clip end MQ: lead wire type(Acc to yoke) MQ= GE Multi-Link 90: 0.9m (wire length) A: AHA code (I= IEC code)



ECG Cable & lead wire	Pictures	REF# examples	Reference to REF#
Din trunk Cable (shielded) Used with Din lead wire Especially for neonate ECG monitoring		B03018-DN-025AI	B: Product category(B= cable yoke type) 05: Leads (03=3-lead; 05=5-lead; 10=10-lead) 018: GE -11pin Plug code DN: Yoke type (DN= DIN type...) 025: Representative length(25cm) A: AHA code (I=IEC code; AI=AHA+IEC)
Din lead wire (Din 1.5mm, Unshielded) Used with Din trunk cable only		C03SDN-90A	C: Product category (C= Reusable Lead wire) 05: Leads (03=3-lead; 05=5-lead; 10=10-lead) S: Snap end C: Clip end N: Neonate clip end DN: lead wire type(Acc to yoke) DN= DIN type 90: 90cm(wire length) A: AHA code (I= IEC code)
Neonate lead wire (Din 1.5mm, Unshielded) Used with Din trunk cable only		C03NDN-90A	C: Product category (C= Reusable Lead wire) 05: Leads (03=3-lead; 05=5-lead; 10=10-lead) N: Neonate clip end DN: lead wire type(Acc to yoke) DN= DIN type 90: 90cm(wire length) A: AHA code (I=IEC code)


ECG Color Code on electrodes



AHA(A)	Alternate AHA(H)	Standard IEC(I)	Alternate IEC(IH)	AHA(A)	Standard IEC(I)
 RA  LL  LA	 RA  RL  LA	 R  F  L	 R  N  L	 RA  RL  V  LL  LA	 R  N  C  F  L
 RA  LL  LA	 RA  RL  LA	 R  F  L	 R  N  L	 RA  RL  V  LL  LA	 R  N  C  F  L
 RA  LL  LA	 RA  RL  LA	 R  F  L	 R  N  L	 RA  RL  V  LL  LA	 R  N  C  F  L




1.1. ECG Cables & Lead wires cross reference




OEM & Models	Monitor Plug	Leads	Trunk Cable	Yoke Conn.	Ld wire Snap	Ld wire Clip	Fixed Snap	Fixed Clip	Din Trunk
AAMI-6PIN									
Without resistance BCI, Mindray Spacelabs Welch Allyn GE		3	/		C03SDN-90A	C03CDN-90A	A03S001-034A	A03C001-034A	B03001-DN-025A
			/		C03SDN-90I	C03CDN-90I	A03S001-034I	A03C001-034I	B03001-DN-025I
		Price	/						A8
		5	/		C05SDN-90A	C05CDN-90A	A05S001-034A	A05C001-034A	B05001-DN-025A
			/		C05SDN-90I	C05CDN-90I	A05S001-034I	A05C001-034I	B05001-DN-025I
		Price	/						A9
With 1K resistance Biolight M8000/9000 BP 308/308E/408		3	/		C03SDN-90A	C03CDN-90A	A03S002-034A	A03C002-034A	B03002-DN-025A
			/		C03SDN-90I	C03CDN-90I	A03S002-034I	A03C002-034I	B03002-DN-025I
		Price	/						A8
		5	/		C05SDN-90A	C05CDN-90A	A05S002-034A	A05C002-034A	B05002-DN-025A
			/		C05SDN-90I	C05CDN-90I	A05S002-034I	A05C002-034I	B05002-DN-025I
		Price	/						A9
Spacelabs LL style		3	B03002-LL-025A		C03SLL-90A	C03CLL-90A	A03S002-034A	A03C002-034A	B03002-DN-025A
			B03002-LL-025I		C03SLL-90I	C03CLL-90I	A03S002-034I	A03C002-034I	B03002-DN-025I
		Price	A0						A9
		5	B05002-LL-025A		C05SLL-90A	C05CLL90A	A05S002-034A	A05C002-034A	B05002-DN-025A
			B05002-LL-025I		C05SLL-90I	C05CLL-90I	A05S002-034I	A05C002-034I	B05002-DN-025I
		Price	B1						A0
Artema/S&W									
8000, Diascope1,-2, Diascope NT, Athena		5	B05039-DX-025A		C05SDX-90A	C05CDX90A	A05S039-034A	A05C039-034A	B05039-DN-025A
			B05039-DX-025I		C05SDX-90I	C05CDX-90I	A05S039-034I	A05C039-034I	B05039-DN-025I
		Price	B3						B2




OEM &Models	Monitor Plug	Leads	Trunk Cable	Lead yokeConn.	Ld wire Snap	Ld wire Clip	Fixed Snap	Fixed Clip	Din Trunk
Biolight									
M7000,M9500, AS Modular		3	B03003-LL-025A		C03SLL-90A	C03CLL-90A	A03S003-034A	A03C003-034A	B03003-DN-025A
			B03003-LL-025I		C03SLL-90I	C03CLL-90I	A03S003-034I	A03C003-034I	B03003-DN-025I
		Price	B2						B1
		5	B05003-LL-025A		C05SLL-90A	C05CLL-90A	A05S001-034A	A05C001-034A	B05003-DN-025A
			B05003-LL-025I		C05SLL-90I	C05CLL-90I	A05S001-034I	A05C001-034I	B05003-DN-025I
		Price	B3						B2
AnyView A2/A3/A5/A6/A8 Q3/Q5/Q7 V6 Vital Sign Monitor		3	B03004-LL-025A		C03SLL-90A	C03CLL-90A	A03S004-034A	A03C004-034A	B03004-DN-025A
			B03004-LL-025I		C03SLL-90I	C03CLL-90I	A03S004-034I	A03C004-034I	B03004-DN-025I
		Price	B3						B2
		5	B05004-LL-025A		C05SLL-90A	C05CLL-90A	A05S004-034A	A05C004-034A	B05004-DN-025A
			B05004-LL-025I		C05SLL-90I	C05CLL-90I	A05S004-034I	A05C004-034I	B05004-DN-025I
		Price	B4						B3
Bionet(Korea)									
BM3, BM3VET		3	B03005-AA-025A		C03SAA-90A	C03CAA-90A	A03S005-034A	A03C005-034A	B03005-DN-025A
			B03005-AA-025I		C03SAA-90I	C03CAA-90I	A03S005-034I	A03C005-034I	B03005-DN-025I
		Price	B1						A0
		5	B05005-AA-025A		C05SAA-90A	C05CAA-90A	A05S005-034A	A05C005-034A	B05005-DN-025A
			B05005-AA-025I		C05SAA-90I	C05CAA-90I	A05S005-034I	A05C005-034I	B05005-DN-025I
		Price	B2						B1
BM5, BM7		3	B03006-AA-025A		C03SAA-90A	C03CAA-90A	A03S006-034A	A03C006-034A	B03006-DN-025A
			B03006-AA-025I		C03SAA-90I	C03CAA-90I	A03S006-034I	A03C006-034I	B03006-DN-025I
		Price	B1						A0
		5	B05006-AA-025A		C05SAA-90A	C05CAA-90A	A05S006-034A	A05C006-034A	B05006-DN-025A
			B05006-AA-025I		C05SAA-90I	C05CAA-90I	A05S006-034I	A05C006-034I	B05006-DN-025I
		Price	B2						B1



OEM & Models	Monitor Plug	Leads	Trunk Cable	Yoke Conn.	Ld wire Snap	Ld wire Clip	Fixed Snap	Fixed Clip	Din Trunk
Biosys(Korea)									
BPM-103		3	B03007-AA-025A		C03SAA-90A	C03CAA-90A	A03S007-034A	A03C007-034A	B03007-DN-025A
			B03007-AA-025I		C03SAA-90I	C03CAA-90I	A03S007-034I	A03C007-034I	B03007-DN-025I
		price	B2						B1
		5	B05007-AA-025A		C05SAA-90A	C05CAA-90A	A05S007-034A	A05C007-034A	B05007-DN-025A
			B05007-AA-025I		C05SAA-90I	C05CAA-90I	A05S007-034I	A05C007-025I	B05007-DN-025I
		price	B3						B2
Bruker									
Bruker (Odam) all, as DC-3002, IH, M3 SM784, SM785		3	B03008-AA-025A		C03SAA-90A	C03CAA-90A	A03S008-034A	A03C008-034A	B03008-DN-025A
			B03008-AA-025I		C03SAA-90I	C03CAA-90I	A03S008-034I	A03C008-034I	B03008-DN-025I
		price	B1						A0
		5	B05008-AA-025A		C05SAA-90A	C05CAA-90A	A05S008-034A	A05C007-034A	B05008-DN-025A
			B05008-AA-025I		C05SAA-90I	C05CAA-90I	A05S008-034I	A05C007-025I	B05008-DN-025I
		price	B2						B1
Comen									
STAR-8000A/B/C/D /E/H, STAR-5000C, C80		3	B03001-AA-025A		C03SAA-90A	C03CAA-90A	A03S001-034A	A03C001-034A	B03001-DN-025A
			B03001-AA-025I		C03SAA-90I	C03CAA-90I	A03S001-034I	A03C001-034I	B03001-DN-025I
		price	A0						A9
		5	B05001-AA-025A		C05SAA-90A	C05CAA-90A	A05S001-034A	A05C001-034A	B05001-DN-025A
			B05001-AA-025I		C05SAA-90I	C05CAA-90I	A05S001-034I	A05C001-034I	B05001-DN-025I
		price	B1						A0
C20, C21, C30, C50, C70, C90, C100, C100A, NC3, NC5, NC8, NC10, NC12, NC19		3	B03009-AA-025A		C03SAA-90A	C03CAA-90A	A03S009-034A	A03C009-034A	B03009-DN-025A
			B03009-AA-025I		C03SAA-90I	C03CAA-90I	A03S009-034I	A03C009-034I	B03009-DN-025I
		price	B2						B1
		5	B05009-AA-025A		C05SAA-90A	C05CAA-90A	A05S009-034A	A05C009-034A	B05009-DN-025A
			B05009-AA-025I		C05SAA-90I	C05CAA-90I	A05S009-034I	A05C009-025I	B05009-DN-025I
		price	B3						B2





OEM & Models	Monitor Plug	Leads	Trunk Cable	Yoke Conn.	Ld wire Snap	Ld wire Clip	Fixed Snap	Fixed Clip	Din Trunk
CAS MEDICAL									
750		3	B03010-AA-025A		C03SAA-90A	C03CAA-90A	A03S010-034A	A03C010-034A	B03010-DN-025A
			B03010-AA-025I		C03SAA-90I	C03CAA-90I	A03S010-034I	A03C010-034I	B03010-DN-025I
		price	B3						B2
		5	B05010-AA-025A		C05SAA-90A	C05CAA-90A	A05S010-034A	A05C010-034A	B05010-DN-025A
			B05010-AA-025I		C05SAA-90I	C05CAA-90I	A05S010-034I	A05C010-025I	B05010-DN-025I
		price	B4						B3
511 Apnea Monitor		3	B03043-AA-025A		C03SAA-90A	C03CAA-90A	A03S043-034A	A03C043-034A	B03043-DN-025A
			B03043-AA-025I		C03SAA-90I	C03CAA-90I	A03S043-034I	A03C043-034I	B03043-DN-025I
		price	B1						A0
		5	B05043-AA-025A		C05SAA-90A	C05CAA-90A	A05S043-034A	A05C043-034A	B05043-DN-025A
			B05043-AA-025I		C05SAA-90I	C05CAA-90I	A05S043-034I	A05C043-025I	B05043-DN-025I
		price	B2						B1
Contec									
CMS 6800, CMS 8000, CMS 9000		3	B03001-AA-025A		C03SAA-90A	C03CAA-90A	A03S001-034A	A03C001-034A	B03001-DN-025A
			B03001-AA-025I		C03SAA-90I	C03CAA-90I	A03S001-034I	A03C001-034I	B03001-DN-025I
		price	B1						A9
		5	B05001-AA-025A		C05SAA-90A	C05CAA-90A	A05S001-034A	A05C001-034A	B05001-DN-025A
			B05001-AA-025I		C05SAA-90I	C05CAA-90I	A05S001-034I	A05C001-034I	B05001-DN-025I
		price	B2						B1
Criticare									
504US, 602A-A, 506, 507,507E, 508, 1100, 1508,8100, Poet 2, Scholar EL,Scholar II 507EL, 8100E		3	B03001-AA-025A		C03SAA-90A	C03CAA-90A	A03S001-034A	A03C001-034A	B03001-DN-025A
			B03001-AA-025I		C03SAA-90I	C03CAA-90I	A03S001-034I	A03C001-034I	B03001-DN-025I
		price	B1						A9
		5	B05001-AA-025A		C05SAA-90A	C05CAA-90A	A05S001-034A	A05C001-034A	B05001-DN-025A
			B05001-AA-025I		C05SAA-90I	C05CAA-90I	A05S001-034I	A05C001-034I	B05001-DN-025I
		price	B2						B1




OEM & Models	Monitor Plug	Leads	Trunk Cable	Yoke Conn.	Ld wire Snap	Ld wire Clip	Fixed Snap	Fixed Clip	Din Trunk
Colin/Omron(Japan)									
BP 308/308E/408		3	B03001-AA-025A		C03SAA-90A	C03CAA-90A	A03S001-034A	A03C001-034A	B03001-DN-025A
			B03001-AA-025I		C03SAA-90I	C03CAA-90I	A03S001-034I	A03C001-034I	B03001-DN-025I
		price	B1						A9
		5	B05001-AA-025A		C05SAA-90A	C05CAA-90A	A05S001-034A	A05C001-034A	B05001-DN-025A
			B05001-AA-025I		C05SAA-90I	C05CAA-90I	A05S001-034I	A05C001-034I	B05001-DN-025I
		price	B2						B1
BP88S		3	B05034-CL-025A		C03SCL-90A	C03CCL-90A	A03C034-034A	A03S034-034A	B03034-DN-025A
			B05034-CL-025I		C03SCL-90I	C03CCL-90I	A03C034-034I	A03S034-034I	B03034-DN-025I
		Price	B5						B4
		5	B05034-CL-025A		C05SCL-90A	C05CCL-90A	A05C034-034A	A05S034-034A	B05034-DN-025A
			B05034-CL-025I		C05SCL-90I	C05CCL-90I	A05C034-034I	A05S034-034I	B05034-DN-025I
		Price	B6						B5
BP88 BP306		3	B03040-CL-025A		C03SCL-90A	C03CCL-90A	A03C040-034A	A03S034-034A	B03034-DN-025A
			B03040-CL-025I		C03SCL-90I	C03CCL-90I	A03C040-034I	A03S034-034I	B03034-DN-025I
		Price	B5						B4
		5	B05040-CL-025A		C05SCL-90A	C05CCL-90A	A05C040-034A	A05S040-034A	B05040-DN-025A
			B05040-CL-025I		C05SCL-90I	C05CCL-90I	A05C040-034I	A05S040-034I	B05040-DN-025I
		Price	B6						B5
Creative									
PC9000S,PC9000, UP6000		3	/		C03SDN-90A	C03CDN-90A	A03C011-034A	A03S011-034A	B03011-DN-025A
			/		C03SDN-90I	C03CDN-90I	A03C011-034I	A03S011-034I	B03011-DN-025I
		Price	/						B1
		5	/		C05SDN-90A	C05CDN-90A	A05C011-034A	A05S011-034A	B05011-DN-025A
			/		C05SDN-90I	C05CDN-90I	A05C011-034I	A05S011-034I	B05011-DN-025I
		Price							B2




OEM & Models	Monitor Plug	Leads	Trunk Cable	Yoke Conn.	Ld wire Snap	Ld wire Clip	Fixed Snap	Fixed Clip	Din Trunk
Dixtal									
BP 308/308E/408		3	B03001-AA-025A		C03SAA-90A	C03CAA-90A	A03S001-034A	A03C001-034A	B03001-DN-025A
			B03001-AA-025I		C03SAA-90I	C03CAA-90I	A03S001-034I	A03C001-034I	B03001-DN-025I
		price	B0						A9
		5	B05001-AA-025A		C05SAA-90A	C05CAA-90A	A05S001-034A	A05C001-034A	B05001-DN-025A
			B05001-AA-025I		C05SAA-90I	C05CAA-90I	A05S001-034I	A05C001-034I	B05001-DN-025I
		price	B1						B0
EP12		3	B03041-AA-025A		C03SAA-90A	C03CAA-90A	A03S041-034A	A03C041-034A	B03041-DN-025A
			B03041-AA-025I		C03SAA-90I	C03CAA-90I	A03S041-034I	A03C041-034I	B03041-DN-025I
		price	B2						B1
		5	B05041-AA-025A		C05SAA-90A	C05CAA-90A	A05S041-034A	A05C041-034A	B05041-DN-025A
			B05041-AA-025I		C05SAA-90I	C05CAA-90I	A05S041-034I	A05C041-034I	B05041-DN-025I
		price	B3						B2
Dräger/Siemens									
SC5000,SC6000,SC7000, SC8000,SC9000 Multi-Parameter Cable		3	B06055-MT-025AI (Dräger 3368391)		C03SDG-90A	C03CDG-90A	/	/	/
					C03SDG-90I	C03CDG-90I	/	/	/
		6			C05SDG-90A	C05CDG-90A	/	/	/
					C05SDG-90I	C05CDG-90I	/	/	/
		Price	E5				/	/	/
		3	B05055-MT-025AI (Dräger MS20093)		C03SSM-90A	C03CSM-90A	/	/	/
					C03SSM-90I	C03CSM-90I	/	/	/
		Price							
		5			C05SSM-90A	C05CSM-90A	/	/	/
					C05SSM-90I	C05CSM-90I	/	/	/
		Price	E5						

OEM &Models	Monitor Plug	Leads	Trunk Cable	Yoke Conn.	Ld wire Snap	Ld wire Clip	Fixed Snap	Fixed Clip	Din Trunk
Dräger/Siemens									
SC 5000, SC 6000, SC 7000, SC 8000, SC 9000		3	B03012-SM-025A		C03SSM-90A	C03CSM-90A	A03S012-034A	A03C012-034A	B03012-DN-025A
			B03012-SM-025I		C03SSM-90I	C03CSM-90I	A03S012-034I	A03C012-034I	B03012-DN-025I
		price	B2						B1
		5	B05012-SM-025A		C05SSM-90A	C05CSM-90A	A05S012-034A	A05C012-034A	B05012-DN-025A
			B05012-SM-025I		C05SSM-90I	C05CSM-90I	A05S012-034I	A05C012-034I	B05012-DN-025I
		price	B3						B2
Sirecust series 400 Sirecust341,610/620/630/7 20/730/732/900/960/961/1 260/1281/1801		3	B03013-SM-025A		C03SSM-90A	C03CSM-90A	A03S013-034A	A03C013-034A	B03013-DN-025A
			B03013-SM-025I		C03SSM-90I	C03CSM-90I	A03S013-034I	A03C013-034I	B03013-DN-025I
		price	B2						B1
		5	B05013-SM-025A		C05SSM-90A	C05CSM-90A	A05S013-034A	A05C013-034A	B05013-DN-025A
			B05013-SM-025I		C05SSM-90I	C05CSM-90I	A05S013-034I	A05C013-034I	B05013-DN-025I
		price	B3						B2
Vitalert 2000, Vitalert 3000, Vista 120 Siemens : Somatom Sensation		3	B03001-SM-025A		C03SSM-90A	C03CSM-90A	A03S001-034A	A03C001-034A	B03001-DN-025A
			B03001- SM -025I		C03SSM-90I	C03CSM-90I	A03S001-034I	A03C001-034I	B03001-DN-025I
		price	B1						A0
		5	B05001- SM -025A		C05SSM-90A	C05CSM-90A	A05S001-034A	A05C001-034A	B05001-DN-025A
			B05001- SM -025I		C05SSM-90I	C05CSM-90I	A05S001-034I	A05C001-034I	B05001-DN-025I
		price	B2						B1
Dräger /Air Shields									
DX 2022, DX 2515		3	B03001-SM-025A		C03SSM-90A	C03CSM-90A	A03S001-034A	A03C001-034A	B03001-DN-025A
			B03001- SM -025I		C03SSM-90I	C03CSM-90I	A03S001-034I	A03C001-034I	B03001-DN-025I
		price	B1						A0
		5	B05001- SM -025A		C05SSM-90A	C05CSM-90A	A05S001-034A	A05C001-034A	B05001-DN-025A
			B05001- SM -025I		C05SSM-90I	C05CSM-90I	A05S001-034I	A05C001-034I	B05001-DN-025I
		price	B2						B1


OEM &Models	Monitor Plug	Leads	Trunk Cable	Yoke Conn.	Ld wire Snap	Ld wire Clip	Fixed Snap	Fixed Clip	Din Trunk
DRE									
Waveline EZ MAX	See AAMI-6PIN								
Waveline Pro Touch-Screen Anesthesia monitor		3	B05014-MQ-025A		C03SMQ-90A	C03CMQ-90A	A03S014-034A	A03C014-034A	B03014-DN-025A
			B05014-MQ-025I		C03SMQ-90I	C03CMQ-90I	A03S014-034I	A03C014-034I	B03014-DN-025I
		price	B3						B2
		5	B05014-MQ-025A		C05SMQ-90A	C05CMQ-90A	A05S014-034A	A05C014-034A	B05014-DN-025A
			B05014-MQ-025I		C05SMQ-90I	C05CMQ-90I	A05S014-034I	A05C014-034I	B05014-DN-025I
		price	B4						B3
ASM 6000, Trax, Vida, Waveline		3	B05022-MQ-025A		C03SMQ-90A	C03CMQ-90A	A03S022-034A	A03C022-034A	B03022-DN-025A
			B05022-MQ-025I		C03SMQ-90I	C03CMQ-90I	A03S022-034I	A03C022-034I	B03022-DN-025I
		price	B2						B1
		5	B05022-MQ-025A		C05SMQ-90A	C05CMQ-90A	A05S022-034A	A05C022-034A	B05022-DN-025A
			B05022-MQ-025I		C05SMQ-90I	C05CMQ-90I	A05S022-034I	A05C022-034I	B05022-DN-025I
		price	B3						B2
		3	B03042-AA-025A		C03SAA-90A	C03CAA-90A	A03S042-034A	A03C042-034A	B03042-DN-025A
			B03042-AA-025I		C03SAA-90I	C03CAA-90I	A03S042-034I	A03C042-034I	B03042-DN-025I
		price	B2						B1
		5	B05042-AA-025A		C05SAA-90A	C05CAA-90A	A05S042-034A	A05C042-034A	B05042-DN-025A
			B05042-AA-025I		C05SAA-90I	C05CAA-90I	A05S042-034I	A05C042-034I	B05042-DN-025I
		price	B3						B2
EDAN									
M3,M3A,M3B, M8, M8A, M8B,M9, M9A		3	B03002-MR-025A		C03SMR-90A	C03CMR-90A	A03S002-034A	A03C002-034A	B03002-DN-025A
			B03002- MR -025I		C03SMR-90I	C03CMR-90I	A03S002-034I	A03C002-034I	B03002-DN-025I
		price	B1						A0
		5	B05002- MR -025A		C05SMR-90A	C05CMR-90A	A05S002-034A	A05C002-034A	B05002-DN-025A
			B05002- MR -025I		C05SMR-90I	C05CMR-90I	A05S002-034I	A05C002-034I	B05002-DN-025I
		price	B2						B1



OEM &Models	Monitor Plug	Leads	Trunk Cable	Yoke Conn.	Ld wire Snap	Ld wire Clip	Fixed Snap	Fixed Clip	Din Trunk
EDAN									
elite V5, elite V6, elite V8, iM50, iM60, iM70, iM80, iM8, iM8A,iM9, iM9A		3	B03026-MR-025A		C03SMR-90A	C03CMR-90A	A03S026-034A	A03C026-034A	B03026-DN-025A
			B03026-MR-025I		C03SMR-90I	C03CMR-90I	A03S026-034I	A03C026-034I	B03026-DN-025I
		price	B2						B1
		5	B05026-MR-025A		C05SMR-90A	C05CMR-90A	A05S026-034A	A05C026-034A	B05026-DN-025A
			B05026-MR-025I		C05SMR-90I	C05CMR-90I	A05S026-034I	A05C026-034I	B05026-DN-025I
		price	B3						B2
Emtel									
FX2000		3	B03044-AA-025A		C03SAA-90A	C03CAA-90A	A03S044-034A	A03C044-034A	B03044-DN-025A
			B03044-AA-025I		C03SAA-90I	C03CAA-90I	A03S044-034I	A03C044-034I	B03044-DN-025I
		price	B3						B2
		5	B05044-AA-025A		C05SAA-90A	C05CAA-90A	A05S044-034A	A05C044-034A	B05044-DN-025A
			B05044-AA-025I		C05SAA-90I	C05CAA-90I	A05S044-034I	A05C044-034I	B05044-DN-025I
		price	B4						B3
FX3000,FX3000C, FX3000MD,E1510		3	B03045-AA-025A		C03SAA-90A	C03CAA-90A	A03S045-034A	A03C045-034A	B03045-DN-025A
			B03045-AA-025I		C03SAA-90I	C03CAA-90I	A03S045-034I	A03C045-034I	B03045-DN-025I
		price	B3						B2
		5	B05045-AA-025A		C05SAA-90A	C05CAA-90A	A05S045-034A	A05C045-034A	B05045-DN-025A
			B05045-AA-025I		C05SAA-90I	C05CAA-90I	A05S045-034I	A05C045-034I	B05045-DN-025I
		price	B4						B3
Fukuda Denshi									
DS-7000,DS-7100, DS-7200, DS-7300		3	B03015-SM-025A		C03SSM-90A	C03CSM-90A	A03S015-034A	A03C015-034A	B03015-DN-025A
			B03015-SM-025I		C03SSM-90I	C03CSM-90I	A03S015-034I	A03C015-034I	B03015-DN-025I
		price	B2						B1
		5	B05015-SM -025A		C05SSM-90A	C05CSM-90A	A05S015-034A	A05C015-034A	B05015-DN-025A
			B05015-SM-025I		C05SSM-90I	C05CSM-90I	A05S015-034I	A05C015-034I	B05015-DN-025I
		price	B3						B2






OEM & Models	Monitor Plug	Leads	Trunk Cable	Yoke Conn.	Ld wire Snap	Ld wire Clip	Fixed Snap	Fixed Clip	Din Trunk
Fukuda Denshi									
DS5100E,DS5300W, DS5400,DS7001 AS/3, CS/3,Cardiicap I, Cardiicap II,Cardiicap/5, Light,S/5, Satliteplus		3	B03016-SM-025A		C03SSM-90A	C03CSM-90A	A03S016-034A	A03C016-034A	B03016-DN-025A
			B03016-SM-025I		C03SSM-90I	C03CSM-90I	A03S016-034I	A03C016-034I	B03016-DN-025I
		price	B3						B2
		5	B05016-SM -025A		C05SSM-90A	C05CSM-90A	A05S016-034A	A05C016-034A	B05016-DN-025A
			B05016-SM-025I		C05SSM-90I	C05CSM-90I	A05S016-034I	A05C016-034I	B05016-DN-025I
		price	B4						B3
DS-7000,DS-7100, DS-7200, DS-7300		10	B10015-MQ-025AI	 MV(V2-V6)	C05SMQ-90A C05SMV-90A	C05CMQ-90A C05CMV-90A	A10S015-034A	A10C015-034A	/
					C05SMQ-90I C05SMV-90I	C05CMQ-90I C05CMV-90I	A10S015-034I	A10C015-034I	/
		price	B8						
DS-7000,DS-7100, DS-7200, DS-7300	 	3	B03017-SM-025A		C03SSM-90A	C03CSM-90A	A03S017-034A	A03C017-034A	B03017-DN-025A
			B03017-SM-025I		C03SSM-90I	C03CSM-90I	A03S017-034I	A03C017-034I	B03017-DN-025I
		price	B7						B6
		5	B05017-SM-025A		C05SSM-90A	C05CSM-90A	A05S017-034A	A05C017-034A	B05017-DN-025A
			B05017-SM-025I		C05SSM-90I	C05CSM-90I	A05S017-034I	A05C017-034I	B05017-DN-025I
		price	B8						B7
		10	B10017-MQ-025A	 MV(V2-V6)	C05SMQ-90A C05SMV-90I	C05CMQ-90A C05CMV-90I	A10S017-034A	A10C017-034A	/
			B10017-MQ-025I		C05SMQ-90A C05SMV-90I	C05CMQ-90A C05CMV-90I	A10S017-034I	A10C017-034I	/
		price	D0						
GE-Critikon									
Dianmap Pro 1000 MPS Select Portable monitors		5	B05002-DX-025A		C05SDX-90A	C05CDX-90A	A05S002-034A	A05C002-034A	B05002-DN-025A
			B05002-DX-025I		C05SDX-90I	C05CDX-90I	A05S002-034I	A05C002-034I	B05002-DN-025I
		price	B2						B1


OEM & Models	Monitor Plug	Leads	Trunk Cable	Yoke Conn.	Ld wire Snap	Ld wire Clip	Fixed Snap	Fixed Clip	Din Trunk
GE-DateX Ohmeda									
AS/3, CS/3, Cardiacap I, Cardiacap II, Cardiacap/5, Light, S/5, Satliteplus		3	B03020-DX-025A		C03SDX-90A	C03CDX-90A	A03S020-034A	A03C020-034A	B03020-DN-025A
			B03020-DX-025I		C03SDX-90I	C03CDX-90I	A03S020-034I	A03C020-034I	B03020-DN-025I
		price	B2						B1
		5	B05020-DX-025A		C05SDX-90A	C05CDX-90A	A05S020-034A	A05C020-034A	B05020-DN-025A
			B05020-DX-025I		C05SDX-90I	C05CDX-90I	A05S020-034I	A05C020-034I	B05020-DN-025I
		price	B3						B2
AS/3, CS/3, Cardiacap II, Satlite Plus S/5		3	B03019-DX-025A		C03SDX-90A	C03CDX-90A	A03S019-034A	A03C019-034A	B03019-DN-025A
			B03019-DX-025I		C03SDX-90I	C03CDX-90I	A03S019-034I	A03C019-034I	B03019-DN-025I
		price	B2						B1
		5	B05019-DX-025A		C05SDX-90A	C05CDX-90A	A05S019-034A	A05C019-034A	B05019-DN-025A
			B05019-DX-025I		C05SDX-90I	C05CDX-90I	A05S019-034I	A05C019-034I	B05019-DN-025I
		price	B3						B2
Dash series, EagleSystems, Solar 3000/8000/9500, Tram 100/ 200/300/ 400/ 451/500/600/800, Caescap B450/650/850 B20/B40/B60 B20i/B40i/B60i		3	B03018-MQ-025A		C03SMQ-90A	C03CMQ-90A	A03S018-034A	A03C018-034A	B03018-DN-025A
			B03018-MQ-025I		C03SMQ-90I	C03CMQ-90I	A03S018-034I	A03C018-034I	B03018-DN-025I
		price	B3						B1
		5	B05018-MQ-025A		C05SMQ-90A	C05CMQ-90A	A05S018-034A	A05C018-034A	B05018-DN-025A
			B05018-MQ-025I		C05SMQ-90I	C05CMQ-90I	A05S018-034I	A05C018-034I	B05018-DN-025I
		price	B3						B2
		6	B06018-MQ-025A		C06SMQ-90A	C06CMQ-90A	/	/	/
			B06018-MQ-025I		C06SMQ-90I	C06CMQ-90I	/	/	/
		Price	B3						
		10	B10018-MQ-025A	 MV(V2-V6)	C05SMQ-90A	C05CMQ-90A	A10S018-034A	A10C018-034A	/
			B10018-MQ-025I		C05SMQ-90A	C05CMQ-90A	A10S018-034I	A10C018-034I	/
		Price	B8		C05SMV-90I	C05CMV-90I			

OEM &Models	Monitor Plug	Leads	Trunk Cable	Yoke Conn.	Ld wire Snap	Ld wire Clip	Fixed Snap	Fixed Clip	Din Trunk
GE-Datex Ohmeda									
Vivid 7, Vivid S6, Vivid i, Vivid q, Vivid e		3	B05047-MQ-025A		C03SMQ-90A	C03CMQ-90A	A03S047-034A	A03C047-034A	B03047-DN-025A
			B05047-MQ-025I		C03SMQ-90I	C03CMQ-90I	A03S047-034I	A03C047-034I	B03047-DN-025I
		Price	B2					B1	
Logiq 5, Logiq 7, Logiq P6, Logiq P6 Pro, Logiq e9		3	B05048-MQ-025A		C03SMQ-90A	C03CMQ-90A	A03S048-034A	A03C048-034A	B03048-DN-025A
			B05048-MQ-025I		C03SMQ-90I	C03CMQ-90I	A03S048-034I	A03C048-034I	B03048-DN-025I
			Price		B3				
GE-Hellige									
TRAM 7200		3	B03046-DX-025A		C03SDX-90A	C03CDX-90A	A03S046-034A	A03C046-034A	B03046-DN-025A
			B03046-DX-025I		C03SDX-90I	C03CDX-90I	A03S046-034I	A03C046-034I	B03046-DN-025I
		price	B6					B5	
		5	B05046-DX-025A		C05SDX-90A	C05CDX-90A	A05S046-034A	A05C046-034A	B05046-DN-025A
			B05046-DX-025I		C05SDX-90I	C05CDX-90I	A05S046-034I	A05C046-034I	B05046-DN-025I
		price	B7					B6	
CardioServ,Generation, SCP 844, SCP 851,SCP 852, SCS 804,Servomed 4,Vicom SM, SMS, SMC, VICOM-SM		3	B03021-DX-025A		C03SDX-90A	C03CDX-90A	A03S021-034A	A03C021-034A	B03021-DN-025A
			B03021-DX-025I		C03SDX-90I	C03CDX-90I	A03S021-034I	A03C019-034I	B03021-DN-025I
		price	B2					B1	
		5	B05021-DX-025A		C05SDX-90A	C05CDX-90A	A05S021-034A	A05C021-034A	B05021-DN-025A
			B05021-DX-025I		C05SDX-90I	C05CDX-90I	A05S021-034I	A05C021-034I	B05021-DN-025I
		price	B3					B2	
Kontron									
7000, 7135, 7135B, 7141, 7210, 7250, 7251, K2000, KSLT, Kolormon,Micromon, Minimon,Supermon, Vitalmon		3	B03022-SM-025A		C03SSM-90A	C03CSM-90A	A03S022-034A	A03C022-034A	B03022-DN-025A
			B03022-SM-025I		C03SSM-90I	C03CSM-90I	A03S022-034I	A03C022-034I	B03022-DN-025I
		price	B2					B1	
		5	B05022-SM-025A		C05SSM-90A	C05CSM-90A	A05S022-034A	A05C022-034A	B05022-DN-025A
			B05022-SM-025I		C05SSM-90I	C05CSM-90I	A05S022-034I	A05C022-034I	B05022-DN-025I
		price	B3					B2	

OEM &Models	Monitor Plug	Leads	Trunk Cable	Yoke Conn.	Ld wire Snap	Ld wire Clip	Fixed Snap	Fixed Clip	Din Trunk
Kontron									
105,107,108,125,12,VSM5 10,Old generation		3	B03049-SM-025A		C03SSM-90A	C03CSM-90A	A03S049-034A	A03C049-034A	B03049-DN-025A
			B03049-SM-025I		C03SSM-90I	C03CSM-90I	A03S049-034I	A03C049-034I	B03049-DN-025I
		price	B2						B1
		5	B05049-SM-025A		C05SSM-90A	C05CSM-90A	A05S049-034A	A05C049-034A	B05049-DN-025A
			B05049-SM-025I		C05SSM-90I	C05CSM-90I	A05S049-034I	A05C049-034I	B05049-DN-025I
		price	B3						B2
Innocare									
IM200, InnoCare-VET		3	B03050-MQ-025A		C03SMQ-90A	C03CMQ-90A	A03S050-034A	A03C050-034A	B03050-DN-025A
			B03050-MQ-025I		C03SMQ-90I	C03CMQ-90I	A03S050-034I	A03C050-034I	B03050-DN-025I
		price	B2						B1
		5	B05050-MQ-025A		C05SMQ-90A	C05CMQ-90A	A05S050-034A	A05C050-034A	B05050-DN-025A
			B05050-MQ-025I		C05SMQ-90I	C05CMQ-90I	A05S050-034I	A05C050-034I	B05050-DN-025I
		price	B3						B2
Infinium									
Mira, Vega		3	B03022-SM-025A		C03SSM-90A	C03CSM-90A	A03S022-034A	A03C022-034A	B03022-DN-025A
			B03022-SM-025I		C03SSM-90I	C03CSM-90I	A03S022-034I	A03C022-034I	B03022-DN-025I
		price	B2						B1
		5	B05022-SM-025A		C05SSM-90A	C05CSM-90A	A05S022-034A	A05C022-034A	B05022-DN-025A
			B05022-SM-025I		C05SSM-90I	C05CSM-90I	A05S022-034I	A05C022-034I	B05022-DN-025I
		price	B3						B2
Lohmerier									
M010,M011/M211 Defi 501		5	B05051-MQ-025A		C05SMQ-90A	C05CMQ-90A	A05S051-034A	A05C051-034A	B05051-DN-025A
			B05051-MQ-025I		C05SMQ-90I	C05CMQ-90I	A05S051-034I	A05C051-034I	B05051-DN-025I
		price	B6						B5

OEM &Models	Monitor Plug	Leads	Trunk Cable	Yoke Conn.	Ld wire Snap	Ld wire Clip	Fixed Snap	Fixed Clip	Din Trunk
MEK									
MP500, MP600, MP1000		5	B05052-AA-025A		C05SAA-90A	C05CAA-90A	A05S052-034A	A05C052-034A	B05052-DN-025A
			B05052-AA-025I		C05SAA-90I	C05CAA-90I	A05S052-034I	A05C052-034I	B05052-DN-025I
		price	B2						B1
360697, M160246, MP1000NT		3	B03053-MQ-025A		C03SAA-90A	C03CAA-90A	A03S053-034A	A03C053-034A	B03053-DN-025A
			B03053-MQ-025I		C03SAA-90I	C03CAA-90I	A03S053-034I	A03C053-034I	B03053-DN-025I
		price	B1						A0
		5	B05053-AA-025A		C05SAA-90A	C05CAA-90A	A05S053-034A	A05C053-034A	B05053-DN-025A
			B05053-AA-025I		C05SAA-90I	C05CAA-90I	A05S053-034I	A05C053-034I	B05053-DN-025I
		price	B2						B1
Mennen									
MP500, MP600, MP1000		3	B03025-AA-025A		C03SAA-90A	C03CAA-90A	A03S025-034A	A03C025-034A	B03025-DN-025A
			B03025-AA-025I		C03SAA-90I	C03CAA-90I	A03S025-034I	A03C025-034I	B03025-DN-025I
		price	B2						B1
		5	B05025-AA-025A		C05SAA-90A	C05CAA-90A	A05S025-034A	A05C025-034A	B05025-DN-025A
			B05025-AA-025I		C05SAA-90I	C05CAA-90I	A05S025-034I	A05C025-034I	B05025-DN-025I
		price	B3						B2
Mindray									
PM7000/8000/9000	SEE AAMI-6PIN								
PM5000,PM6000 BeneViewT5/T6/T8 BeneVision,BeneView, iPM,iMEC, uMEC series BeneHeart,defibrillator		3	B03026-MR-025A		C03SMR-90A	C03CMR-90A	A03S026-034A	A03C026-034A	B03026-DN-025A
			B03026-MR-025I		C03SMR-90I	C03CMR-90I	A03S026-034I	A03C026-034I	B03026-DN-025I
		price	B2						B1
		5	B05026-MR-025A		C05SMR-90A	C05CMR-90A	A05S026-034A	A05C026-034A	B05026-DN-025A
			B05026-MR-025I		C05SMR-90I	C05CMR-90I	A05S026-034I	A05C026-034I	B05026-DN-025I
		price	B3						B2

OEM &Models	Monitor Plug	Leads	Trunk Cable	Yoke Conn.	Ld wire Snap	Ld wire Clip	Fixed Snap	Fixed Clip	Din Trunk
Mindray-Datascope									
Datascope -Mindray V12, V21		5	B05027-DT-025A		C05SDT-90A	C05CDT-90A	A05S027-034A	A05C027-034A	B05027-DN-025A
			B05027-DT-025I		C05SDT-90I	C05CDT-90I	A05S027-034I	A05C027-034I	B05027-DN-025I
		price	B3						
NihonKohden									
Def. Cardiolife Lifescope serieos 6-8-9-11,6102-6205-6120		5	B05029-NA-025A		C05SNA-90A	C05CNA-90A	A05S029-034A	A05C029-034A	B05029-DN-025A
			B05029-NA-025I		C05SNA-90I	C05CNA-90I	A05S029-034I	A05C029-034I	B05029-DN-025I
		price	B3						
		5	B05030-NA-025A		C05SNA-90A	C05CNA-90A	A05S030-034A	A05C030-034A	B05030-DN-025A
			B05030-NA-025I		C05SNA-90I	C05CNA-90I	A05S030-034I	A05C030-034I	B05030-DN-025I
		Price	B3						
TEC-7621C/TEC-7631C TEC-7721C/TEC-7731C 7700		3	B03031-NK-025A		C03SNK-90A	C03CNK-90A	A03S031-034A	A03C031-034A	B03031-DN-025A
			B03031- NK -025I		C03SNK-90I	C03CNK-90I	A03S031-034I	A03C031-034I	B03031-DN-025I
		Price	B2						
		5	B05031-NK-025A		C05SNK-90A	C05CNK-90A	A05S031-034A	A05C031-034A	B05031-DN-025A
			B05031-NK-025I		C05SNK-90I	C05CNK-90I	A05S031-034I	A05C031-034I	B05031-DN-025I
		Price	B3						
BSM-2300 lifescope I BSM-4100 lifescope P BSM-5100 lifescope A BSM-9510 lifescope M BSM-9800 lifescope S		3	B03028-NK-025A		C03SNK-90A	C03CNK-90A	A03S028-034A	A03C028-034A	B03028-DN-025A
			B03028- NK -025I		C03SNK-90I	C03CNK-90I	A03S028-034I	A03C028-034I	B03028-DN-025I
		Price	B3						
		5	B05028-NK-025A		C05SNK-90A	C05CNK-90A	A05S028-034A	A05C028-034A	B05028-DN-025A
			B05028-NK-025I		C05SNK-90I	C05CNK-90I	A05S028-034I	A05C028-034I	B05028-DN-025I
		Price	B4						

OEM & Models	Monitor Plug	Leads	Trunk Cable	Yoke Conn.	Ld wire Snap	Ld wire Clip	Fixed Snap	Fixed Clip	Din Trunk
Nihon Kohden									
Life Scope:OPV-1500		3	B03054-NK-025A (JC-103TA)		C03SNK-90A	C03CNK-90A	A03S054-034A	A03C054-034A	B03054-DN-025A
			B03054-NK-025I (JC-103T)		C03SNK-90I	C03CNK-90I	A03S054-034I	A03C054-034I	B03054-DN-025I
		Price	B4						B3
BSM-7102A, BSM-7300A, BSM-8108A, BSM-8301A, BSM-8302A, BSM-8500A, BSM-8502A, BSM-8800A, Life Scope 12, Life Scope 6, Life Scope 7, Life Scope 8, Life Scope 9, MU-831RA, OEC-6102A, OEC-6105A, OEC-7102A, OEC-8108A, TEC 5531E, TEC-5200A, TEC-7100A, TEC-7200A, TEC-7300A		3	B03002-NK-025A		C03SNK-90A	C03CNK-90A	A03S002-034A	A03C002-034A	B03002-DN-025A
			B03002-NK-025I		C03SNK-90I	C03CNK-90I	A03S002-034I	A03C002-034I	B03002-DN-025I
		Price	B1						A0
		5	B05002-NK-025A		C05SNK-90A	C05CNK-90A	A05S002-034A	A05C002-034A	B05002-DN-025A
			B05002-NK-025I		C05SNK-90I	C05CNK-90I	A05S002-034I	A05C002-034I	B05002-DN-025I
		Price	B2						B1
		3	B03002-NA-025A		C03SNA-90A	C03CNA-90A	A03S002-034A	A03C002-034A	B03002-DN-025A
			B03002-NA-025I		C03SNA-90I	C03CNA-90I	A03S002-034I	A03C002-034I	B03002-DN-025I
		Price	B1						A0
		5	B05002-NA-025A		C05SNA-90A	C05CNA-90A	A05S002-034A	A05C002-034A	B05002-DN-025A
			B05002-NA-025I		C05SNA-90I	C05CNA-90I	A05S002-034I	A05C002-034I	B05002-DN-025I
		Price	B2						B1
Philips									
Merlin78352C/78354C 78834C, ML65A,66A M1001A/M1002A, Virdia, Ominicare CMA, Heartstream Codemaster XL		3	B03035-AA-025A		C03SAA-90A	C03CAA-90A	A03S035-034A	A03C035-034A	B03035-DN-025A
			B03035-AA-025I		C03SAA-90I	C03CAA-90I	A03S035-034I	A03C035-034I	B03035-DN-025I
		price	B2						B1
		5	B05035-AA-025A		C05SAA-90A	C05CAA-90A	A05S035-034A	A05C035-034A	B05035-DN-025A
			B05035-AA-025I		C05SAA-90I	C05CAA-90I	A05S035-034I	A05C035-034I	B05035-DN-025I
		price	B3						B2

OEM &Models	Monitor Plug	Leads	Trunk Cable	Yoke Conn.	Ld wire Snap	Ld wire Clip	Fixed Snap	Fixed Clip	Din Trunk
Philips									
MP5, MP 20, MP 30, MP 50, MP 60,MP 70, Merlin, SureSigns VM6, SureSigns VM8, V24, V24E, V26		3	B03035-HP-025A		C03SHP-90A	C03CHP-90A	A03S035-034A	A03C035-034A	B03035-DN-025A
			B03035-HP-025I		C03SHP-90I	C03CHP-90I	A03S035-034I	A03C035-034I	B03035-DN-025I
		Price	B1						A0
		5	B05035-HP-025A		C05SHP-90A	C05CHP-90A	A05S035-034A	A05C035-034A	B05035-DN-025A
			B05035-HP-025I		C05SHP-90I	C05CHP-90I	A05S035-034I	A05C035-034I	B05035-DN-025I
		Price	B2						B1
Multi-Measurement Modules: M3000A, M3001A M3002A Patient Monitors: MP5,MP2,ECG Modules: M1001B, M1002B, 860335, 860338		10	B10035-HP-025A		C05SHC-90A Chest V1-V6	C05CHC-90A Chest V1-V6	A10S035-034A	A10C035-034A	/
					C05SHP-90A	C05CHP-90A			
					C05SHC-90I Chest C1-C6	C05CHC-90I Chest C1-C6	A10S035-034I	A10C035-034I	/
					C05SHP-90I	C05CHP-90I			
		Price	B5						
Physio Control									
Lifepak 11, Lifepak 12, Lifepak 15, Lifepak 20, Lifepak 20E		3	/	/	/	/	A03S024-034A	A03C024-034A	B03024-DN-025A
			/		/	/	A03S024-034I	A03C024-034I	B03024-DN-025I
		Price							
		5	/	/	/	/	A05S024-034A	A05C024-034A	B05024-DN-025A
			/		/	/	A05S024-034I	A05C024-034I	B05024-DN-025I
		Price							
		10			A04S024-034A	/	A10S024-034A	/	/
					A04S024-034I	/	A10S024-034I	/	/
		Price	A04S024-034A C05SPB-90A (V1-V6)						

OEM & Models	Monitor Plug	Leads	Trunk Cable	Yoke Conn.	Ld wire Snap	Ld wire Clip	Fixed Snap	Fixed Clip	Din Trunk
Spacelabs									
90701,90303,90305, 90308,90309,90385	SEE AAMI-6PIN								
Ultraview 90492, 91220, 90496, 91496		3	B06037-LB-025A		C03SLB-90A	C03CLB-90A	A03S037-034A	A03C037-034A	B03037-DN-025A
			B06037-LB-025I		C03SLB-90I	C03CLB-90I	A03S037-034I	A03C037-034I	B03037-DN-025I
		Price	B3						B2
		5	B06037-LB-025A		C05SLB-90A	C05CLB-90A	A05S037-034A	A05C037-034A	B05037-DN-025A
			B06037-LB-025I		C05SLB-90I	C05CLB-90I	A05S037-034I	A05C037-034I	B05037-DN-025I
		Price	B4						B3
		3	B03037-SL-025A		C03SSL-90A	C03CSL-90A	A03S037-034A	A03C037-034A	B03037-DN-025A
			B03037-SL-025I		C03SSL-90I	C03CSL-90I	A03S037-034I	A03C037-034I	B03037-DN-025I
		Price	B2						B1
		5	B05037-SL-025A		C05SSL-90A	C05CSL-90A	A05S037-034A	A05C037-034A	B05037-DN-025A
			B05037-SL-025I		C05SSL-90I	C05CSL-90I	A05S037-034I	A05C037-034I	B05037-DN-025I
		Price	B3						B2
		10	B10037-MQ-025A	 MV(V2-V6)	C05SMQ-90A	C05CMQ-90A	A10S018-034A	A10C018-034A	/
			B10037-MQ-025I		C05SMV-90I	C05CMV-90I			
					C05SMQ-90A	C05CMQ-90A	A10S018-034I	A10C018-034I	/
		Price	B8						
Welch Allyn	See AAMI								
Protocol	See AAMI-6PIN								
Micropaq Propaq LT		5	/		C05SDN-90A	C05CDN-90A	A05S038-034A	A05C038-034A	B05038-DN-025A
			/		C05SDN-90I	C05CDN-90I	A05S038-034I	A05C038-034I	B05038-DN-025I
		Price							B5

OEM &Models	Monitor Plug	Leads	Trunk Cable	Yoke Conn.	Ld wire Snap	Ld wire Clip	Fixed Snap	Fixed Clip	Din Trunk
ZOLL									
R SERIES		10	/	/	/	/	A10S039-034A	A10C039-034A	/
E SERIES			/	/	/	/	A10S039-034I	A10C039-034I	/
M SERIES		Price							
M CCT SERIES									
X SERIES		10	/	/	/	/	A10S040-034A	A10C040-034A	
			/	/	/	/	A10S040-034I	A10C040-034I	
		Price							
		4	/	/	/	/	A04S040-034A	A04C040-034A	
			/	/	/	/	A04S040-034I	A04C040-034I	
		Price							
		6	/	/	C06SZL-90A	C06CZL-90A	/	/	
			/	/	C06SZL-90I	C06CZL-90I	/	/	
		Price							



A04S040-034A

+



C06SZL-90A (Chest)

=

A10S040-034A

One Piece cables with neonate electrode end for most brands are available upon request. For neonate application, Din Style neonate leads are also the alternatives.

1.2. Major types or brands ECG lead wires



Specification

Cable material	TPU Jacket, Light gray
Connector to Yoke	Compatible with ECG cables from the world's leading brands of monitors
Lead wire length	0.9m

Compatibility

Major brands ECG or multi-parameter monitoring system

Types of Major ECG lead wires	Leads	Snap	Descriptions	Price code	Clip	Descriptions	Price code
Din1.5mm Leads							
	3	C03SDN-90A	3-Ld, Snap,AHA,35in.		C03CDN-90A	3-Ld, Clip,AHA,35in.	
		C03SDN-90I	3-Ld, Snap,IEC,35in.		C03CDN-90I	3-Ld, Clip,IEC,35in.	
	5	C05SDN-90A	5-Ld, Snap,AHA,35in.		C05CDN-90A	5-Ld, Clip,AHA,35in.	
		C05SDN-90I	5-Ld, Snap,IEC,35in.		C05CDN-90I	5-Ld, Clip,IEC,35in.	
	7	C05SDN-90A	7-Ld, Snap,AHA,35in.		C05CDN-90A	7-Ld, Clip,AHA,35in.	
		C05SDN-90I	7-Ld, Snap,IEC,35in.		C05CDN-90I	7-Ld, Clip,IEC,35in.	
	10	C10SDN-90A	10-Ld, Snap,AHA,35in.		C10CDN-90A	10-Ld, Clip,AHA,35in.	
		C10SDN-90I	10-Ld, Snap,IEC,35in.		C10CDN-90I	10-Ld, Clip,IEC,35in.	

Types of Major ECG lead wires	Leads	Snap	Descriptions	Price code	Clip	Descriptions	Price code
Din1.5mm Leads For Neonatal							
	3	/	/		C03NDN-90A	3-Ld, Neo clip,AHA,35in.	
		/	/		C03NDN-90I	3-Ld, Neo Clip,IEC,35in.	
LL style Leads							
	3	C03SLL-90A	3-Ld, Snap,AHA,35in.		C03CLL-90A	3-Ld, Clip,AHA,35in.	
		C03SLL-90I	3-Ld, Snap,IEC,35in.		C03CL-90I	3-Ld, Clip,IEC,35in.	
	5	C05SLL-90A	5-Ld, Snap,AHA,35in.		C05CLL-90A	5-Ld, Clip,AHA,35in.	
		C05SL-90I	5-Ld, Snap,IEC,35in.		C05CLL-90I	5-Ld, Clip,IEC,35in.	
Mindray Leads							
	3	C03SMR-90A	3-Ld, Snap,AHA,35in.		C03CMR-90A	3-Ld, Clip,AHA,35in.	
		C03SMR-90I	3-Ld, Snap,IEC,35in.		C03CMR-90I	3-Ld, Clip,IEC,35in.	
	5	C05SMR-90A	5-Ld, Snap,AHA,35in.		C05CMR-90A	5-Ld, Clip,AHA,35in.	
		C05SMR-90I	5-Ld, Snap,IEC,35in.		C05CMR-90I	5-Ld, Clip,IEC,35in.	

Types of Major ECG lead wires	Leads	Snap	Descriptions	Price code	Clip	Descriptions	Price code
Philips Leads							
	3	C03SAA-90A	3-Ld, Snap,AHA,35in.		C03CAA-90A	3-Ld, Clip,AHA,35in.	
		C03SAA-90I	3-Ld, Snap,IEC,35in.		C03CAA-90I	3-Ld, Clip,IEC,35in.	
	5	C05SAA-90A	5-Ld, Snap,AHA,35in.		C05CAA-90A	5-Ld, Clip,AHA,35in.	
		C05SAA-90I	5-Ld, Snap,IEC,35in.		C05CAA-90I	5-Ld, Clip,IEC,35in.	
	3	C03SHP-90A	3-Ld, Snap,AHA,35in.		C03CHP-90A	3-Ld, Clip,AHA,35in.	
		C03SHP-90I	3-Ld, Snap,IEC,35in.		C03CHP-90I	3-Ld, Clip,IEC,35in.	
	5	C05SHP-90A	5-Ld, Snap,AHA,35in.		C05CHP-90A	5-Ld, Clip,AHA,35in.	
		C05SHP-90I	5-Ld, Snap,IEC,35in.		C05CHP-90I	5-Ld, Clip,IEC,35in.	

Types of Major ECG lead wires	Leads	Snap	Descriptions	Price code	Clip	Descriptions	Price code
GE-Multilink Leads							
	3	C03SMQ-90A	3-Ld, Snap,AHA,35in.		C03CMQ-90A	3-Ld, Clip,AHA,35in.	
		C03SMQ-90I	3-Ld, Snap,IEC,35in.		C03CMQ-90I	3-Ld, Clip,IEC,35in.	
	5	C05SMQ-90A	5-Ld, Snap,AHA,35in.		C05CMQ-90A	5-Ld, Clip,AHA,35in.	
		C05SMQ-90I	5-Ld, Snap,IEC,35in.		C05CMQ-90I	5-Ld, Clip,IEC,35in.	
GE Chest V2-V6 C2-C6	5	C05SMV-90A	5-Ld, Snap,AHA,35in.		C05CMV-90A	5-Ld, Clip,AHA,35in.	
		C05SMV-90I	5-Ld, Snap,IEC,35in.		C05CMV-90I	5-Ld, Clip,IEC,35in.	
	10	C10SMQ-90A	10-Ld, Snap,AHA,35in.		C10CMQ-90A	10-Ld, Clip,AHA,35in.	
		C10SMQ-90I	10-Ld, Snap,IEC,35in.		C10CMQ-90I	10-Ld, Clip,IEC,35in.	

Types of Major ECG lead wires	Leads	Snap	Descriptions	Price code	Clip	Descriptions	Price code
NihonKohden Leads							
	3	C03SNK-90A	3-Ld, Snap,AHA,35in.		C03CNK-90A	3-Ld, Clip,AHA,35in.	
		C03SNK-90I	3-Ld, Snap,IEC,35in.		C03CNK-90I	3-Ld, Clip,IEC,35in.	
	6	C06SNK-90A	6-Ld, Snap,AHA,35in.		C06CNK-90A	5-Ld, Clip,AHA,35in.	
		C06SNK-90I	6-Ld, Snap,IEC,35in.		C06CNK-90I	5-Ld, Clip,IEC,35in.	
Dräger-Siemens Leads							
	3	C03SDG-90A	3-Ld, Snap,AHA,35in.		C03CDG-90A	3-Ld, Clip,AHA,35in.	
		C03SDG-90I	3-Ld, Snap,IEC,35in.		C03CDG-90I	3-Ld, Clip,IEC,35in.	
	5	C05SDG-90A	5-Ld, Snap,AHA,35in.		C05CDG-90A	5-Ld, Clip,AHA,35in.	
		C05SDG-90I	5-Ld, Snap,IEC,35in.		C05CDG-90I	5-Ld, Clip,IEC,35in.	
	6	C06SDG-90A	6-Ld, Snap,AHA,35in.		C06CDG-90A	6-Ld, Clip,AHA,35in.	
		C06SDG-90I	6-Ld, Snap,IEC,35in.		C06CDG-90I	6-Ld, Clip,IEC,35in.	
Spacelabs Leads							
	5	C05SSL-90A	5-Ld, Snap,AHA,35in.		C05CSL90A	5-Ld, Clip,AHA,35in.	
		C05SSL-90I	5-Ld, Snap,IEC,35in.		C05CSL-90I	5-Ld, Clip,IEC,35in.	

2. Telemetry ECG Cables And Holter ECG Cables



Compatibility

Major brands ECG transmitter & Holter system

Specification

Cable material	TPU Jacket, Light gray
Lead wire length	0.9m
Connector & Yoke	Golden plated pins

OEM & Models	Photo	Lead#	Parts-Med REF#	Descriptions	OEM PN	Price code
GE Healthcare Apexpro FH Telemetry Patientnet wireless telemetry systems		3	C03SAP-90A	3-Ld, Snap,AHA,35in.	394111-002	
			C03SAP-90I	3-Ld, Snap,IEC,35in.		
			C03CAP-90A	3-Ld, Clip,AHA,35in.		
			C03CAP-90I	3-Ld, Clip,IEC,35in.		
		5	C05SAP-90A	5-Ld, Snap,AHA,35in.	394111-006	
			C05SAP-90I	5-Ld, Snap,IEC,35in.		
			C05CAP-90A	5-Ld, Clip,AHA,35in.		
			C05CAP-90I	5-Ld, Clip,IEC,35in.		
		6	C06SAP-90A	6-Ld, Snap,AHA,35in.	394111-008	
			C06SAP-90I	6-Ld, Snap,IEC,35in.		
			C06CAP-90A	6-Ld, Clip,AHA,35in.		
			C06CAP-90I	6-Ld, Clip,IEC,35in.		
GE Hralthcare ApexPro Telemetry		6	C06SMQ-90A	6-Ld, Snap,AHA,35in.	421930-001	
			C06SMQ-90I	6-Ld, Snap,IEC,35in.	421932-001	
		6	C06CMQ-90A	6-Ld, Clip,AHA,35in.	421930-002	
			C06CMQ-90I	6-Ld, Clip,IEC,35in.	421930-002	
Mindray TMS-6012,TEL-100		5	C05STM-90A	5-Ld, Snap,AHA,35in.	115-004869-00	
			C05STM-90I	5-Ld, Snap,IEC,35in.		
			C05CTM-90A	5-Ld, Clip,AHA,35in.	115-004873-00	
			C05CTM-90I	5-Ld, Clip,IEC,35in.		

OEM & Models	Photo	Lead#	Parts-Med REF#	Descriptions	OEM PN	Price code
Philips Legacy M2600A/ M2601A, M2636A Transmitters 		3	C03SAT-90A	3-Ld, Snap,AHA,35in.	M2590A	
			C03SAT-90I	3-Ld, Snap,IEC,35in.	M2594A	
			C03CAT-90A	3-Ld, Clip,AHA,35in.	M2591A	
			C03CAT-90I	3-Ld, Clip,IEC,35in.	M2595A	
		5	C05SAT-90A	5-Ld, Snap,AHA,35in.	M2592A	
			C05SAT-90I	5-Ld, Snap,IEC,35in.	M2596A	
			C05CAT-90A	5-Ld, Clip,AHA,35in.	M2593A	
			C05CAT-90I	5-Ld, Clip,IEC,35in.	M2597A	
Philips IntelliVue Telemetry M4841A M2601B 		3	C03SHT-90A	3-Ld, Snap,AHA,35in.	989803151991	
			C03SHT-90I	3-Ld, Snap,IEC,35in.	989803152001	
			C03CHT-90A	3-Ld, Clip,AHA,35in.	989803173141	
			C03CHT-90I	3-Ld, Clip,IEC,35in.	989803151981	
		5	C05SHT-90A	5-Ld, Snap,AHA,35in.	989803173151	
			C05SHT-90I	5-Ld, Snap,IEC,35in.	989803152081	
			C05CHT-90A	5-Ld, Clip,AHA,35in.	989803152071	
			C05CHT-90I	5-Ld, Clip,IEC,35in.	989803152061	
Philips Medical MX40 Without SPO2 		3	C03SMX-90A	3-Ld, Snap,AHA,35in.		
			C03SMX-90I	3-Ld, Snap,IEC,35in.	989803171801	
			C03CMX-90A	3-Ld, Clip,AHA,35in.		
			C03CMX-90I	3-Ld, Clip,IEC,35in.	989803171901	
		5	C05SMX-90A	5-Ld, Snap,AHA,35in.	989803171821	
			C05SMX-90I	5-Ld, Snap,IEC,35in.		
			C05CMX-90A	5-Ld, Clip,AHA,35in.	989803171831	
			C05CMX-90I	5-Ld, Clip,IEC,35in.	989803171931	
		6	C06SMX-90A	6-Ld, Snap,AHA,35in.		
			C06SMX-90I	6-Ld, Snap,IEC,35in.		
			C06CMX-90A	6-Ld, Clip,AHA,35in.	989803171861	
			C06CMX-90I	6-Ld, Clip,IEC,35in.	989803171961	

OEM & Models	Photo	Lead#	Parts-Med REF#	Descriptions	OEM PN	Price code
Philips Medical MX40 With SPO2 Senear 		3	C03SMS-90A	3-Ld, Snap,AHA,35in.		
			C03SMS-90I	3-Ld, Snap,IEC,35in.		
			C03CMS-90A	3-Ld, Clip,AHA,35in.	989803171811	
			C03CMS-90I	3-Ld, Clip,IEC,35in.	989803171911	
		5	C05SMS-90A	5-Ld, Snap,AHA,35in.		
			C05SMS-90I	5-Ld, Snap,IEC,35in.		
			C05CMS-90A	5-Ld, Clip,AHA,35in.	989803171851	
			C05CMS-90I	5-Ld, Clip,IEC,35in.	989803171951	
		6	C06SMS-90A	6-Ld, Snap,AHA,35in.		
			C06SMS-90I	6-Ld, Snap,IEC,35in.		
			C06CMS-90A	6-Ld, Clip,AHA,35in.	989803171871	
			C06CMS-90I	6-Ld, Clip,IEC,35in.	989803171971	
Mortara H12+™ Digital Holter Recorder Mortara X12+, Mortara 		10	C10SMT1-90A	10-Ld, Snap,AHA,35in.	9293-034-50	
			C10SMT1-90I	10-Ld, Snap,IEC,35in.	9293-017-50	
Mortara H3+ 		5	C05SMT2-90A	5-Ld, Snap,AHA,35in.		
			C05SMT2-90I	5-Ld, Snap,IEC,35in.		

OEM & Models	Photo	Lead#	Parts-Med REF#	Descriptions	OEM PN	Price code
Philips Medical DigiTrak XT Zymed DigiTrak-Plus		5	C05SPM-90A	5-Ld, Snap,IEC,35in.	989803157491	
			C05SPM-90I	5-Ld, Snap,IEC,35in.		
Jincomed MIC-12H Series MIC-3H Series		7	C07SJC-90A	7-Ld, Snap,IEC,35in.		
			C07SJC-90I	7-Ld, Snap,IEC,35in.		
		10	C10SJC-90A	10-Ld, Snap,IEC,35in.		
			C10SJC-90I	10-Ld, Snap,IEC,35in.		
GE Healthcare SEER Light		5	C05SSL-90A	5-Ld, Snap,IEC,35in.	2008594-001	
			C05SSL-90I	5-Ld, Snap,IEC,35in.	2008594-003	
		7	C07SSL-90A	7-Ld, Snap,IEC,35in.	2008594-002	
			C07SSL-90I	7-Ld, Snap,IEC,35in.	2008594-004	
Bi-biomed BI6603,BI6612,BI6803, BI6812,BI9000,BI9800		3	C03SBI-90A	3-Ld, Snap,IEC,35in.		
			C03SBI-90I	3-Ld, Snap,IEC,35in.		
		5	C05SBI-90A	5-Ld, Snap,IEC,35in.		
			C05SBI-90I	5-Ld, Snap,IEC,35in.		
		7	C07SBI-90A	7-Ld, Snap,IEC,35in.		
			C07SBI-90I	7-Ld, Snap,IEC,35in.		
		10	C10SBI-90A	10-Ld, Snap,IEC,35in.	H2111	
			C10SBI-90I	10-Ld, Snap,IEC,35in.		

3. Disposable ECG Monitoring Lead Wires

3.1 Disposable ECG lead wires

Single patient use (disposable) from patient to monitor Reduces the risk of cross contamination from reusable cables and leads, improving patient care and cost control









- Full compatibility of the lead wires enable to be used with original ECG trunk Cables or devices, no extra adapter required.
- Built-in ribbon cable for patient comfort and ease of use
- Easy to read labeling facilities easy connection and lead placement
- Unique light blue color indicated the single use



• Patient Safety






• Patient Comfort

• Electrical Performance

• Ease of Use

OEM & Models	Connector	Descriptions	Lead #	REF# Snap(AHA/IEC)	REF# Clip(AHA/IEC)	Packaging	Price code
Mindray -Datascope		Disposable 3-lead shielded, 90cm	3	D03SDT-90A D03SDT-90I	D03CDT-90A D03CDT-90I	10PCS/BOX	
		Disposable 5-lead shielded, 90cm	5	D05SDT-90A D05SDT-90I	D05CDT-90A D05CDT-90I	10PCS/BOX	
Mindray		Disposable 3-lead shielded, 90cm	3	D03SMR-90A D03SMR-90I	D03CMR-90A D03CMR-90I	10PCS/BOX	
		Disposable 5-lead shielded, 90cm	5	D05SMR-90A D05SMR-90I	D05CMR-90A D05CMR-90I	10PCS/BOX	
Philips		Disposable 3-lead shielded, 90cm	3	D03SHP-90A D03SHP-90I	D03CHP-90A D03CHP-90I	10PCS/BOX	
		Disposable 5-lead shielded, 90cm	5	D05SHP-90A D05SHP-90I	D05CHP-90A D05CHP-90I	10PCS/BOX	
		Disposable 5-lead shielded, Chest, 90cm C2-C6;V2-V6	5	D05SHC-90A D05SHC-90I	D05CHC-90A D05CHC-90I	10PCS/BOX	
		Disposable 3-lead shielded, 90cm	3	D03SHT-90A D03SHT-90I	D03CHT-90A D03CHT-90I	10PCS/BOX	
		Disposable 5-lead shielded, 90cm	5	D05SHT-90A D05SHT-90I	D05CHT-90A D05CHT-90I	10PCS/BOX	

OEM & Models	Connector	Descriptions	Lead #	REF# Snap(AHA/IEC)	REF# Clip(AHA/IEC)	Packaging	Price code
Dräger		Disposable 3-lead shielded, 90cm	3	D03SDG-90A D03SDG-90I	D03CDG-90A D03CDG-90I	10PCS/BOX	
		Disposable 5-lead shielded, 90cm	5	D05SDG-90A D05SDG-90I	D05CDG-90A D05CDG-90I	10PCS/BOX	
		Disposable 6-lead shielded, 90cm	6	D06SDG-90A D06SDG-90I	D06CDG-90A D06CDG-90I	10PCS/BOX	
Dräger -Siemens		Disposable 3-lead shielded, 90cm	3	D03SSM-90A D03SSM-90I	D03CSM-90A D03CSM-90I	10PCS/BOX	
		Disposable 5-lead shielded, 90cm	5	D05SSM-90A D05SSM-90I	D05CSM-90A D05CSM-90I	10PCS/BOX	
GE Multi-link		Disposable 3-lead shielded, 90cm	3	D03SMQ-90A D03SMQ-90I	D03CMQ-90A D03CMQ-90I	10PCS/BOX	
		Disposable 5-lead shielded, 90cm	5	D05SMQ-90A D05SMQ-90I	D05CMQ-90A D05CMQ-90I	10PCS/BOX	
		Disposable 5-lead shielded, Chest, 90cm C2-C6;V2-V6	5	D05SMV-90A D05SMV-90I	D05CMV-90A D05CMV-90I	10PCS/BOX	

OEM & Models	Connector	Descriptions	Lead #	REF# Snap(AHA/IEC)	REF# Clip(AHA/IEC)	Packaging	Price code
GE Apex Pro Telemetry		Disposable 5-lead shielded for Apex, 90cm	6	D06SMQ-90A D06SMQ-90I	D06CMQ-90A D06CMQ-90I	10PCS/BOX	
NihonKohden		Disposable 3-lead shielded, 90cm	3	D03SNK-90A D03SNK-90I	D03CNK-90A D03CNK-90I	10PCS/BOX	
		Disposable 6-lead shielded, 90cm	6	D06SNK-90A D06SNK-90I	D06CNK-90A D06CNK-90I	10PCS/BOX	
Spacelabs		Disposable 3-lead shielded, 90cm	3	D03SSL-90A D03SSL-90I	D03CMQ-90A D03CMQ-90I	10PCS/BOX	
		Disposable 5-lead shielded, 90cm	5	D05SSL-90A D05SSL-90I	D05CSL-90A D05CSL-90I	10PCS/BOX	
Philips Legacy M2600A/ M2601A, M2636A Transmitters		Disposable 3-lead shielded, 90cm	3	D03SAT-90A D03SAT-90I	D03CAT-90A D03CAT-90I	10PCS/BOX	
		Disposable 3-lead shielded, 90cm	5	D05SAT-90A D05SAT-90I	D05CAT-90A D05CAT-90I	10PCS/BOX	

4.Diagnosis ECG Cables & Accessories



Specification







Cable material	TPU Jacket, Light gray
Main Cable length	Fixed cable 2.0m;Trunk 2.5m
Connector	plug &Golden plated pins
Lead wire length	1.2m(limb lead);0.8m(Chest lead)







Compatibility

Major brands ECG diagnosis system.







Product name	Photo	Part # Example	Reference to part number
One piece EKG fixed cable (A)		EA10B-MQ1-020A EA10B-MQ1-020I (Example)	E: EKG diagnostic cable or lead wire A: One piece EKG fixed cable 10: 10-lead(12=12 lead) B: banana end (N= needle, C=Clip, S=snap) MQ: GE Marquette factory/model code 1:Connector configure; B: banana end 020: cable length 2.0m I: IEC code (A= AHA code)
EKG Trunk Cable (B)		EB10-MQ1-025A EB10-MQ1-025I (Example)	E: EKG diagnostic cable or lead wire B: Trunk cable 10: 10-lead(12=12 lead) MQ: GE Marquette, factory/model code 1:Connector configure 025: cable length 2.5m I: IEC code (A= AHA code)
EKG Lead wire set (C)		EC10B-MQ-012A EC10B-MQ-012I (Example)	E: EKG diagnostic cable or lead wire C: EKG Lead wire set 10: 10-lead(12=12 lead) B: banana end (N= needle, C=Clip, S=snap) MQ: GE Marquette VS (Connector configure) 012: cable length 1.2m I: IEC code (A= AHA code)





All one piece cable or lead wires are available with banana end, needle end, Clip end and snap end





Code 代码	AHA 美标	IEC 欧标	IH 替代欧标
4mm Banana end (B)			
3mm Needle end (N)			

Code 代码	AHA 美标	IEC 欧标	IH 替代欧标
Clip End (C)			
Snap End (S)			






4.1. One piece EKG fixed cable






OEM& Models	EKG machine Plug	Lead#	4.0mm banana end	3.0mm needle end	Clip end	Snap end	Price code
Mortara/Burdick/Cardiac Science							
EK-10 and Elite, E350, E350i, E550, E560, Eclipse™ LE/LE.		10	EA10B-BK1-020A EA10B-BK1-020I (012-0700-00)	EA10N-BK1-020A EA10N-BK1-020I	EA10C-BK1-020A EA10C-BK1-020I	EA10S-BK1-020A EA10S-BK1-020I	
Burdick Quinton Eclipse 4,8,400,800, 850,850i,plus,Eclipse LE II		10	EA10B-BK2-020A EA10B-BK2-020I	EA10N-BK2-020A EA10N-BK2-020I	EA10C-BK2-020A EA10C-BK2-020I	EA10S-BK2-020A EA10S-BK2-020I	
Mortara > Burdick ELI-200、ELI-50		10	EA10B-MT1-020A EA10B-MT1-020I	EA10N-MT1-020A EA10N-MT1-020I	EA10C-MT1-020A EA10C-MT1-020I	EA10S-MT1-020A EA10S-MT1-020I	
Fukuda Denshi							
FX101, FX102F, FX328U, FCP2201		10	EA10B-FD1-020A EA10B-FD1-020I	EA10N-FD1-020A EA10N-FD1-020I	EA10C-FD1-020A EA10C-FD1-020I	EA10S-FD1-020A EA10S-FD1-020I	
Fukuda ME							
KP-500,KP-500		10	EA10B-FD2-020A EA10B-FD2-020I	EA10N-FD2-020A EA10N-FD2-020I	EA10C-FD2-020A EA10C-FD2-020I	EA10S-FD2-020A EA10S-FD2-020I	
GE-Marquette							
MAC 500,MAC 1100,MAC 1200		10	EA10B-MQ1-020A EA10B-MQ1-020I	EA10N-MQ1-020A EA10N-MQ1-020I	EA10C-MQ1-020A EA10C-MQ1-020I	EA10S-MQ1-020A EA10S-MQ1-020I	



OEM& Models	EKG machine Plug	Lead#	4.0mm banana end	3.0mm needle end	Clip end	Snap end	Price code
Kenz							
CARDIOLINE 103,106 (Italy)		10	EA10B-KZ1-020A EA10B-KZ1-020I	EA10N-KZ1-020A EA10N-KZ1-020I	EA10C-KZ1-020A EA10C-KZ1-020I	EA10S-KZ1-020A EA10S-KZ1-020I	
ECG 108/110 Cardico 302/601/1210/1211		10	EA10B-KZ2-020A EA10B-KZ2-020I	EA10N-KZ2-020A EA10N-KZ2-020I	EA10C-KZ2-020A EA10C-KZ2-020I	EA10S-KZ2-020A EA10S-KZ2-020I	
Hellige / Siemens Hörmann / Bosch							
Hellige:Cardiosys, EK36, EK403, EK41, EK413, EK43,EK512, EK53 Siemens Hörmann:Cardiostat 1 / 3 T / 703 / 710,Cardina C,Cardirex 3T/4/6T/62, Cathcor, MegacardElemaA/B/R/E Bosch:Bioset3500/3600/3700/6 000/8000		10	EA10B-HL1-020A EA10B-HL1-020I	EA10N-HL1-020A EA10N-HL1-020I	EA10C-HL1-020A EA10C-HL1-020I	EA10S-HL1-020A EA10S-HL1-020I	
Mortara > Burdick ELI-200、ELI-50		10	EA10B-MT1-020A EA10B-MT1-020I	EA10N-MT1-020A EA10N-MT1-020I	EA10C-MT1-020A EA10C-MT1-020I	EA10S-MT1-020A EA10S-MT1-020I	
NihonKohden							
CaduoFax 6151,6851,6353 CadioFax 8110		10	EA10B-NK1-020A EA10B-NK1-020I	EA10N-NK1-020A EA10N-NK1-020I	EA10C-NK1-020A EA10C-NK1-020I	EA10S-NK1-020A EA10S-NK1-020I	

OEM& Models	EKG machine Plug	Lead#	4.0mm banana end	3.0mm needle end	Clip end	Snap end	Price code
NihonKohden							
CaduoFax Q,ECG-9110K, ECG-9130K,Doing Juang Cardioline (New)		10	EA10B-NK2-020A EA10B-NK2-020I	EA10N-NK2-020A EA10N-NK2-020I	EA10C-NK2-020A EA10C-NK2-020I	EA10S-NK2-020A EA10S-NK2-020I	
Philips/HP							
M1770A,M1771A,M1772A Pagewrite100/200/ 200i		10	EA10B-HP1-020A EA10B-HP1-020I	EA10N-HP1-020A EA10N-HP1-020I	EA10C-HP1-020A EA10C-HP1-020I	EA10S-HP1-020A EA10S-HP1-020I	
Schiller							
AT1, AT2, AT2 Plus,AT4 Resistor 10K Long Screw		10	EA10B-SH1-020A EA10B-SH1-020I	EA10N-SH1-020A EA10N-SH1-020I	EA10C-SH1-020A EA10C-SH1-020I	EA10S-SH1-020A EA10S-SH1-020I	
AT3,AT6,CS6,AT5, AT10,AT60 Short Screw		10	EA10B-SH2-020A EA10B-SH2-020I	EA10N-SH2-020A EA10N-SH2-020I	EA10C-SH2-020A EA10C-SH2-020I	EA10S-SH2-020A EA10S-SH2-020I	

4.2. EKG Trunk Cable

OEM& Models Mortara/Burdick	EKG machine Plug	Lead#	REF#	Used with lead wire set/GE Multi-link set	Price code
EK10-E350		10	EB10-BK1-025AI	10-lead, chest /0.9m;Limb/1.2m EC10B-MQ-012A, EC10N-MQ-012A,EC10C-MQ-012A, EC10S-MQ-012A,EC10B-MQ-012I, EC10N-MQ-012I, EC10C-MQ-012I,EC10S-MQ-012I (See the next page of 10-lead wire sets for detail)	
Burdick Quinton Eclipse 4,8,400,800, 850,850i,plus,Eclipse LE II and ATRIA PLUS		10	EB10-BK1-025AI	Same as above	
ELI-200、ELI-50		10	EB10-MT1-025AI	Same as above	
Fukuda Denshi					
FX101, FX102F, FX328U, FCP2201		10	EB10-FD1-025AI	Same as above	
Fukuda ME					
KP-500、KP-500D		10	EB10-FD2-025AI	Same as above	
GE-Marquette					
MAC 500,MAC 1100, MAC 1200		10	EB10-MQ1-025AI	Same as above	

Hellige/Siemens					
Hellige:Cardiosys, EK36, EK403, EK41, EK413, EK43,EK512, EK53		10	EB10-HL1-025AI	Same as above	
Kenz					
CARDIOLINE 103,106 (Italy)		10	EB10-KZ1-025AI	Same as above	
ECG 108/110 Cardico 302/601/1210/1211		10	EB10-KZ2-025AI	Same as above	
NihonKohden					
CaduoFax 6151,6851,6353 CadioFax 8110		10	EB10-NK1-025AI	Same as above	
CaduoFax Q,ECG-9110K ECG-9130K,Doing Juang Cardioline (New)		10	EB10-NK2-025AI	Same as above	
Philips					
M1770A,M1771A,M1772A Pagewrite100/200/200i		10	EB10-HP1-025AI	Same as above	

OEM& Models	EKG machine Plug	Lead#	REF#	Used with lead wire set/GE Multi-link set	Price code
Schiller					
AT1, AT2, AT2 Plus,AT4 Resistor 10K Long Screw		10	EB10-SH1-025AI	Same as above	
AT3,AT6,CS6,AT5, AT10,AT60 Short Screw		10	EB10-SH2-025AI	Same as above	
GE Healthcare-Resting sys.					
CAM14 coiled patient Cable Used with GE MAC 5000 resting system with CAM14 Acquisition Module			EBCAM-GE14 (2016560-001)4.5M+0.65M EBCAM-GE14-S (2016560-001) 1.4M+0.3M	Used with CAM 14 patient lead wires on next page.	
CASE 12, CASE 15, MAC PC, MAC 6, MAC 12, MAC 15, CENTRA Systems with AM4 or AM5 Acquisition Module			EBAM-GE4-S (700044-202)0.25+4.6M EBAM-GE4-S (700044-203) 0.25+1.3M	Used with CAM 14patient lead wires on next page	
GE MAX-1 Exercise Testing System with AM4 or AM5 Acquisition Module			EBMAX1-GE4 (700044-201)		

4.2 EKG Lead wire

GE Multi-Link 10-Lead-integrated leads	OEM & Models	Electrode end	REF# (AHA)	REF# (IEC)	Pieces/Set	Price code
	Banana End 900112-001(AHA) 900112-002(IEC)	4.0mm banana	EC10B-MQ-120A	EC10B-MQ-120I	10	
		3.0mm needle	EC10N-MQ-120A	EC10N-MQ-120I	10	
		Clip	EC10C-MQ-120A	EC10C-MQ-120I	10	
		Snap	EC10S-MQ-120A	EC10S-MQ-120I	10	
GE CAM 14 patient lead wire	OEM & Models	Electrode end	REF# (AHA)	REF# (IEC)	Pieces/Set	Price code
		4mm banana	EC10B-VS-120A	EC10B-VS-120I	10	
		3mm needle	EC10N-VS-120A	EC10N-VS-120I	10	
		Clip	EC10C-VS-120A	EC10C-VS-120I	10	
		Snap	EC10S-VS-120A	EC10S-VS-120I	10	
	420101-002	CAM/Universal	EC10D-VS-120A	EC10D-VS-120I	10	
	420101-001 Used with 900178-001/900178-002 or Unimed Adapter as electrodes end (see below)	CAM/Universal	EC14D-VS-120A	EC14D-VS-120I	14	
	2001925-005(66cm)	CAM/Universal	EC10D-VS-066	66cm,no color coded	10	
	2001925-002(74cm)	CAM/Universal	EC10D-VS-074	80cm,no color coded	10	
	2001925-006(90cm)	CAM/Universal	EC10D-VS-090	90cm,no color coded	10	
	2001925-003(102cm)	CAM/Universal	EC10D-VS-102	100cm,no color coded	10	
	2001925-004(130cm)	CAM/Universal	EC10D-VS-130	120cm,no color coded	10	

GE AM4/AM5 patient lead wires	OEM & Models	Electrode end	REF# (AHA)	REF# (IEC)	Pieces/Set	Price code
		4mm banana	EC10B-AM-120A	EC10B-AM-120I	10	
		3mm needle	EC10N-AM-120A	EC10N-AM-120I	10	
		Clip	EC10C-AM-120A	EC10C-AM-120I	10	
		Snap	EC10S-AM-120A	EC10S-AM-120I	10	
	900177-001 Used with 900178-001/900178-002 or Unimed Adapter as electrodes end (see below)	CAM/Universal	EC12D-AM-120A	EC12D-AM-120I	12	
	406581-004(90cm)	CAM/Universal	EC12D-AM-090	80cm,no color coded	10	
	406581-006(130cm)	CAM/Universal	EC12D-AM-130	120cm,no color coded	10	
GE CAM14/AM4/AM5 Electrode Adapter	OEM & Models	Electrode end	REF# (AHA)	REF# (IEC)	Pieces/Set	Price code
	900178-001,AHA 900178-101,IEC 900179-201,AHA 900178-201,IEC For the above leads set Or GE original CAM14/AM4/AM5	4mm banana,10 pcs/set	EAD10-DB-A	EAD10-DB-I	10	
		4mm banana,14 pcs/set	EAD14-DB-A	EAD14-DB-I	14	
		4mm banana,12 pcs/set	EAD12-DB-A	EAD12-DB-I	12	

Electrode Adapter	OEM & Models	Electrode end	REF# (AHA)	REF# (IEC)	Pieces/Set	Price code
	4mm banana -Snap	Snap,10 pcs/set	EAD10-BS	EAD10-DB	10	
	4mm banana -Clip	Clip,10 pcs/set	EAD10-BC	EAD10-DB	10	
	Adapter-Alligator clip	Alligator clip, 10 pcs/set	EAD10-DA	EAD10-DA	10	
Philips Lead wire	OEM & Models	Electrode end	REF# (AHA)	REF# (IEC)	Pieces/Set	Price code
	PageWriter Touch: 860284, PageWriter Trim Cardiographs, 860286, 860287, 860288, 860290, 860297	4mm banana	EC10B-PW1-012A (989803129161)	EC10B-PW1-012I (989803129191)	10	
		3mm needle	EC10N-PW1-012A	EC10N-PW1-012I	10	
		Clip	EC10C-PW1-012A	EC10C-PW1-012I	10	
		Snap	EC10S-PW1-012A	EC10S-PW1-012I	10	

Philips Lead wire	OEM & Models	Electrode end	REF# (AHA)	REF# (IEC)	Pieces/Set	Price code
	989803151631/989803151641	4mm banana	EC10B-PW2-012A	EC10B-PW2-012I	10	
		3mm needle	EC10N-PW2-012A	EC10N-PW2-012I	10	
		Clip	EC10C-PW2-012A	EC10C-PW2-012I	10	
		Snap	EC10S-PW2-012A	EC10S-PW-012I	10	
Mortara Lead wire						
	9293-046-60	4mm banana	EC10B-MT-012A	EC10B-MT-012I	10	
		3mm needle	EC10N-MT-012A	EC10N-MT-012I	10	
		Clip	EC10C-MT-012A	EC10C-MT-012I	10	
		Snap	EC10S-MT-012A	EC10S-PW-012I	10	
AA style Lead wire						
		4mm banana	EC10B-AA-012A	EC10B-AA-012I	10	
		3mm needle	EC10N-AA-012A	EC10N-AA-012I	10	
		Clip	EC10C-AA-012A	EC10C-AA-012I	10	
		Snap	EC10S-AA-012A	EC10S-AA-012I	10	
Welch Allyn Lead wire						
	Monitor:CP100/CP200 Used for Welch Allyn 400294	4mm banana	EC10B-WA-012A	EC10B-WA-012I	10	
		3mm needle	EC10N-WA-012A	EC10N-WA-012I	10	
		Clip	EC10C-WA-012A	EC10C-WA-012I	10	
		Snap	EC10S-WA-012A	EC10S-WA-012I	10	

4.3 EKG Accessories-1

REF#	ECG Electrodes	Package	Price code	REF to Figures
ECL01-A-I	Adult Limb Clamps, IEC, universal	4 pcs/box		A
ECL01-P-I	Pediatric Limb Clamps, IEC, universal	4 pcs/box		B
ESE01-A-Green	Adult suction electrodes, Universal	6 pcs/box		C
ESE01-P-Blue	Pediatric suction electrodes	6 pcs/box		D
ELC01-A-Green	Adult Limb Clamps/AHA, all clamps in green	4 pcs/box		E
ECL01-P-Green	Pediatric Limb Clamps/AHA, all clamps in green	4 pcs/box		F
SEA-Blue	Adult blue bulb	6 pcs/box		G
SEP-Green	Pediatric green bulb	6 pcs/box		H
ESE02-A-Blue	Adult suction electrodes, Universal	6 pcs/box		I
ESE02-P-Gray	Pediatric suction electrodes	6 pcs/box		J
ECL02-A-I	Adult Limb Clamps, IEC, 4.0mm banana,3.0mm needle	4 pcs/box		K



4.3 EKG Accessories-2

REF#	ECG Electrodes	Package	Price code	REF to Figures
ECL03-A-I	Adult Limb Clamps/IEC (snap, grabber, 3mm needle, 4mm banana)	4 pcs/box		A
ECL03-A-I	Adult limb clamps/IEC(3mm needle, 4mm banana)	4 pcs/box		B
ECL04-A-I	Adult limb clamps/IEC/smooth surface(3mm needle, 4mm banana)	4 pcs/box		C
ECL04-A-I	Adult limb clamps/IEC/Agcl plating (multifunctional)	4 pcs/box		D
ECL05-A-I	Adult limb clamps/IEC/Nickel plating (multifunctional)	4 pcs/box		E
ECL06-A-Blue	Adult limb clamps/blue (Grabber, snap, 3mm needle, 4mm banana)	4 pcs/box		F
ESE02-P-Gray	Adult suction electrodes/gray (universal)	6 pcs/box		G
ECL02-P-I	Pediatric limb clamps/IEC/Agcl plating (Snap, grabber, 3mm needle and 4mm banana)	4 pcs/box		H



5. EEG&EMG Lead Wires



1.Cup 1-10 Ld (1082PD-10A)



2.Discoid Electrodes 1-10 Ld(1082PD-10B)



3. Clip 1-10 Ld(1082GD-10A)



4. Spring Clip 1-10 Ld(1082GD-04B)



5.Alligators Clips 1-10 Ld(1082GD-06C)



6.Alligators Clip 1-10 Ld(1082GD-07D)



7.Alligators Clip1-10 Ld(Red)(1082GD-10E)



8.Alligators Clip 1-10 Ld(Black) (1082GD-10F)



9.Hat 1-10 Ld(1082PD-10C)



10.Clip Electrodes 2 Ld(1082GD-02G)



11.Ear Clip Electrodes 2 Ld(1082PD-02A)



12. Discoid Electrodes 1Ld(Yellow) (1082PD-01B)



13. Black Cup Electrodes 1Ld(Red) (1082PD-01D)



14. Sterile Disposable Subdermal
Needle Electrodes 2Ld(UTP wires) (1082ND-02B)



15. Sterile Disposable Subdermal
Needle Electrodes 1Ld(1082ND-01A)