

# 2025-2026

Customers First... Friends for Life

# **TERMS AND CONDITIONS**

Payment Terms.....Net 30 days

Credit Cards: Visa, Master Card, Discover & American Express \*\*fees apply\*\*

#### **SHIPPING**

All shipments are F.O.B. shipping point. Shipping charges will be charged on all invoices when applicable. All shipping UPS ground service. Shipments of size and weight, that do not meet UPS guidelines, will be handled by freight line transportation. Handling fee will be applied when receiver is 3rd party.

Please Note: Carriers are responsible for losses and all damages to shipped products. Acceptance of a shipment as delivered is your acknowledgment that the shipment was received in good condition. In case of visible damage, do not sign for shipment or accept shipment unless the delivering agent notes the nature and extent of the damages. Then notify the delivering carrier and file claim for damages at once. If you find concealed damages, follow the same procedure of contacting the delivering carrier. Claims for damage visible or concealed cannot be filed by us. Any shortages in shipment must be reported by Purchasing and/or Receiving Department within 10 days of delivery.

#### RETURNED ITEMS

No material will be accepted for credit without prior approval. Please call for an RGA (Return Goods Authorization) number. Items must be returned in good condition, (i.e., not used, damaged or altered). For any item returned and, upon inspection, determined to have been used, altered or damaged, the customer will be contacted before settlement. Returns are subject to the manufacturer's restocking charge where applicable. Custom-made positioners are not returnable. All returns should be made within 30 days after receipt of order. No C.O.D.'s accepted. No cash refunds. Credits are issued against future orders. Credits are valid for 1 year. Any and all credits will be voided if not used within one year of invoice date.

#### **TAXES**

Sales tax charged: AL, CA, FL, KY, IL, LA, MD, NC, NV, NY, & TX.

Customer is responsible for supplying tax exemption form via mail, fax, or email.

PRICES SUBJECT TO CHANGE WITHOUT NOTICE

#### MINIMUM ORDER

A \$15.00 handling fee will be added to ordered \$30.00 and below.

#### LIMITED WARRANTY

VIP, Inc. assumes no liability for consequential damages of any kind. Upon acceptance of products, customer assumes liability for all direct or indirect results and consequences of the use or misuse of any product by anyone. This warranty is expressly in lieu of all other warranties, guarantees, obligations, or liabilities expressed or implied by VIP, Inc. or representatives, with all statutory/implied warranties, other than title, hereby negated and excluded.


P.O. Box 10853 Pensacola, FL32524

1-800-326-4761 • 850-476-1078 • 24 Hour Fax 850-476-8426 vanarsdale.ar@gmail.com (Accounting) vanarsdale.cs@gmail.com (Customer Service)

# **TABLE OF CONTENTS**

A		D		Н	
Abduce Knee Positioner Set Accessory Holders/Caddy Adhesive Remover Wipes Air Freshener Alcohol Prep Pads Alginate Impression Material Alloy Alloy Dispenser Aloe Gel Aluminum Cooling Tray Aluwax Anatomical Drawings Anatomical Rubber Stamps Anti-Static Spray Apron-PVC-Nylon Aquaplast Pellets	29 16 14 9 16 41 3 3 9 4 42 51 & 52 52 53 2 42	Digital Level, pocket/w magner DisCide D-Wipes Towels Disinfectant/Sanitizer Disposable Polyethylene Overboots Disposable Tissues Disposable U-Frames Drill Dust Mask  E Elasto-Gel Bolus	t 21 8 7 8 &9 44 18 47-50 5 2	Hammers Hand Sanitizer (Aloe Gel) Hand Soap (D-Lead) Head Rest (Plastic/Foam) Heat Gun Heating Pans HEPA Vacuum Cleaner Hot Wire/Cutting Wire  I Ink (Tattoo) Ink Cups Immobilization System Impression Materials Incontinence Clamps	5 9 7 26 5 & 46 46 7 1 16 49 &, 50 42 38
Thermoplastic Arm Protector Sleeves Arm Shields Arm & Shoulder Strap Attenuate Arm Positioner Aviation Snips  B	47&48 2 45 27 28 5	Electron Cone (Inserts) Electron Cone (Templates) Electron Foam Cutter Electron Styrofoam Elevated Arm Positioner Exact Enhance Bar Exact Lok-Bar (Varian) Eye Shields Eye Shield Sterilization Tray	15 15 1 1 19 35 35 17	Infection Control I.V. Stand IZI Port Skin Markers - Call Fo	8 55
BB's and Wire Badge Board Belly Boards Bite Pieces Blocks (Oncology Styrofoam) Bolus Material Bolus (Solid) BOLX	10-14 19 30 43 1 40 - 42 40 41	F Face Protector-Shield/ Visor Feet Immobilizer File (8" half round) Flexible Curve Rulers Foam Blocks Foam Cutter Foam Cutter Wire Floor Mats	56 2 44 5 25 1 1	Lead BB's Lead Letters, Numbers Lead Markers Lead Simulator Set Lead Letters/Storage Board Lead Sheets Lead Test Kit Lead Weight Lead Wire Leg Immobilizer Levels	12 & 24 24 12 24 24 5 7 4 12 33 21
Calipers-Curved, Flat Blades	00	Foam Positioners/Soft	6 32, 33, 34	Lubricating Jelly	22
Carfusion Dye,Holder,Sticks Carts (Heavy Duty) Cassette Grips Chin Straps China Markers Clamps (penile) Contour Wire Cooling Trays	20 16 4 24 39 14 38 25	Glass Bead Warmer Gloves (Neoprene) Gloves (Latex) Powder Free Glue Dots	38 2 41 36	Jawlocks <b>M</b>	43
CT Markers CT Rectal Markers CT Vaginal Depth Scale	4 10 & 11 22 22	Goggles Goniometers Grip Rings	2 25 31	Magnification Rulers Marking Pens Mask (Respirator) Masking Tape Mats (Safety/Comfort) Tacky Mayo Stand Medi-Air Metal Alloy	25 14 2 5 6 55 9 3

Miscellaneous Mold Room Items Mold Release Mole Stick	9 3 50	Snips Soft Sponge Soldering Iron Stools (Step)	5 50 5 55
N		Styrofoam Blocks Superflab Bolus	1 40
Needles (for Tattoo)	16	Superstuff Bolus Suremark BB's & Wire 10, 11,	40
0		Carcinant BB 3 a Wile 10, 11,	12, 0 10
Odor Eliminator Overhead Arms Positioner	9 28	T Table Pads (Custom Made)	54
Patient Transfer Board Penlight (Disposable)	53 14	Table Straps Tacky Mats Tape (Radiation Warning) Tattoo Accessories Technician 20cm Rulers	45 6 37 16 25
Pens & Pencils Pink Dental Wax Pitcher - Insulated Pliers	14 42 3 5	Testicle Shield Stand	15 23 Available 23
Pocket Rulers Positioners, Abduce Knee Positioners, Prone Positioners (Soft) Prone Pillo 32,	25 29 27 33, 34 30	Thermometer (Digital/Glass) Thermoplastic Materials Timo Head Supports Tongs Tongue Depressors	3 47,48 26 3 43
Prone Face Positioner	27	Tools Towels, Disposable (D-Wipe) Transfer Boards, Straps & Wall Rack	5 7 53
RAD Guide Templates RADTAT - Suremark Rectal Markers Reprosil Rubber Mallet Rubber Stamps Rulers (Stainless Steel & Plastic) Rust Cleaner (Water Bath)	15 13 22, 23 42 5 52 25 46 13	U-frame Wall Hanger U-frame Immobilization Systems U-frames Unlimited Scissors Utility Wax	48 49, 50 47, 48 46
S			42
Safety Floor Mats Sanizide Plus Disinfectant Scissors Screw Driver Set Securline Pens Shadowform Rectal Markers Shadowform Vaginal Markers Sharpies Shields (Testicular) Shoulder Retractor, Straps	6 8 46 5 14 23 23 14 23 27	V Vacuum HEPA Vaginal Depth Scale & Holders Vaginal Marker, Shadowform Velcro One-Wrap	7 22 23 44
Signs (Safety/Warning) Silicone Sealant Silicone Mold Spray Silky Rock Silverman Supports Simulation Lead Letters Skin Markers Skimmer Spoon	37 50 3 41 26 24 10-13 3	Wall Rack - Transfer Board Warning Signs, Tape Waterbaths Waxes, (Pink, Utility, Aluwax) Wing Board	53 37 46 42 19





Remit to address P.O. Box 10853 Pensacola, FL 32524




Accounting: vanarsdale.ar@gmail.com

Customer Service: vanarsdale.cs@gmail.com



Phone: 1-800-326-4761 1-850-476-1078





# "Every small gesture combines into a flood of caring and support"

## Personalized Sympathy Cards

All cards are the same size (5" x 7") portrait or landscape All cards are the same price.


Blank envelopes included with order -

Printed envelopes (on flap) (3 lines)

Your choice of verse, font, and ink colors (shown) for inside of cards











Contact us for a copy of our complete Sympathy Card collection or visit our website:


www.vanarsdaleproducts.com 800-326-4761 vanarsdale.cs@gmail.com

# Certificates to Recognize Your Latients' Achievements

- Border Colors: Royal Blue, Burgundy, Grey, and Gold
  - \*Teal is available with additional fees\*
- Borders and text are inter-changeable
- Any text may be modified or an original text may be submitted for printing
- Please provide the following information when ordering:
  - 1) Border number, 2) Border color, 3) Text choice,
  - 4) Name of facility as it should read, 5) Quantity

Catalog # CON-16





Text A

# Certificate of Merit Appreciation

Be It Declared That

HAS COMPLETED THE PRESCRIBED COURSE OF RADIATION THERAPY WITH THE HIGHEST DEGREE OF COURAGE, DETERMINATION AND GOOD NATURE. WE APPRECIATE THE CONFIDENCE PLACED IN US AND THE OPPORTUNITY TO SERVE YOU.

#### REGIONAL SOUTHWEST MEDICAL CENTER

Radiation Oncology Department Physicians and Staff

Text B

# REGIONAL SOUTHWEST MEDICAL CENTER Radiation Oncology Staff

hresents

# Award of Achievement

It's nice to share in your success To see what you've achieved You made it happen all because You worked hard and believed.

Text C

#### Regional West Medical Center DEPARTMENT OF RADIATION THERAPY

Be It Declared to All Present That

HAS COMPLETED THE PRESCRIBED COURSE OF RADIATION THERAPY WITH A HIGH ORDER OF PROFICIENCY IN THE SCIENCE AND ART OF BEING CHEERFUL, OUTSTANDING IN HIGH COURAGE, TOLERANT AND DETERMINED IN ALL ORDERS GIVEN.

#### **CERTIFICATE OF MERIT**

By Order of The Diploma Committee

# PREMIUM PRECISION SURFACED FOAM

<b>BLOCK SIZE</b>	<b>BOX COUNT</b>	CATALOG #
8"×8"×3"	13	870-3008
10" × 10" × 3"	13	870-3010
12" × 12" × 3"	13	870-3012

\*\*Additional sizes avaible. Contact us for a quote\*\*

# Adjustable Heat Foam Cutter


The Adjustable Heat Foam Cutter has a slippery work surface that is Alu Corbond compound makes it easier to push Styrofoam around the large table. The solid aluminum over arm has a clearance of 5 1/2' for cutting Styrofoam blocks from 3/4" thick to 3" thick. The holder and wire spool may be shifted adjusted along the over arm to enable miter cutting.

#### Features:

- \* An LED indicates operation and thus reminds you not to touch the hot wire to prevent burnt fingers (the wire heats to the maximum temperature in less than 1 second)
- \* Printed grid and protractor on worktable assist with division and cutting.
- \* On/Off switch
- \* Electronic variable temperature control knob.
- \* Lockable feed bar/fence for angle ciuts
- \* One spool of NiChrome wire 0.008" diameter x 98 ft.
- \* One year warranty

#### **NiChrome Wire**

Catalog # VIP-869-9101





**Adjustable Heat Foam** Cutter

Catalog # 869-910



# **Electron Cone Sizes**

BLOCK SIZE	<b>BOX COUNT</b>	CATALOG #
8" × 8" × ¾"	49	870-0708
10" × 10" × ¾"	49	870-0710
11 ½" × 11 ½" × ¾"	49	870-0711
8" × 8" × 1"	38	871-1008
10" × 10" × 1"	38	871-1010
12" × 12" × 1"	38	871-1012

# **MOLD ROOM SAFETY PRODUCTS**

# **Protection at its Best!**

#### **FACE PROTECTION**

Protecto-Shield® Headgear style V5, plastic brow guard ratchet adjustment and vinyl sweatband. Constructed of tough thermoplastic materials, our molded Protecto-Shield® headgear is designed with a special brim to prevent fluids from splashing down inside the visor. Protecto-Shield® visor M86 is constructed of propionate material to provide the required style and tint for almost any application. Light weight and attaches easily with our combination of positive lock and snap fasteners.

Protecto-Shield® and Visor

Size: 8" x 13" Color: medium green Catalog # VIP-176

#### 3-M DUST AND FUME/MIST RESPIRATOR

TC-84A-1298 approved

Used in applications such as soldering, metal cutting or metal pouring where a metal fume hazard exists. Approved for lead dust and fumes. (2) Easily fits under face protector. Exhalation valve improves worker comfort and acceptance.

Catalog # VIP-177

#### **ARM PROTECTOR/SLEEVES**

Disposable — Have elastic at wrist and elbow. Color: white Catalog # VIP-194 24/pkg

#### **PVC/NYLON APRON**

Full front body protection. Rugged PVC apron with nylon backing. Resists most acids, caustics and punctures. Construction features: overlapped and sewn edges with reinforced. heat-sealed vinyl grommets and cotton ties.

Size: 35" x 45" Color: green

Catalog # VIP-179

#### **CLEAR GOGGLES**

Fit over glasses Catalog # VIP-SEE

NIOSH . . . National Institute for Occupational Safety and Health

OSHA . . . Occupational Safety and Health Administration MSAH . . . Mine Safety and Health Administration

PEL . . . . Permissible Exposure Limit

(1) NIOSH/MSHA approved for dusts and mists with PEL not less than 0.05 mg/m3

(2) NIOSH/MSHA approved for dust fumes and mists with a PEL not less than 0.05 mg/m3 and for radon daughters attached to dusts, fumes and mists.



#### **NEOPRENE BLEND GLOVES**

Neoprene blend over natural latex provides excellent protection against acids, alcohols, animal fats, caustics, salts, petroleum-based solvents and detergents. Gloves have 100% cotton flocking to absorb perspiration. Diamond embossed for excellent grip on wet or dry surfaces. Twenty-eight mil thick for extra protection against chemicals, abrasions and punctures.

Length: 13"

Catalog # VIP-178LG (size 10) Catalog # VIP-178MED (size 9) Catalog # VIP-178SM (size 8)

# **ALLOY DISPENSER BUILT TO LAST!**

# Catalog # VIP-878-060

1.5 gallon capacity 120 VAC, 8 AMPS



#### Features Include:

- Alloy temperature stability ± 4°
- Microprocessor PID temperature control
- "J" thermocouple
- Open thermocouple sensor will shut unit down
- Temperature is displayed in fahrenheit or celsius
- Easy view temperature controls
- Set point reading is green and alloy temperature is red
- Insulated container preserves power consumption
- Exterior stays cool even after hours of operation
- Easy to clean round container
- Stainless steel corrosion-resistant interior and exterior
- One year manufacturer's warranty

# METAL ALLOY FOR RADIATION SHIELDING





•	
Bismuth	50%
Lead	27%
Tin	13%
Cadmium	10%
Melting Point 158°F-70	°C
Density: .3458 lb./cubic inch	

Catalog # VIP-879-158 50 lb. box

Composition:



Catalog # VIP-879-202 50 lb. box

Density: .3500 lb./cubic inch

Due to fluctuation of the metal market, prices subject to change. Call for quote.

# **MOLD ROOM ACCESSORIES**

CADFREE



Catalog #

VIP-184 Digital Thermometer (pocket size)

-40°-450°F

VIP-878-152 Thermometer (glass)

0°-230°F Red Spirit Filled

VIP-878152C Thermometer (glass)

-20°-50°C Blue Spirit Filled

VIP-185 Ladle, stainless steel

VIP-186 Skimmer Spoon, stainless steel
VIP-4TM Metal Tongs w/plastic tip

VIP-4TM Metal Tongs w/plastic tip VIP-188 Wire Cleaning Brush



STAINLESS STEEL PITCHER

Excellent for pouring metal alloy.

Easy, safe handle.

Three different sizes.

Catalog # VIP-99-20 oz. Catalog # VIP-99-33 oz. Catalog # VIP-99-48 oz.



Food grade, non-staining, non-CFC, odorless. Simply spray on molds to prevent parts from sticking. Can be used on 158 and 203 metal alloy.

Net Wt. 13 oz.

Catalog # VIP-110



# **COOLING PLATES & STYROFOAM WEIGHT**



#### **BLOCK COOLING TRAY**

Constructed of aluminum, uses an external water supply to rapidly cool alloy blocks. Unit can be attached to a faucet and sink or be plumbed directly into current plumbing. 1" x 12" x 24"

Catalog # VIP-878-271



#### STYROFOAM LEAD WEIGHT

Weight to assist in block stabilization.
16 lbs
Catalog # VIP-878-274

# **UTILITY CARTS**



# NEXEL STAINLESS STEEL UTILITY CART

Noncorrosive cart with 20-gauge stainless steel shelves and 16-gauge uprights. Three shelves, 12" clearance between shelves. Top two shelves have lips on three sides for secure transport; bottom shelf is flush for oversized items.

400lb. capacity

Catalog # VIP-SSC3



# LAKESIDE STANDARD DUTY PLASTIC UTILITY CART

Stain and odor-resistant material maintains sanitation. Durable impact-resistant polyethylene construction.
Three Shelves.

Self Size: 25-1/2"W x 16"D Shelf Clearance: 12-1/8" 300lb. weight capacity

Dimensions: 31-3/4"W x 38"H x 17-1/2"D 4' all swivel casters - **Assembly required** 

Catalog # VIP-BPC3

# **LEAD SHEETS**



1/16" (1.6mm) x 12" x 24" Catalog # VIP- 062L

1/8" (3.2mm) x 12" x 24' Catalog # VIP-125L

1/16" x 12" x 12" Catalog # VIP-28L

1/8" x 12" x 12" Catalog # VIP-29L

1/24" (1.1mm) x 12" x 24" Catalog # VIP-1224124

# MOLD ROOM MISCELLANEOUS ITEMS

#### **MASKING TAPE**

High Temperature, 3" wide x 60 yards long Catalog # 7627A27

High Temperature, 1" wide x 60 yards long **Catalog # 7627A23** 

#### **PISTOL HEAT GUN**

650 degree F Temperature Setting. 120 Volts, 475 W, three prong. \*\*Adjustable Stand Included\*\*

Catalog # 3414K11

Replacement Element for Silent Heat Gun Catalog # 341412

#### **ELECTRIC SOLDERING IRON**

Heavy Duty, 120 Vac, 200/260 Watts, 900/1000° F **Catalog # 7728A1** 

## **LIGHT-DUTY SOLDERING IRON**

80 Watts, 900° F, Chisel Tip Design Catalog # 7662A697

**SLIP JOINT PLIERS** — 6 3/4"

Catalog # 45378

**EIGHT INCH HALF ROUND FILE w/handle** 

Catalog # 931311

**AVIATION SNIPS** — Strait Cut

Catalog # 942773

**ELECTRIC DRILL** — Variable Speed, 3/8"

Catalog # 10114

RUBBER MALLET

Catalog # M45789

NAIL HAMMER — 16 oz.

Catalog # 38126

**BALL PEEN HAMMER** — 4 oz.

Catalog # 38467

**SCREWDRIVER SET** — Six Slotted & four Phillips

Catalog # 941809

# **Safety/Comfort Mats**



# **ERGO MAT**™

- 1/8" non-skid corrugated surface adds safety to work station areas.
- <sup>3</sup>/8" buoyant sponge underside provides maximum anti-fatigue relief.
- Permanent lamination of materials for extended service life.
- <sup>1</sup>/<sub>2</sub>" overall thickness provides deluxe comfort for tough industrial applications.
- Beveled edges on all four sides reduce tripping hazards.
- \* Not guaranteed against spiked heels or stool legs.

Catalog # EM-474-3-5 3' x 5'

# COMFORT-EZE™

Offers fatigue relief for the person who stands all day

- Molded single piece construction is tough and durable; withstands abuse.
- Resilient rubber compound provides superb cushioning to relieve fatigue.
- Uniquely engineered 5/16" knob pattern allows air to circulate beneath feet.
- Double sided design extends the service life of the mat.
- Safety edges are beveled on all four sides for extra safety.
- 3/8" thick; Color is black.

Catalog # 447-30-10 30" x 10'





# FRAMELESS TACKY MATS

Blue tear-off sheets

These low-profile mats are made up of layers of adhesive-coated sheets. Each sheet has a nonleeching, non toxic, FDA compliant antimicrobial agent that inhibits the growth of bacteria, mold, and mildew for the most effective contamination control offered. When the top sheet gets dirty, simply pull the numbered corner tab to remove it and expose the next clean sheet. For indoor use only.

#### 30 layers per mat

Catalog	# TM-1836	18" x 36" x .05"
Catalog	# TM-1845	18" x 45" x .05"
Catalog	# TM-2545	25" x 45" x .05"
Catalog	# TM-3645	36" x 45" x .05"
Catalog	# TM-3660	36" x 60" x .05"

#### 60 layers per mat

Catalog # TM-1836-60	18" x 36" x .09"
Catalog # TM-1845-60	18" x 45" x .09"
Catalog # TM-2545-60	25" x 45" x .09"
Catalog # TM-3645-60	36" x 45" x .09"
Catalog # TM-3660-60	36" x 60" x .09"

# REUSABLE STICKY-SURFACE MATS

When these sticky mats get dirty, don't throw them away—simply use the cleaner (sold below) to restore the surface's original dirt- and dust-grabbing characteristics. Made of cut- and scratchresistant cushioned vinyl, mats are not affected by heat, cold, or humidity. An adhesive backing keeps them in place. For indoor use. Color is blue.



Catalog # TM-2424R 24" x 24" x 1/8" Catalog # TM-2448R 24" x 48" x 1/8" Cleaner for Reusable Mats

1 qt. Container Catalog # TM-RC

# Four Gallon HEPA Dry Vacuum



#### Catalog # UZ-934

Note: For dry applications only.

#### **FEATURES**

- Four gallon capacity means fewer bag changes
- Powerful suction for more thorough cleaning
- Lightweight, only 13 ½ lbs. for easy portability
- Microstat filtration prevents exhausting particulates back into the air
- · Will not tip over due to a low center of gravity

#### **SPECIFICATIONS**

Two stage fan and easy-access HEPA filter. Includes: 35 foot long power cord, 6 foot long hose, floor tool, and a two-piece aluminum wand.

Horsepower: 1.3 Amps: 8.7 cfm: 74

Voltage: 115-120, 1000 watts Dimensions: 15" x 15" x 22"; 0.16oz

Available Accessories: VIP-48316 Replacement HEPA filter

VIP-126913 Replacement Dust Bags: Pkg of 5

VIP-121698 Exhaust Filter



## **D-Lead Test Kit**

#### Detects the Presence of Lead on Skin in 3 Seconds

D-Lead Test Kit is an easy and inexpensive personal hygiene training tool. The economical, reliable, and easy-to-use kit safely tests your hands, skin, and all wettable surfaces for fast detection of lead from metal alloys, oxides, salts, and dusts. A brilliant yellow color appears immediately when lead is present. Sensitive detection limit can recognize concentrations as low as 20mg total lead. No chemicals to mix and no unpleasant odors.

Catalog # VIP-72368 D-Lead Test Kit

# **D-Wipe Towels**

Quickly and safely removes multiple layers of metals and contamination from the skin and all wipeable surfaces. 6" x 8".





# **Hand Soap**

Unique cleaners quickly and completely remove arsenic, cadmium, chromium, lead, mercury, silver, and zinc from workers' skin and other washable surfaces. Will not irritate skin or interfere with waterwaste treatment systems.

Catalog # VIP-23394

Hand soap, 8 oz. bottle with flip cap

Catlagg # VIP-84151

Hand soap, 1 gallon refill bottle



# SaniZide Plus™ Germicidal Solution

SaniZide Plus is a Quaternary ammonium compound formulated for use as a broad spectrum disinfectant deodorizer on environmental surfaces. Premixed and easy to use, let SaniZide Plus help you comply with the OSHA Bloodborne Pathogens Standard which requires the use of an "appropriate disinfectant" that is tuberculocidal and virucidal against HIV-1.

- Ready-to-use hospital grade disinfectant
- Kills HIV and TB
- · Virucidal, Tuberculocidal, Fungicidal, Bactericidal
- Controls mold/mildew
- · Neutralizes odors and freshens in one step
- Helps prevent cross contamination
- EPA registered

Catalog #422-3691
SaniZide Plus Disinfectant • 32 oz. Trigger Sprayer

# DisCide ULTRA — Palmero Surface Disinfectant

Effectively kills HBV, MRSA, VRE, HIV-1 (AIDS virus), herpes, adenovirus, staphylococcus aureus, pseudomonas, salmonella, E. coli and is fungicidal in one minute! DisCide ULTRA is also tuberculocidal. DisCide ULTRA is the solution to all your infection control problems, combining a powerful formula that kills effectively and quickly.

	Catalog#	<u>Unit</u>
Towelettes 6" x 6¼" (160 count)	VIP-8515	1 each
10" x 10" (60 count)	VIP-8523	1 each
Gallon Refill	VIP-8580	1 each
Quart w/Sprayer	VIP-8572	1 each
Canister Holder	VIP-1394	1 each

# - U it's safe enough for two-year-old Brandon, . .









# Medi-Aire® Rejuvenates the Air

A spray or two of Medi-aire's concentrated formula terminates unpleasant odors present in hospitals, nursing homes, or clinics.

Patients and staff will appreciate Medi-aire.

Available in easy to carry 1 oz. spray bottles and economical 8 oz. re ill bottles.

Catalog # VIP-140F Fresh Air 1 oz. Spray Bottle

Catalog # VIP-148F 8 oz. Refill Bottle Catalog # VIP-148L 8 oz. Refill Bottle







# **Bye Bye Odor Eliminator**

- Reduces airborne bacteria for cleaner air
- Neutralizes offensive odors and freshens the air
- For use in industrial, insititutional, or house-hold settings
- Non-toxic for household safety
- Mild and unobtrusive scent
- 7.5 oz. spay bottle

# Catalog # VIP-142



# **Fresh Odor Eliminator**

MEDLINE, Fresh Odor Eliminator, 2oz spray

Catalog # VIP-141F





# **Precision Skin** k Marking and Labeling **Systems**

# **CLEARLINE**™ Semi-lucent Scar Markers

Clearline™ low-density, non-metallic, semi-lucent scar markers provide distinct outlining of scars and other areas of concern without obscuring micro-calcifications in mammography.















#### Clearline™ Labels Available in Six Sizes

CL-8—0.8mm semi-lucent line—packaged with 336 cm per box CL-10—1.0mm semi-lucent line—packaged with 336 cm per box CL-15—1.5mm semi-lucent line—packaged with 320 cm per box CL-20—2.0mm semi-lucent line—packaged with 320 cm per box CL-25—2.5mm semi-lucent line—packaged with 203 cm per box CL-30—3.0mm semi-lucent line—packaged with 165 cm per box





ensure easy and accurate placement without contaminating the adhesive. Packaged with 110 labels per box. \*PM-400 Series uses an alternate material and adhesion for more humid environments and particular skin surfaces.

The newest members of the PortalMark™ label family offer additional options for delineating radiation treatment areas. Their premium medical grade adhesive is hypo-allergenic and latex-free, resulting in less irritation for the

patient. Constructed of waterproof, yet breathable material for maximum performance and comfort, PortalMarks hold strong for up to two weeks. Our unique Suregrip™ applicator discs

Catalog # PM-400

Catalog # PM-425

Catalog # PM-450 Black 110 labels per box Black 110 labels per box Blue 110 labels per box Catalog # PM-475 Blue 110 labels per box



Vision Line™ products offer the perfect combination of tack surfaces and easy-peel corners, keeping the markers correctly positioned with minimal discomfort to the patient upon removal.

These specially formulated NON-METALLIC markers are clearly detectable on film, reducing read time.

Visionmark™ CT markers are clear and consistent with virtually no scatter.

Our patent pending Vision Line™ products are safer for you, your patients and the environment.

# CT-Mark™ and CT-Line™


The CT-Mark<sup>™</sup> and CT-Line<sup>™</sup> have been created specifically for CT imaging. This specially formulated NON-METALLIC material is the result of years of research and development, resulting in a virtually scatter-free opaque marker that is perfect for nearly all CT marking applications.





CT-Mark Available in 2.3mm and 4.0mm ball sizes. CT-23 Packaged with 110 markers per box. CT-40 Packaged with 110 markers per box.

# **CT-Line**™




CT-W10 - 1.0mm line - packaged with 300 cm per box. CT-W20 - 2.0mm line - packaged with 300 cm per box.



# Visionmark™

The Visionmark<sup>™</sup> is the next generation in general use skin markers. These specially formulated NON-METALLIC markers perform consistently and accurately in a variety of applications. An excellent tool for distinguishing between nipple shadow and lesion, the Visionmark<sup>™</sup> reduces the need for repeat examinations. Easy Peel Tabs.

Available in **2.0mm** and **2.5mm** ball sizes. Packaged with 110 markers per box.



# Visionmark™CT

These specially formulated NON-METALLIC markers reduce read time on CT scans with virtually no scatter. These sizes and ball densities are perfect for all CT skin marking applications. Easy Peel Tabs.

Available in **3.0mm**, **4.0mm**, and **5.0mm** ball sizes. Packaged with 50 markers per box. Packaged with 50 markers per box.

## Catalog #

V-20, 25, CT-23 V-30 V-40, 50, CT-40, CC-15, CC-25 CT-W10 CT-W20



# **Markers and Accessories**

## Relief Tabs™

sure mari





sure

**Double Sided Foam Tape** 



SL-15RT

SI -20RT

SL-25RT

Give your patients some relief with Suremark's new Relief Tabs™. The unique adhesive-free center won't stick to sensitive areas of the skin. Available in 1.5mm and 2.0mm and 2.5mm ball sizes. Packaged with 110 markers per box.



# **Storage Box**

This convenient box holder provides a permanent visible storage area for your Suremark X-ray Skin Marker boxes. Holds five or ten standard boxes or two R and L boxes. Place on a desktop or mount to a wall or imaging machine.

BH-5 BH-1

Holds five boxes

BH-10 Holds ten boxes

# These reusable ½" x 1" double sided tape strips are ideal for fastening film markers to cas-

settes. Packaged with 110 strips per box. SL-DT1



## Suremark®













The most common style for general use,
Suremark 1.0mm, 1.2mm, 1.5mm, 1.8mm,
2.0mm and 2.5mm labels are perfect for
distinguishing between a nipple shadow and a lesion,
thereby reducing costly repeat examinations.
Packaged with 110 markers per box.

# **Powermark**™







This design maintains all of the advantages of our regular Suremark labels, with larger lead ball sizes for visibility in radiation therapy. Available in 3.0mm, 4.0mm, and 5.0mm ball sizes. Packaged with 50 markers per box.

**SL-30** 

SL-40

SL-50

# Suremark® Tabs™

This design maintains all of the advantages of our regular Suremark\* labels, with the added benefit of a "tabbed" feature that enables natural and conforming placement on more prominent nipples or where more skin contact is preferred.

Available in both 1.5mm, 2.0mm and 2.5mm ball sizes. Packaged with 110 markers per box.







#### Catalog #

SL-15RT, 20RT, 25RT SL-DT1 SL-10, 12, 15, 18, 20, 25 SL-30, 40, 50 SL-15T, 20T, 25T

# Suremark® Wire™


This valuable tool is used for marking scar tissue, outlining palpable nodules, moles, entry points, and in CT scans.

Available in .3mm, .4mm


Packaged with 336 cm per box

Available in .8mm, 1.2mm, and 1.5mm wire sizes.

Packaged with 336cm per box.



(.3mm wire) (.4mm wire) (.8mm wire) (1.2mm wire) (1.5mm wire)



The SL-Wire can be placed in any orientation to mark scar tissue or placed vertically (as shown above) to provide an accurate positioning reference on CT scans.

#### Catalog #

SL-W3, SL-W4, SL-W8, SL-W12, SL-W15, TM-25

# **Thermoplastic Markers**







Theramark TM markers stay securely where you put them throughout the entire treatment process. They are repositionable during the initial set up, giving you time to get the markers in just the right spot. The adhesive bonds stronger as time goes on. Theramarks adhere directly to immobillization masks without having to use tape.

Only for use on immobilization masks.

Not medical grade adhesive.

Adhesive made with natural rubber latex.

Packaged with 110 markers per box.

# **RadTat Skin Marking**

Catalog # RTC - 100





**RED** 

Catalog # VIP-33RD

**BLACK** Catalog # VIP-33BK **GREEN** 

Catalog # VIP-33GR

**BLUE** 

Catalog # VIP-33BL



# **CHINA MARKERS**

**BLACK** Catalog # VIP-34BK

**GREEN** 

Catalog # VIP-34GR Catalog # VIP-34YW Catalog # VIP-34BL

**RED CRIMSON RED** Catalog # VIP-34WH Catalog # VIP-34RD Catalog # VIP-34CR



# **SECURLINE**

Fine tip, non-toxic, gentian violet ink. High visibility. Prep resistant.

Catalog # VIP-440



# **DISPOSABLE PEN LIGHT**

Catalog # VIP-PL Pen Light

Catalog # VIP-PKT **Pocket Protector** 



# PAINT MARKING PEN







Fine line, oil based marking pen. Writes on anything. Long-lasting . . . ink will not clog or thicken. A must for every department.

Catalog # VIP-45BK (Black)

Catalog # VIP-45WH (White)

Catalog # VIP-45GR (Green)

Catalog # VIP-45RB (Royal Blue)

Catalog # VIP-45AQ (Aqua)

Catalog # VIP-45OR (Orange)

Catalog # VIP-45VI (Violet)

BLACK, WHITE AND COLORS: **GREEN** 

ROYAL BLUE AQUA ORANGE VIOLET

# **CardinalHealth™ ADHESIVE REMOVER**

CardinalHealth™Adhesive Remover Wipes are excellent for removing paint marks from Deco Pens. In the event that lines need to be changed, the lines can be removed quite easily (without any unnecessary scrubbing). Adhesive Remover Wipes do leave a slightly oily residue after removing the preexisting lines. We strongly recommend that you clean the area gently with a little alcohol before remarking.



**Adhesive Remover Wipes** 75 individual packets per box



# **ELECTRON CONE INSERTS**

Inserts fabricated from 1/4" thick acrylic. Provides a flat surface for marking on aperature templates.


**Siemens** 


# Varian Catalog # Description VIP-125-6 6 cm (7.6 x 7.6 cm) VIP-125-10 10 cm (11.3 x 11.3 cm) VIP-125-15 15 cm (16.1 x 16.1 cm) VIP-125-20 20 cm (20.9 x 20.9 cm) VIP-125-25 25 cm (25.7 x 25.7 cm)

Catalog #	Description
VIP-126-6	(6 x 5 cm)
VIP-126-10	(11.5 x 11.5 cm)
VIP-126-15	(15.5 x 15.5 cm)
VIP-126-20	(20 x 20 cm)
VIP-126-25	(24.6 x 24.6 cm)

Varian 2	2100-C Inserts
Catalog	# Description
11-400	(8.9 x 8.9 cm)
11-402	(11.8 x 11.8 cm)
11-403	(16.5 x 16.5 cm)
11-404	(21.2 x 21.2 cm)
11-405	(25.9 x 25.9 cm)







# **ELECTRON CONE TEMPLATES**

- Templates fabricated from clear mylar
- Simplifies electron cone aperature block fabrication
  - -Recording field outlines is easier
  - CAX clearly marked on templates
  - Sold separately (in sets of 25 each)



<u>Varian</u>		
Catalog #	Description	
VIP-127-6	6 cm	
VIP-127-10	10 cm	
VIP-127-15	15 cm	
VIP-127-20	20 cm	
VIP-127-25	25 cm	


<u>Siemens</u>
Catalog #
VIP-128-6
VIP-128-10
VIP-128-15
VIP-128-20
VIP-128-25

Varian 2100-C		
Catalog #	Description	
11-400T	6 cm	
11-402T	10 cm	
11-403T	15 cm	
11-404T	20 cm	
11-405T	25 cm	

# **RAD GUIDE TEMPLATES**

Clear mylar templates for recording blocked fields. Uniformly cut and printed for convenient use.

Catalog #	<u>Description</u>
VIP-101	100 ea. 8" x 6" + Target
VIP-109	100 ea. 8" x 11" + Target



# **CARFUSION DYE**

Catalog # VIP-464 Dye 4 oz. bottle

Catalog # VIP-464C (Case of 6)

Catalog # VIP-463 Dye Holder (anti-spill)

Catalog # VIP-468 Wood Sticks (100/pkg)



# TATTOO/MARKING INK

Catalog # VIP-1PI Black Tattoo Ink ½ oz. bottle

Catalog # VIP-1EB Blue Tattoo Ink (only color available) ½ oz. bottle



# **ACRYLIC STAND**

Convenient accessory used in conjunction with ink cups. Storage for ink cups and shelf to hold cups when in use.

Catalog # VIP-1A



Functions as an ink well, neat, easy, & sanitary. 100/pkg

Catalog # VIP-1C





## ALCOHOL PREP PADS

200 individually wrapped. Size 1 ½" x 1 ½"
70% Isopropyl alcohol for optimum anti-bacterial action.

Catalog # VIP-1APP



Catalog #VIP-168LG 10"H x 9 5/8"L x 5 3/8"W

# Acrylic Accessory Holders


These organizers have separate compartments that conveniently hold several items routinely used daily. Made of clear plastic acrylic.



Catalog #VIP-168SM 81/4"H x 6 1/4"L x 41/4"W

# **TUNGSTEN EYE SHIELDS**

#### 2 and 3 mm Thick



#### Example


The surface dose on a particular accelerator is 72% for a 6 MeV with 80% at the inner surface of the 3mm thick eyelid which is a 20% dose deficit. Using an eye shield consisting of 3mm of tungsten and 1mm of aluminum, there is a 29% increase in dose at the surface of the shields aluminum cap due to backscatter. This gives a total dose of 109% at the inner surface of the eyelid and 101% at the outer surface of a 3mm thick eyelid.



Specifications

Tungsten Density: 17 g/cm<sup>3</sup>
Aluminum Density: 2.718 g/cm<sup>3</sup>

	TRANSMISSION Using XV-2 Film 6 MeV 9 MeV	
2mm Tungsten	3.4% 5.6%	108% 111%
2mm Tungsten + 0.5 Aluminum	3.0% 4.8%	103% 106%
2mm Tungsten + 1mm Aluminum	3.0% 4.4%	95% 102%
3mm Tungsten	2.5% 3.3%	112% 113%
3mm Tungsten + 0.5 Aluminum	2.4% 2.9%	102% 105%
3mm Tungsten + 1mm Aluminum	2.5% 2.8%	97% 106%



	Model #	Description	Thickness	Dimension A	Dimension B	Dimension (	Dimension D	Price
S	936-583	Tungsten Eye Shield	2mm	11.6mm	15.7mm	17.8mm	18.8mm	Quote
	936-585	Tungsten Eye Shield	2mm	13.3mm	17.4mm	19.5mm	20.5mm	Quote
	936-587	Tungsten Eye Shield	2mm	15.0mm	19.1mm	21.1mm	22.1mm	Quote
	936-589	Tungsten Eye Shield	2mm	16.7mm	20.8mm	22.9mm	23.9mm	Quote
	936-591	Tungsten Eye Shield	2mm	18.4mm	22.5mm	24.6mm	25.6mm	Quote
М	936-596	Tungsten Eye Shield	3mm	11.0mm	17.4mm	19.5mm	20.5mm	Quote
	936-598	Tungsten Eye Shield	3mm	12.7mm	19.1mm	21.2mm	22.2mm	Quote
L	936-601	Tungsten Eye Shield	3mm	14.4mm	20.8mm	22.9mm	23.9mm	Quote
	936-623	Tungsten Eye Shield	3mm	16.1mm	22.5mm	24.6mm	25.6mm	Quote
	936-627	Tungsten Eye Shield	3mm	17.8mm	24.2mm	26.3mm	27.3mm	Quote

# **Eyeshield Sterilization Tray**

Consists of Base, Lid, and Mat

Trays can be sterilized by all standard methods of sterilization (autoclave, steam, dry heat and chemical). Trays are molded from General Electric's ULTEM® resin, using mold flow analysis to guarantee product strength, structural integrity and extended life cycle. Flexibility of layout with finger-tip mat inserts to protect instruments. Trays are stackable. Integral feet in bases and corner ridges on lids allow for easier handling in autoclaves and storage. All bases have pebble-type surface to keep mats from sticking. 2.5" x 6.0" x 1.25"

Catalog # VIP-936-500

**Call for Quote** 



# **Prone Head Supports**



Catalog# VC307SF

\*\*More information on page 27\*\*
CALL FOR QUOTE

# Disposable Paper Covers

# **DISPOSABLE PAPER—DUNCAN**

Catalog # VIP-84-B 100 sheets/box

Catalog # VIP-84-C 6 boxes/case 100 sheets/box



For best long-term use of MG head supports and to protect patients we recommend the use of disposable paper covers. This absorbent tissue has been pre-cut to it all headrests, and is available by the box (100 sheets/box) or case (six boxes/case).

Precut paper is available for all Duncan supports.



# **ELEVATED ARM POSITIONER**

Constructed of Lexan material for the base and unbreakable delrin plastic for the hand grips. Two built-in pegs for securing your Silverman plastic head and neck supports (template included for drilling ¼" holes). Timo head and neck supports may be used when taped/glued to the template(s).

Alpha Cradle foam mixture set may be used to construct a custom form positioner for the head and neck and secured by the two base pegs. The handgrips, pegs and 24" x 8" base complete the unit.





# **EXTENDED WING BOARD**


The Extended Wing Board is constructed of durable, lightweight ABS and can be used with a T- or U-Grip handle. Its extended length allows for additional arm positioning flexibility. It is CT compatible and indexable for enhanced reproducibility from imaging through treatment. It easily attaches to the MT-350 and MT-250 breastboards.

**Catalog # MT-WBO8T** Extended Wing Board w/T-Grip - 22.5" wide

Catalog # MT-WBO8U Extended Wing Board w/U-Grip - 22.5" wide

**CALL FOR QUOTE** 





EASY TO READ ACCURATE
MEASUREMENTS FROM THE
TABLE TOP TO PATIENT'S SURFACE







# **BODY CALIPERS**

#### MANUFACTURING DATA:

- Precision machined anodized aluminum.
- Components hand tapped, counter bored and assembled with steel screws for strength, stability and durability. Unlike rivets used in some products steel screws may be re-tightened.
- Easy-to-read screen printed scale is protected from wear by a recessed sliding scale indicator.
- Knurled thumbscrew positioner firmly grips the smooth operating slide at any site on scale.

#### **APPLICATION FEATURES:**

- Can be utilized for table top distances to determine central-axis and off-axis measurement . . . wide base guarantees perpendicularity.
- Reversing caliper orientation assists in determining separations such as laterals, obliques and off-table measurements.
- Scale color coding on both sides of Standard Body Calipers and Extended Body Calipers eliminates confusion with scale orientation.
- Scale resolution defined into one millimeter increments.
- Eliminates dosage variations of up to 4 percent.
- Pointer(s) easily reach around extended body contours.

Catalog#	DESCRIPTION
VIP-95	Standard Body Caliper (45 cm)
VIP-96	Extended Body Caliper (55 cm)
VIP-97	Head & Neck Caliper (25cm)

# **CURVED BLADE BODY CALIPERS**

This lighweight precision Curved Blade Body Caliper has one scale and a non-removable slider which prevents errors caused by reading the wrong scale on a dual scale caliper. The slider can be locked at any position. The rigid leg that slips under the patient is 1/8" thick and 11" long.

Catalog #	Description	Scale
705-546	Curved Blade Body Caliper	45 cm
705-551	Curved Blade Body Caliper	55 cm

# **BODY CALIPER with Flat Blade**

This lighweight Body Caliper with flat blade has one scale on a non-removable slider that prevents errors caused by reading the wrong scale on a dual scale caliper. The slider can be locked at any position. The rigid leg that slips under the patient is 1/8" thick and 11" long.

Catalog #	Description	Scale
705-545	Body Caliper	45 cm
705-550	<b>Body Caliper</b>	55 cm

## **DIGITAL LEVEL**

This digital level is an ideal product to accurately measure the 0.0° to 90.0° angles. This level can assist in set-ups and treatment of patients that require such items as: tilt boards, cassette holders, breast squeeze, treatment tables, gantry angles, and head and neck immobilizers. It is lightweight and only 9" long.

#### Features include:

#### - Angle-Reading Capability

Reads all angles — not just level and plumb.

#### - Recalibration

There are two ways to recalibrate Smart Tool.

#### - Display

Liquid crystal display can be read from the side and up to 10 feet away. Display even flips over with the level so it always reads "right side up."

#### - Accuracy

Smart Tool allows you to choose between two levels of accuracy: 0.1° and 0.2°.

#### Module

The "brain" that powers Smart Tool is a rugged, solid-state sensor sealed inside a weather-resistant polycarbonate module. It has no moving parts or vial to break.

#### - Battery

Smart Tool uses a standard 9-volt battery with a life of 300-400 hours. It will give a "low battery" message on the display when it's time for a replacement, and it shuts itself off after five minutes if not used.

#### - Maintenance

Clean with water by rinsing, spraying or wiping with a wet cloth. Do not clean plastic parts with solvents.



#### Catalog # VIP-192 Smart Tool

One-year limited warranty. See owner's manual for details.

# **POCKET LEVEL**



The pen-size, 12 cm Pocket Level comes in a hexagonal aluminum case. The bottom surface is perfectly flat for accurate surface leveling. A large easy to read glass vial is mounted in an anodized hexagonal tube with a black plastic clip.

Catalog #	Description

352-112

12 cm Pocket Level with Clip

# **MAGNETIC GANTRY LEVEL WITH GLOW-VIEW**



The easy viewing Magnetic Gantry Level is 9" long with three lined vials (plumb, level and 45°) with "Glow-View." It is made of tough shock resistant material. The base utilizes a strong magnetic strip that secures to iron and steel surfaces.

Catalog #	Description
VIP-A67	Magnetic Gantry Level

#### **VAGINAL - RECTAL MARKER**

Vaginal - Rectal - External: Position and magnification can be determined accurately by counting the lead balls which are spaced 1 cm apart. This marker can be placed on top or inside of the patient and used with adjustable external locator during simulation or port filming on treatment machines. Sterilization by gas.

#### **SPECIFICATIONS:**

#### Catalog # VIP-132

- 5/16" diameter x 30 cm flexible tubing 3/16" diameter lead balls spaced every cm
- · Adjustable stainless steel external locator



#### CT RECTAL MARKER

## Rectal • Vaginal • External

The NON-METALLIC CT Rectal Marker is a flexible tube packed with teflon balls spaced at 1 cm intervals from center to center. The Rectal Marker is used to accurately obtain both the rectum position and magnification by counting the balls. This determines the rectum location relative to the radiation field. An adjustable anus marker can be utilized during CT simulation. The Rectal Marker can be used in the vagina or as an external marker.

#### **SPECIFICATIONS:**

#### Catalog # 460-010

Sterilization: Gas; Tubing: 5/16" (8 mm) dia. x 30 cm flexible

Balls: 3/16" (5 mm) dia. spaced every cm; Anus Locator: Adjustable delrin



#### VAGINAL DEPTH SCALE w/INTROITUS MARKER

The Vaginal Depth Scale can be used during simulation or port filming on treatment machines. Position and magnification can be determined accurately by counting the lead balls which are spaced 1 cm apart. An Introitus Marker is included which can be positioned and then locked to the scale with a nylon thumb screw. The Introitus Marker (Model 707-010) is also sold separately.

#### **SPECIFICATIONS:**

Catalog #707- 000

Material: Acrylic Rod; 16 mm dia. x 26 cm; Lead Balls: 3/16" spaced every cm; Scale Length:15 cm



#### CT VAGINAL DEPTH SCALE

This NON-METALLIC CT Vaginal Depth Scale is 1.6 cm in diameter and 26 cm long and has teflon balls which are spaced 1 cm apart and are 5 mm diameter. An introitus marker is included with the 1.6 cm diameter vaginal scale. A thumb screw holds the introitus marker in position.

#### **SPECIFICATIONS:**

#### Catalog # VIP-707-145

Material: Acrylic; Density: 1.1859 g/cm3, Scale Length: 15 cm



#### VAGINAL DEPTH SCALE HOLDER

This holder has a swivel/angling joint that slides up/down on the base plate post. This allows the Vaginal Depth Scale to be adjusted and held at the proper angle. The heavy steel base has a nickel chrome finish. The Vaginal Depth Scale must be bought separately.



Catalog # 707-020

## SHADOWFORM Radiopaque rectal and vaginal markers

ShadowForm rectal and vaginal markers are a sanitary and convenient means of localizing and documenting the exact location of anatomical structures when planning radiation treatment of the pelvic area. They assist the radiation oncologist in developing a specific treatment program that accounts for the location of tissues in the rectum, sigmoid and cervix.

Suitable for use with both port films and the low-energy beams found in simulation rooms. ShadowForm markers are made from a soft, smooth and flexible barium-impregnated polymer that provides the greatest possible patient comfort. They also are individually packaged and disposable—eliminating the potential for cross-contamination.

#### SHADOWFORM RECTAL MARKERS

Available in 10 cm and 35 cm lengths, these tubes are ribbed at 1 cm intervals to provide a convenient reference when compensation for film magnification. A detachable T-bar handle provided as an external anchor, and an open distal end permits the marker to function as pof a barium delivery system.

CATALOG #	DESCRIPTION	UNITS PER CASE
VIP-1710	10 cm Rectal Marker with T-Bar Handle	10
VIP-1740	35 cm Rectal Marker with T-Bar Handle	10

Provided as an 18 cm tube, these markers may be trimmed to any desired length.



## SHADOWFORM VAGINAL MARKER



The detachable T-bar handle acts as an external anchor, and can be used to position a 10 cm, ribbed rectal marker inside of the vaginal marker whenever a reference for compensating film magnification is required.

CATALOG #	DESCRIPTION	<b>UNITS PER CASE</b>
VIP-1735	18 cm Vaginal Marker	10
	with T-Bar Handle	

# **TESTICLE SHIELDS**

These lead shields, available in three sizes, are designed to eliminate or reduce scatter radiation to the testicles. The sections of the shield are secured with two silicone rubber straps supplied with each unit. The halves may be secured with standard tape.

#### SPECIFICATIONS:

VIP-926-100

\* Each unit is cast to provide a lead wall thickness of  $\frac{1}{2}$ " (1.27 cm).

CT Adjustable Testicle Shield Stand

- \* A 90° open sector is provided in each shield to allow narrow but comfortable attachment to the patient.
- \* Lugs cast into each half of the shield unit provide anchor points for the two silicone rubber straps used to hold the unit together.
- \* All surfaces of the shields are painted with a blue-colored epoxy to reduce lead exposure to both the patient and therapist.

CATALOG # VIP-928-200	DESCRIPTION Round Testicle Shield (small)	I.D. 2" (5.08 cm), O.D. 3" (7.62 cm) Weight: 4 lbs.
VIP-928-250	Round Testicle Shield (medium)	I.D. 2 ½" (6.35 cm), O.D. 3 ½" (8.89 cm) 5.7 lbs.
VIP-928-300	Round Testicle Shield (large)	I.D. 3" (7.63 cm), O.D. 4" (10.16 cm) Weight: 7.25 lbs.

Replacement Straps Available - Please call






# **LEAD LETTERS/NUMBERS**

Large ¼" thick x ¾" high lead letters/numbers are mounted on plastic and sealed with epoxy sealant. Each piece has a pre-drilled hole for choice of mounting to shoulder hooks, OR affixed to Velcro® for ease of mounting and remounting to your storage board and/ or film cassette. When ordering, PLEASE specify the shoulder hook style or Velcro style mounting.

These lead pieces are *virtually unbreakable*. However, if one should break, return the item and it will be replaced.

#### SPECIAL ORDERS AVAILABLE:

- \* Choice of letters and/or numbers, i.e., SETUP, CLINAC-20, etc.
- \* Personalized with initials, i.e., L TC, R TC, L MKV, R MKV, etc.



## SIMULATOR - ACCELERATOR ROOM

Catalog #

VIP-79-1	One (1) Letter
VIP-79-2	Two (2) Letters
VIP-79-3	Three (3) Letters
VIP-79-4	Four (4) Letters
VIP-79-5	Complete Set of 1

Set includes the following abbreviated letters: L (left), R (right), AR (anterior right), PR (posterior right), AL (anterior left), PL (posterior left), LAO (left anterior oblique), LPO (left posterior oblique), RAO (right anterior oblique), RPO (right posterior oblique), MED (medial), LAT (lateral), TANG (tangential).

# **LEAD MARKERS**



These easy to use lead markers are composed of rounded 2mm lead, encapsulated in plastic with a soft, self-sticking, peel-off backing. Excellent for identifying, localizing anatomical structures, or used as a magnification factor with the non-stick side next to the skin.

Catalog # LNA-20 100/pkg

# **STORAGE BOARD**



Store your lead letters/numbers on this attractive clear acrylic board that can be wall mounted near your work area. Each board has two strips of Velcro or 16 shoulder hooks with space enough to mount a full set of lead pieces.

Size: 7"h x 16"w x  $\frac{1}{4}$ " with predrilled  $\frac{1}{4}$ " holes for mounting.

Catalog # VIP-79-6V w/velcro Catalog # VIP-79-6H w/hooks

# **CONTOUR WIRE**

- \* Lead-free (97% tin, 3% copper)
- \* 1/16" diameter x 25.5'
- \* Weight: Qne pound

1/16" Catalog #VIP-172-16



# **GONIOMETER**

Angle ruler, plastic • 0-13cm/0-5" • Opaque with black lettering

Catalog #VIP-718



# **FLEXIBLE CURVE**

Used for those irregular shapes needed for contouring patient skin surfaces before transferring to graph paper. Graduations show inches and centimeters.

Catalog #VIP-165-20 20" curve Catalog #VIP-165-32 32" curve



# **RULERS**



#### Catalog #

Stainless Steel: 30cm scale VIP-730-30

60cm scale VIP-731-60

Transparent Plastic: 12" 30cm VIP-405-12

18" 45cm VIP-405-18

VIP-720 12/pkg White Plastic: 20cm

45cm VIP-722 3/pkg

# **POCKET RULER**



Designed especially for radiation therapists. Four centimeter scales on each 3" x 6" mylar card. Will enhance field and ODI lights. Rounded corners on one end to prevent user injury. Opposite end square corner starting at .0mm.

#### Catalog #VIP-121


# **MAGNIFICATION RULER**

#### CT RADIOGRAPHIC MRI PORTAL

Choose from three different long-lasting variations of the laminated magnification rulers. Each set comes with 10 clear mylar rulers for easy reading of radiographs.

- -The R-Mag consists of 15 scales and is graduated in magnification from 1.0 cm to 1.75 cm.
- -The D-Mag consists of 11 scales, graduated in demagnification from 1.0 cm to .50 cm (ideal for demagnified treatment plans).
- -The CT-Mag consists of 22 point scales and has been expanded to 40 cm for use in CT and MRI readings.

Catalog #	Description
VIP-250	R-Mag(Magnifiction)
VIP-251	D-Mag (Demagnification)
VIP-252	CT-Mag (CT & MRI)



## **NEW IMPROVED HEAD REST-For Rounded End Style**

Catalog # VIP-48 set A-F

Separately:

Catalog # VIP-48 (then add letter you need)





#### **SILVERMAN HEAD-NECK SUPPORTS**

Catalog # VIP-50 (set of six A-F)

Separately: add desired letter (A-F) to

VIP-50, i.e.: VIP50A



#### FOR CIVCO/MED-TEC'S ROUNDED END STYLE

Separately add desired letter to MT-VIP-50: i.e., MT-VIP-50-A

Set:MT-VIP-50 (set A-F)



#### **TIMO HEAD-NECK SUPPORTS**

Catalog # VIP-49 (A-F Set)

Colors: Blue, Red, Green, Brown, Yellow, Gray or Assorted (Colors Limited)



### CIVCO/MED-TEC STYLE WITH ACRYLIC PLATES



Add MT before our catalog number, i.e.: Catalog # MT-VIP-49 (set of six A-F)

Ordered separately, i.e.: MT-VIP-49A

Color choice not available

#### SILVERMAN HEAD-NECK SUPPORTS FOR VANARSDALE'S SYSTEMS

Catalog # VIP-50 APHS (Set)



Ordered Separately, i.e: VIP-50-APHS-B

Each support is lettered A-F for easy identification.

These supports are used with our APHS I & II systems. Head rests are punched with two holes on each side to line up pegs on units.

See page 52&53 for our units.

VIP-50 AAP (A-F Set)
Two holes punched on ends for the
Attenuate Arm Positioner



See page 28 for AAP System

### **ARM & SHOULDER EXTENSION STRAPS**

Used for extending the patient's shoulders down and out of the treatment field for head and neck technique. Place straps around patient's wrists and pull each strap length down and under his/her feet. You can control the amount of tension needed for best results. Fabricated with blue velvet nylon strapping and velcro-like fastening system. Measures 50" long with 12" wrist straps. Now availabe in black vinyl for easy cleaning!

Catalog # VIP- 43 (Blue Velour) Catalog # VIP- 43-V (Black Vinyl)



## **SHOULDER RETRACTOR**

Used for extending the patient's shoulders down and out of the treatment field for head and neck technique. For those patients who need the extra sturdiness, this right angle polycarbonate frame is the perfect item. Straps with Velcro' fasteners are placed around patient's wrists, then a nylon rope keeps tension on the shoulders with a self-locking bracket. A tension spring in the rope line keeps the tension on the shoulders without "over strain." Releasing the tension is done by simply pulling ropes laterally away from the patient.



Catalog # VIP-47SV (Spring & Velcro Replacement)



### PRONE FACE POSITIONER II

This set may be used similar to the PRONE PILLO POSITIONER. The complete set includes a total of three (3) positioners. One (1) straight face positioner (9" square x 4.5" high), one (1) 20° cervical/flexion face piece (12" square by 4.25" and 2.75" high) and one (1) chest angle support positioner (8" w x 18" L x 2.5" and 1" H). This system is very helpful for patients that need prone positioning. Excelent for treatment of medulloblastoma, occipital fields, pa spines, pelvic and hips. The face pieces are constructed of semi-rigid foam material. May be used as a disposable/permanent item(s). For face protection, the Duncan disposable sheets or disposable plastic shower caps may be used. Shower caps are used by placing them over the foam piece(s) and making a small split in the center for patient breathing. The angled sponge is covered with microvent material. (See page 14C for information). Can be used with or without the foam pieces depending on the circumstances of patient positioning.

-Easy to position and re-position -Comfortable contour face rest -May be used with or without chest angle support Immediate visual check -Cut-out for "easy-to-see" anatomical structure -Adequate breathing space.

Prone Face Positioner II can be used six (6) ways:

- Straight face piece with or without chest support piece.
- 20° face piece in cervical flexion with or without chest support piece.
- 20° face piece in cervical extension with or without chest support piece.



Prone Face Positioner II - Set



Straight Face Positioning



Flexion Positioning



**Extension Positioning** 

#### **VINYL COATED**

Covered Chest Support

Duncan Disposable Paper

Disposable Shower Caps

Catalog # 306A

Catalog # VIP-84-B

Catalog # VIP-198

Prone Face Positioner II - Set Straight Face Positioner Angle Face Positioner Catalog # VC307 Catalog # VC307SF Catalog # VC307AF

# **ATTENUATE ARMS POSITIONER™**

Patent Pending

For patients required to lay in a supine position with their arms positioned above their head, feeling comfortable is almost an impossibility. That is—until now.

This arm(s) positioner, called the "ATTENUATE ARMS POSITIONER," offers improvements over similar products. It's constructed of high density foam—soft, yet firm. A polyethylene material supports the base and maintains rigidity (with handle). Head and neck supports attach to the base for reproducible positioning.

Patients may be positioned with proper head support, arms extended above their head and placed onto the foam positioner. A slightly elevated angle allows the arms to rest, with an additional angle for the wrist area. Patients have a choice: to hold onto a soft "Grip Ring" or just relax.

This positioner is excellent for esophagus, chest, pancreas, ribs, and breast treatment technique. The patient can be set-up and treated in a fixed position throughout the treatment course. Excellent for oblique and lateral technique. It is also helpful for any treatment technique which requires the patient to have only one arm in the extended position.



Dimensions: 20" wide and 9.5" high
Can easily be placed into the opening of the
Computed Tomography diagnostic unit. Dosimetry
can benefit by confirming their patient treatment plan(s),
and there is no metal component.



Extended position.

Looks Good Feels Good Works Good



Washable vinyl covering.



This accessory provides added support to the humerous (upper limb) portion of the patient's arm(s). May be attached to the Positioner via Velcro fasteners. Available in blue vinyl spray protective coating.



A handle is constructed in the polyethylene base for ease of handling and to store on counter or wall hanging.



A set of six (6) head and neck positions attach easily onto the two (2) pegs.

## VIP-63VC Attenuate Arms Positioner, Blue Vinyl Coated Set

Includes:(One set of Silverman Head Supports (A-F) and one Grip Ring(VIP-182))

VIP-63AC-V AAP Arm Support, Blue Vinyl Coated

Replacements can be ordered separately

VIP50AAP Head & Neck Support (A-F)
VIP-182 Medium Grip Ring

# **Abduce Knee Positioner Set**


# Includes two Foam Positioners and one Grip Ring

The purpose of this product is to arrange patients in an abduction pelvic position with maximum comfort.

This blue vinyl coated semi-firm sponge can be very useful for a treatment setting of the entire pelvis, single hip, and vulva/scrotum area.

Place the sponge(s) under the femur/knee area with the lower side of the sponge, a 10° slant toward the patient's trunk. The knee(s) are seated into the sponge's concave area.

The patient's feet can be held together with the support of the vinyl coated grip ring.



Catalog # VIP-300 Abduce Knee Positioner
Includes 2 Foam Positioners and 1 Lg Grip Ring (VIP-183)

Catalog # VIP-183 Large Grip Ring

# PRON-PILLO®



A molded soft foam positioning device designed to provide maximum comfort for the patient when lying in the prone position. The unique feature of Pron-Pillo is its ability to support both head and shoulders to assure a comfortable position for the neck.

Eliminates use of conventional pillows and additional time involved in placing your patient in a comfortable prone position. Raised base allows patient to breathe freely while enjoying a relaxed prone position.

**PRON-PILLO** 

**CATALOG #FFF306** 





A Belly Board is used to allow the intestinal tract to drop out of the treatment field. Treatment through the foam material has a low absorption of the beam. This unit is made of lightweight, no-breakable high density foam and is coated with a durable blue finish for easy cleaning. The light weight allows for easy handling of the Belly Board. The soft flexible surface provides patient comfort and helps retain excellent body positioning. The bottom surface is anti-skid on treatment and simulator couch-es. The Belly Board includes one opening reducer that when placed in the opening makes the size of the opening smaller. When the EXACT ENHANCE BAR is attached to the Belly Board with THE GLUE DOTS, the belly board can then be indexed using a Lock-Bar or a similar bar. This will allow positioning the belly board in the same position daily.

Cleaning Instructions: Soap and water, 50/50 Bleach/Water or an ammonia solution can be used.

Additional Specification: Opening: 12 1/4" W x 13 3/8"L (31 x 34cm)

Blue Reducer: 11.8"W x 3.23"L x 7"T (30 x 8.2 x 17.8cm)

Catalog # VIP-256-035 Belly Boards

## GRIP RINGS™ — A Comfortable Idea







Now you can give your patients that "extra something" to help them feel more comfortable. By using Grip Rings™ you offer patients something to make them feel more secure in an already stressful situation. Another benefit of using the Grips is that patients are more relaxed, and original set-up and treatment position techniques are easier to verify.

#### **GUARANTEED FOR 90 DAYS**

#### How is it Designed and Constructed?

The circular design of Grip Rings allows pressure to be evenly distributed, preventing harsh marks on the skin. The attractive and smooth outer skin becomes an integral part of the cushion through a high strength bond to the foam core. Only high tech, non-absorbing, closed cell cushing mediums are encapsulated to ensure the highest measure of quality. The outer skin of tough blue vinyl is smooth and easy to wipe clean with most anti-bacterial solutions, and free of surface textures that trap and hold bacteria from soil and body fluids. If the soft, supple outer skin is cut or torn, the inner cushion will not "sponge-up" liquids/fluids as some products that utilize conventional upholstery materials. Because of the smooth vinyl surface, the ring(s) can be placed over the patient's wrist (and removed) easily and quickly without causing additional trauma to the patient (such as cuts, burns, abrasions, or hair removal as when tape is used).



Five Sizes Offered





#### **SETS AVAILABLE:**

**Catalog # VIP-166** Grip Rings™(3) Contains 3 rings: Small, Medium, and Large

**Catalog # VIP-166-4** Grip Rings<sup>™</sup>(4) Contains 4 rings: Small, Medium, Large, and X-Large

#### SOLD SEPARATELY:

Catalog # VIP-181 SMALL: children & small adult hands (I.D. 3")

Catalog # VIP-182 MEDIUM: average size adult hands (I.D. 4/1/4")

Catalog # VIP-183 LARGE: large adult hands and small feet (I.D. 6")

Catalog # VIP-083 X-LARGE: large feet or with shoes (I.D. 8")

XX-LARGE: extra large feet or with boots (I.D. 10")

## **POSITIONERS**

## Covered Positioners

All covered foam blocks are enclosed with non-seam velcro strips. Blue Microvent Soft™ Pillow Ticking Fabric covers these polyfoam positioning blocks. Wipe them clean, no laundering required.

PROTECTION • ECONOMY PERFORMANCE • COMFORT


## "Microvent Soft" Pillow Ticking Fabric . . .

- Antibacterial
- Antifungal
- Noiseless
- · Non-allergenic
- Stain resistant
- · Flame retardent
- · Self-deodorizing
- Easy-to-clean
- · Long lasting
- Unaffected by discharges from eyes, ears, nose and mouth
- Helps reduce risks of cross infection
- Lowers housekeeping, laundering, replacement, inventory and disposal costs.
- Unaffected by typical hospital disinfectants when used in those dilutions recommended by the manufacturer.

Custom positioners made to your specifications ... most any size, any shape, and are covered in black vinyl material.

Call us for more information.


#### Standard Set



Standard Set Covered Catalog # FFF-702-C Set consists of: FFF101-C, FFF102-C, FFF103-C, FFF104-C, FFF105-C, FFF106-C

#### **Individual Pieces**

35°/55° Angle 9¾" x 24" x 6¾"




20° Wedge 11" x 7" x 4" 4" 11" 20°


Covered Catalog # FFF-101-C

Covered Catalog # FFF-104-C

45° Triangle 7" x 24" x 7"



Rectangle 11" x 7" x 2"



Covered Catalog # FFF-102-C

Covered Catalog # FFF-105-C

15° Wedge
11" x 7" x 3"

3"

11"


7"

Covered Catalog # FFF-103-C


Circle 6" x 11/2"

Covered Catalog # FFF-106-C


## POSITIONERS continued . . .




Angular Bolsters FF301-C covered 21" FFF301-CM 24"




Angular Bolsters
Catalog #FFF302-C covered



Bolster, Knee and Back
Catalog # FFF-303-C covered 15"
FFF-303-CM 21"




Bolster, Myelogram
Catalog # FFF309-C covered 18"
FFF309-CM 24"




Leg Immobilizer - Covered In Black Vinyl Catalog # FFF-704-C

Catalog # FFF-704-CE 14"W x 16"L x 6" x 11" H






45° Triangle 9" x 9" x 24" x 12" **Catalog#FFF-107-C covered** 11.5 x 11.5 x 21 x 16 **FFF-107-CM** 



Torso Blocks
Catalog # FFF401-C covered

In addition to the positioners presented in our catalog, we would be glad to try to fill any custom order you may have.


If you need a positioner size not listed, please give us a call so we can discuss meeting your needs.




gles, covered
Description
1" H x 10" L x 7" W 2" H x 10" L x 7" W 3" H x 10" L x 7" W 4" H x 10" L x 7" W 4" H x 20" L x 10" W

## POSITIONERS continued . . .


#### Extra Long Wedges



Catalog # FFF405-C covered




Catalog # FFF407-C covered




Catalog # FFF124-C covered

#### Wedges




Catalog # FFF203-C covered




Catalog # FFF205-C covered 10" x 10" x 10" Catalog # FFF205-CM 7" x 7" x 7"

#### **Chest Positioners**




Chest Flat

Catalog # FFF402-C covered



Chest Positioner
Catalog # FFF403-C covered

#### Chest/Leg Elevation Wedges



15" x 24" x 24"

Catalog # FFF502-C covered

In addition to the positioners presented in our catalog, we would be gl ad to try to fill any custom order you may have.


If you need a positioner size not listed, please give us a call so we can discuss meeting your needs.

## Exact Enhance Bar

assists in the setup of

## VAC FIX, VAC-LOK<sup>™</sup> ALPHA CRADLE<sup>®</sup> SUPPORTS

This product was designed to "ENHANCE" the EXACT BAR used with the Varian Accelerator Exact Treatment Couch. The EXACT BAR can be adapted to the vacuum bags and Alpha Cradle type of patient immobilization forms with the use of three self-adhesive Glue Dots.



Glue Dots Exact Enhance Bar Varian Exact Lok-Bar



May be used with any plastic Vac Bag system again and again because of easy removal.


#### Exact Enhance Bar Catalog # EXC-EB-1

**Varian Exact** 

Catalog # MT-VAR-20

The Enhance Bar is made of high-density polyethylene material. Measurements are 1"W x 16"L and  $\frac{1}{2}$ " in thickness. There are holes drilled. You have the choice of centering or shifting  $2\frac{3}{4}$ " in either direction, laterally.

- · Consistent in accurate positioning
- · Easier location of "Exact Bar" and vacuum bags/Alpha Cradle
- · Replacement guarantee against chipping, breaking or warping



Seat your Varian Exact Lock-Bar or similar bar in a "notch" along the length of the couch/table top.



Place the Exact Enhance Bar onto the Loc-Bar, matching the "two pegs" or similar design unit.



Now use "Glue Dots" to adhere the Exact Enhance Bar to your immobilization form. Take two or three "Glue Dots" and press the open side into the Exact Enhance Bar (as shown). Remove the wax paper and repeat for each Glue Dot. (Do not touch with fingers.)



Place the patient immobilization form on top of the Bar/Glue Dots.



Position patient onto immobilization vacuum/mold to construct form.



Remove patient from completed form. The Exact Enhance Bar will be attached to the form and can be easily removed and reseated onto the Lok-Bar/couch when needed.

The Exact Enhance Bar can be used repeatedly.

The Glue Dots are easily replaced by simply "rolling" them off the bar or form and discarding them.

#### Catalog # EXC-EB-1 Exact Enhance Bar



Easy to use, simply touch item to glue dot and lift. The Glue Dot will transfer to the item. Do not remove the Glue Dot with fingers.

Glue Dots are 3/4" round. Dots cannot be reused. Samples are available.

Catalog#XD43-601 Glue Dots 200 count/10sheets

Page 36

1-800-326-4761 • 24 Hour FAX 850-476-8426

## **WARNING SIGNS & TAPE**

#### SILK SCREENED SIGNS

Signs in this section are silk screen quality, 1/8" plexiglas (acrylic) backboard with double-faced tape attached for mounting.

For a combined sign Catalog # PS04ES (Call for Quote)

#### **PREGNANT**

IF YOU ARE PREGNANT, OR THINK YOU MAY BE, TELL THE RADIATION THERAPIST BEFORE HAVING RADIOGRAPHS TAKEN.

#### **EMBARAZO?**

SI USTED ESTA EMBARAZADA O CREE QUE LO ESTA NOTIFIQUE-SELO AL TECNOLOGO ANTES DE QUE LE TOMEN LOS RAYOS-X. 10" x 4"


8" x 10"

#### Catalog # PS-04E

Size: 10" x 8"

#### Catalog # PS-04S

Size: 10" x 8"





VIP-657P



VIP-7710P



VIP-204P



VIP-445P

#### **Peel & Stick Caution Labels**

Label size: 21/4" x 21/4" 10 per sheet





Catalog # VIP-CL215



VIP-202P



CRA-17 10"x4"

# CAUTION THIS MACHINE PRODUCES RADIATION WHEN ENERGIZED

CRE-18 10"x4"



VIP-203P



VIP-205P



CBH-14 10"x4"

## **INCONTINENCE CLAMPS/MEATUS CLAMPS**



Whitver Clamp Catalog # VIP-U410



Black Clamp Catalog # VIP-U411



Hyams Clamp Catalog # VIP-U412



Zipser Clamp Catalog # VIP-U413





Cunningham Incontinence Clamp Catalog # VIP-U414R (Regular) Catalog # VIP-U414L (Large)





Baumrucker Incontinence Clamp Catalog # VIP-U415


(No metal parts)

#### PATIENT CHIN STRAPS

#### Catalog # VIP-85


Help immobilize your patients when head and neck areas are under treatment with this easy-to-use chin strap. a contoured chin cup that fits comfortably over the patient's chin; while the velcro hook and loop fasteners allow you to pull desired tension needed for patient position-alignment.

An easy-to-use Velcro® restraining system provide secure immobilization time after time without gumming up the treatment couch. Our padded straps are comfortable, washable and reusable.



ADULT CHIN STRAP

Catalog # CS-A



PEDIATRIC CHIN STRAP

Catalog # CS-P



## LATEX FREE TUBULAR ELASTIC NET

For the effective fixation of dressing without restricting circulation, irritating the patient's skin, or binding body movement.

- -Tubular material(s): Latex free elastic tubular net
- 25 yards per roll (22.9 meters)

A very suitable undercover bandage when a protective barrier is needed for covering head and neck area. Simply cut to size indicated; pull over area to be covered. To cut holes for appendages, pull on bandage, cut strands where appendages will pass through.



Catalog #:	<u>Size</u>	<u>Width</u>	<u>Application</u>
VIP-91-1	1	5/8 inch	Adult fingers and toes
VIP-91-2	2	3/4 inch	Infant feet or arms
VIP-91-3	3	7/8 inch	Small hands and arms
VIP-91-4	4	1 inch	Medium arms, small ankles
VIP-91-5	5	1 1/4 inch	Large arms, small knees, medium knees
VIP-91-6	6	1 5/8 inch	Large ankles, medium knees, small thighs
VIP-91-7	7	2 inch	Large knees, medium thighs
VIP-91-8	8	2 1/4 inch	Large thighs
VIP-91-9	9	2 1/2 inch	Small trunks
VIP-91-10	10	2 7/8 inch	Medium trunks
VIP-91-11	11	3 inch	Large trunks

#### **ELASTO-GEL**

#### ... The stick-to Gel on purpose ...

Materials are frequently used in high energy radiation therapy in order to deliver the prescribed dose to the patient's skin surface. Elasto-gel is easy to work with, reusable on the same patient, is mildly adhesive, and conforms exceptionally well to body contours. It may be easily cut with scissors to the desired size and shape. Standard backing is stretch cloth on one side and clear removable plastic on the other. The stretch backing may be removed by first moistening with a damp cloth and then peeling the cloth from the gel. Because of the self-adhesive nature of the gel pieces, they may be layered together for various thicknesses. Density 1.20 gm. Sample upon request.

Catalog #	Size (inches)	Size (cm)	Thickness
BL9705	12" x 12"	30 x 30	0.5cm
BL9710	12" x 12"	30 x 30	1.0cm
BL9803	8" x 16"	20 x 40	0.3cm
BL9805	8" x 16"	20 x 40	0.5cm
BL9810	8" x 16"	20 x 40	1.0cm



#### SUPERFLAB BOLUS MATERIAL

- For maximum dose built-up to skin
- Conforms well to patient's contour
- Very elastic and quite "FLABBY"
- Specific gravity: 1.02
- Maintains uniformity to thickness
- · Can be washed
- · Can be cut
- Approved for human contact
- · Made of synthetic oil gel
- · Will not dry out



The SUPERFLAB BOLUS MATERIAL is offered in thicknesses which provide maximum dose build-up for relevant photon energies. Since the material does not suffer inelastic strain from normal stresses, it does not have to be bagged or wrapped in plastic film to maintain its shape. At the option of the user, however, SUPERFLAB may be wrapped in disposable plastic film for cleanliness and use; or it may be washed with soap and water as needed, followed by an application of talcum powder or corn starch.

Thickness (cm)	Standard Sizes <u>Size (cm)</u>	Catalog#
0.2	30 x 30	8117-0.2
0.3	30 x 30	8117-0.3
0.5	30 x 30	8117-0.5
1.0	30 x 30	8117-1.0
1.5	30 x 30	8117-1.5
2.0	30 x 30	8117-2.0
2.5	30 x 30	8117-2.5
3.0	30 x 30	8117-3.0
4.0	30 x 30	8117-4.0

Thickness (cm)	Large Sizes Size (cm)	Catalog#
0.3	40 x 40	8117-03-4040
0.5	40 x 40	8117-05-4040
0.5	40 x 60	8117-05-4060
1.0	40 x 40	8117-10-4040
1.0	40 x 60	8117-10-4060

#### SUPER STUFF BOLUS MATERIALS

Super Stuff is one of the best materials to be used for bolus to date. Super Stuff comes in pink, blue, green and white powder form. One ounce of powder is packaged in a plastic bag with a black line on the bag to indicate the water fill line. (One cup of water is added to one ounce of powder). After mixing, by squeezing on the plastic bag, Super Stuff takes on the consistency of Jello. It can then be shaped to fit the contours of the patient. Long shelf life. Density: 1.02g/cm<sup>3</sup>

Catalog # VIP-169 PK (pink)	50/1 oz. pkgs
Catalog # VIP-169 BL (blue)	50/1 oz. pkgs
Catalog # VIP-169 GR (green)	50/1 oz. pkgs
Catalog # VIP-169 WH (white)	50/1 oz. pkgs



#### **BOLX**

with sealed skin/film to prevent sticking — easy to clean

A superior contour bolus is now available. Uniquely self adhering and tissue equivalent. Available in 0.5cm-4.0cm thickness. Effective density relative to water is 1.02g/cm<sup>3</sup>. Bolx without skin is also available at the same price—use same catalog number and add w/o.



Catalog #	<b>Thickness</b>	<u>Size</u>
VIP-58-131	0.5cm	30cm x 30cm
VIP-58-132	1.0cm	30cm x 30cm
VIP-58-133	1.5cm	30cm x 30cm
VIP-58-134	2.0cm	30cm x 30cm
VIP-58-135	3.0cm	30cm x 30cm

Catalog #	<b>Thickness</b>	<u>Size</u>
VIP-58-131S	0.5cm	15cm x 15cm
VIP-58-132S	1.0cm	15cm x 15cm
VIP-58-133S	1.5cm	15cm x 15cm
VIP-58-134S	2.0cm	15cm x 15cm

#### **ALGINATE IMPRESSION MATERIAL**

A <u>fast setting</u> alginate impression material that is used to produce facial molds, dental and gum impressions. This fast setting material uses water as a base and when applied to a facial area or bite block, will produce an impression within 1–2 minutes.

Catalog #0485-001 1 lb. can

Catalog #VIP-79420 Spatula, plaster #11R
Catalog #VIP-H-00565 Flexibole, large 60<sup>cc</sup>





#### **SILKY ROCK**



Catalog # 6834-01996

Yellow Microstone 33 lbs.

This material is used to produce solid patient facial molds post alginate impression material. Silky Rock powder, when mixed with regular tap water, forms a "rock-like" solid mold.

After hardening of the mold, thin sheets of lead can be placed onto the section/ treatment interest with a small hammer, and contoured into the shape needed.

Small rounded end instruments(s) can shape the smaller contoured areas such as the nose/eye. Small wooden dowels are sometimes used for this purpose by simply sanding the ends to a convex shape.

#### "POWDER FREE" LATEX GLOVES

Non-chlorinated Powder Free and Hypoallergenic

Catalog # VIP-70100 Extra Small Catalog # VIP-74100 Medium

Catalog # VIP-72100 Small Catalog # VIP-76100 Large 100/box

VanArsdale Innovative Products, Inc.

1-800-326-4761

24 Hour FAX 850-476-8426

#### **AQUAPLAST PELLETS**



Soften these pellets quickly in hot water for easy shaping! Fabricate custom bite blocks and bolus with these easy to mold pellets. Simply soften in hot water for less than 1 minute and contour—either directly against the skin or a plaster or dental stone mold. Produces a rigid, unbreakable mold, which can be repeatedly modified by reheating.

Catalog #RT-7001 1 lb. Catalog #RT-7003 3 lbs.

When forming bite blocks, place the softened pellets in a small plastic bag prior to inserting in the mouth. (DO NOT PLACE SOFTENED PELLETS DIRECTLY INTO THE MOUTH). Once the mold is completed and cooled, remove the plastic bag.

#### **UTILITY WAX**



This paraffin wax material is soft and pliable and can be molded into the smallest detailed configuration. It comes packaged in a box of 55 square rope strips (11" x 3/16" x 3/16") and is ready to use. Starting with one rope strip, and by adding rope to rope, you can construct the shape you need without having to mix or melt a material/chemical. The wax will maintain its shape as long as it isn't pressured out of shape.

Catalog # VIP-171 55 ropes/box white Catalog # VIP-171R 55 ropes/box red



#### PINK DENTAL WAX

This wax may be used as bolus material, or used for coating eye shields. 1.6 mm thick.

Catalog #VIP-158 1 pound box - Approx. 25 sheets per box
Catalog #VIP-158-5 5 pound box - Approx. 125 sheets per box



#### **ALUWAX - IMPRESSIONAL WAX**

MAKE A LASTING IMPRESSION . . . A simple, yet ef-fective method to use with your BITE BLOCK HEAD AND NECK IMMOBILIZATION SYSTEM. This dental wax (gray in color) is easy to form by heating in a water bath (or your room faucet water) 90° to 155°F. The wax may be used in sheet strip form or rolled into a mass. After softened, place around the bite piece and into the patient's mouth. Setting time, 1–2 minutes. Sets hard—guaranteed you'll like it!

Catalog # VIP-163 25 sheets/box








#### **BITE PIECES**

An accessory for head and neck immobilization systems. We recommend using Aluwax as your molding compound. We can supply bite pieces for the Hughes-Davidson, Radiation Products Design, and others.



A

Catalog #VIP-274-100 (12pk) - Adult Disposable Jaw Lock

B

Catalog #VIP-164Hd (12pk) - Plexiglass Bite Pieces

C

Catalog #VIP-164RPD(12pk)- Plexiglass Bite Pieces

D

Catalog #VIP-164QG (12pk) - Quick Grip

#### PLEXIGLAS BITE PIECES/TONGUE DEPRESSOR

This item can be used with or without a bite block head and neck immobilization system. The length of this bite block is 4cm longer than the average 9cm, making the total length 12.8cm.


Aluwax or other molding material can be molded around the bite block and positioned in the usual manner. As the patient closes his/her mouth, the extended piece will force the tongue to stay in a lower position. Outline (parallel) of the treatment field can easily be identified by the bite block stem end, reassuring tongue depression.

Features: A 4.5cm length with a .7cm i.d. plastic tube permitting free airway of patient breathing.


Another use: By wrapping gauze or other material around the bite block, it can be used as tongue depressor or an airway in an emergency situation.

Can be cleaned and reused if gas sterilized.

VIP-52 (minimum order 6)







Using regular bite block

Using extended bite block

### VELCRO<sup>®</sup> ONE-WRAP™

The Ultimate in Patient Comfort Restraining

- Velcro One-Wrap is extremely soft with no sharp edges.
- Strong as tape without the adhesive tape pain to the patient.
- Easy and quick to use for your strapping and positioning needs.
- Numerous positioning uses Economical Reusable

A Velcro strap that has the *loop side* and *hook side* back-to-back. You'll never have to look for "matching pieces" of Velcro again; just cut the length you need. Black and 1½"w. Minimum 5 yards.

Catalog #VIP-180



#### DISPOSABLE POLYETHYLENE OVERBOOTS



Slip-on, disposable, 0.003"-gauge polyethylene overboots fit over any shoe size. Easy-tie closure—just rip along the perforations and tie around your leg. Overboots come in convenient dispenser packages of 50 (25 pair). Durable thermoplastic dispenser is clear and features four keyholes so it can be mounted on the wall.

Measures  $10\frac{4}{4}$  H x  $12\frac{1}{5}$  W x  $2\frac{1}{2}$  D.


Catalog #	Description
8BB-150866	Overboots, XL
8BB-150865	Overboots, XXL

#### TABLE RAIL SAFETY STRAP

VIP's safety straps help to immobilize the patient during radiation treatment and prevent falling from the table. The straps are constructed of soft, yet strong blue velvet cloth material which is completely washable.

To prevent the straps from falling completely to the floor while not in use, fasten each strap to itself as illustrated. This enables you to leave the straps attached to the side rails at all times.

To secure the patient in position, simply pull the strap with Velcro loop side facing up first and then the hook side strap second. Easy, fast and comfortable for all.



The straps come in a two piece set. Each piece is 6" wide and 30" long. Slip the proper end of each piece through each side rail and fasten with the Velcro closures. Position each strap as illustrated.

#### Catalog # VIP-170

#### TILT TABLE STRAP



Specially designed tilt table strap is made of Velcro\* brand loop material, double-faced and cushioned for patient comfort, with Velcro\* brand hook material at each end of strap for greater control and adjustability. Can be washed and sterilized.

- Catalog #VIP-281045
- Color: Beige
- Size: 4" x 72"

#### **IMPROVED TABLE STRAP**



Extra strong table strap with more Velcro® to accommodate very small and extra large patients. Made of soft blue velour with 18″ hook and 24″ loop Velcro® restraint sewn on. Can be washed and sterilized. Also available in black vinyl for easy cleaning!


- Catalog # VIP-156
- Color: Blue
- Size: 6" x 82"
- Material: Blue Velour
- Catalog # VIP-156V
- Color: Black
- Size: 6" x82'
- Material: Black Vinyl

#### **TABLE STRAP**

Strong, reliable, cushioned, and soft against the skin. For the patient, this strap is designed to ensure patient safety and comfort. For the therapist, it provides the greatest degree of adjustability and control. Made of soft blue velour with Velcro\* brand loop. Can be washed and sterilized. Also available in black vinyl for easy cleaning!

- Catalog # VIP-173
- Color: Blue
- Size: 4" x 72"
- Material : Blue Velour
- Catalog # VIP-173V
- Color: Black
- Size: 4" x 76"
- Material: Black Vinyl





#### PATIENT ARM SHIELDS

Patients may feel more comfortable by placing this acrylic plastic shield by their side while being transported from place to place.

- Prevents personnel from leaning on patient's arm or I.V. lines and tubes.
- Slides beneath table pad.
- Lightweight clear plastic that is RADIOLUCENT.
- Easily cleaned with hospital disinfectants.
- Supplied with carrying/hanging ports.

Size: 18" long x 6.5 high x 9" deep x 1/4" thick.

Catalog # VIP-134 (Set of Two)



## **Heating Pans**



Heating Pan, electric, temperature controlled.
I.D. dimensions 11.5"w x 15.5"L x 2.5"D
\*Specific brand will depend on availability
Catalog # VIP-EHP



Heating Pan, electric, temperature controlled, W-3V, 5 gallon capacity. 700 watts.
I.D. dimensions 12"w x 20"L x 6 3/4"D
O.D. dimensions 14"w x 22"L x 9 ½"D

#### Catalog # VIP-HWS-1



#### **Water Bath**

- Overall: 31.25 x 24 x 13" (79.37 x 60.96 x 33.02 cm)
- Internal: 25 x 19 x 5" (63.5 x 48.3 x 12.7 cm)
- 120 Volts
- Digital Readout Level

The Water Bath accomodates all thermoplastic sizes including Type-S head, neck and shoulders, HipFix, and ThermFix thermoplastics.

#### Catalog # MTRH-1500

\*Dimensions do not include shut off valve mounted on side.

## **Mold Room Accessories**

#### **Unlimited Scissors**



Cuts through a completed immobilizer and unsoftened Aquaplast™ with powerful leverage.
Flexible, non-slip handles for a comfortable grip.
Unlimited Scissors (common model for right and left hands)
Catalog # RT4932



Complete system, 650°F heating element, adjustable stand & heat shrink attachment.

Catalog # VIP-341411 Heat Gun
Catalog # VIP-341412 Replacement heat element

#### Waterbath Rust Prevention C.R.E.S



Many customers believe stainlesss steel to be impervious to rusting, but that is not the case. The unit can rust from being exposed to harsh cleaning chemicals, foreign debris, or lack of proper and regular cleaning.

C.R.E.S. is recommended by the manufacturer to remove rust from the unit.

#### Catalog # VIP-CRES-K


Kit includes sponge, scratch pad, gloves, and 16 oz. cleaner

## Aquaplast RT™ Miscellaneous Items

Catalog #	Description
RT-4454	12" Opti-Handles, pair, for use with breast boards
RT-4454-14	14" Opti-Handles, pair, for use with AccuFix Quest breast board
RT-4932	Unlimited Scissors, universal left or right handed
RT-7001	Adapt It pellets, 1 lb. resealable package
RT-7003	Adapt It pellets, 3 lb. jar/container
RT-2068	Stainless Steel Aquapan for S-Frames (7.5 gallon capacity 3" deep
	water level) with hospital grade plug, 110 volts
	(22" x 26" x 5" ID - inside diameter)







## Aquaplast RT™ Bolus Material

Catalog #	Description
RT-1908-4	Ivory, 4.8mm Aquaplast RT Bolus Material, 7.6cm x 22cm (3" x 9"), solid
RT-1910-4	Ivory, 4.8mm Aquaplast RT Bolus Material, 20cm x 22cm (8" x 9"), solid
RT-1913-2	Ivory, 2.4mm Aquaplast RT Bolus Material, 30cm x 30cm (12" x 12"), solid
RT-1913-4	Ivory, 4.8mm Aquaplast RT Bolus Material, 30cm x 30cm (12" x 12"), solid
RT-1922	Ivory, 4.8mm Aquaplast RT Bolus Material, 45cm x 60cm (18" x 24"), solid
RT-1930-4	Ivory, 4.8mm Aquaplast RT Bolus Material, 43cm x 43cm (17" x 17"), solid
RT-1931-4	Ivory, 4.8mm Aquaplast RT Bolus Material, 30cm x 45cm (12" x 18"), solid

## Aquaplast RT™ Thermoplastic - 2.4mm Thick

RT-1779O Original Ivory, disposable U-Frame, perforate

RT-1779SO Original Ivory, Slimline, disposable U-Frame, perforated

RT-1782SD Ivory, disposable head & shoulders "S" Frame, perforated

RT-1783P3 Ivory, 12" x 18" (holes punched in 12" side for Quest Board), perforated

RT-1788 Ivory, 9" x 12" sheet, perforated

RT-1788P Ivory, 9" x 12" sheet (prepunched for U-Frame), perforated



## Aquaplast RT™ Thermoplastic - 3.2mm Thick

RT-1822 Ivory, 18" x 24" sheet, solid RT-1824 Ivory, 24" x 36" sheet, solid RT-1882 Ivory, less stretch, 18" x 24" sheet, perforated RT-1889O Original Ivory, disposable U-Frame, perforated RT-1889SO Original Ivory, Slimline disposable U-Frame (wider opening), perforated RT-1889SXO Original Ivory, Slimline disposable U-Frame extended 2", perforated RT-1889SD Ivory, disposable head only S-Frame, perforated RT-1889SLO Ivory, Slimline U-Frame 9" x 14", long RT-1891P Ivory, 9" x 10" sheet (prepunched for U-Frame), perforate RT-1892SD Ivory, disposable head & shoulders S-Frame, perforated



Disposable Head Frame



Slimline Frame



Extended Frame

## **U-Frame Wall Hanger**



A product that is made to help keep your patient molded U-Frames off of counter tops, etc. and mounted on your treatment room wall in an organized manner.

The unit has five  $1\frac{1}{2}$ " pegs for U-Frame hangings, spaced 11" apart. Three drilled wall mounting holes. Overall dimensions: 3" H x 48" W x  $^3$ /8". Clear acrylic plastic.

Catalog # VIP-94 U-Frame Wall Hanger

## **U-Frame Immobilization System**

#### Product Features:


Minimal time needed for immobilization patient set-up technique, and can be utilized by one person n Reproduction patient set-up positioning for day-to-day treatments n Patient multiple positioning, 0° (flat) to 38.5° (flexion) using the 0 ring on the angle rod Two (2) openings, 5.25" x 9". One in the base and one in the angle to accommodate through the table technique n Base plate has non-skid tabs n U-Frames available as resusable or disposable n Uses standard Silverman plastic head and neck supports

Can accommodate the use of Timo head and neck supports by simply glueing to additional templates

Entire unit can be utilized in CT and MRI diagnostic units for verification of treatment set-ups due to all acrylic/plastic material

This is a dedicated system for patients undergoing external radiation treatments of the head and neck region. The system uses a perforated thermoplastic molding material. It's safe and easy to use when softened in temperature controlled water, and can be silicone sealed using polycarbonate plastic reusable frames. With no preparation to be done to the patient, the U-Frame can easily be draped over the patient's face and secured by four (4) quick release clamps. Alignment pegs are seated for precision repositioning, and the thermoplastic molding material can be shaped to conform to treatment adaptation and patient comfort.

Added feature: The APHS-1 and the APHS-2 have a milled-out area in the superior portion of each unit. This milled-out area can be seated onto the Exact Bar of the Varian Accelerator Exact Couch



APHS- 2 Unit

If you prefer to use disposable U-Frames, you'll like the VIP 1764s. Although the size states 8" x 9", it refers to only the back side of the frame. The patient coverage remains the same. The look is different because of the rounded upper corners, made special for VIP. We have drilled two additional peg seating holes 2.7 cms above the regular peg seating hole. This "extra" hole will allow positioning the mold close to the apex of the head plus additional coverage below the chin. See the photo below.

#### Package #2 (disposable)

System includes: one base unit, one template (for head rest alignment and thermoplastic shrinkage if needed), one angle rod with "0" ring for flexion angle setting, three disposable U-Frames with perforated thermoplastic sheets.




#### Package #1 (reusable)


System includes: one base unit, one template (for head rest alignment and thermoplastic shrinkage if needed), one angle rod with "0" ring for flexion angle setting, two reusable U-Frames with thermoplastic sheets, one scraper and one tube of silicone sealant.

Outside dimensions: 12" x 18" x <sup>3</sup>/<sub>4</sub>".

#### **Prone Positioning Accessory**



Now available:
A soft foam vinyl covered prone positioner made to assist patients in a prone head position.
Pad size: 8 ½" x 6 ½" x ½"



The pad is attached with Velcro to a 3/s" thick polycarbonate plastic frame with a cutout for the patient's face and two 1/4" holes for securing onto the base plate. Safe and easily cleaned. Constructed for the APHS II Immobilization Unit, but may possibly be modified for other units. (CIVCO/MEDTEC, Inc.)


#### APHS-1 Immobilization System • Same as APHS-2, without the angle

Catalog #	<u>Description</u>
APHS-2	U-Frame Immobilization System with angling feature
APHS-1	U-Frame Immobilization System
APHS-3	Prone Positioner Pad
APF-1	Reusable U-Frame, Lexan
VIP-1764-D2	Disposable U-Frame (TP) 2.4mm
VIP-1764-D3	Disposable U-Frame (TP) 3.2mm
VIP-200	Silicone Sealant—3 oz.
AP-711	Templates 7" x 11" x 1/8"


Vac-Lok Bar for Exact Couch

Head & Neck supports—see page 28 Heating Pan—see page 52

MT-VAR-20



APHS-1





APHS-1 with Exact Vac-Loc Bar

	Catalog #	<u>Description</u>	
A	VIP-1764 D2 VIP-1764 D3	Ivory, 2.4mm, 9" x 10" sheet Ivory, 3.2mm, 9" x 10" sheet	
В	VIP-1764-D2-S VIP-1764-D3-S	Ivory, Slimline frame 2.4mm, 9" x 10" sheet Ivory, Slimline frame 3.2mm, 9" x 10" sheet	
C	VIP-1764-D3X APHS-RS	Ivory, Extended Sheet 3.2mm, 9" x 12" Ivory 1" Replacement Screws w/Black Knob (5/package)	





### **Soft Sponge and Molestick**

Now your patients' rough spots can be soft spots when they have thermoplastic material (TM) formed over their body parts to produce immobilization. A patient immobilized under any TM can become very uncomfortable. This can be due to swelling, pressure points, sores or lacerations.

By placing a small piece of an adhesive sponge padding to your cured TM over these areas, a much softer point of pressure is created, giving you a much happier patient.


**Molestick**, an open-cell foam with an adhesive back, is available in a roll size 6" x 6 yds, 1/16" plus thickness that compresses under weight. May also be used on the plastic head supports to create "a feeling" of softness. Beige color.

Catalog #VIP-T4700 Call for sample


**Soft Sponge**, a closed-cell foam with an adhesive back which can be cleaned with soap and water is available in sheet size 18" x 23.5" x 1/8", 4 sheets/case. Beige color. Provides a more consistent "spacer," with your choice of either removing or not removing it. Helpful when you want to add to your head support due to shrinkage of the tumor, etc.

Catalog #VIP-T4445 Call for sample









# ANATOMICAL DRAWING LABELS & RUBBER STAMPS

With these anatomical drawings you can record your treatment fields in the chart of patients undergoing radiation therapy treatment. Each drawing comes on a printed sheet with eight easy-to-peel, pressure sensitive, self-adhesive labels, or may be ordered on a rubber stamp.

■ Easy to use ■ Easy to draw/sketch on


■ No gluing necessary ■ Less expensive than photos

Illustration shown actual size.

Order by illustration number(s).

Each sheet contains eight identical drawings (actual size 10.2cm x 6.3cm).


Go to **www.vanarsdaleproducts.com** for a complete list of all anatomical drawings available for labels and rubber stamps.







#### Drawings 70 and 71 shown actual size.





\*Each sheet contains **six** identical whole body drawings.





\*Each sheet contains  $\underline{\textbf{six}}$  identical whole body drawings.

71. Whole Body Posterior Anterior

## **Anatomical Drawings Ordering Information**

#### **Label Sheets**


Catalog #AD followed by the number of drawing being ordered. Example: Catalog # AD-17

#### Order by sheet number

ie: three each of #3, 4, 5 and four each of #16, 17, 18, 19, etc.

#### **AD-17 Example**

Left Tangential Breast



#### PATIENT TRANSFER BOARD

Designed to ease patient transfer from bed or stretcher to treatment table with minimum effort. Made of long lasting polyethylene material.

Catalog # VIP-40

Size: 72" x 22" x 1/4"

Straps for Transfer Board

Catalog # VIP-40-S

#### **WALL RACK**

An excellent way to store your patient transfer board.

Catalog # VIP-40-A 16"w x 4"h x 3/4"d peg holes 12-13" apart







#### WHEELCHAIR PATIENT TRANSFER MOVER

If your wheelchair has arms that are removable, there is a way to help save your back from strain when moving a patient from wheelchair onto the treatment table and back to the wheelchair. The "mover" is made of polyethylene material and just flexible enough to bend over table corners and under the patient, yet strong enough to withstand heavy abuse. A nice item to have around!









#### **ANTI-STATIC CLEANER**

Use Staticide for maximum static electricity control. Great for your patient transfer boards, as well as computers, glass, fabric, or plastic.
Non-toxic, non-flammable, biodegradable.
Long lasting, safe to use.



32 oz. Trigger Spray Bottle Catalog # VIP-42


## RADIATION THERAPY TABLE PADS

#### What Makes Our Pads Different?

- Minimizes beam contamination with the open Rx area
- Minimizes SSD drift as associated with solid pads
- Velcro-like bonding holds pads "in place"
- Easy to remove and return to desired position
- Easy to clean—washable vinyl
- Standard size 72"L x 18"W x 1" thick

Patients receiving radiation treatments have reason to complain of discomfort due to the hard surface of the treatment table. In severe cases, this could even cause treatment error due to patient motion during the course of treatment.


We recommend that the simulator and treatment table pads be identical due to set-up and treatment duplication, in addition to patient comfort.



\*Additional charge for pads thicker than 1", and longer than 82", and/or wider than 24".

## The Electron Block System

EBS is designed to make electron block pouring easy and more precise. The aluminum base plate promotes faster cooling and will not allow bulging of the frames enabling the blocks to fit into the accessory with ease and without breaking code reader pins.



## No more:

- C clamps
- angle bars
- cooling plates



- bulged frames
- broken pins
- embarrassing treatment delays

Catalog #	<b>Dimensions</b>	
EBS-6	(6X6X)	
EBS-10	(10X10)	
EBS-15	(15X15)	
EBS-20	(20X20)	
EBS-25	(25X25)	
EBS-5	All 5 sizes - S	Set

# STEP STOOL Catalog #BSS-60



- Durable stool with flared legs for stability.9" overall height.
- 18 gauge, chrome plated tubular steel constructionSteel top with skid resistant black mat.
  - Shipping weight: 8 lbs., UPS.

STEP STOOL with handrail Catalog #BSH-65



- Same stool as the BSS-60
- 34" overall assembled height.
- Handrail is 7/8" dia. chrome plated tubular steel.
  - Comes with hardware for easy assembly.
  - Shipping weight: 1 box@ 4 lbs. ea., UPS.

# VanArsdale Innovative Products, Inc. P.O. Box 10853 Pensacola, FL 32524 1-800-326-4761 Fax 850-476-8426

## **ORDER FORM**

'		lale.cs@gmail.com	DA	TE	
FROM:			P.O. #		
SHIPPING			BILLING		
ADDRESS:			ADDRESS:		
- I		1			
QTY ITEM# DESC			RIPTION	UNIT PRICE	
1		If Paying By Credit	Card Pleas	se Fill Out	
Name	on Credit Card				
Card #		Ex	p. Date		
Securi	ty Code on Back _				
	in LIDS Craused	Callagt #			
	ip UPS Ground ip Pre-pay &  Ad			Fax number ————————————————————————————————————	



#### 10% Discount



when you mention this coupon
Compliments from

VanArsdale Innovative Products, Inc.

#### January—February Sale

- Tattoo Ink
- Ink Cups
- Tattoo Stand
- Alcohol Pads

Items found on Page 16



#### 10% Discount



when you mention this coupon

Compliments from

VanArsdale Innovative Products, Inc.

#### March—April Sale

All Rulers

Items found on Page 25



#### 10% Discount



when you mention this coupon

Compliments from

VanArsdale Innovative Products, Inc.

#### May—June Sale

All Positioners

Items found on Pages 32-34



#### 10% Discount



when you mention this coupon

Compliments from

VanArsdale Innovative Products, Inc.

#### July—August Sale

Grip Rings All sizes

Items found on Page 31



#### 10% Discount



when you mention this coupon

Compliments from

VanArsdale Innovative Products, Inc.

#### <u>September—October Sale</u>

Pen Marking Products

Items found on Page 14



#### 10% Discount



when you mention this coupon

Compliments from

VanArsdale Innovative Products, Inc.

#### **November—December Sale**

 Mold Room (Accessories Only)
 Items found on Page 3 and 5 VanArsdale Innovative Products Inc. P.O. Box 10853 • Pensacola, FL 32524

> 1-800-326-4761 • 850-476-1078 24 Hour Fax 850-476-8426

vanarsdale.cs@gmail.com (Customer Service)
vanarsdale.ar@gmail.com (Accounts Receivable/Payable)

## Customers First... Friends for Life