

AVRASYA MEDICAL CONSUMBLES CATALOG



www.avrasyaexport.com

OpSet Adult Blood Pressure Cuffs

It is Latex-free, single patient use a cuff intended for adults

The cuffs can reduce mutual infection catching and hospitalacquired infections. Disposable, for single patient.

- Does not contain LaTeX Comes in 5 sizes
- Customized adapters for ensuring equality





OpSet Blood Pressure Bag

Blood Pressure Bag is safe, high-quality and reliable product for the treatment of intravenous or intrarterial infusion.

It is designed for single patient use and is for the control and prevention of mutual infection and the formation of hospital infections.

The safety valve prevents excessive inflation (330 mmHg). - Does not contain LaTeX Maintains the desired infusion pressure.

- Colour-coded indicator provides simple pressure monitoring (0-300 mmHg).
- 3-way stopcock provides precision pressure control.
- Hook slide provides one-hand use. 100% tested.
- 500 ml, 1000 ml, 3000 ml





OpMask Single Blood Bag

The capacity is 450 ml.

Donor needle gauge size: 16G

Each blood bag kit is in a PE and PET vacuum packed bag, and 5-set packs of aluminum foil.

Anticoagulant solution should contain 63 ml of CPDA-1 (citrate phospate dextrose -adenin 1) anticoagulant solution. Packaging must be in packaging packed in a sterile way and individually.





OpSet I.V. CANNULA

Catheter manufactured high quality FEP Teflon
Flexible wings for a better placement
Safe placement on the skin, balanced
and secure grip
Catheter durable against twisting with silicone end
14G Orange/16G Grey/18G Green/20G
Pink/22G Blue/24G Yellow/26G Purple
BOX of 100 pcs /CARTON of 2000 pcs.





OpSet 3-Way Stopcock

Easy Stopcock Rotation.

High Pressure Leakproofing Test.

Understandable surface format regulates maximum fluid flow with minimum turbulence.

Thanks to 6% luer-cone, compatibility with any standard product.

On Lipid Resistance Level.

Box of 50 pcs./ Pack of 1000 pcs.



OpMask Blood Transfusion Kit

Air inlet drilling end 20 drops = 1 ml 200-µm filter Roller drop tuner

Luer lock connector with Injection port, 18G needle

Hose length: 150 cm - in single sterile packaging

Packaging: 400 pcs/pack





Ethox Blood Pressure Device

Surgi-Press is a durable and effective device for appropriate and correct infusion.

A piston ensuring a circulation movement has been designed for consistent and accurate infusion. Surgi-Press, eliminating the need for continuous monitoring, creates a constant pressure without periodic pumping.

The bulb can provide automatic infusion without meter and pumping. Easy to install, easy to use.

No need for periodical pumping and setting. Minimum pressure change/difference





OpMask Infusion Set with Burette

100 ml-150 ml application

Suitable for administering liquid drugs into the body when baby transfusion where infusion dose needs to be administered in an accurate and controlled and timely way.

The main structure and Properties;

Burette has a standard scale. Liquid can be seen clearly on the tubing wall;

Drip Chamber is transparent, easy to observe;





OpMask Drop Set Kit

Y-injection port, with tube cap, male/female luer lock with protective cap

Y-injection port LaTeX-containing/latex-free

Usage: used to set the exact flow rate

sterilized with E.O.

Single package: Blister

BOX of 50 pcs./PACK of 400 pcs.





OpMask Serum Set

Air-wayed drilling bit for infusion bottles and bags Airline filter with cap - 20 drops = 1 ml, 15 μ m fluid filter

Roller clamp - with injection port - Luer slip connector 21 G needle- Hose length: 150 cm

Single sterile packaging Packing: 500 pcs/pack





OpMask Adult Urine Bag with Valve

One-way valve - Capped tapered input connection Soft and kink resistant PVC tubing

ML type scale on the bag - Sterile and disposable

Tube length: 90 cm, 100 cm

Packaging: PE bundle or paper-plastic packaging

2000 ml, anti-flow device and output

With white Floors or transparent bag

sterilized with E.O.

Stopcocked/stopcockless (hanger hole)

Single package: PE/Blister

Package: 250 PCs





OpMask Speculum

Small, Medium, Large.

Made from polystyrene (PS), suitable for medical use

non-toxic, does not irritate.
Sterile or Clean Room.
Each product is packaged as single.
Only for single use, and should be discarded after using

Types: Binding, spinned, Gear Speculum Pack of 100 or 200 PCs.



OpSet T Spiral 380 A

Contraceptive!


It contains 176 mg. of copper wire on the vertical arm of each t-shaped plastic frame.

Additionally, there is a cuff containing 66.5 mg of approximately copper placed on each horizontal arm.

The total surface of the device is 380 mm2.

In boxes of 25, in ready-to-use single packaging.

Sterile





OpSet 2-3 Way Latex Foley Probes

Traumatic intubation is provided with super sof tapered eyes which have non-distal-ends, and that are proportional to each other.

This is achieved by the ultra-thin high elasticity of the balloon.

100% silicone elastomer.

26 Dimensions between 26 Fr 16Fr, Depending on the customer needs, small balloon capacity is also available.

FG Color Codes:

8 Black, 10 Grey, 12 White, 14 Green 16 Orange, 18 Red, 20 Yellow, 22 Cream, 24 Blue, 26 Pink. Sizes



OpSet Intubation Tube- Cuffed

It providing easily locating the radio opaque line and tube position.

The tructure to prevent King

At E. T.O., a standard proximal end of 15 mm is available.

It is sterile and non-progenic.

Dimensions: 2.5mm, 3.0mm, 3.5mm, 4.0mm,

4.5mm, 5.0mm, 5.5mm, 6.0mm, 6.5mm, 7.0mm,

7.5mm, 8.0mm, 8.5mm, 9.0mm, 9.5mm, 10mm





OpSet Tracheostomy Tube - cuffed and uncuffed

Used in patients whose airway is through a tracheostomy and whose ventilation is mechanical.

Made from Non-toxic, irritant-free medical grade PVC.

Thin-walled, heat-sensitive tube softens with the body heat and is resistant to bending.

Superior ease of application provided with lockable obturator in guide wire technique, Seldinger guide wire is widely used for comforting the intubation.

Distal ends having rounded corners in atraumatic intubation and intubation is intended for increasing the comfort of the patient with softer inputs and outputs.

It can determine the positon of the tube along the radio opaque line and the marked part.

SIZES

6.0 mm, 6.5 mm, 7.0 mm, 7.5 mm, 8.0 mm, 8.5 mm, 9.0 mm





OpMask Cautery Pen

Package: 100 PCs

170 mm-ABS pen view, flow design 7 cm Overall Stainless Steel Tip 2.36 mm- tip diameter- 3 m- PVC medical cable Maximum output current: 1A - Frequency region: 0-1.0 MHZ Flow resistance: 200Ω - disposable, sterile package Each of them is packed with an internal PE plastic bag, and an outer sterile bag.



OpMask Cautery Plate (Single-Double)

Properties:

Self adhesive, soft
Superconductivity and adherence
With electricity protection
preventive against infections
Compatible with the any brand device
Box:50 pcs./Pack:400 pcs.



OpMask Stapler Remover

DISPOSABLE/STERILIZED
Can be sterilized with ETO
In single sterile blister packaging
Packaging: 200 PCs/pack



OpMask Vinyl Examination Gloves (Powder-Free)

Double layer of PVC/PU film provides a double protection. It has a special formulation, is odorless.

Thanks to superb slippery internal surfaces, it can be worn easily.

Small, Medium, Large. Box 100 pcs/Pack: 2000 pcs.





CUTANPLAST® Standard Spongostan

Bloodstopper Pad In the boxes of 20 pieces Sterile.

Areas of use;

Parenchymal surgery and Breast Surgery

Dimensions:70x50x10 mm

Bloodstopper Pad

In boxes of 10, Sterile.

It is used in all types of anal surgery.

Dimensions: 80x30 mm





SURUWAX ™ Sterile Bone Wax

SURUWAX ™, is a sterile mixture of beeswax, paraffin wax and isoproply palmitate;

Indications

SURUWAX [™] is used for the mechanical haemostasis of the bone in injuries/surgeries.

Used in applications to patch the holes in the skull.

Open heart and thoracic surgery, brain surgery,

Orthopedics Surgery, Contra Indications,

On infected areas where bone regeneration and/or fusion are requested

Advantages: easy to shape and implement easy-to-open-cover, disposable, sterile packaging does not effect the biological structure of the body Minimal absorbable, sterilization,

Gamma Radiation

Packing: 12 PCs. 2.5 gms





OpMask Skin Stapler

DISPOSABLE/STERILIZED
Stainless steel
35 pieces of wide (w) staples
In single sterile blister packaging
Packaging: 5 PCs/box, 100 PCs/Pack





OpMask Automatic Tourniquet

Applied to single bones.

(Applied to Biceps bone in the arm, to the thigh bone in the leg)

Tourniquet should be applied in serious and life-threatening bleedings that cannot be prevented otherwise.

Tourniquet is only needed in case of a large vein (artery) cut.

Tourniquet should be positioned a little above the wound, in a way not to contact the joint and the wound.



