EMERGENCY EYEWASH & SHOWER EQUIPMENT

SUPPLEMENTAL EYEWASHES include bottled eyewash and other devices that do not meet 15-minute drench criteria.

These supplemental solutions are designed to treat nuisance hazards and/or as an initial first aid to injured workers during the path to certified eyewash, but they do not meet the 15 minute drench criteria of the ANSI standard Z358.1. Bottled eyewash is designed to be positioned close to the point of need.

Acceptable for immediate first-aid, bottled eyewashes should support, but not substitute, a primary eyewash or shower unit.


SALINE EYEWASH BOTTLES

- Contains a buffered saline solution that is sterile and superior to tap water for emergency eye care
- Blow-Fill-Seal (BFS) bottle is tamper resistant, and maintains the sterility and integrity of the eyewash
- · Easy to remove and open in emergency situations
- 1, 4, and 8 oz. bottles have user controlled delivery, delivering drops or an irrigating stream of solution, depending on the squeeze pressure applied to the bottle
- Trilingual instructions (English, French and Spanish)
- Does not meet ANSI standard Z358.1 15-minute flushing requirement



Honeywell

Model No.	Mfg. No.	Cap. Oz.	Price /Each
SEC467	32-000451-0000	1	
SEC468	32-000452-0000	4	
SEC469	32-000445-0000	8	
SEC470	32-000454-0000	16	
SEC471	32-000455-0000	32	



SALINE EYEWASH WALL STATIONS

- Holds bottle firmly in place while allowing easy access in an emergency
- Mounts to wall with provided screws
- Integral overhand protects bottles from dust and debris
- · Durable plastic is easy to clean
- Simple, highly visible operating instructions are provided on each unit
- Bottle contains a buffered saline solution that is sterile and superior to tap water for emergency eye care
- Blow-Fill-Seal (BFS) bottle is tamper resistant, and maintains the sterility and integrity of the eyewash
- Delivers drops or an irrigating stream of solution, depending on the squeeze pressure applied to the bottle
- Trilingual instructions
- (English, French and Spanish)

 Does not meet ANSI standard Z358.1
- 15-minute flushing requirement





Honeywell

Model	Mfg.	Station	Bottle Cap.	Price
No.	No.	Type	Oz.	/Each
SEC472	32-000460-0000	Single	16	
SEC473	32-000465-0000	Double	16	
SEC474	32-000461-0000	Single	32	
SEC475	32-000462-0000	Double	32	

NEUTRALISER™ EYEWASH & FACE-RINSE

- Provides immediate first aid
- and neutralized chemical burns

 Ingredients neutralize acids and alkalis
- better than standard saline solution
- Responds to all accidents with one "all risks" formula, which includes antiseptics to prevent infection
- Quick rinse design 1/4 twist to open
- Bottle pods allow you to locate the solution close to the potential risk







Model	Mfg.		Price /Each	
No.	No.	Description		
SEI502	80-1600-1	1L Neutraliser™ Solution		
SEI503	80-1600-2	474 ml Neutraliser™ Solution		
SEI504	80-1600-1BN	1L BottlePod™ w/ Neutraliser™ Sign		
SEI505	80-1600-2BN	474 ml BottlePod™ w/ Neutraliser™ Sign		