



ALTUS HEALTHCARE SOLUTIONS

Hygiene is Our Foundation – Quality You Can Trust

SURGEONS GOWN & SMMS FABRIC

Altus Healthcare Solutions specializes in manufacturing a variety of single-use surgeons gowns made from SMMS fabric, available in multiple GSM options. Our products are designed for maximum safety, comfort, and durability in various surgical environments.



Licence Number:
MFG/MD/2025/000432

TYPES OF SURGEONS GOWNS

Type of Gown	Type of Gown Description
Surgeons Gown	Standard gown designed for general surgical procedures.
Wrap Around Gown	Features a wrap design for ease of wearing and comfort.
Reinforced Wraparound Gown	Enhanced protection with reinforced areas, suitable for high-risk procedures.
Special Protection Gown	Special Protection Gown Designed for maximum barrier protection in critical environments.

WHY CHOOSE OUR SURGEONS GOWNS?

- **Tailored Variety:** A range of gown types catering to different procedural needs, ensuring safety and comfort across various surgical environments.
- **High-Quality Materials:** Manufactured from SMMS fabric for reliable protection against contaminants.
- **Custom Packaging Options:** Offered in both regular and book folding designs for easy handling and storage.
- **Single-Use Design:** Each gown is designed for one


Contact Us:

Mail: Marketing@altus-healthcare.com
Visit US: www.altus-healthcare.com

+91 9497701278
+91 7034433566



INSIDE



BREATHABLE NON-WOVEN FABRIC




VELCRO ADJUSTABLE NECKLINE



KNITTED CUFF



REINFORCED FROM INSIDE



WATER REPELLENT



ULTRASONIC SEALING

PRODUCT SPECIFICATIONS

Attribute	35 GSM	40 GSM	43 GSM
Material	SMMS (Spunbond, Meltblown, Spunbond)	SMMS (Spunbond, Meltblown, Spunbond)	SMMS (Spunbond, Meltblown, Spunbond)
Weight	35 grams per square meter	40 grams per square meter	43 grams per square meter
Barrier Protection	Moderate liquid and particulate barrier	Good liquid and particulate barrier	Superior liquid and particulate barrier
Breathability	High	Moderate	Lower than 35 and 40 GSM
Comfort	Very comfortable for short procedures	Comfortable for extended wear	Ideal for lengthy surgical sessions
Applications	Minor surgeries, outpatient procedures	General surgeries, procedural use	Critical surgeries, high-risk environments
Single-use	Yes	Yes	Yes
Size	130 x 145cm	130 x 145 cm	130 x 145 cm
Color	Medical Blue	Medical Blue	Medical Blue
Packaging	Individually wrapped in sterile	Individually wrapped in sterile	Individually wrapped in sterile
Attribute	35 GSM	40 GSM	43 GSM
	Options with regular and book folding	Options with regular and book folding	Options with regular and book folding
Sterilization	Options Available in sterile options	Available in sterile options	Available in sterile options
Folding Options	Regular Folding, Book Folding	Regular Folding, Book Folding	Regular Folding, Book Folding
Price Range	₹43 - ₹100	₹47 - ₹110	₹50 - ₹115

SPECIAL ORDERS

We also accommodate special orders for additional GSM options such as 50 GSM (₹52 - ₹120) and 60 GSM (₹57 - ₹125) SMMS, tailored to meet your specific requirements.



Manufactured By :

Altus Health Care Solutions

No -60 G/2- Trichy Road, Pappampatti Piruvu, Coimbatore, Tamilnadu - 641 103

Web : www.altus-healthcare.com | Email : altushealthcaresolutions@gmail.com

Customer Care : +91 95663 95077/ 70344 33566