ENGINEERING SPECIFICATIONS: TOC Systems ON-LINE ULTRA-PURE TOC

Model/Name	StarTOC ON-LINE Ultra-Pure Analyzer	Notes
Amalyzaia Mashad	Hitas Divas (Non Diagrapsias Infrased (NDID) CO. Datastics	
Analysis Method	Ultra-Pure/Non-Dispersive Infrared (NDIR) CO ₂ Detection	
Possible Analytes Measured	Total Organic Carbon; Total Carbon; Total Inorganic Carbon; Non-Purgeable Organic Carbon; TOC-True (including purgeable Organics)	
Detector Type (CO ₂)	NDIR (solid state; no moving parts; computer-controlled; non-reflective sample cell-impervious to corrosion and guaranteed for 5 years)	
Control/Data Handling	Industrial Microsoft Windows CE Computer, Touch Screen, Paperless Chart Recorder	
Sample Introduction	Digital Injection	
Measurement Specifications		
Measurement Ranges (ppb)	1000 10,000	!
Accuracy/Repeatability (%)	3% of Full Scale or +/-40ppb whichever is greater	
Carrier Gas Flow (mL/min.)	300 mL/max - Computer Controlled Mass Flow Controller CO ₂ & HC - FREE AIR, OR O ₂ 15 +/- 2 PSI	Bottled Gas Supply Preferred
Average Analysis/Response Time (minutes) Including TIC	8-12	
Outputs	RS-232 RS-485(option) 4-20 mA TOC Output Additional 4-20 mA TOC Output available as option Alarm Relays (3) 2 Level Alarms 1 Malfunction Alarm Loss of Carrier	
	Out of Service Signal (when unit in calibration or validation check)	

Flat Panel Color Paperless Chart Recorder Display Data Storage Internal Memory Sample Requirements TOC Systems offers Temperature Range 4° - 90°C cooler for higher 40° - 194° F temperatures Fast by-pass loop Sample flow rate 20 mL/min suggested Inlet pressure 0-0.5 psi Sample Drain gravity/air break UTILITIES Required Power: 100/240 VAC 50/60 HZ. 15 Amp Service Carrier Gas CO₂ & HC - free air or O₂ (300 mL/minute-max.); 15 +/- 2 PSI Reagents Phosphoric Acid Calibration Standards D.I. Water Sodium Persulfate **Environment** Operating Temperature:10° - 50°C 50° - 122°F Construction Dual Compartment FRP Cabinet **Enclosure** Dimensions (HXWXD) 36 x 24 x 12 (in)

Mounting

Rack or Wall Mount

Weight	34 Kg 75 Lbs	
Area Classification Option	NEMA - 4, IP 65 Enclosure	
Compliance/Certification	EPA Methods	
	Complying with all International Standards, such as: DIN-EN 1484, DIN-ENV 12260, DIN 38409-H3, ISO 8245, Standard Method 5310B, Standard Method 5310C, Standard Method 5310D, USEPA 415, USEPA 9060, ASTM D5173, EN 13137	
General Features	Historical Data/Time and Date Stamped	
	Multi-Stream Analysis Auto-Calibration Auto-Validation Auto-Cleaning	User to specify number of streams
	TOC Systems also offers custom options.	

(All performance specifications have been verified in a controlled laboratory environment. Actual field performance may vary with application measuring range and detection limits depend on the method, injection volume, vessel purity, chemicals and gases used, and the qualification of the operators.)