

**AJL**  
vet



## FELINE INTRAOCULAR LENS

### PRODUCT DESCRIPTION AND INDICATION

Ophthalmic medical device designed for implantation in feline crystalline extraction procedures.

The product is made of HEMA (UV), a biocompatible acrylic copolymer that offers excellent geometric and folding properties allowing the lens implantation through a small incision.

Indicated for implantation in the posterior chamber, replacing the feline crystalline extracted for the correction of aphakia.

Veterinary use only.

### HOW SUPPLIED

Single and sterile, in a glass vial, immersed in high purified water.

A single-use injector set may be provided with each intraocular lens.



### PRODUCT SPECIFICATIONS AND MODELS

#### Optic Design

Biconvex  
7 mm diameter  
+53.5 D

#### Available Lengths and References

12.00 mm [ref. CAT701200]  
13.00 mm [ref. CAT701300]  
14.00 mm [ref. CAT701400]