


AJL
vet

RIGID CANINE IOL

DESCRIPTION AND INDICATIONS

AJL Rigid Canine Intraocular Lens, made of polymethacrylate (PMMA) with UV filter, is an ophthalmic device designed for implantation in canine crystalline extraction procedures.

Indicated for implantation in the posterior chamber, replacing the canine crystalline extracted for the correction of aphakia.



DIFFERENT MODELS

Ref: P651400

14mm (length), 6.5mm (optic), 41D (power)

Ref: P701600

16mm (length), 7mm (optic), 41D (power)

HOW SUPPLIED

Single and sterile. Ready to use.