



Knotless Design Implant


Screw (PLGA/ $\beta$ -TCP)

Insertion Tip (PEEK)

Hybrid Dual Thread Pattern

# AlternatiV+<sup>TM</sup> Max

Knotless Anchor for Use with BioBrace<sup>®</sup>



BioBrace<sup>®</sup> Implants Sold Separately

PRODUCT OVERVIEW

 CONMED

# AlternatiV+™ Max

Knotless Anchor for Use with BioBrace®



**The AlternatiV+™ Max Knotless Anchor** is intended for fixation of suture (soft tissue) to bone in the shoulder, foot/ankle, knee, and elbow in the following procedures:

» **Shoulder:** Rotator Cuff Repair, Bankart Repair, SLAP Lesion Repair, Biceps Tenodesis, Acromio-Clavicular Separation Repair, Deltoid Repair, Capsular Shift or Capsulolabral Reconstruction

» **Foot/Ankle:** Lateral Stabilization, Medial Stabilization, Achilles Tendon Repair

» **Knee:** Medial Collateral Ligament Repair, Lateral Collateral Ligament Repair, Posterior Oblique Ligament Repair, Iliotibial Band Tenodesis

» **Elbow:** Biceps Tendon Reattachment, Ulnar or Radial Collateral Ligament Reconstruction

## Ordering Information

Description	Catalog Number
AlternatiV+™ Max Knotless 4.75mm Anchor	20KPN4751
BioBrace® FA 5mm x 250mm	BB5X250-FA
BioBrace® FA 23mm x 30mm	BB23X30-FA



Distributed by:



6000 Poplar Ave., Suite 115  
Memphis, TN 38119  
844 602 6637  
www.conmed.com/footandankle

28, chemin du Petit Bois  
69130 Ecully  
France  
+33 (0)4 72 29 26 26

CONMED Corporation  
1311 Concept Blvd.  
Largo, FL 33773  
USA

AlternatiV+™ is a trademark of, and manufactured by Genesis Biologics, Inc. 17791 Mitchell North, Suite D, Irvine, CA 92614. All other ™ and ® denote Trademarks and Registered Trademarks of CONMED or its affiliates. ©2024 CONMED Corporation. All Rights Reserved.

Some products not available in all markets. This document provides information regarding medical devices and instruments in surgical procedures. It is not medical advice and healthcare professionals should use their own judgment before using to treat a particular patient. Healthcare professionals should be trained in the use of such devices before surgery and should always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any product.

PB-GB001-01 Rev A 092024