


# HYALURONIC ACID


We introduce our innovative brand of hyaluronic acid-based fillers, with different concentrations to satisfy all needs and guarantee an immediate, safe and non-invasive lifting effect.

From correction of expression lines to the treatment of large imperfections, our facial and body fillers meet the needs of an increasingly demanding clientele and professionals who are constantly concerned about the volume-enhancing and moisturizing properties of the products.

Our medical devices stand out for their innovative formula that uses latest generation cross-linked hyaluronic acids, extremely thin needles of the best Japanese steel and CE-certified formulas that guarantee total resorption.



ACNE SCARS



PERIOcular  
WRINKLES

BAR  
CODE





EVERY 4 MONTHS

Packaging

1x1ml

Volumizing effect



PH

7

Needle

30G



EVERY 4 MONTHS

20MG/ML. ts special formula is indicated to treat

the **SMILE LINE**, especially to support the

problems of medium-young skins. It is

recommended to use a **30G needle.**




SUPERFICIAL MIDDLE  
FOLDS (Glabella)

NOSE

DEEP MIDDLE FOLDS

LIPS





EVERY 4-6 MONTHS

Packaging

1x1ml

Volumizing effect




PH

7

Needle

27G




EVERY 4-6 MONTHS

22MG/ML. Its special formula is indicated to treat the **T ZONE** especially to support the problems of medium-young skins. It is recommended to use a

**27G needle.**









EVERY 6 MONTHS

Packaging

1x1ml

Volumizing effect



PH

7

Needle

27G



EVERY 6 MONTHS

24MG/ML. Opera III is based on an high concentration of hyaluronic acid that make the product suitable to support the problems of **MATURE SKINS**, especially the areas that

outline the smile. It is recommended to use a **27G**

**needle.**



DEEP  
FURROWES

CHEEKS

MALAR  
AREA

CHIN







EVERY 6-12 MONTHS

Packaging

1x1ml

Volumizing effect



PH

7

Needle

25G



EVERY 6-12 MONTHS


25MG/ML. Its high concentration make the product

especially suitable to support problems of

**MATURE AND VERY MATURE SKINS.** It is

recommended to use a **25G needle.**






ANTI AGING



DAMAGES  
CAUSED BY  
CHRONO AND  
PHOTO AGING





1 TREATMENT PER WEEK  
4 TREATMENT IN TOTAL

Packaging

1x1ml

Volumizing effect


○

PH

7

Needle

30G



1 TREATMENT PER WEEK  
4 TREATMENT IN TOTAL


20MG/ML. Opera Rev is the ideal product to correct

**MILD CUTANEOUS DEPRESSIONS** of the

face. It is recommended to use a **30G** with the BAP

technique





## BODY FILLER

20 mg/ml. 10 ml doses.

Revolutionary, advanced macromolecular hyaluronic acid body filler, not of animal origin, designed to reshape the body in a completely safe way, with no need for surgery.

The body filler makes it possible to harmonise large volumes such as buttocks, calves, forearms and areas requiring natural reshaping.

