

A Perfect Smile Starts with Great Data!

Pre-Op extraoral data –
Photo or face scan

- Pre-op and/or Mock Up Full Face Natural Smile
- Pre-op and/or Mock Up Full Face Retracted Smile

Corresponding intraoral data

- Pre-op and/or Mock Up upper, lower, and occlusal scans

AND/OR

- 360 degree scan of any existing prosthetics with wash impression and bite registration scan

All photos must correspond to intraoral data and patient cannot be numb. The intraoral and extraoral data work as a pair.

WWW.BBDL.DESIGN



Day of surgery data

BEFORE AND AFTER implant placement:

- Jaw scan of TADs, breadcrumb, or other stationary surgical marker scan.

AFTER IMPLANT PLACEMENT

- Photogrammetry (preferred), scan ladder, or verification bar
- Implant and scanbody information, pontic sites preferences regarding shape and pressure, and desired screw type (recommend Vortex, SIN smart angle, Rosen, or Dess)

WhatsApp: 661-803-9695