

Parts of ECG Machine

Electrode Cable Connector

- Leads are color-coded:
 - RA (Red)**: Right Arm
 - LA (Yellow)**: Left Arm
 - RL (Green)**: Right Leg
 - LL (Black)**: Left Leg
 - V1-V6 (White)**: Chest Leads

Lead Wires


- Connect the electrodes to the machine.

ECG Paper

- Prints the heart's electrical activity as a waveform.

Control Buttons & Dials

- Used to start, adjust, & control the ECG machine.



ECG Paper

- Small Square = 0.04 sec
- Large Square = 0.20 sec

Power Button

- Turns the ECG machine ON and OFF.

Power Cable

- Supplies power to the ECG machine.

Monitor/Display Screen

- Shows the waveform of the heart's electrical activity and the heart rate.

Leads are color-coded:

- RA (Red)**: Right Arm
- LA (Yellow)**: Left Arm
- RL (Green)**: Right Leg
- LL (Black)**: Left Leg
- V1-V6 (White)**: Chest Leads

