





People-oriented" technology



SPINAL SYSTEM PRODUCT CATALOG



www.hof-scientific.com





COMPANY PROFILE



FAQ

HOF Scientific (UK) & Bensberg Sciences (USA) – A Global Holding Group
The main Base of the Head Office in United Kingdom's
HOF Scientific, in partnership with Bensberg Sciences, proudly leads a dynamic international consortium of companies headquartered across the United Kingdom, the United States, the United Arab Emirates, Germany.
Specializing in medical devices, food supply, and energy solutions (Oil and Gas), our group brings together deep industry knowledge, regulatory expertise, and a commitment to innovation across critical sectors.
With over 15 years of experience, we serve as a strategic bridge between global manufacturers and high-demand markets—particularly across the Middle East ensuring seamless supply chains, market access, and regulatory compliance.
As a group, we are united by a vision of long-term partnerships, quality excellence, and market leadership, The group also leads specialized consultancy services in project management, engineering design, construction oversight, and geotechnical and surveying solutions. From planning to execution, we help clients deliver high-quality, well-engineered infrastructure with confidence and efficiency.



PRODUCT CONTENT

ZION TM 5.5 Posterior Thoracolumbar Internal Fixation System	01-03
ATLASH 5.5 Posterior Thoracolumbar Internal Fixation Syste	04-05
ETNA TM 5.5/6.0 Posterior Thoracolumbar Internal Fixation System	06-13
EVEREST™ 5.5 MIS Spine System	14-15
CERVINO TM Anterior Cervical Plate	16
ALPST™ Posterior Cervical Internal Fixation System	17-19
Titanium Mesh	20
Cage System	21
Instrument Set	22-30



ZION™ 5.5 Posterior Thoracolumbar Internal Fixation System



ZION™ Multi Axial Pedicle Screw



ZION™ Reduction Multi Axial
Pedicle Screw



Indications

- Disc degenerative disease, Spondylolisthesis, Stenosis, Spinal fracture, Spinal tumor, Spinal deformity

Features & Benefits

- Low profile, anatomical contoured implants and smooth surface allow for increased versatility and reduce potential of irritation
- Single axial and multi axial screw options for optimal screw insertion trajectory
- Top opening and tightening screws and hooks
- Multi axial screw motion range: ±18°
- Increase flexibility in screw and rod placement
- Reduce screw head impact

P/N	Material	Spec.
010013021025	TC4	Φ4.0×25
010013021030	TC4	Φ4.0×30
010013021035	TC4	Φ4.0×35
010013021040	TC4	Φ4.0×40
010013022025	TC4	Φ4.5×25
010013022030	TC4	Φ4.5×30
010013022035	TC4	Φ4.5×35
010013022040	TC4	Φ4.5×40
010013023025	TC4	Φ5.0×25
010013023030	TC4	Φ5.0×30
010013023035	TC4	Φ5.0×35
010013023040	TC4	Φ5.0×40
010013023045	TC4	Φ5.0×45
010013023050	TC4	Φ5.0×50
010013024025	TC4	Φ5.5×25
010013024030	TC4	Φ5.5×30
010013024035	TC4	Φ5.5×35
010013024040	TC4	Φ5.5×40
010013024045	TC4	Φ5.5×45
010013024050	TC4	Φ5.5×50
010013025025	TC4	Φ6.0×25

P/N	Material	Spec.
010013025030	TC4	Φ6.0×30
010013025035	TC4	Φ6.0×35
010013025040	TC4	Φ6.0×40
010013025045	TC4	Φ6.0×45
010013025050	TC4	Φ6.0×50
010013025055	TC4	Φ6.0×55
010013026025	TC4	Φ6.5×25
010013026030	TC4	Φ6.5×30
010013026035	TC4	Φ6.5×35
010013026040	TC4	Φ6.5×40
010013026045	TC4	Φ6.5×45
010013026050	TC4	Φ6.5×50
010013026055	TC4	Φ6.5×55
010013027025	TC4	Φ7.0×25
010013027030	TC4	Φ7.0×30
010013027035	TC4	Φ7.0×35
010013027040	TC4	Φ7.0×40
010013027045	TC4	Φ7.0×45
010013027050	TC4	Φ7.0×50
010013027055	TC4	Φ7.0×55


P/N	Material	Spec.
010013640025	TC4	Φ4.0×25
010013640030	TC4	Φ4.0×30
010013640035	TC4	Φ4.0×35
010013640040	TC4	Φ4.0×40
010013640250	TC4	Φ4.5×25
010013640300	TC4	Φ4.5×30
010013640350	TC4	Φ4.5×35
010013640400	TC4	Φ4.5×40
010013650250	TC4	Φ5.0×25
010013650300	TC4	Φ5.0×30
010013650350	TC4	Φ5.0×35
010013650400	TC4	Φ5.0×40
010013650450	TC4	Φ5.0×45
010013650500	TC4	Φ5.0×50
010013650250	TC4	Φ5.5×25
010013650300	TC4	Φ5.5×30
010013650350	TC4	Φ5.5×35
010013650400	TC4	Φ5.5×40
010013650450	TC4	Φ5.5×45
010013650500	TC4	Φ5.5×50
010013660250	TC4	Φ6.0×25

P/N	Material	Spec.
010013660300	TC4	Φ6.0×30
010013660350	TC4	Φ6.0×35
010013660400	TC4	Φ6.0×40
010013660450	TC4	Φ6.0×45
010013660500	TC4	Φ6.0×50
010013660550	TC4	Φ6.0×55
010013665250	TC4	Φ6.5×25
010013665300	TC4	Φ6.5×30
010013665350	TC4	Φ6.5×35
010013665400	TC4	Φ6.5×40
010013665450	TC4	Φ6.5×45
010013665500	TC4	Φ6.5×50
010013665550	TC4	Φ6.5×55
010013670250	TC4	Φ7.0×25
010013670300	TC4	Φ7.0×30
010013670350	TC4	Φ7.0×35
010013670400	TC4	Φ7.0×40
010013670450	TC4	Φ7.0×45
010013670500	TC4	Φ7.0×50
010013670550	TC4	Φ7.0×55


ZION™ Single Axial Pedicle Screw



ZION™ Reduction Single Axial Pedicle Screw



ZION™ 5.5mm Rod



P/N	Material	Spec.
010013028025	TC4	Φ4.0×25
010013028030	TC4	Φ4.0×30
010013028035	TC4	Φ4.0×35
010013028040	TC4	Φ4.0×40
010013029025	TC4	Φ4.5×25
010013029030	TC4	Φ4.5×30
010013029035	TC4	Φ4.5×35
010013029040	TC4	Φ4.5×40
010013030025	TC4	Φ5.0×25
010013030030	TC4	Φ5.0×30
010013030035	TC4	Φ5.0×35
010013030040	TC4	Φ5.0×40
010013030045	TC4	Φ5.0×45
010013030050	TC4	Φ5.0×50
010013031025	TC4	Φ5.5×25
010013031030	TC4	Φ5.5×30
010013031035	TC4	Φ5.5×35
010013031040	TC4	Φ5.5×40
010013031045	TC4	Φ5.5×45
010013031050	TC4	Φ5.5×50
010013032025	TC4	Φ6.0×25

P/N	Material	Spec.
010013032030	TC4	Φ6.0×30
010013032035	TC4	Φ6.0×35
010013032040	TC4	Φ6.0×40
010013032045	TC4	Φ6.0×45
010013032050	TC4	Φ6.0×50
010013032055	TC4	Φ6.0×55
010013033025	TC4	Φ6.5×25
010013033030	TC4	Φ6.5×30
010013033035	TC4	Φ6.5×35
010013033040	TC4	Φ6.5×40
010013033045	TC4	Φ6.5×45
010013033050	TC4	Φ6.5×50
010013033055	TC4	Φ6.5×55
010013034025	TC4	Φ7.0×25
010013034030	TC4	Φ7.0×30
010013034035	TC4	Φ7.0×35
010013034040	TC4	Φ7.0×40
010013034045	TC4	Φ7.0×45
010013034050	TC4	Φ7.0×50
010013034055	TC4	Φ7.0×55

P/N	Material	Spec.
010013640025	TC4	Φ4.0×25
010013640030	TC4	Φ4.0×30
010013640035	TC4	Φ4.0×35
010013640040	TC4	Φ4.0×40
010013640045	TC4	Φ4.0×45
010013640050	TC4	Φ4.0×50
010013640055	TC4	Φ4.0×55
010013640125	TC4	Φ4.5×25
010013640130	TC4	Φ4.5×30
010013640135	TC4	Φ4.5×35
010013640140	TC4	Φ4.5×40
010013640145	TC4	Φ4.5×45
010013640150	TC4	Φ4.5×50
010013640155	TC4	Φ4.5×55
010013640250	TC4	Φ5.0×25
010013640300	TC4	Φ5.0×30
010013640350	TC4	Φ5.0×35
010013640400	TC4	Φ5.0×40
010013640450	TC4	Φ5.0×45
010013640500	TC4	Φ5.0×50
010013640550	TC4	Φ5.0×55
010013645250	TC4	Φ5.5×25
010013645300	TC4	Φ5.5×30
010013645350	TC4	Φ5.5×35
010013645400	TC4	Φ5.5×40
010013645450	TC4	Φ5.5×45
010013645500	TC4	Φ5.5×50
010013645550	TC4	Φ5.5×55
010013650250	TC4	Φ6.0×25


P/N	Material	Spec.
010013650300	TC4	Φ6.0×30
010013650350	TC4	Φ6.0×35
010013650400	TC4	Φ6.0×40
010013650450	TC4	Φ6.0×45
010013650500	TC4	Φ6.0×50
010013650550	TC4	Φ6.0×55
010013655250	TC4	Φ6.5×25
010013655300	TC4	Φ6.5×30
010013655350	TC4	Φ6.5×35
010013655400	TC4	Φ6.5×40
010013655450	TC4	Φ6.5×45
010013655500	TC4	Φ6.5×50
010013655550	TC4	Φ6.5×55
010013660250	TC4	Φ7.0×25
010013660300	TC4	Φ7.0×30
010013660350	TC4	Φ7.0×35
010013660400	TC4	Φ7.0×40
010013660450	TC4	Φ7.0×45
010013660500	TC4	Φ7.0×50
010013660550	TC4	Φ7.0×55

P/N	Material	Spec.
010013660050	TC4	L=50
010013660051	TC4	L=60
010013660052	TC4	L=70
010013660053	TC4	L=80
010013660054	TC4	L=90
010013660055	TC4	L=100
010013660056	TC4	L=110
010013660057	TC4	L=120
010013660058	TC4	L=130
010013660059	TC4	L=140
010013660060	TC4	L=150
010013660061	TC4	L=160

P/N	Material	Spec.
010013660062	TC4	L=180
010013660063	TC4	L=200
010013660064	TC4	L=210
010013660065	TC4	L=220
010013660066	TC4	L=240
010013660067	TC4	L=250
010013660068	TC4	L=260
010013660069	TC4	L=300
010013660070	TC4	L=350
010013660071	TC4	L=400
010013660072	TC4	L=500



ZION™ Crosslink, I



P/N	Material	Spec.
03001 3600073	TC4	L=50
03001 3600074	TC4	L=55
03001 3600075	TC4	L=60

ZION™ Crosslink, II



P/N	Material	Spec.
03001 3600015	TC4	L=50 + 65

ZION™ Pedicle Hook



P/N	Material	Spec.
03001 3600076	TC4	L
03001 3600077	TC4	M
03001 3600078	TC4	S

ZION™ Lumbar Laminar Hook




P/N	Material	Spec.
03001 3600079	TC4	L
03001 3600080	TC4	M
03001 3600081	TC4	S

ZION™ Thoracic Laminar Hook




P/N	Material	Spec.
03001 3600082	TC4	L
03001 3600083	TC4	M
03001 3600084	TC4	S

ZION™ Left Offset Hook



P/N	Material	Spec.
03001 3600085	TC4	L

ZION™ Right Offset Hook




P/N	Material	Spec.
03001 3600086	TC4	R

ZION™ Set Screw



P/N	Material	Spec.
03001 3600087	TC4	Ø8×1

ATLAS™ 5.5 Posterior Thoracolumbar Internal Fixation System



Indications

- Disc degenerative disease, Spondylolisthesis, Spondylosis, Spinal fracture, Spinal tumor, Spinal deformity

Features & Benefits

- Low profile screw head makes the base well-stacked to reduce potential of irritation
- StarDrive recess provides an improved torque transmission to the set screw
- Extended tabs with notches on both sides, facilitate the tab removal and provides break-off surface
- Polyaxial screw motion range: ±22°, to facilitate the placement of pedicle screws and rods.
- Dual lead thread design for optimal bone purchase
- The screw tip facilitates the insertion into the pedicle and realizes the bicortical fixation

ATLAS™ Polyaxial Pedicle Screw



P/N	Material	Spec.
03001 3037030	TC4	Φ4.5×30
03001 3037035	TC4	Φ4.5×35
03001 3037040	TC4	Φ4.5×40
03001 3037045	TC4	Φ4.5×45
03001 3038030	TC4	Φ5.0×30
03001 3038035	TC4	Φ5.0×35
03001 3038040	TC4	Φ5.0×40
03001 3038045	TC4	Φ5.0×45
03001 3038050	TC4	Φ5.0×50
03001 3038055	TC4	Φ5.0×55
03001 3039030	TC4	Φ5.5×30
03001 3039035	TC4	Φ5.5×35
03001 3039040	TC4	Φ5.5×40
03001 3039045	TC4	Φ5.5×45
03001 3039050	TC4	Φ5.5×50

P/N	Material	Spec.
03001 3039055	TC4	Φ5.5×55
03001 3041035	TC4	Φ6.0×35
03001 3041040	TC4	Φ6.0×40
03001 3041045	TC4	Φ6.0×45
03001 3041050	TC4	Φ6.0×50
03001 3041055	TC4	Φ6.0×55
03001 3042035	TC4	Φ6.5×35
03001 3042040	TC4	Φ6.5×40
03001 3042045	TC4	Φ6.5×45
03001 3042050	TC4	Φ6.5×50
03001 3042055	TC4	Φ6.5×55
03001 3043035	TC4	Φ7.0×35
03001 3043040	TC4	Φ7.0×40
03001 3043045	TC4	Φ7.0×45
03001 3043050	TC4	Φ7.0×50

ATLAS™ Reduction Polyaxial Pedicle Screw



P/N	Material	Spec.
03001 3450030	TC4	Φ5.0×30
03001 3450035	TC4	Φ5.0×35
03001 3450040	TC4	Φ5.0×40
03001 3450050	TC4	Φ5.0×45
03001 3450055	TC4	Φ5.0×55
03001 3450060	TC4	Φ5.0×60
03001 3450070	TC4	Φ5.0×70
03001 3450080	TC4	Φ5.0×80
03001 3450090	TC4	Φ5.0×90
03001 3450100	TC4	Φ5.0×100

P/N	Material	Spec.
03001 3460045	TC4	Φ6.0×45
03001 3460050	TC4	Φ6.0×50
03001 3460055	TC4	Φ6.0×55
03001 3460060	TC4	Φ6.0×60
03001 3460065	TC4	Φ6.0×65
03001 3460070	TC4	Φ6.0×70
03001 3460075	TC4	Φ6.0×75
03001 3460080	TC4	Φ6.0×80
03001 3460085	TC4	Φ6.0×85
03001 3460090	TC4	Φ6.0×90



ETNA™ 5.5 Posterior Thoracolumbar Internal Fixation System

05

ETNA™ 5.5 Posterior Thoracolumbar Internal Fixation System

06

ATLAS™ Monoaxial Pedicle Screw



P/N	Material	Spec.
030013044030	TC4	Φ5.0x30
030013044035	TC4	Φ5.0x35
030013044040	TC4	Φ5.0x40
030013044045	TC4	Φ5.0x45
030013045035	TC4	Φ5.5x35
030013045040	TC4	Φ5.5x40
030013045045	TC4	Φ5.5x45
030013045050	TC4	Φ5.5x50
030013046035	TC4	Φ6.0x35
030013046040	TC4	Φ6.0x40

P/N	Material	Spec.
030013046045	TC4	Φ6.0x45
030013046050	TC4	Φ6.0x50
030013047035	TC4	Φ6.5x35
030013047040	TC4	Φ6.5x40
030013047045	TC4	Φ6.5x45
030013047050	TC4	Φ6.5x50
030013048035	TC4	Φ7.0x35
030013048040	TC4	Φ7.0x40
030013048045	TC4	Φ7.0x45
030013048050	TC4	Φ7.0x50

ATLAS™ Reduction Monoaxial Pedicle



P/N	Material	Spec.
030013250030	TC4	Φ5.0x30
030013250035	TC4	Φ5.0x35
030013250040	TC4	Φ5.0x40
030013250045	TC4	Φ5.0x45
030013255035	TC4	Φ5.5x35
030013255040	TC4	Φ5.5x40
030013255045	TC4	Φ5.5x45
030013255050	TC4	Φ5.5x50
030013260035	TC4	Φ6.0x35
030013260040	TC4	Φ6.0x40

P/N	Material	Spec.
030013260045	TC4	Φ6.0x45
030013260050	TC4	Φ6.0x50
030013265035	TC4	Φ6.5x35
030013265040	TC4	Φ6.5x40
030013265045	TC4	Φ6.5x45
030013265050	TC4	Φ6.5x50
030013270035	TC4	Φ7.0x35
030013270040	TC4	Φ7.0x40
030013270045	TC4	Φ7.0x45
030013270050	TC4	Φ7.0x50

ATLAS™ 5.5mm Rod



P/N	Material	Spec.
030013055050	TC4	L=50
030013055060	TC4	L=60
030013055070	TC4	L=70
030013055080	TC4	L=80
030013055090	TC4	L=90
030013055100	TC4	L=100
030013055110	TC4	L=110
030013055120	TC4	L=120
030013055130	TC4	L=130
030013055140	TC4	L=140
030013055150	TC4	L=150
030013055160	TC4	L=160

P/N	Material	Spec.
030013055180	TC4	L=180
030013055200	TC4	L=200
030013055210	TC4	L=210
030013055220	TC4	L=220
030013055240	TC4	L=240
030013055250	TC4	L=250
030013055260	TC4	L=260
030013055300	TC4	L=300
030013055350	TC4	L=350
030013055400	TC4	L=400
030013055500	TC4	L=500

ATLAS™ 5.5 Crosslink II



ATLAS™ 5.5 Set Screw



P/N	Material	Spec.
030013000060	TC4	L=60
030013000070	TC4	L=70
030013000080	TC4	L=80

P/N	Material	Spec.
030013000001	TC4	/

ETNA™ 5.5 Posterior Thoracolumbar Internal Fixation System



Indications

- Intervertebral disc degenerative disease, Fracture or dislocation, Spinal stenosis, Scoliosis, Kyphosis and lordosis, Spinal tumor, Fetal joint formation and previous fusion failure, Degenerative scoliosis, Lumbar degeneration instability, Intervertebral disc herniation or spine stenosis, Discogenic back pain, Lumbar spondylolisthesis requires simultaneous decompression and reduction of spinal canal

Features & Benefits

- The bigger diameter with upper shaft segment of screws, the lower broken rate of screws.
- Opened connector, easy to connect beak rod reducer to pedicle screw
- Minimized design of pedicle screw head, minimum tissue irritation
- Slot design at the top of shaft
- Easy to connect beak rod reducer to head of screws, fixed securely. The reduce of joint axis, Reducing the loss of clench force, Adding the construct of stopping loosening.
- Low Buttons, the proximal is corresponded to pedicle screws, The distal is corresponded to reduction screws

ETNA™ 5.5 Polyaxial Pedicle Screw



P/N	Material	Spec.
030013049025	TC4	Φ4.0x25
030013049030	TC4	Φ4.0x30
030013049035	TC4	Φ4.0x35
030013049040	TC4	Φ4.0x40
030013049045	TC4	Φ4.0x45
030013049050	TC4	Φ4.0x50
030013049055	TC4	Φ4.0x55
030013049060	TC4	Φ4.0x60
030013051025	TC4	Φ4.5x25
030013051030	TC4	Φ4.5x30
030013051035	TC4	Φ4.5x35
030013051040	TC4	Φ4.5x40
030013051045	TC4	Φ4.5x45
030013051050	TC4	Φ4.5x50
030013051055	TC4	Φ4.5x55
030013051060	TC4	Φ4.5x60
030013052025	TC4	Φ5.0x25
030013052030	TC4	Φ5.0x30
030013052035	TC4	Φ5.0x35
030013052040	TC4	Φ5.0x40
030013052045	TC4	Φ5.0x45
030013052050	TC4	Φ5.0x50
030013052055	TC4	Φ5.0x55
030013052060	TC4	Φ5.0x60
030013057025	TC4	Φ7.0x25
030013057030	TC4	Φ7.0x30
030013057035	TC4	Φ7.0x35
030013057040	TC4	Φ7.0x40
030013057045	TC4	Φ7.0x45
030013057050	TC4	Φ7.0x50
030013057055	TC4	Φ7.0x55
030013057060	TC4	Φ7.0x60

P/N	Material	Spec.
030013053045	TC4	Φ5.5x45
030013053050	TC4	Φ5.5x50
030013053055	TC4	Φ5.5x55
030013053060	TC4	Φ5.5x60
030013054025	TC4	Φ6.0x25
030013054030	TC4	Φ6.0x30
030013054035	TC4	Φ6.0x35
030013054040	TC4	Φ6.0x40
030013054045	TC4	Φ6.0x45
030013054050	TC4	Φ6.0x50
030013054055	TC4	Φ6.0x55
030013054060	TC4	Φ6.0x60
030013056025	TC4	Φ6.5x25
030013056030	TC4	Φ6.5x30
030013056035	TC4	Φ6.5x35
030013056040	TC4	Φ6.5x40
030013056045	TC4	Φ6.5x45
030013056050	TC4	Φ6.5x50
030013056055	TC4	Φ6.5x55
030013056060	TC4	Φ6.5x60
030013057025	TC4	Φ7.0x25
030013057030	TC4	Φ7.0x30
030013057035	TC4	Φ7.0x35
030013057040	TC4	Φ7.0x40
030013057045	TC4	Φ7.0x45
030013057050	TC4	Φ7.0x50
030013057055	TC4	Φ7.0x55
030013057060	TC4	Φ7.0x60



ETNA™ 5.5S Pedicle Threaded
Internal Fixation System

07

ETNA™ 5.5 Long Arm Polyaxial Pedicle Screw



P/N	Material	Spec.
000013058025	TC4	Ø4.0x25
000013058030	TC4	Ø4.0x30
000013058035	TC4	Ø4.0x35
000013058040	TC4	Ø4.0x40
000013058045	TC4	Ø4.0x45
000013058050	TC4	Ø4.0x50
000013058055	TC4	Ø4.0x55
000013058060	TC4	Ø4.0x60
000013059025	TC4	Ø4.5x25
000013059030	TC4	Ø4.5x30
000013059035	TC4	Ø4.5x35
000013059040	TC4	Ø4.5x40
000013059045	TC4	Ø4.5x45
000013059050	TC4	Ø4.5x50
000013059055	TC4	Ø4.5x55
000013059060	TC4	Ø4.5x60
000013061025	TC4	Ø5.0x25
000013061030	TC4	Ø5.0x30
000013061035	TC4	Ø5.0x35
000013061040	TC4	Ø5.0x40
000013061045	TC4	Ø5.0x45
000013061050	TC4	Ø5.0x50
000013061055	TC4	Ø5.0x55
000013061060	TC4	Ø5.0x60
000013062025	TC4	Ø5.5x25
000013062030	TC4	Ø5.5x30
000013062035	TC4	Ø5.5x35
000013062040	TC4	Ø5.5x40

P/N	Material	Spec.
000013062045	TC4	Ø5.5x45
000013062050	TC4	Ø5.5x50
000013062055	TC4	Ø5.5x55
000013062060	TC4	Ø5.5x60
000013063025	TC4	Ø6.0x25
000013063030	TC4	Ø6.0x30
000013063035	TC4	Ø6.0x35
000013063040	TC4	Ø6.0x40
000013063045	TC4	Ø6.0x45
000013063050	TC4	Ø6.0x50
000013063055	TC4	Ø6.0x55
000013063060	TC4	Ø6.0x60
000013064025	TC4	Ø6.5x25
000013064030	TC4	Ø6.5x30
000013064035	TC4	Ø6.5x35
000013064040	TC4	Ø6.5x40
000013064045	TC4	Ø6.5x45
000013064050	TC4	Ø6.5x50
000013064055	TC4	Ø6.5x55
000013064060	TC4	Ø6.5x60
000013066025	TC4	Ø7.0x25
000013066030	TC4	Ø7.0x30
000013066035	TC4	Ø7.0x35
000013066040	TC4	Ø7.0x40
000013066045	TC4	Ø7.0x45
000013066050	TC4	Ø7.0x50
000013066055	TC4	Ø7.0x55
000013066060	TC4	Ø7.0x60

ETNA™ 5.5 Long Arm Monoaxial Pedicle Screw




P/N	Material	Spec.
000013075025	TC4	Ø4.0x25
000013075030	TC4	Ø4.0x30
000013075035	TC4	Ø4.0x35
000013075040	TC4	Ø4.0x40
000013075045	TC4	Ø4.0x45
000013075050	TC4	Ø4.0x50
000013075055	TC4	Ø4.0x55
000013075060	TC4	Ø4.0x60
000013076025	TC4	Ø4.5x25
000013076030	TC4	Ø4.5x30
000013076035	TC4	Ø4.5x35
000013076040	TC4	Ø4.5x40
000013076045	TC4	Ø4.5x45
000013076050	TC4	Ø4.5x50
000013076055	TC4	Ø4.5x55
000013076060	TC4	Ø4.5x60
000013077025	TC4	Ø5.0x25
000013077030	TC4	Ø5.0x30
000013077035	TC4	Ø5.0x35
000013077040	TC4	Ø5.0x40
000013077045	TC4	Ø5.0x45
000013077050	TC4	Ø5.0x50
000013077055	TC4	Ø5.0x55
000013077060	TC4	Ø5.0x60
000013078025	TC4	Ø5.5x25
000013078030	TC4	Ø5.5x30
000013078035	TC4	Ø5.5x35
000013078040	TC4	Ø5.5x40

P/N	Material	Spec.
000013078045	TC4	Ø5.5x45
000013078050	TC4	Ø5.5x50
000013078055	TC4	Ø5.5x55
000013078060	TC4	Ø5.5x60
000013079025	TC4	Ø6.0x25
000013079030	TC4	Ø6.0x30
000013079035	TC4	Ø6.0x35
000013079040	TC4	Ø6.0x40
000013079045	TC4	Ø6.0x45
000013079050	TC4	Ø6.0x50
000013079055	TC4	Ø6.0x55
000013079060	TC4	Ø6.0x60
000013080025	TC4	Ø6.5x25
000013080030	TC4	Ø6.5x30
000013080035	TC4	Ø6.5x35
000013080040	TC4	Ø6.5x40
000013080045	TC4	Ø6.5x45
000013080050	TC4	Ø6.5x50
000013080055	TC4	Ø6.5x55
000013080060	TC4	Ø6.5x60
000013081025	TC4	Ø7.0x25
000013081030	TC4	Ø7.0x30
000013081035	TC4	Ø7.0x35
000013081040	TC4	Ø7.0x40
000013081045	TC4	Ø7.0x45
000013081050	TC4	Ø7.0x50
000013081055	TC4	Ø7.0x55
000013081060	TC4	Ø7.0x60

ETNA™ 5.5S Pedicle Threaded
Internal Fixation System

08

ETNA™ 5.5 Monoaxial Pedicle Screw



P/N	Material	Spec.
000013067025	TC4	Ø4.0x25
000013067030	TC4	Ø4.0x30
000013067035	TC4	Ø4.0x35
000013067040	TC4	Ø4.0x40
000013067045	TC4	Ø4.0x45
000013067050	TC4	Ø4.0x50
000013067055	TC4	Ø4.0x55
000013067060	TC4	Ø4.0x60
000013068025	TC4	Ø4.5x25
000013068030	TC4	Ø4.5x30
000013068035	TC4	Ø4.5x35
000013068040	TC4	Ø4.5x40
000013068045	TC4	Ø4.5x45
000013068050	TC4	Ø4.5x50
000013068055	TC4	Ø4.5x55
000013068060	TC4	Ø4.5x60
000013069025	TC4	Ø5.0x25
000013069030	TC4	Ø5.0x30
000013069035	TC4	Ø5.0x35
000013069040	TC4	Ø5.0x40
000013069045	TC4	Ø5.0x45
000013069050	TC4	Ø5.0x50
000013069055	TC4	Ø5.0x55
000013069060	TC4	Ø5.0x60
000013071025	TC4	Ø5.5x25
000013071030	TC4	Ø5.5x30
000013071035	TC4	Ø5.5x35
000013071040	TC4	Ø5.5x40

P/N	Material	Spec.
000013071045	TC4	Ø5.5x45
000013071050	TC4	Ø5.5x50
000013071055	TC4	Ø5.5x55
000013071060	TC4	Ø5.5x60
000013072025	TC4	Ø6.0x25
000013072030	TC4	Ø6.0x30
000013072035	TC4	Ø6.0x35
000013072040	TC4	Ø6.0x40
000013072045	TC4	Ø6.0x45
000013072050	TC4	Ø6.0x50
000013072055	TC4	Ø6.0x55
000013072060	TC4	Ø6.0x60
000013073025	TC4	Ø6.5x25
000013073030	TC4	Ø6.5x30
000013073035	TC4	Ø6.5x35
000013073040	TC4	Ø6.5x40
000013073045	TC4	Ø6.5x45
000013073050	TC4	Ø6.5x50
000013073055	TC4	Ø6.5x55
000013073060	TC4	Ø6.5x60
000013074025	TC4	Ø7.0x25
000013074030	TC4	Ø7.0x30
000013074035	TC4	Ø7.0x35
000013074040	TC4	Ø7.0x40
000013074045	TC4	Ø7.0x45
000013074050	TC4	Ø7.0x50
000013074055	TC4	Ø7.0x55
000013074060	TC4	Ø7.0x60

ETNA™ 5.5 Crosslink, III



ETNA™ 5.5 Pedicle Hook



ETNA™ 5.5 Extended Body Hook



P/N	Material	Spec.
000013100050	TC4	L=50
000013100055	TC4	L=55
000013100060	TC4	L=60
000013100065	TC4	L=65
000013100070	TC4	L=70
000013100075	TC4	L=75
000013100080	TC4	L=80

P/N	Material	Spec.
000013100085	TC4	L=85
000013100090	TC4	L=90
000013100095	TC4	L=95
000013100100	TC4	L=100
000013100105	TC4	L=105
000013100110	TC4	L=110
000013100115	TC4	L=115

P/N	Material	Spec.
000013100004	TC4	Large
000013100005	TC4	Medium
000013100006	TC4	Small

P/N	Material	Spec.
000013100007	TC4	Large
000013100008	TC4	Medium
000013100009	TC4	Small



ETNA™ 5.5 Ramped Blade Hook



P/N	Material	Spec.
030013100010	TC4	Wide
030013100011	TC4	Narrow

ETNA™ 5.5 Lumbar Laminar Hook



P/N	Material	Spec.
030013100012	TC4	Supra, Medium
030013100013	TC4	Infra, Large
030013100014	TC4	Infra, Small

ETNA™ 5.5 Thoracic Angled Hook



P/N	Material	Spec.
030013100015	TC4	Left
030013100016	TC4	Right

ETNA™ 5.5 Offset Hook



P/N	Material	Spec.
030013100017	TC4	Left
030013100018	TC4	Right

ETNA™ 5.5 Domino Connector



P/N	Material	Spec.
030013100019	TC4	Straight
030013100020	TC4	Square

ETNA™ 5.5mm Rod



P/N	Material	Spec.	P/N	Material	Spec.
030013600050	TC4	L=50	030013600062	TC4	L=180
030013600051	TC4	L=60	030013600063	TC4	L=200
030013600052	TC4	L=70	030013600064	TC4	L=210
030013600053	TC4	L=80	030013600065	TC4	L=220
030013600054	TC4	L=90	030013600066	TC4	L=240
030013600055	TC4	L=100	030013600067	TC4	L=250
030013600056	TC4	L=110	030013600068	TC4	L=260
030013600057	TC4	L=120	030013600069	TC4	L=300
030013600058	TC4	L=130	030013600070	TC4	L=350
030013600059	TC4	L=140	030013600071	TC4	L=400
030013600060	TC4	L=150	030013600072	TC4	L=500
030013600061	TC4	L=160			

ETNA™ 5.5 Set Screw



P/N	Material	Spec.
030013100021	TC4	Bx1

ETNA™ 6.0 Posterior Thoracolumbar Internal Fixation System

ETNA™ 6.0 Polyaxial Pedicle Screw



P/N	Material	Spec.	P/N	Material	Spec.
030013201025	TC4	Φ4.0×25	030013204045	TC4	Φ5.5×45
030013201030	TC4	Φ4.0×30	030013204050	TC4	Φ5.5×50
030013201035	TC4	Φ4.0×35	030013204055	TC4	Φ5.5×55
030013201040	TC4	Φ4.0×40	030013204060	TC4	Φ5.5×60
030013201045	TC4	Φ4.0×45	030013205025	TC4	Φ6.0×25
030013201050	TC4	Φ4.0×50	030013205030	TC4	Φ6.0×30
030013201055	TC4	Φ4.0×55	030013205035	TC4	Φ6.0×35
030013201060	TC4	Φ4.0×60	030013205040	TC4	Φ6.0×40
030013202025	TC4	Φ4.5×25	030013205045	TC4	Φ6.0×45
030013202030	TC4	Φ4.5×30	030013205050	TC4	Φ6.0×50
030013202035	TC4	Φ4.5×35	030013205055	TC4	Φ6.0×55
030013202040	TC4	Φ4.5×40	030013205060	TC4	Φ6.0×60
030013202045	TC4	Φ4.5×45	030013206025	TC4	Φ6.5×25
030013202050	TC4	Φ4.5×50	030013206030	TC4	Φ6.5×30
030013202055	TC4	Φ4.5×55	030013206035	TC4	Φ6.5×35
030013202060	TC4	Φ4.5×60	030013206040	TC4	Φ6.5×40
030013203025	TC4	Φ5.0×25	030013206045	TC4	Φ6.5×45
030013203030	TC4	Φ5.0×30	030013206050	TC4	Φ6.5×50
030013203035	TC4	Φ5.0×35	030013206055	TC4	Φ6.5×55
030013203040	TC4	Φ5.0×40	030013206060	TC4	Φ6.5×60
030013203045	TC4	Φ5.0×45	030013207025	TC4	Φ7.0×25
030013203050	TC4	Φ5.0×50	030013207030	TC4	Φ7.0×30
030013203055	TC4	Φ5.0×55	030013207035	TC4	Φ7.0×35
030013203060	TC4	Φ5.0×60	030013207040	TC4	Φ7.0×40
030013204025	TC4	Φ5.5×25	030013207045	TC4	Φ7.0×45
030013204030	TC4	Φ5.5×30	030013207050	TC4	Φ7.0×50
030013204035	TC4	Φ5.5×35	030013207055	TC4	Φ7.0×55
030013204040	TC4	Φ5.5×40	030013207060	TC4	Φ7.0×60


ETNA™ 6.0 Long Arm Polyaxial Pedicle Screw



P/N	Material	Spec.	P/N	Material	Spec.
030013208025	TC4	Φ4.0×25	030013211045	TC4	Φ5.5×45
030013208030	TC4	Φ4.0×30	030013211050	TC4	Φ5.5×50
030013208035	TC4	Φ4.0×35	030013211055	TC4	Φ5.5×55
030013208040	TC4	Φ4.0×40	030013211060	TC4	Φ5.5×60
030013208045	TC4	Φ4.0×45	030013212025	TC4	Φ6.0×25
030013208050	TC4	Φ4.0×50	030013212030	TC4	Φ6.0×30
030013208055	TC4	Φ4.0×55	030013212035	TC4	Φ6.0×35
030013208060	TC4	Φ4.0×60	030013212040	TC4	Φ6.0×40
030013209025	TC4	Φ4.5×25	030013212045	TC4	Φ6.0×45
030013209030	TC4	Φ4.5×30	030013212050	TC4	Φ6.0×50
030013209035	TC4	Φ4.5×35	030013212055	TC4	Φ6.0×55
030013209040	TC4	Φ4.5×40	030013212060	TC4	Φ6.0×60
030013209045	TC4	Φ4.5×45	030013213025	TC4	Φ6.5×25
030013209050	TC4	Φ4.5×50	030013213030	TC4	Φ6.5×30
030013209055	TC4	Φ4.5×55	030013213035	TC4	Φ6.5×35
030013209060	TC4	Φ4.5×60	030013213040	TC4	Φ6.5×40
030013210025	TC4	Φ5.0×25	030013213045	TC4	Φ6.5×45
030013210030	TC4	Φ5.0×30	030013213050	TC4	Φ6.5×50
030013210035	TC4	Φ5.0×35	030013213055	TC4	Φ6.5×55
030013210040	TC4	Φ5.0×40	030013213060	TC4	Φ6.5×60
030013210045	TC4	Φ5.0×45	030013214025	TC4	Φ7.0×25
030013210050	TC4	Φ5.0×50	030013214030	TC4	Φ7.0×30
030013210055	TC4	Φ5.0×55	030013214035	TC4	Φ7.0×35
030013210060	TC4	Φ5.0×60	030013214040	TC4	Φ7.0×40
030013211025	TC4	Φ5.5×25	030013214045	TC4	Φ7.0×45
030013211030	TC4	Φ5.5×30	030013214050	TC4	Φ7.0×50
030013211035	TC4	Φ5.5×35	030013214055	TC4	Φ7.0×55
030013211040	TC4	Φ5.5×40	030013214060	TC4	Φ7.0×60



ETNA™ 6.0 Monoaxial Pedicle Screw



P/N	Material	Spec.
000013215025	TC4	Ø4.0x25
000013215030	TC4	Ø4.0x30
000013215035	TC4	Ø4.0x35
000013215040	TC4	Ø4.0x40
000013215045	TC4	Ø4.0x45
000013215050	TC4	Ø4.0x50
000013215055	TC4	Ø4.0x55
000013215060	TC4	Ø4.0x60
000013216025	TC4	Ø4.5x25
000013216030	TC4	Ø4.5x30
000013216035	TC4	Ø4.5x35
000013216040	TC4	Ø4.5x40
000013216045	TC4	Ø4.5x45
000013216050	TC4	Ø4.5x50
000013216055	TC4	Ø4.5x55
000013216060	TC4	Ø4.5x60
000013217025	TC4	Ø5.0x25
000013217030	TC4	Ø5.0x30
000013217035	TC4	Ø5.0x35
000013217040	TC4	Ø5.0x40
000013217045	TC4	Ø5.0x45
000013217050	TC4	Ø5.0x50
000013217055	TC4	Ø5.0x55
000013217060	TC4	Ø5.0x60
000013218025	TC4	Ø5.5x25
000013218030	TC4	Ø5.5x30
000013218035	TC4	Ø5.5x35
000013218040	TC4	Ø5.5x40

P/N	Material	Spec.
000013218045	TC4	Ø5.5x45
000013218050	TC4	Ø5.5x50
000013218055	TC4	Ø5.5x55
000013218060	TC4	Ø5.5x60
000013219025	TC4	Ø6.0x25
000013219030	TC4	Ø6.0x30
000013219035	TC4	Ø6.0x35
000013219040	TC4	Ø6.0x40
000013219045	TC4	Ø6.0x45
000013219050	TC4	Ø6.0x50
000013219055	TC4	Ø6.0x55
000013219060	TC4	Ø6.0x60
000013220025	TC4	Ø6.5x25
000013220030	TC4	Ø6.5x30
000013220035	TC4	Ø6.5x35
000013220040	TC4	Ø6.5x40
000013220045	TC4	Ø6.5x45
000013220050	TC4	Ø6.5x50
000013220055	TC4	Ø6.5x55
000013220060	TC4	Ø6.5x60
000013221025	TC4	Ø7.0x25
000013221030	TC4	Ø7.0x30
000013221035	TC4	Ø7.0x35
000013221040	TC4	Ø7.0x40
000013221045	TC4	Ø7.0x45
000013221050	TC4	Ø7.0x50
000013221055	TC4	Ø7.0x55
000013221060	TC4	Ø7.0x60

ETNA™ 6.0 Long Arm Monoaxial Pedicle Screw



P/N	Material	Spec.
000013222025	TC4	Ø4.0x25
000013222030	TC4	Ø4.0x30
000013222035	TC4	Ø4.0x35
000013222040	TC4	Ø4.0x40
000013222045	TC4	Ø4.0x45
000013222050	TC4	Ø4.0x50
000013222055	TC4	Ø4.0x55
000013222060	TC4	Ø4.0x60
000013223025	TC4	Ø4.5x25
000013223030	TC4	Ø4.5x30
000013223035	TC4	Ø4.5x35
000013223040	TC4	Ø4.5x40
000013223045	TC4	Ø4.5x45
000013223050	TC4	Ø4.5x50
000013223055	TC4	Ø4.5x55
000013223060	TC4	Ø4.5x60
000013224025	TC4	Ø5.0x25
000013224030	TC4	Ø5.0x30
000013224035	TC4	Ø5.0x35
000013224040	TC4	Ø5.0x40
000013224045	TC4	Ø5.0x45
000013224050	TC4	Ø5.0x50
000013224055	TC4	Ø5.0x55
000013224060	TC4	Ø5.0x60
000013225025	TC4	Ø5.5x25
000013225030	TC4	Ø5.5x30
000013225035	TC4	Ø5.5x35
000013225040	TC4	Ø5.5x40

P/N	Material	Spec.
000013225045	TC4	Ø5.5x45
000013225050	TC4	Ø5.5x50
000013225055	TC4	Ø5.5x55
000013225060	TC4	Ø5.5x60
000013226025	TC4	Ø6.0x25
000013226030	TC4	Ø6.0x30
000013226035	TC4	Ø6.0x35
000013226040	TC4	Ø6.0x40
000013226045	TC4	Ø6.0x45
000013226050	TC4	Ø6.0x50
000013226055	TC4	Ø6.0x55
000013226060	TC4	Ø6.0x60
000013227025	TC4	Ø6.5x25
000013227030	TC4	Ø6.5x30
000013227035	TC4	Ø6.5x35
000013227040	TC4	Ø6.5x40
000013227045	TC4	Ø6.5x45
000013227050	TC4	Ø6.5x50
000013227055	TC4	Ø6.5x55
000013227060	TC4	Ø6.5x60
000013228025	TC4	Ø7.0x25
000013228030	TC4	Ø7.0x30
000013228035	TC4	Ø7.0x35
000013228040	TC4	Ø7.0x40
000013228045	TC4	Ø7.0x45
000013228050	TC4	Ø7.0x50
000013228055	TC4	Ø7.0x55
000013228060	TC4	Ø7.0x60

ETNA™ 6.0 Crosslink, III



ETNA™ 6.0 Pedicle Hook



ETNA™ 6.0 Extended Body Hook



ETNA™ 6.0 Ramped Blade Hook



ETNA™ 6.0 Lumbar Laminar Hook I



P/N	Material	Spec.
000013200005	TC4	L=50
000013200006	TC4	L=55
000013200007	TC4	L=60
000013200008	TC4	L=65
000013200009	TC4	L=70
000013200010	TC4	L=75
000013200011	TC4	L=80

P/N	Material	Spec.
000013200012	TC4	L=85
000013200013	TC4	L=90
000013200014	TC4	L=95
000013200015	TC4	L=100
000013200016	TC4	L=105
000013200017	TC4	L=110
000013200018	TC4	L=115

P/N	Material	Spec.
000013200019	TC4	Large
000013200020	TC4	Medium
000013200021	TC4	Small

P/N	Material	Spec.
000013200022	TC4	Large
000013200023	TC4	Medium
000013200024	TC4	Small

P/N	Material	Spec.
000013200025	TC4	Wide
000013200026	TC4	Narrow

P/N	Material	Spec.
000013200027	TC4	Supra, Medium
000013200028	TC4	Infra, Large
000013200029	TC4	Infra, Small



ETNA™ 6.0 Thoracic Angled Hook



ETNA™ 6.0 Offset Hook



ETNA™ 6.0 Domino Connector



ETNA™ 6.0mm Rod



ETNA™ 6.0 Set Screw



P/N	Material	Spec.
030013200010	TC4	Left
030013200011	TC4	Right

P/N	Material	Spec.
030013200012	TC4	Left
030013200013	TC4	Right

P/N	Material	Spec.
030013200014	TC4	Straight
030013200015	TC4	Square


P/N	Material	Spec.	P/N	Material	Spec.
030013200016	TC4	L=50	030013200048	TC4	L=180
030013200017	TC4	L=60	030013200049	TC4	L=200
030013200018	TC4	L=70	030013200050	TC4	L=210
030013200019	TC4	L=80	030013200051	TC4	L=220
030013200040	TC4	L=90	030013200052	TC4	L=240
030013200041	TC4	L=100	030013200053	TC4	L=250
030013200042	TC4	L=110	030013200054	TC4	L=260
030013200043	TC4	L=120	030013200055	TC4	L=300
030013200044	TC4	L=130	030013200056	TC4	L=350
030013200045	TC4	L=140	030013200057	TC4	L=400
030013200046	TC4	L=150	030013200058	TC4	L=500
030013200047	TC4	L=160			

P/N	Material	Spec.
030013200059	TC4	Bx1


EVEREST™ 5.5 MIS Spine System



EVEREST™ 5.5 MIS Monoaxial Pedicle Screw



EVEREST™ 5.5 MIS Polyaxial Pedicle Screw



Indications

- Disc degenerative disease, Spondylolisthesis, Stenosis, Spinal fracture, Spinal tumor, Spinal deformity

Features & Benefits

- Percutaneous pedicle screw fixation, Small incision, less stripping of paraspinal muscles, minimal blood loss and rapid postoperative recovery
- Apply for multi-level constructs
- Anatomically precontoured connecting rod, also allows additional rod bending
- Simplified and unique rod placement techniques allow arbitrarily curved rod insertion
- Double-load threaded and alternatively cement injectable cannulated
- Pedicle screws ease screw insertion while adding stability, which is especially important for patients with osteoporosis

P/N	Material	Spec.	P/N	Material	Spec.
030013750025	TC4	Φ5.0x25	030013760045	TC4	Φ6.0x45
030013750030	TC4	Φ5.0x30	030013760050	TC4	Φ6.0x50
030013750035	TC4	Φ5.0x35	030013760055	TC4	Φ6.0x55
030013750040	TC4	Φ5.0x40	030013760060	TC4	Φ6.0x60
030013750045	TC4	Φ5.0x45	030013765025	TC4	Φ6.5x25
030013750050	TC4	Φ5.0x50	030013765030	TC4	Φ6.5x30
030013750055	TC4	Φ5.0x55	030013765035	TC4	Φ6.5x35
030013750060	TC4	Φ5.0x60	030013765040	TC4	Φ6.5x40
030013755025	TC4	Φ5.5x25	030013765045	TC4	Φ6.5x45
030013755030	TC4	Φ5.5x30	030013765050	TC4	Φ6.5x50
030013755035	TC4	Φ5.5x35	030013765055	TC4	Φ6.5x55
030013755040	TC4	Φ5.5x40	030013765060	TC4	Φ6.5x60
030013755045	TC4	Φ5.5x45	030013770025	TC4	Φ7.0x25
030013755050	TC4	Φ5.5x50	030013770030	TC4	Φ7.0x30
030013755055	TC4	Φ5.5x55	030013770035	TC4	Φ7.0x35
030013755060	TC4	Φ5.5x60	030013770040	TC4	Φ7.0x40
030013760025	TC4	Φ6.0x25	030013770045	TC4	Φ7.0x45
030013760030	TC4	Φ6.0x30	030013770050	TC4	Φ7.0x50
030013760035	TC4	Φ6.0x35	030013770055	TC4	Φ7.0x55
030013760040	TC4	Φ6.0x40	030013770060	TC4	Φ7.0x60

P/N	Material	Spec.	P/N	Material	Spec.
030013850025	TC4	Φ5.0x25	030013860045	TC4	Φ6.0x45
030013850030	TC4	Φ5.0x30	030013860050	TC4	Φ6.0x50
030013850035	TC4	Φ5.0x35	030013860055	TC4	Φ6.0x55
030013850040	TC4	Φ5.0x40	030013860060	TC4	Φ6.0x60
030013850045	TC4	Φ5.0x45	030013865025	TC4	Φ6.5x25
030013850050	TC4	Φ5.0x50	030013865030	TC4	Φ6.5x30
030013850055	TC4	Φ5.0x55	030013865035	TC4	Φ6.5x35
030013850060	TC4	Φ5.0x60	030013865040	TC4	Φ6.5x40
030013855025	TC4	Φ5.5x25	030013865045	TC4	Φ6.5x45
030013855030	TC4	Φ5.5x30	030013865050	TC4	Φ6.5x50
030013855035	TC4	Φ5.5x35	030013865055	TC4	Φ6.5x55
030013855040	TC4	Φ5.5x40	030013865060	TC4	Φ6.5x60
030013855045	TC4	Φ5.5x45	030013870025	TC4	Φ7.0x25
030013855050	TC4	Φ5.5x50	030013870030	TC4	Φ7.0x30
030013855055	TC4	Φ5.5x55	030013870035	TC4	Φ7.0x35
030013855060	TC4	Φ5.5x60	030013870040	TC4	Φ7.0x40
030013860025	TC4	Φ6.0x25	030013870045	TC4	Φ7.0x45
030013860030	TC4	Φ6.0x30	030013870050	TC4	Φ7.0x50
030013860035	TC4	Φ6.0x35	030013870055	TC4	Φ7.0x55
030013860040	TC4	Φ6.0x40	030013870060	TC4	Φ7.0x60



EVEREST™ 5.5mm Rod, Straight



P/N	Material	Spec.
00001395030	TC4	Ø5.5x30
00001395035	TC4	Ø5.5x35
00001395040	TC4	Ø5.5x40
00001395045	TC4	Ø5.5x45
00001395050	TC4	Ø5.5x50
00001395055	TC4	Ø5.5x55
00001395060	TC4	Ø5.5x60
00001395065	TC4	Ø5.5x65
00001395070	TC4	Ø5.5x70
00001395075	TC4	Ø5.5x75
00001395080	TC4	Ø5.5x80
00001395085	TC4	Ø5.5x85
00001395090	TC4	Ø5.5x90
00001395095	TC4	Ø5.5x95
00001395100	TC4	Ø5.5x100
00001395105	TC4	Ø5.5x105
00001395110	TC4	Ø5.5x110
00001395115	TC4	Ø5.5x115

P/N	Material	Spec.
00001395120	TC4	Ø5.5x120
00001395125	TC4	Ø5.5x125
00001395130	TC4	Ø5.5x130
00001395135	TC4	Ø5.5x135
00001395140	TC4	Ø5.5x140
00001395145	TC4	Ø5.5x145
00001395150	TC4	Ø5.5x150
00001395155	TC4	Ø5.5x155
00001395160	TC4	Ø5.5x160
00001395165	TC4	Ø5.5x165
00001395170	TC4	Ø5.5x170
00001395175	TC4	Ø5.5x175
00001395180	TC4	Ø5.5x180
00001395185	TC4	Ø5.5x185
00001395190	TC4	Ø5.5x190
00001395195	TC4	Ø5.5x195
00001395200	TC4	Ø5.5x200

EVEREST™ 5.5mm Rod, Curved



P/N	Material	Spec.
000013002025	TC4	Ø5.5x25
000013002030	TC4	Ø5.5x30
000013002035	TC4	Ø5.5x35
000013002040	TC4	Ø5.5x40
000013002045	TC4	Ø5.5x45
000013002050	TC4	Ø5.5x50
000013002055	TC4	Ø5.5x55
000013002060	TC4	Ø5.5x60
000013002065	TC4	Ø5.5x65
000013002070	TC4	Ø5.5x70

P/N	Material	Spec.
000013002075	TC4	Ø5.5x75
000013002080	TC4	Ø5.5x80
000013002085	TC4	Ø5.5x85
000013002090	TC4	Ø5.5x90
000013002100	TC4	Ø5.5x100
000013002110	TC4	Ø5.5x110
000013002120	TC4	Ø5.5x120
000013002130	TC4	Ø5.5x130
000013002140	TC4	Ø5.5x140
000013002150	TC4	Ø5.5x150

EVEREST™ 5.5 Set Screw



P/N	Material	Spec.
000013002	TC4	Ø6.6

CERVINO™ Anterior Cervical Plate



Indications

- Cervical spine fixation of 1-3 segments in C2-C7, Cervical degenerative disease, Cervical spine trauma, Cervical vertebrae deformity, Cervical spine tumor

Features & Benefits

- 1.8mm thick, low profile and high intensity
- Reversible and adjustable locking design for easy adjustment and revision
- Self-locking mechanism with one step rotation of 90° for convenient performance
- Anatomically precontoured plate and smooth surface reduce potential for irritation
- Self-tapping and self-drilling screw design help to ease surgical operation

CERVINO™ Anterior Cervical Plate



P/N	Material	Spec.
000013104022	TC4	Ø6, L=22mm
000013104024	TC4	Ø6, L=24mm
000013104026	TC4	Ø6, L=26mm
000013104028	TC4	Ø6, L=28mm
000013106033	TC4	Ø6, L=33mm
000013106036	TC4	Ø6, L=36mm
000013106038	TC4	Ø6, L=38mm
000013106041	TC4	Ø6, L=41mm
000013106043	TC4	Ø6, L=43mm
000013106046	TC4	Ø6, L=46mm
000013106051	TC4	Ø6, L=51mm
000013108053	TC4	Ø6, L=53mm
000013108056	TC4	Ø6, L=56mm
000013108059	TC4	Ø6, L=59mm
000013108061	TC4	Ø6, L=61mm
000013108064	TC4	Ø6, L=64mm
000013108069	TC4	Ø6, L=69mm
000013108071	TC4	Ø6, L=71mm
000013108075	TC4	Ø6, L=75mm

CERVINO™ Anterior Cervical Screw, Self-tapping



P/N	Material	Spec.
000013118013	TC4	L=13mm
000013118014	TC4	L=14mm
000013118015	TC4	L=15mm
000013118016	TC4	L=16mm
000013118017	TC4	L=17mm


CERVINO™ Anterior Cervical Screw, Self-drilling



P/N	Material	Spec.
000013128013	TC4	L=13mm
000013128014	TC4	L=14mm
000013128015	TC4	L=15mm
000013128016	TC4	L=16mm
000013128017	TC4	L=17mm



ALPS™ Posterior Cervical Internal Fixation System



Indications

Instability in occipitocervical segment or upper cervical spine

- Rheumatoid Arthritis
- Deformity
- Trauma
- Tumor

Instability in low cervical spine

- Trauma
- Tumor
- Latrogenic instability following laminectomy
- Degeneration of cervical spine

Features & Benefits

- The implants provide the flexibility required to accommodate variations in patient anatomy
- Variable angle self-tapping cortical and cancellous screws
- System includes a complete set of implants and instruments designed to optimize fixation to the occiput

ALPS™ 4.0mm Multi Axial Cancellous Screw



P/N	Material	Spec.
030013540010	TC4	4.0x10
030013540012	TC4	4.0x12
030013540014	TC4	4.0x14
030013540016	TC4	4.0x16
030013540018	TC4	4.0x18
030013540020	TC4	4.0x20
030013540022	TC4	4.0x22
030013540024	TC4	4.0x24

P/N	Material	Spec.
030013540026	TC4	4.0x26
030013540028	TC4	4.0x28
030013540030	TC4	4.0x30
030013540032	TC4	4.0x32
030013540034	TC4	4.0x34
030013540036	TC4	4.0x36
030013540038	TC4	4.0x38
030013540040	TC4	4.0x40

ALPS™ 4.0mm Occipital Screw



P/N	Material	Spec.
030013640006	TC4	4.0x6
030013640008	TC4	4.0x8
030013640010	TC4	4.0x10
030013640012	TC4	4.0x12
030013640014	TC4	4.0x14
030013640016	TC4	4.0x16
030013640018	TC4	4.0x18
030013640020	TC4	4.0x20


ALPS™ 3.5mm Multi Axial Cortical Screw



P/N	Material	Spec.
030013035010	TC4	3.5x10
030013035012	TC4	3.5x12
030013035014	TC4	3.5x14
030013035016	TC4	3.5x16
030013035018	TC4	3.5x18
030013035020	TC4	3.5x20
030013035022	TC4	3.5x22
030013035024	TC4	3.5x24
030013035026	TC4	3.5x26
030013035028	TC4	3.5x28
030013035030	TC4	3.5x30

P/N	Material	Spec.
030013035032	TC4	3.5x32
030013035034	TC4	3.5x34
030013035036	TC4	3.5x36
030013035038	TC4	3.5x38
030013035040	TC4	3.5x40
030013035042	TC4	3.5x42
030013035044	TC4	3.5x44
030013035046	TC4	3.5x46
030013035048	TC4	3.5x48
030013035050	TC4	3.5x50
030013035052	TC4	3.5x52

ALPS™ 4.0mm Multi Axial Cortical Screw



P/N	Material	Spec.
030013040010	TC4	4.0x10
030013040012	TC4	4.0x12
030013040014	TC4	4.0x14
030013040016	TC4	4.0x16
030013040018	TC4	4.0x18
030013040020	TC4	4.0x20
030013040022	TC4	4.0x22
030013040024	TC4	4.0x24
030013040026	TC4	4.0x26
030013040028	TC4	4.0x28
030013040030	TC4	4.0x30

P/N	Material	Spec.
030013040032	TC4	4.0x32
030013040034	TC4	4.0x34
030013040036	TC4	4.0x36
030013040038	TC4	4.0x38
030013040040	TC4	4.0x40
030013040042	TC4	4.0x42
030013040044	TC4	4.0x44
030013040046	TC4	4.0x46
030013040048	TC4	4.0x48
030013040050	TC4	4.0x50
030013040052	TC4	4.0x52

ALPS™ 3.5mm Multi-Axial Cortical Bone Shaft Screw




P/N	Material	Spec.
030013135022	TC4	3.5x22
030013135024	TC4	3.5x24
030013135026	TC4	3.5x26
030013135028	TC4	3.5x28
030013135030	TC4	3.5x30
030013135032	TC4	3.5x32
030013135034	TC4	3.5x34
030013135036	TC4	3.5x36

ALPS™ 4.5mm Occipital Screw



P/N	Material	Spec.
030013145006	TC4	4.5x6
030013145008	TC4	4.5x8
030013145010	TC4	4.5x10
030013145012	TC4	4.5x12
030013145014	TC4	4.5x14
030013145016	TC4	4.5x16
030013145018	TC4	4.5x18
030013145020	TC4	4.5x20

ALPS™ 3.2 Lateral Connector (Closed)



P/N	Material	Spec.
030013032010	TC4	3.2x10

ALPS™ 3.2 Lateral Connector (Open)




P/N	Material	Spec.
030013032011	TC4	3.2x11

ALPS™ Crosslink Rod



P/N	Material	Spec.
030013001040	TC4	40mm
030013001045	TC4	45mm
030013001050	TC4	50mm
030013001055	TC4	55mm
030013001060	TC4	60mm
030013001070	TC4	70mm



ALPS™ Cross Link Clamp



P/N	Material	Spec.
030013001001	TC4	/

ALPS™ Domino Connector



P/N	Material	Spec.
030013001002	TC4	/

ALPS™ Occipital Screw Connector



P/N	Material	Spec.
030013001005	TC4	/

ALPS™ 4.5mm Laminar Hook



P/N	Material	Spec.
030013001006	TC4	/

ALPS™ 6.0mm Laminar Hook



P/N	Material	Spec.
030013001007	TC4	/

ALPS™ Rod



P/N	Material	Spec.
030013001120	TC4	120mm
030013001240	TC4	240mm
030013001300	TC4	300mm

ALPS™ Occipito Cervical Plate




P/N	Material	Spec.
030013001100	TC4	L=100
030013001200	TC4	L=200

ALPS™ Midline Occipital Plate



P/N	Material	Spec.
030013001003	TC4	30°
030013001004	TC4	40°

Titanium Mesh




Indications

- Reconstruct the vertebral column after corpectomy due to tumor or trauma
- Infectious disease, e.g. restore spinal stability after evacuation of vertebral tuberculosis

Features & Benefits

- Titanium implants in various dimensions provide option to choose the configuration that is best suitable to the patient's individual pathology and anatomy
- Mesh implant allows for custom fit and for bone in-growth
- Mesh implants can be inserted anteriorly, laterally or antero-laterally

Titanium Mesh




P/N	Material	Spec.
030073010010	TC4	Φ10×10
030073010030	TC4	Φ10×30
030073010050	TC4	Φ10×50
030073010100	TC4	Φ10×100
030073012010	TC4	Φ12×10
030073012030	TC4	Φ12×30
030073012050	TC4	Φ12×50
030073012100	TC4	Φ12×100
030073014010	TC4	Φ14×10
030073014030	TC4	Φ14×30
030073014050	TC4	Φ14×50
030073014100	TC4	Φ14×100
030073016010	TC4	Φ16×10
030073016030	TC4	Φ16×30
030073016050	TC4	Φ16×50
030073016100	TC4	Φ16×100
030073018010	TC4	Φ18×10
030073018030	TC4	Φ18×30
030073018050	TC4	Φ18×50
030073018100	TC4	Φ18×100
030073020010	TC4	Φ20×10
030073020030	TC4	Φ20×30
030073020050	TC4	Φ20×50
030073020100	TC4	Φ20×100



Cage System



Cervical Cage



PLIF Cage



TLIF Cage



Cervical Cage Indications

- Cervical intervertebral fusion after discectomy, Cervical syndrome, Cervical disc herniation, Anterior cervical instability, Revision surgery for failed disc operation, post-operative instability

Features & Benefits

- Radiolucent PEEK material allows the growth of bone through the center hole of the cage to be visualized-built-in X-rays markers and cage visualization
- Shape teeth on implant surface promote primary stability and prevent cage migration
- Curved and wedge-shaped design accommodates endplate shape variance

Lumbar Cage Indications

- Degenerative discopathy and spinal instability, Degenerative or isthmic spondylolisthesis, Revision surgery for failed disc operation, Lumbar Pseudarthrosis or failed interbody fusion

Features & Benefits




- Radiolucent PEEK material allows the growth of bone through the center hole of the cage to be visualized-built-in X-rays markers and cage visualization
- Shape teeth on implant surface promote primary stability and prevent cage migration
- Curved and wedge-shaped design accommodates endplate shape variance

P/N	Material	Spec.
030103011005	PEEK	5x11
030103011006	PEEK	6x11
030103011007	PEEK	7x11
030103011008	PEEK	8x11
030103011009	PEEK	9x11
030103011010	PEEK	10x11
030103011011	PEEK	11x11
030103011012	PEEK	12x11

P/N	Material	Spec.	P/N	Material	Spec.
030103008022	PEEK	8x22	030103012022	PEEK	12x22
030103008026	PEEK	8x26	030103012026	PEEK	12x26
030103008032	PEEK	8x32	030103012032	PEEK	12x32
030103008036	PEEK	8x36	030103012036	PEEK	12x36
030103010022	PEEK	10x22	030103014022	PEEK	14x22
030103010026	PEEK	10x26	030103014026	PEEK	14x26
030103010032	PEEK	10x32	030103014032	PEEK	14x32
030103010036	PEEK	10x36	030103014036	PEEK	14x36

P/N	Material	Spec.
030103107030	PEEK	7x30
030103109030	PEEK	9x30
030103111030	PEEK	11x30
030103112030	PEEK	12x30
030103113030	PEEK	13x30
030103115030	PEEK	15x30
030103107034	PEEK	7x34
030103109034	PEEK	9x34
030103111034	PEEK	11x34
030103112034	PEEK	12x34
030103113034	PEEK	13x34
030103115034	PEEK	15x34




ATLAS 5.5 Posterior Thoracolumbar Internal Fixation System Instrument Set



No.	REF	Description	Qty
1	22001301	Double-Head Probe, Straight & Curved	1
2	22001302	Single-Head Probe, Straight	1
3	22001303	Tap-Q4.5	1
4	22001304	Tap-Q5.5	1
5	22001305	Tap-Q6.5	1
6	22001306	Tap-Q7.5	1
7	22001307	Trail Rod 200mm	1
8	22001308	Crosslink Screwdriver SW3.5	1
9	22001309	Screw Cutter for Long Arm Screw	1
10	22001310	Crosslink Nut Holder SW3.5	1
11	22001311	Screw Nut Holder T27	1
12	22001312	In-situ Bending Iron L	1
13	22001313	In-situ Bending Iron L	1
14	22001314	Pedicle Probe Straight	1
15	22001315	Pedicle Probe Bent	1
16	22001316	Screwdriver T27	1
17	22001317	Fixation Pin Ball-type	1
18	22001318	Fixation Pin Ball-type	1
19	22001319	Fixation Pin Pillar-type	1
20	22001320	Fixation Pin Pillar-type	1
21	22001321	Crosslink Measurement	1
22	22001322	Screw Measurement	1
23	22001323	Rod Holding Forceps	1
24	22001324	Rod Pushing Forceps	1
25	22001325	5.5mm System Spinal Persuader Sleeve	1
26	22001326	5.5mm System Spinal Persuader	1
27	22001327	Compressor	1
28	22001328	Spreader	1
29	22001329	Rod Twister	1
30	22001330	Counter Torque for Screw Cutter	1
31	22001331	12 Nm Torque Wrench	1
32	22001332	T-handle Quick Coupling	1
33	22001333	Straight Rechat Screwdriver	1
34	22001334	AWL	1
35	22001335	Rod Bender	1
36	22001336	Rod Pusher/Retractor	1
37	22001337	Insert Device for Fixation Pin	1
38	22001338	Rod Holder	1
39	22001339	Screwdriver for Polysial Screw	1
40	22001340	Screwdriver for Monosial Screw	1
41	22001341	Aluminium Box	1







ETNA 5.5 Posterior Thoracolumbar Internal Fixation System Instrument Set






No.	REF	Description	Size	Qty
1	321001	Thoracic Awl		1
2	321002	Curved Probe, for thoracic	Curved	1
3	321003	Straight Probe, for thoracic	Straight	1
4	321004	Curved Probe, for lumbar	Curved	1
5	321005	Straight Probe, for lumbar	Straight	1
6	321006	Double-head Feeler Probe		1
7	321007	Enlarged Feeler Probe		1
8	321008	Pedicle Marker	Cylindrical	4
9	321009	Pedicle Marker	Sphere	4
10	321010	Pedicle Marker Insert		1
11	321011	Tap, Ø4.5	4.5mm	1
12	321012	Tap, Ø5.0	5.0mm	1
13	321013	Tap, Ø5.5	5.5mm	1
14	321014	Tap, Ø6.0	6.0mm	1
15	321015	Tap, Ø6.5	6.5mm	1
16	321016	Tap, Ø7.0	7.0mm	1
17	321017	Monoaxial Screwdriver	Mono	1
18	321018	Polyaxial Screwdriver	Poly	1
19	321019	Provisional Plug Driver	T25	1
20	321020	Final Plug Driver	T25	2
21	321021	Screw Wrench		1
22	321022	Ring Countertorque		1
23	321023	French Bender		1
24	321024	Trail Rod	300mm	1
25	321025	Rod Holder I		1
26	321026	Rod Holder II		2
27	321027	Hex Rod Rotation Wrench	SW4.5	1
28	321028	Rod Pusher	Curved	1
29	321029	In-situ Bender	Left	1
30	321030	In-situ Bender	Right	1
31	321031	Rocker	Short	1
32	321032	Rocker	Long	1
33	321033	Compressor		1
34	321034	Spreader		1
35	321035	Vertebra Distractor	L	1
36	321036	Vertebra Distractor	R	1
37	321037	Beak Rod Reducer		1
38	321038	Reduction Sleeve	Long	1
39	321039	Ratchet Handle		2
40	321040	Quick Connect Handle	T-Shape	1
41	321041	Torque Limiting Handle	12Nm	1
42	321042	Crosslink Nut Driver	SW1.5	1
43	321043	Crosslink Nut Driver	SW8	1
44	321044	Depth Gauge		1
45	321045	Reduction Tab Breaker		1
46	321046	Rod Holder		1
47	321099	Aluminum Box		1

EVEREST 5.5 MIS Spine System






No.	REF	Description	Size	Qty
1	20505000	Reamer	Ø8	1
2	20505105	Guide Pin	Ø11.5	4
3	20505201	Tap, Ø5	Ø5	1
4	20505202	Tap, Ø6	Ø6	1
5	20505203	Tap, Ø7	Ø7	1
6	20505300	Polyaxial Screwdriver	M9	1
7	20505400	Monoaxial Screwdriver	M9	1
8	20505502	Anti-torque Wrench	5.5	1
9	20505600	Multilayer Expander	Ø18.8	1
10	20505700	Quick Connect Handle	Cannulated	2
11	20505800	Reduction Tab Breaker		1
12	20505900	Guide Rod		1
13	20506000	Provisional Tighter	SW4.0	1
14	20506700	Set Screw Driver with T-handle	SW4.0	1
15	20504100	Puncture Needle	Ø9	2
16	20505403	Rod Holder, Long	Ø5.5	1
17	20505404	Rod Holder, Short	Ø5.5	1
18	20507000	Set Screw Driver	SW4.0	1
19	20506900	Reduction Pole A		2
20	20506800	Reduction Pole B		2
21	20507000	Reduction Sleeve		1
22	20507200	Reduction Bar		1
23	20507100	Distraction & Compression Device		1
24	20506500	Compressor		1
25	20506600	Distractor		1
26	23801100	Rod Bender		1
27	90503801	Guide Rack		1
28	20500002	Aluminum Box		1

CERVINO Anterior Cervical Plate Instrument Set



No.	REF	Description	Qty
1	22000201	Plate and Screw Box	1
2	22000202	Drill Bit Ø2.5	1
3	22000203	Drill Bit Ø2.5	1
4	22000204	Screwdriver Sleeve	1
5	22000205	Screwdriver Sleeve	1
6	22000206	Curette	1
7	22000207	Screwdriver Hex SW2.5	1
8	22000208	Tap with Stopper	1
9	22000209	AWL	1
10	22000210	Nut Screwdriver	1
11	22000211	Brace Screwdriver	1
12	22000212	Screwdriver Hex SW2.5	1
13	22000213	Hand Drill	1
14	22000214	Plate Holder Forcep	1
15	22000215	Cervical Hook	1
16	22000216	Cervical Hook	1
17	22000217	Guide	1
18	22000218	Distractor	1
19	22000219	Location Screw Screwdriver	1
20	22000220	Brace Screw	1
21	22000221	Plate Bender	1
22	22000222	Location Screw	4
23	22000223	Aluminium Box	1


ALPS Posterior Cervical Fixation System Instrument Set



No.	REF	Description	Qty
1	22000301	In-situ Bending Iron Left	1
2	22000302	In-situ Bending Iron Right	1
3	22000303	Distractor	1
4	22000304	Feeler for Screw Channel Bent	1
5	22000305	Feeler for Screw Channel Straight	1
6	22000306	Hex Nut Holder SW3.0	1
7	22000307	Hex Screwdriver SW3.0 Long	1
8	22000308	Drill Bit Ø2.4	1
9	22000309	Drill Bit Ø2.7	1
10	22000310	Tap Ø3.5	1
11	22000311	Tap Ø4.0	1
12	22000312	Mould Rod Ø3.5	1
13	22000313	Depth Gauge 0-40mm	1
14	22000314	Crosslink Hex Screwdriver SW2.5 Short	1
15	22000315	Quick Coupling T-handle	1
16	22000316	Rod Pusher	1
17	22000317	Hole Open Forcep	1
18	22000318	AWL	1
19	22000319	Screw/Hook Holder Forcep	1
20	22000320	Rod Holder Forcep	1
21	22000321	Counter Torque	1
22	22000322	Pedicle Screw Screwdriver	1
23	22000323	Insert Device for Fixation Pin	1
24	22000324	Protection Sleeve	1
25	22000325	Drill Guide	1
26	22000326	Rod Cutter	1
27	22000327	Parallel-Compression Forcep	1
28	22000328	Rod Twist	1
29	22000329	Quick Coupling Straight Handle	1
30	22000330	Distractor Forcep	1
31	22000331	Rod Bender	1
32	22000332	Aluminium Box	1




Titanium Mesh Instrument Set



No.	REF	Description	Qty
1	22000501	Cage Stand	1
2	22000502	Pressure 6mm	1
3	22000503	Pressure 18mm	1
4	22000504	Pusher Straight	1
5	22000505	Osteotome	1
6	22000506	Pressure 12mm	1
7	22000507	Pusher Curved	1
8	22000508	Cage Cutter	1
9	22000509	Cage Holding Forcep	1
10	22000510	Implant Measure 10/12mm	1
11	22000511	Implant Measure 16/18mm	1
12	22000512	Implant Measure 22/25mm	1
13	22000513	Aluminium Box	1



Cervical Cage Instrument Set



No.	REF	Description	Qty
1	22000601	Bone Graft Inserter	1
2	22000602	Trail Holder	1
3	22000603	Peek Cage Holder	1
4	22000604	Bone Graft Cage	1
5	22000605	Cage Box	1
6	22000606	Trail 10mm	1
7	22000607	Trail 8mm	1
8	22000608	Trail 6mm	1
9	22000609	Trail 7mm	1
10	22000610	Trail 5mm	1
11	22000611	Trail 5mm	1
12	22000612	Trail 4mm	1
13	22000613	Aluminium Box	1





PLIF Cage Instrument Set



No.	REF	Description	Qty
1	22000701	Trial 20x8mm	1
2	22000702	Trial 20x12mm	1
3	22000703	Trial 22x8mm	1
4	22000704	Trial 12x22mm	1
5	22000705	Trial 26x8mm	1
6	22000706	Trial 26x12mm	1
7	22000707	Trial 20x10mm	1
8	22000708	Trial 20x14mm	1
9	22000709	Trial 22x10mm	1
10	22000710	Trial 22x14mm	1
11	22000711	Trial 26x10mm	1
12	22000712	Reamer 7mm	1
13	22000713	Reamer 8mm	1
14	22000714	Reamer 9mm	1
15	22000715	Bone Graft Insertion	1
16	22000716	Hammer	1
17	22000717	Trial 26x14mm	1
18	22000718	Trial 32x8mm	1
19	22000719	Trial 32x10mm	1
20	22000720	Trial 32x12mm	1
21	22000721	Trial 32x14mm	1
22	22000722	Distractor	1
23	22000723	Bone Graft Cage	1
24	22000724	Cage Holder	1
25	22000725	Curette	1
26	22000726	Nerve Dissector	1
27	22000727	Reamer 8x6.5mm	1
28	22000728	Reamer 8.5x7.5mm	1
29	22000729	Reamer 9.5x8.5mm	1
30	22000730	Reamer 10x10.5mm	1
31	22000731	Reamer 10.5x11.5mm	1
32	22000732	Implant Cage Box	1
33	22000733	Aluminium Box	1

TLIF Cage Instrument Set



No.	REF	Description	Qty
1	22000801	T-handle Quick Coupling	1
2	22000802	Insert Guide	1
3	22000803	Flat Reamer 13mm	1
4	22000804	Flat Reamer 11mm	1
5	22000805	Trial 7x25mm	1
6	22000806	Flat Reamer 9mm	1
7	22000807	Trial 9x25mm	1
8	22000808	Flat Reamer 9mm	1
9	22000809	Trial 11x25mm	1
10	22000810	Trial 13x25mm	1
11	22000811	Trial 15x25mm	1
12	22000812	Cage Holding Forcep	1
13	22000813	Bone Graft Insertion	1
14	22000814	Bone Graft Cage	1
15	22000815	Curved Impactor	1
16	22000816	Curette 8mm Right	1
17	22000817	Side Curette 8mm	1
18	22000818	Straight Impactor	1
19	22000819	Sliding Hammer	1
20	22000820	Ring Curette 7mm	1
21	22000821	Curette 8mm Right	1
22	22000822	Side Curette 6mm	1
23	22000823	Curved Impactor	1
24	22000824	Implant Cage Box	1
25	22000825	Aluminium Box	1

HOF SCIENTIFIC LTD

Legal Manufacturer & EU Authorized Representative

Designed and Manufactured by HOF Scientific, Ltd

128 City Road, EC1V 2NX, London, United Kingdom

+44 7771 033343 - business@hof-scientific.com