





www.hof-scientific.com



TRAUMA SYSTEM PRODUCT CATALOG

People-oriented"technology





COMPANY PROFILE



FAQ

HOF Scientific (UK) & Bensberg Sciences (USA) – A Global Holding Group
The main Base of the Head Office in United Kingdom's
HOF Scientific, in partnership with Bensberg Sciences, proudly leads a dynamic international consortium of companies headquartered across the United Kingdom, the United States, the United Arab Emirates, Germany.
Specializing in medical devices, food supply, and energy solutions (Oil and Gas), our group brings together deep industry knowledge, regulatory expertise, and a commitment to innovation across critical sectors.
With over 15 years of experience, we serve as a strategic bridge between global manufacturers and high-demand markets—particularly across the Middle East ensuring seamless supply chains, market access, and regulatory compliance. As a group, we are united by a vision of long-term partnerships, quality excellence, and market leadership, The group also leads specialized consultancy services in project management, engineering design, construction oversight, and geotechnical and surveying solutions. From planning to execution, we help clients deliver high-quality, well-engineered infrastructure with confidence and efficiency.

PRODUCT CONTENT

Locking Plate System	01-14
Mini-Plate System	15-21
Mini Locking Plate System	22-24
Cannulated Screw System	25-29
Intramedullary Nail System	30-39



Locking Plate System



Indication

The locking plate system is indicated for internal fixation of fractures of the clavicle, humeral, ulna&radius, tibia, femur, calcaneus, metacarpals, pelvis, etc. It may also be used in revision surgery of pseudoarthrosis, non-union and mal-union, osteotomies.

Features & Benefits

- Anatomically pre-contoured plate.
- Reduced plate thickness in the periarticular region.
- Rounded and tapered plate ends.
- K-wire or suture holes for temporary plate fixation.
- Cortical screw can be used.
- Limited contacts design on the plate's under-surface.

Femur LC-LCP



P/N	Material	Spec.	Length
030071001006	TA3	6H	133
030071001007	TA3	7H	151
030071001008	TA3	8H	169
030071001009	TA3	9H	187
030071001010	TA3	10H	205
030071001011	TA3	11H	223
030071001012	TA3	12H	241
030071001014	TA3	14H	277
030071001016	TA3	16H	313
030071001018	TA3	18H	349

Using 5.0mm locking screws, 4.5mm cortical screws.

Proximal Lateral Femur Locking Plate III



P/N	Material	Spec.	Length
030071002005	TA3	5H, L	118
030071002007	TA3	7H, L	154
030071002009	TA3	9H, L	190
030071002011	TA3	11H, L	226
030071002013	TA3	13H, L	262
030071002105	TA3	5H, R	118
030071002107	TA3	7H, R	154
030071002109	TA3	9H, R	190
030071002111	TA3	11H, R	226
030071002113	TA3	13H, R	262

Using 5.0mm/6.5mm locking screws, 4.5mm cortical screws.

Proximal Femur Condylar Locking Plate



P/N	Material	Spec.	Length
030071074104	TA3	4H, L	149
030071074106	TA3	6H, L	187
030071074108	TA3	8H, L	225
030071074110	TA3	10H, L	263
030071074112	TA3	12H, L	301
030071074114	TA3	14H, L	339
030071074116	TA3	16H, L	377
030071074004	TA3	4H, R	149
030071074006	TA3	6H, R	187
030071074008	TA3	8H, R	225
030071074010	TA3	10H, R	263
030071074012	TA3	12H, R	301
030071074014	TA3	14H, R	339
030071074016	TA3	16H, R	377

Using 5.0mm/6.5mm locking screws, 4.5mm cortical screws.

Distal Lateral Femoral Supracondylar Locking Plate



P/N	Material	Spec.	Length	P/N	Material	Spec.	Length
030071075104	TA3	4H, L	110	030071075004	TA3	4H, R	110
030071075105	TA3	5H, L	128	030071075005	TA3	5H, R	128
030071075106	TA3	6H, L	146	030071075006	TA3	6H, R	146
030071075107	TA3	7H, L	163.5	030071075007	TA3	7H, R	163.5
030071075108	TA3	8H, L	182	030071075008	TA3	8H, R	182
030071075109	TA3	9H, L	200	030071075009	TA3	9H, R	200
030071075110	TA3	10H, L	218	030071075010	TA3	10H, R	218
030071075111	TA3	11H, L	236	030071075011	TA3	11H, R	236
030071075112	TA3	12H, L	254	030071075012	TA3	12H, R	254
030071075113	TA3	13H, L	270	030071075013	TA3	13H, R	270
030071075114	TA3	14H, L	286	030071075014	TA3	14H, R	286

Using 5.0mm/6.5mm locking screws, 4.5mm cortical screws.

Distal Femur Locking Plate II



P/N	Material	Spec.	Length	P/N	Material	Spec.	Length
030071032005	TA3	5H, L	122	030071032105	TA3	5H, R	122
030071032006	TA3	6H, L	140	030071032106	TA3	6H, R	140
030071032007	TA3	7H, L	158	030071032107	TA3	7H, R	158
030071032008	TA3	8H, L	176	030071032108	TA3	8H, R	176
030071032009	TA3	9H, L	194	030071032109	TA3	9H, R	194
030071032010	TA3	10H, L	212	030071032110	TA3	10H, R	212
030071032011	TA3	11H, L	230	030071032111	TA3	11H, R	230
030071032012	TA3	12H, L	248	030071032112	TA3	12H, R	248
030071032013	TA3	13H, L	266	030071032113	TA3	13H, R	266
030071032014	TA3	14H, L	284	030071032114	TA3	14H, R	284

Using 5.0mm locking screws, 4.5mm cortical screws.

Distal Lateral Femur Locking Plate IV



P/N	Material	Spec.	Length
030071003005	TA3	5H, L	127
030071003007	TA3	7H, L	159
030071003009	TA3	9H, L	191
030071003011	TA3	11H, L	223
030071003013	TA3	13H, L	255
030071003015	TA3	15H, L	287
030071003017	TA3	17H, L	319
030071003105	TA3	5H, R	127
030071003107	TA3	7H, R	159
030071003109	TA3	9H, R	191
030071003111	TA3	11H, R	223
030071003113	TA3	13H, R	255
030071003115	TA3	15H, R	287
030071003117	TA3	17H, R	319

Using 5.0mm locking screws, 4.5mm cortical screws.



Tibia LC-LCP



P/N	Material	Spec.	Length
030071004005	TA3	5H	103
030071004006	TA3	6H	120.5
030071004007	TA3	7H	138
030071004008	TA3	8H	155.5
030071004009	TA3	9H	173
030071004010	TA3	10H	190.5
030071004011	TA3	11H	208
030071004012	TA3	12H	225.5
030071004013	TA3	13H	243
030071004015	TA3	15H	278
030071004017	TA3	17H	313

Using 5.0mm locking screws, 4.5mm cortical screws.

Proximal Medial Tibia Locking Plate



P/N	Material	Spec.	Length
030071076105	TA3	5H, L	125
030071076107	TA3	7H, L	161
030071076109	TA3	9H, L	197
030071076111	TA3	11H, L	233
030071076113	TA3	13H, L	269
030071076005	TA3	5H, R	125
030071076007	TA3	7H, R	161
030071076009	TA3	9H, R	197
030071076011	TA3	11H, R	233
030071076013	TA3	13H, R	269

Using 5.0mm locking screws, 4.5mm cortical screws.

4.5 Proximal Medial Tibia Locking Plate II



P/N	Material	Spec.	Length
030071008004	TA3	4H, L	81
030071008005	TA3	5H, L	97
030071008006	TA3	6H, L	113
030071008007	TA3	7H, L	129
030071008008	TA3	8H, L	145
030071008104	TA3	4H, R	81
030071008105	TA3	5H, R	97
030071008106	TA3	6H, R	113
030071008107	TA3	7H, R	129
030071008108	TA3	8H, R	145

Using 5.0mm locking screws, 4.5mm cortical screws.

Proximal Lateral Tibia Locking Plate IV



P/N	Material	Spec.	Length
030071005005	TA3	5H, L	137
030071005007	TA3	7H, L	169
030071005009	TA3	9H, L	201
030071005011	TA3	11H, L	233
030071005105	TA3	5H, R	137
030071005107	TA3	7H, R	169
030071005109	TA3	9H, R	201
030071005111	TA3	11H, R	233

Using 5.0mm locking screws, 4.5mm cortical screws.

Proximal Lateral Tibia Locking Plate VI



P/N	Material	Spec.	Length
030071024004	TA3	4H, L	83
030071024005	TA3	5H, L	99
030071024006	TA3	6H, L	115
030071024007	TA3	7H, L	131
030071024008	TA3	8H, L	147
030071024104	TA3	4H, R	83
030071024105	TA3	5H, R	99
030071024106	TA3	6H, R	115
030071024107	TA3	7H, R	131
030071024108	TA3	8H, R	147

Using 3.5mm locking screws, 3.5mm cortical screws.

Proximal Lateral Tibia Locking Plate VII



P/N	Material	Spec.	Length	P/N	Material	Spec.	Length
030071031002	TA3	2H, L	90	030071031002	TA3	2H, R	90
030071031003	TA3	3H, L	108	030071031003	TA3	3H, R	108
030071031104	TA3	4H, L	126	030071031104	TA3	4H, R	126
030071031105	TA3	5H, L	144	030071031105	TA3	5H, R	144
030071031106	TA3	6H, L	162	030071031106	TA3	6H, R	162
030071031107	TA3	7H, L	180	030071031107	TA3	7H, R	180
030071031108	TA3	8H, L	198	030071031108	TA3	8H, R	198
030071031109	TA3	9H, L	216	030071031109	TA3	9H, R	216
030071031110	TA3	10H, L	234	030071031110	TA3	10H, R	234
030071031111	TA3	11H, L	252	030071031111	TA3	11H, R	252
030071031112	TA3	12H, L	270	030071031112	TA3	12H, R	270
030071031113	TA3	13H, L	288	030071031113	TA3	13H, R	288
030071031114	TA3	14H, L	306	030071031114	TA3	14H, R	306

Using 5.0mm locking screws, 4.5mm cortical screws.

3.5 Distal Medial Tibia Locking Plate II



P/N	Material	Spec.	Length
030071009004	TA3	4H, L	109
030071009006	TA3	6H, L	133
030071009008	TA3	8H, L	157
030071009010	TA3	10H, L	181
030071009012	TA3	12H, L	205
030071009014	TA3	14H, L	229
030071009104	TA3	4H, R	109
030071009106	TA3	6H, R	133
030071009108	TA3	8H, R	157
030071009110	TA3	10H, R	181
030071009112	TA3	12H, R	205
030071009114	TA3	14H, R	229

Using 3.5mm locking screws, 3.5mm cortical screws.

Distal Medial Tibia Locking Plate II



P/N	Material	Spec.	Length
030071007005	TA3	5H, L	108
030071007007	TA3	7H, L	140
030071007009	TA3	9H, L	172
030071007011	TA3	11H, L	204
030071007013	TA3	13H, L	236
030071007105	TA3	5H, R	108
030071007107	TA3	7H, R	140
030071007109	TA3	9H, R	172
030071007111	TA3	11H, R	204
030071007113	TA3	13H, R	236

Using 5.0mm locking screws, 4.5mm cortical screws.



3.5 Distal Lateral Tibia Locking Plate II



P/N	Material	Spec.	Length
030071010004	TA3	4H, L	58.5
030071010006	TA3	6H, L	82.3
030071010008	TA3	8H, L	106.5
030071010010	TA3	10H, L	130.5
030071010012	TA3	12H, L	154.5
030071010104	TA3	4H, R	58.5
030071010106	TA3	6H, R	82.3
030071010108	TA3	8H, R	106.5
030071010110	TA3	10H, R	130.5
030071010112	TA3	12H, R	154.5

Using 3.5mm locking screws, 3.5mm cortical screws.

Distal Lateral Tibia Locking Plate II



P/N	Material	Spec.	Length
030071006005	TA3	5H, L	108
030071006007	TA3	7H, L	140
030071006009	TA3	9H, L	172
030071006011	TA3	11H, L	204
030071006013	TA3	13H, L	236
030071006105	TA3	5H, R	108
030071006107	TA3	7H, R	140
030071006109	TA3	9H, R	172
030071006111	TA3	11H, R	204
030071006113	TA3	13H, R	236

Using 5.0mm locking screws, 4.5mm cortical screws.

Distal Lateral Fibular Locking Plate VI



P/N	Material	Spec.	Length	P/N	Material	Spec.	Length
030071025003	TA3	3H, L	73	030071026103	TA3	3H, R	73
030071025004	TA3	4H, L	86	030071026104	TA3	4H, R	86
030071025005	TA3	5H, L	99	030071026105	TA3	5H, R	99
030071025006	TA3	6H, L	112	030071026106	TA3	6H, R	112
030071025007	TA3	7H, L	125	030071026107	TA3	7H, R	125
030071025008	TA3	8H, L	138	030071026108	TA3	8H, R	138
030071025009	TA3	9H, L	151	030071026109	TA3	9H, R	151
030071025010	TA3	10H, L	164	030071026110	TA3	10H, R	164
030071025011	TA3	11H, L	177	030071026111	TA3	11H, R	177
030071025012	TA3	12H, L	190	030071026112	TA3	12H, R	190

Using 3.5mm locking screws, 3.5mm cortical screws.

Calcaneal Locking Plate VIII



P/N	Material	Spec.	Length
030071080051	TA3	L	51
030071080060	TA3	L	60
030071080068	TA3	L	68
030071080151	TA3	R	51
030071080160	TA3	R	60
030071080168	TA3	R	68

Using 3.5mm locking screws, 3.5mm cortical screws.

Straight Reconstruction Locking Plate II



3.5 Straight Reconstruction Locking Plate II



4.5 Straight Reconstruction Locking Plate II



Oblique Radius T-Locking Plates II



P/N	Material	Spec.	Length
030071011004	TA3	4H	48
030071011005	TA3	5H	60
030071011006	TA3	6H	72
030071011007	TA3	7H	84
030071011008	TA3	8H	96
030071011009	TA3	9H	108
030071011010	TA3	10H	120
030071011011	TA3	11H	132
030071011012	TA3	12H	144
030071011014	TA3	14H	168
030071011016	TA3	16H	192
030071011018	TA3	18H	216
030071011020	TA3	20H	240
030071011022	TA3	22H	264

Using 3.5mm locking screws, 3.5mm cortical screws.

P/N	Material	Spec.	Length
030071022006	TA3	6H	86
030071022007	TA3	7H	99
030071022008	TA3	8H	112
030071022009	TA3	9H	125
030071022010	TA3	10H	138
030071022011	TA3	11H	151
030071022012	TA3	12H	164
030071022014	TA3	14H	190
030071022016	TA3	16H	216
030071022018	TA3	18H	242

Using 3.5mm locking screws, 3.5mm cortical screws.

P/N	Material	Spec.	Length
030071023003	TA3	3H	119
030071023004	TA3	4H	136
030071023005	TA3	5H	153
030071023006	TA3	6H	170
030071023007	TA3	7H	187
030071023008	TA3	8H	204
030071023009	TA3	9H	221
030071023011	TA3	11H	255
030071023013	TA3	13H	289
030071023015	TA3	15H	323

Using 3.5mm locking screws, 3.5mm cortical screws.

P/N	Material	Spec.	Length
030071012003	TA3	3H, L	50.2
030071012004	TA3	4H, L	61.6
030071012005	TA3	5H, L	73
030071012006	TA3	6H, L	84.4
030071012007	TA3	7H, L	95.8
030071012008	TA3	8H, L	107.2
030071012103	TA3	3H, R	50.2
030071012104	TA3	4H, R	61.6
030071012105	TA3	5H, R	73
030071012106	TA3	6H, R	84.4
030071012107	TA3	7H, R	95.8
030071012108	TA3	8H, R	107.2

Using 3.5mm locking screws, 3.5mm cortical screws.



Small T-Locking Plate (head holes 4)



P/N	Material	Spec.	Length
030071013003	TA3	3H	46.1
030071013004	TA3	4H	57.5
030071013005	TA3	5H	68.9
030071013006	TA3	6H	80.3
030071013007	TA3	7H	91.7
030071013008	TA3	8H	103.1

Using 3.5mm locking screws, 3.5mm cortical screws.

Small T-Locking Plate (head holes 3)



P/N	Material	Spec.	Length
030071014003	TA3	3H	50
030071014004	TA3	4H	59.5
030071014005	TA3	5H	71
030071014006	TA3	6H	82.5
030071014007	TA3	7H	94
030071014008	TA3	8H	105.5

Using 3.5mm locking screws, 3.5mm cortical screws.

Distal Radius Locking Plate I , 9 Holes Head



P/N	Material	Spec.	Length
030071033003	TA3	3H, L	50
030071033004	TA3	4H, L	59
030071033005	TA3	5H, L	68
030071033006	TA3	6H, L	77
030071033007	TA3	7H, L	86
030071033008	TA3	8H, L	95
030071033103	TA3	3H, R	50
030071033104	TA3	4H, R	59
030071033105	TA3	5H, R	68
030071033106	TA3	6H, R	77
030071033107	TA3	7H, R	86
030071033108	TA3	8H, R	95

Using 2.7mm locking screws, 2.7mm cortical screws.

Variable Angle Two-Column Volar Distal Radius Locking Plate, Narrow



P/N	Material	Spec.	Length
030071071002	TA3	2H, L	42
030071071003	TA3	3H, L	51
030071071004	TA3	4H, L	63
030071071005	TA3	5H, L	72
030071071102	TA3	2H, R	42
030071071103	TA3	3H, R	51
030071071104	TA3	4H, R	63
030071071105	TA3	5H, R	72

Using 2.7mm VA locking screws, 2.7mmlocking screws, 2.7mm cortical screws.

Variable Angle Two-Column Volar Distal Radius Locking Plate, Head Holes 6



Variable Angle Two-Column Volar Distal Radius Locking Plate, Head Holes 7



Ulna Radius LC-DCP



Humerus LC-DCP



P/N	Material	Spec.	Length
030071072002	TA3	2H, L	45
030071072003	TA3	3H, L	54
030071072004	TA3	4H, L	66
030071072005	TA3	5H, L	75
030071072102	TA3	2H, R	45
030071072103	TA3	3H, R	54
030071072104	TA3	4H, R	66
030071072105	TA3	5H, R	75

Using 2.7mm VA locking screws, 2.7mmlocking screws, 2.7mm cortical screws.

P/N	Material	Spec.	Length
030071073002	TA3	2H, L	47
030071073003	TA3	3H, L	55
030071073004	TA3	4H, L	66
030071073005	TA3	5H, L	75
030071073102	TA3	2H, R	47
030071073103	TA3	3H, R	55
030071073104	TA3	4H, R	66
030071073105	TA3	5H, R	75

Using 2.7mm VA locking screws, 2.7mmlocking screws, 2.7mm cortical screws.

P/N	Material	Spec.	Length
030071016004	TA3	4H	67
030071016005	TA3	5H	80
030071016006	TA3	6H	93
030071016007	TA3	7H	106
030071016008	TA3	8H	119
030071016009	TA3	9H	132
030071016010	TA3	10H	145
030071016011	TA3	11H	158
030071016012	TA3	12H	171

Using 3.5mm locking screws, 3.5mm cortical screws.

P/N	Material	Spec.	Length
030071015004	TA3	4H	72.5
030071015005	TA3	5H	87.2
030071015006	TA3	6H	101.9
030071015007	TA3	7H	116.6
030071015008	TA3	8H	131.3
030071015009	TA3	9H	146
030071015010	TA3	10H	160.7
030071015011	TA3	11H	175.4
030071015012	TA3	12H	190.1

Using 3.5mm locking screws, 3.5mm cortical screws.



Proximal Dorsal Ulna Locking Plate II



P/N	Material	Spec.	Length
030071020004	TA3	4H, L	86
030071020006	TA3	6H, L	112
030071020008	TA3	8H, L	138
030071020010	TA3	10H, L	164
030071020012	TA3	12H, L	190
030071020014	TA3	14H, L	216
030071020104	TA3	4H, R	86
030071020106	TA3	6H, R	112
030071020108	TA3	8H, R	138
030071020110	TA3	10H, R	164
030071020112	TA3	12H, R	190
030071020114	TA3	14H, R	216

Using 3.5mm locking screws, 3.5mm cortical screws.

Proximal Lateral Humerus Locking Plate IV



P/N	Material	Spec.	Length
030071017003	TA3	3H, L	80
030071017004	TA3	4H, L	92
030071017005	TA3	5H, L	104
030071017006	TA3	6H, L	116
030071017007	TA3	7H, L	128
030071017008	TA3	8H, L	140
030071017103	TA3	3H, R	80
030071017104	TA3	4H, R	92
030071017105	TA3	5H, R	104
030071017106	TA3	6H, R	116
030071017107	TA3	7H, R	128
030071017108	TA3	8H, R	140

Using 3.5mm locking screws, 3.5mm cortical screws.

3.5 Proximal Lateral Humeral Locking Plate II



P/N	Material	Spec.	Length
030071019002	TA3	2H, Standard	90
030071019003	TA3	3H, Standard	108
030071019004	TA3	4H, Standard	126
030071019005	TA3	5H, Standard	144
030071019006	TA3	6H, Long	162
030071019007	TA3	7H, Long	180
030071019008	TA3	8H, Long	198
030071019009	TA3	9H, Long	216
030071019010	TA3	10H, Long	234
030071019011	TA3	11H, Long	252
030071019012	TA3	12H, Long	270

Using 3.5mm locking screws, 3.5mm cortical screws.

Proximal Humeral Locking Plate



P/N	Material	Spec.	Length	P/N	Material	Spec.	Length
030071077102	TA3	2H, L	91	030071077002	TA3	2H, R	91
030071077103	TA3	3H, L	107	030071077003	TA3	3H, R	107
030071077104	TA3	4H, L	123	030071077004	TA3	4H, R	123
030071077105	TA3	5H, L	139	030071077005	TA3	5H, R	139
030071077106	TA3	6H, L	155	030071077006	TA3	6H, R	155
030071077107	TA3	7H, L	171	030071077007	TA3	7H, R	171
030071077108	TA3	8H, L	187	030071077008	TA3	8H, R	187
030071077109	TA3	9H, L	203	030071077009	TA3	9H, R	203
030071077110	TA3	10H, L	219	030071077010	TA3	10H, R	219
030071077111	TA3	11H, L	235	030071077011	TA3	11H, R	235
030071077112	TA3	12H, L	251	030071077012	TA3	12H, R	251
030071077113	TA3	13H, L	267	030071077013	TA3	13H, R	267
030071077114	TA3	14H, L	283	030071077014	TA3	14H, R	283

Using 3.5mm locking screws, 3.5mm cortical screws.

Distal Medial Humeral Locking Plate II



P/N	Material	Spec.	Length
030071027003	TA3	3H, L	58
030071027005	TA3	5H, L	84
030071027007	TA3	7H, L	110
030071027009	TA3	9H, L	136
030071027011	TA3	11H, L	162
030071027013	TA3	13H, L	188
030071027015	TA3	15H, L	214
030071027103	TA3	3H, R	58
030071027105	TA3	5H, R	84
030071027107	TA3	7H, R	110
030071027109	TA3	9H, R	136
030071027111	TA3	11H, R	162
030071027113	TA3	13H, R	188
030071027115	TA3	15H, R	214

Using 2.7mm/3.5mm locking screws, 3.5mm cortical screws.

Distal Posterior Lateral Humeral Locking Plate II



P/N	Material	Spec.	Length
030071028004	TA3	4H, L	70
030071028006	TA3	6H, L	94
030071028008	TA3	8H, L	120
030071028010	TA3	10H, L	146
030071028012	TA3	12H, L	172
030071028104	TA3	4H, R	70
030071028106	TA3	6H, R	94
030071028108	TA3	8H, R	120
030071028110	TA3	10H, R	146
030071028112	TA3	12H, R	172

Using 2.7mm/3.5mm locking screws, 3.5mm cortical screws.

Distal Humeral Y-Locking Plate II



P/N	Material	Spec.	Length
030071078103	TA3	3H, L	88.5
030071078104	TA3	4H, L	104.5
030071078105	TA3	5H, L	120.5
030071078003	TA3	3H, R	88.5
030071078004	TA3	4H, R	104.5
030071078005	TA3	5H, R	120.5

Using 3.5mm locking screws, 3.5mm cortical screws.

Clavicle Diaphysis Locking Plate II



P/N	Material	Spec.	Length
030071018006	TA3	6H, L	74.1
030071018008	TA3	8H, L	98.1
030071018010	TA3	10H, L	121.4
030071018012	TA3	12H, L	143.8
030071018106	TA3	6H, R	74.1
030071018108	TA3	8H, R	98.1
030071018110	TA3	10H, R	121.4
030071018112	TA3	12H, R	143.8

Using 3.5mm locking screws, 3.5mm cortical screws.



Clavicle Hook Locking Plate IV



P/N	Material	Spec.	Hook Length
030071081005	TA3	5H, L	14
030071081006	TA3	6H, L	14
030071081007	TA3	7H, L	14
030071082005	TA3	5H, L	17
030071082006	TA3	6H, L	17
030071082007	TA3	7H, L	17
030071081105	TA3	5H, R	14
030071081106	TA3	6H, R	14
030071081107	TA3	7H, R	14
030071082105	TA3	5H, R	17
030071082106	TA3	6H, R	17
030071082107	TA3	7H, R	17

Using 3.5mm locking screws, 3.5mm cortical screws.

Clavicle Hook Locking Plate V



P/N	Material	Spec.	Length
030071021004	TA3	4H, L	73.5
030071021104	TA3	4H, R	73.5

Using 3.5mm locking screws, 3.5mm cortical screws.

Clavicle Locking plate II



P/N	Material	Spec.	Length
030071029003	TA3	3H, L	76
030071029004	TA3	4H, L	88
030071029005	TA3	5H, L	100
030071029006	TA3	6H, L	112
030071029007	TA3	7H, L	124
030071029008	TA3	8H, L	135
030071029103	TA3	3H, R	76
030071029104	TA3	4H, R	88
030071029105	TA3	5H, R	100
030071029106	TA3	6H, R	112
030071029107	TA3	7H, R	124
030071029108	TA3	8H, R	135

Using 3.5mm locking screws, 3.5mm cortical screws.

S-Clavicle Locking Plate



P/N	Material	Spec.	Length
030071030006	TA3	6H, L	74
030071030007	TA3	7H, L	86
030071030008	TA3	8H, L	98
030071030010	TA3	10H, L	122
030071030106	TA3	6H, R	74
030071030107	TA3	7H, R	86
030071030108	TA3	8H, R	98
030071030110	TA3	10H, R	122

Using 3.5mm locking screws, 3.5mm cortical screws.

2.7mm Variable Angle Locking Screw



P/N	Material	Spec.
030012227008	TC4	8mm
030012227010	TC4	10mm
030012227012	TC4	12mm
030012227014	TC4	14mm
030012227016	TC4	16mm
030012227018	TC4	18mm
030012227020	TC4	20mm
030012227022	TC4	22mm
030012227024	TC4	24mm
030012227026	TC4	26mm
030012227028	TC4	28mm
030012227030	TC4	30mm

2.7mm Locking Screw



P/N	Material	Spec.
030012327006	TC4	6mm
030012327008	TC4	8mm
030012027010	TC4	10mm
030012027012	TC4	12mm
030012027014	TC4	14mm
030012027016	TC4	16mm
030012027018	TC4	18mm
030012027020	TC4	20mm
030012027022	TC4	22mm
030012027024	TC4	24mm
030012027026	TC4	26mm
030012027028	TC4	28mm
030012027030	TC4	30mm

P/N	Material	Spec.
030012027032	TC4	32mm
030012027034	TC4	34mm
030012027036	TC4	36mm
030012027038	TC4	38mm
030012027040	TC4	40mm
030012027042	TC4	42mm
030012027044	TC4	44mm
030012027046	TC4	46mm
030012027048	TC4	48mm
030012027050	TC4	50mm
030012027055	TC4	55mm
030012027060	TC4	60mm

3.5mm Locking Screw



P/N	Material	Spec.
030012035010	TC4	10mm
030012035012	TC4	12mm
030012035014	TC4	14mm
030012035016	TC4	16mm
030012035018	TC4	18mm
030012035020	TC4	20mm
030012035022	TC4	22mm
030012035024	TC4	24mm
030012035026	TC4	26mm
030012035028	TC4	28mm
030012035030	TC4	30mm
030012035032	TC4	32mm
030012035034	TC4	34mm
030012035036	TC4	36mm

P/N	Material	Spec.
030012035038	TC4	38mm
030012035040	TC4	40mm
030012035042	TC4	42mm
030012035044	TC4	44mm
030012035046	TC4	46mm
030012035048	TC4	48mm
030012035050	TC4	50mm
030012035055	TC4	55mm
030012035060	TC4	60mm
030012035065	TC4	65mm
030012035070	TC4	70mm
030012035075	TC4	75mm
030012035080	TC4	80mm
030012035085	TC4	85mm

5.0mm Locking Screw



P/N	Material	Spec.
030012050018	TC4	18mm
030012050020	TC4	20mm
030012050022	TC4	22mm
030012050024	TC4	24mm
030012050026	TC4	26mm
030012050028	TC4	28mm
030012050030	TC4	30mm
030012050032	TC4	32mm
030012050034	TC4	34mm
030012050036	TC4	36mm
030012050038	TC4	38mm
030012050040	TC4	40mm
030012050042	TC4	42mm

P/N	Material	Spec.
030012050044	TC4	44mm
030012050046	TC4	46mm
030012050048	TC4	48mm
030012050050	TC4	50mm
030012050055	TC4	55mm
030012050060	TC4	60mm
030012050065	TC4	65mm
030012050070	TC4	70mm
030012050075	TC4	75mm
030012050080	TC4	80mm
030012050085	TC4	85mm
030012050090	TC4	90mm



6.5mm Locking Screw



P/N	Material	Spec.
030012065035	TC4	35mm
030012065040	TC4	40mm
030012065045	TC4	45mm
030012065050	TC4	50mm
030012065055	TC4	55mm
030012065060	TC4	60mm
030012065065	TC4	65mm
030012065070	TC4	70mm
030012065075	TC4	75mm

P/N	Material	Spec.
030012065080	TC4	80mm
030012065085	TC4	85mm
030012065090	TC4	90mm
030012065095	TC4	95mm
030012065100	TC4	100mm
030012065105	TC4	105mm
030012065110	TC4	110mm
030012065115	TC4	115mm
030012065120	TC4	120mm

2.7mm Cortical Screws, Full Threaded



P/N	Material	Spec.
030012427006	TC4	6mm
030012427008	TC4	8mm
030012427010	TC4	10mm
030012427012	TC4	12mm
030012427014	TC4	14mm
030012427016	TC4	16mm
030012427018	TC4	18mm
030012427020	TC4	20mm
030012427022	TC4	22mm
030012427024	TC4	24mm
030012427026	TC4	26mm
030012427028	TC4	28mm
030012427030	TC4	30mm

3.5mm Cortical Screws, Full Threaded



P/N	Material	Spec.
030012135010	TC4	10mm
030012135012	TC4	12mm
030012135014	TC4	14mm
030012135016	TC4	16mm
030012135018	TC4	18mm
030012135020	TC4	20mm
030012135022	TC4	22mm
030012135024	TC4	24mm
030012135026	TC4	26mm
030012135028	TC4	28mm
030012135030	TC4	30mm

P/N	Material	Spec.
030012135032	TC4	32mm
030012135034	TC4	34mm
030012135036	TC4	36mm
030012135038	TC4	38mm
030012135040	TC4	40mm
030012135042	TC4	42mm
030012135044	TC4	44mm
030012135046	TC4	46mm
030012135048	TC4	48mm
030012135050	TC4	50mm

4.0mm Cancellous Screws, Partial Threaded



P/N	Material	Spec.
030012140012	TC4	12mm
030012140014	TC4	14mm
030012140016	TC4	16mm
030012140018	TC4	18mm
030012140020	TC4	20mm
030012140022	TC4	22mm
030012140024	TC4	24mm
030012140026	TC4	26mm
030012140028	TC4	28mm
030012140030	TC4	30mm
030012140032	TC4	32mm

P/N	Material	Spec.
030012140034	TC4	34mm
030012140036	TC4	36mm
030012140038	TC4	38mm
030012140040	TC4	40mm
030012140042	TC4	42mm
030012140044	TC4	44mm
030012140046	TC4	46mm
030012140048	TC4	48mm
030012140050	TC4	50mm
030012140052	TC4	52mm
030012140054	TC4	54mm

4.0mm Cancellous Screws, Full Threaded



4.5mm Cortical Screw



6.5mm Cancellous screw, full threaded



6.5mm Cancellous screw, partial threaded



P/N	Material	Spec.
030012240014	TC4	14mm
030012240016	TC4	16mm
030012240018	TC4	18mm
030012240020	TC4	20mm
030012240022	TC4	22mm
030012240024	TC4	24mm
030012240026	TC4	26mm
030012240028	TC4	28mm
030012240030	TC4	30mm
030012240032	TC4	32mm
030012240034	TC4	34mm

P/N	Material	Spec.
030012240036	TC4	36mm
030012240038	TC4	38mm
030012240040	TC4	40mm
030012240042	TC4	42mm
030012240044	TC4	44mm
030012240046	TC4	46mm
030012240048	TC4	48mm
030012240050	TC4	50mm
030012240052	TC4	52mm
030012240054	TC4	54mm

P/N	Material	Spec.
030012045016	TC4	16mm
030012045018	TC4	18mm
030012045020	TC4	20mm
030012045022	TC4	22mm
030012045024	TC4	24mm
030012045026	TC4	26mm
030012045028	TC4	28mm
030012045030	TC4	30mm
030012045032	TC4	32mm
030012045034	TC4	34mm
030012045036	TC4	36mm

P/N	Material	Spec.
030012045038	TC4	38mm
030012045040	TC4	40mm
030012045042	TC4	42mm
030012045044	TC4	44mm
030012045046	TC4	46mm
030012045048	TC4	48mm
030012045050	TC4	50mm
030012045052	TC4	52mm
030012045054	TC4	54mm
030012045056	TC4	56mm
030012045058	TC4	58mm

P/N	Material	Spec.
030012165030	TC4	30mm
030012165035	TC4	35mm
030012165040	TC4	40mm
030012165045	TC4	45mm
030012165050	TC4	50mm
030012165055	TC4	55mm
030012165060	TC4	60mm
030012165065	TC4	65mm
030012165070	TC4	70mm
030012165075	TC4	75mm

P/N	Material	Spec.
030012165080	TC4	80mm
030012165085	TC4	85mm
030012165090	TC4	90mm
030012165095	TC4	95mm
030012165100	TC4	100mm
030012165105	TC4	105mm
030012165110	TC4	110mm
030012165115	TC4	115mm
030012165120	TC4	120mm

P/N	Material	Spec.
030012265030	TC4	30mm
030012265035	TC4	35mm
030012265040	TC4	40mm
030012265045	TC4	45mm
030012265055	TC4	55mm
030012265060	TC4	60mm
030012265065	TC4	65mm
030012265070	TC4	70mm
030012265075	TC4	75mm

P/N	Material	Spec.
030012265080	TC4	80mm
030012265085	TC4	85mm
030012265090	TC4	90mm
030012265095	TC4	95mm
030012265100	TC4	100mm
030012265105	TC4	105mm
030012265110	TC4	110mm
030012265115	TC4	115mm
030012265120	TC4	120mm



Mini-Plate System



Indications

The Mini-Plate System is indicated for treatment of small bone fracture of hand and foot (except calcaneus and talus). It may also be used in revision surgery of non-union and mal-union.

Features & Benefits

- Available 1.5mm, 2.0mm, 2.5mm and 2.7mm implants
- StarDrive recess mates with self-retaining screwdriver improving torque transmission
- Low profile plates reduce soft tissue irritation

Phalangeal Mini-Plate, Straight



P/N	Material	Spec.	Length
030071034003	TA3	3H	19
030071034004	TA3	4H	26
030071034005	TA3	5H	33
030071034006	TA3	6H	40
030071034007	TA3	7H	47

Using 1.5mm Cortical Screw

Metacarpal Mini-Plate, T-type



P/N	Material	Spec.	Length
030071035004	TA3	4H	17
030071035005	TA3	5H	29.5
030071035006	TA3	6H	37
030071035007	TA3	7H	44.5
030071035008	TA3	8H	52

Using 1.5mm Cortical Screw

Phalangeal Mini-Plate, L-type



P/N	Material	Spec.	Length
030071036004	TA3	4H, L	17
030071036104	TA3	4H, R	17

Using 1.5mm Cortical Screw

Metacarpal Mini-Plate, L-type



1.5 Straight Plate



1.5 T Condylar Plate(3 holes head)



1.5 T Condylar Plate (4 holes head)



1.5 H Condylar Plate (4 holes)



P/N	Material	Spec.	Length
030071036005	TA3	5H, L	29.5
030071036105	TA3	5H, R	29.5

Using 1.5mm Cortical Screw

P/N	Material	Spec.	Length
030071037003	TA3	3H	13.8
030071037004	TA3	4H	18.8
030071037005	TA3	5H	23.8
030071037006	TA3	6H	28.8
030071037007	TA3	7H	33.8
030071037008	TA3	8H	38.8
030071037009	TA3	9H	43.8
030071037010	TA3	10H	48.8

Using 1.5mm Cortical Screw

P/N	Material	Spec.	Length
030071038003	TA3	3H	18.8
030071038004	TA3	4H	23.8
030071038005	TA3	5H	28.8
030071038006	TA3	6H	33.8
030071038007	TA3	7H	38.8

Using 1.5mm Cortical Screw

P/N	Material	Spec.	Length
030071039003	TA3	3H	19
030071039004	TA3	4H	24
030071039005	TA3	5H	29
030071039006	TA3	6H	34
030071039007	TA3	7H	39

Using 1.5mm Cortical Screw

P/N	Material	Spec.	Length
030071040004	TA3	4H	9

Using 1.5mm Cortical Screw



1.5 H Condylar Plate (8 holes)



P/N	Material	Spec.	Length
030071040008	TA3	8H, L	22
030071040108	TA3	8H, R	22

Using 1.5mm Cortical Screw

1.5 Condylar Mini Plate



P/N	Material	Spec.	Length
030071041003	TA3	3H, L	23
030071041004	TA3	4H, L	28
030071041005	TA3	5H, L	33
030071041103	TA3	3H, R	23
030071041104	TA3	4H, R	28
030071041105	TA3	5H, R	33

Using 1.5mm Cortical Screw

1.5 Y Condylar Plate



P/N	Material	Spec.	Length
030071042003	TA3	3H	20.9
030071042004	TA3	4H	25.9
030071042005	TA3	5H	30.9
030071042006	TA3	6H	35.9
030071042007	TA3	7H	40.9

Using 1.5mm Cortical Screw

1.5 T Condylar Plate



P/N	Material	Spec.	Length
030071043003	TA3	3H	18.8
030071043004	TA3	4H	23.5
030071043005	TA3	5H	28.8
030071043006	TA3	6H	33.8
030071043007	TA3	7H	38.8

Using 1.5mm Cortical Screw

2.0 H Condylar Plate (4 holes)



P/N	Material	Spec.	Length
030071044004	TA3	4H	12

Using 2.0mm Cortical Screw

2.0 Condylar Mini Plate



2.0 T Condylar Plate



2.0 Straight Plate



2.0 T Condylar Plate (2 holes head)



2.0 T Condylar Plate (3 holes head)



P/N	Material	Spec.	Length
030071045003	TA3	3H, L	24
030071045004	TA3	4H, L	29
030071045005	TA3	5H, L	34
030071045103	TA3	3H, R	24
030071045104	TA3	4H, R	29
030071045105	TA3	5H, R	34

Using 2.0mm Cortical Screw

P/N	Material	Spec.	Length
030071046003	TA3	3H	20
030071046004	TA3	4H	25
030071046005	TA3	5H	30
030071046006	TA3	6H	35
030071046007	TA3	7H	40
030071046008	TA3	8H	45
030071046009	TA3	9H	50

Using 2.0mm Cortical Screw

P/N	Material	Spec.	Length
030071047003	TA3	3H	15
030071047004	TA3	4H	20
030071047005	TA3	5H	25
030071047006	TA3	6H	30
030071047007	TA3	7H	35
030071047008	TA3	8H	40
030071047009	TA3	9H	45
030071047010	TA3	10H	50

Using 2.0mm Cortical Screw

P/N	Material	Spec.	Length
030071048003	TA3	3H	20
030071048004	TA3	4H	25
030071048005	TA3	5H	30
030071048006	TA3	6H	35
030071048007	TA3	7H	40
030071048008	TA3	8H	45

Using 2.0mm Cortical Screw

P/N	Material	Spec.	Length
030071049003	TA3	3H	20
030071049004	TA3	4H	25
030071049005	TA3	5H	30
030071049006	TA3	6H	35
030071049007	TA3	7H	40
030071049008	TA3	8H	45

Using 2.0mm Cortical Screw



2.0 T Condylar Plate (4 holes head)



P/N	Material	Spec.	Length
030071050003	TA3	3H	20
030071050004	TA3	4H	25
030071050005	TA3	5H	30
030071050006	TA3	6H	35
030071050007	TA3	7H	40
030071050008	TA3	8H	45
Using 2.0mm Cortical Screw			

2.0 LC-DCP Plate



P/N	Material	Spec.	Length
030071051003	TA3	3H	22.4
030071051004	TA3	4H	29.6
030071051005	TA3	5H	36.8
030071051006	TA3	6H	44
030071051007	TA3	7H	51.2
030071051008	TA3	8H	58.4
Using 2.0mm Cortical Screw			

2.5 Straight Plate



P/N	Material	Spec.	Length
030071052003	TA3	3H	18
030071052004	TA3	4H	24
030071052005	TA3	5H	30
030071052006	TA3	6H	36
030071052007	TA3	7H	42
030071052008	TA3	8H	48
Using 2.5mm Cortical Screw			

2.5 T Condylar Plate (2 holes head)



P/N	Material	Spec.	Length
030071053003	TA3	3H	24
030071053004	TA3	4H	30
030071053005	TA3	5H	36
030071053006	TA3	6H	42
030071053007	TA3	7H	48
Using 2.5mm Cortical Screw			

2.5 T Condylar Plate (3 holes head)



P/N	Material	Spec.	Length
030071054003	TA3	3H	24
030071054004	TA3	4H	30
030071054005	TA3	5H	36
030071054006	TA3	6H	42
030071054007	TA3	7H	48
Using 2.5mm Cortical Screw			

2.5 Y Condylar Plate



2.5 Condylod Plate



2.5 LC-DCP Plate



2.7 One-fourth Tubular Plate



2.7 Condylod Mini Plate



P/N	Material	Spec.	Length
030071055003	TA3	3H	26
030071055004	TA3	4H	32
030071055005	TA3	5H	38
030071055006	TA3	6H	44
030071055007	TA3	7H	50
Using 2.5mm Cortical Screw			

P/N	Material	Spec.	Length
030071056003	TA3	3H, L	24.6
030071056004	TA3	4H, L	32.6
030071056005	TA3	5H, L	40.6
030071056006	TA3	6H, L	48.6
030071056007	TA3	7H, L	56.6
030071056103	TA3	3H, R	24.6
030071056104	TA3	4H, R	32.6
030071056105	TA3	5H, R	40.6
030071056106	TA3	6H, R	48.6
030071056107	TA3	7H, R	56.6
Using 2.5mm Cortical Screw			

P/N	Material	Spec.	Length
030071057003	TA3	3H	24.6
030071057004	TA3	4H	32.6
030071057005	TA3	5H	40.6
030071057006	TA3	6H	48.6
030071057007	TA3	7H	56.6
030071057008	TA3	8H	64.6
Using 2.5mm Cortical Screw			

P/N	Material	Spec.	Length
030071058003	TA3	3H	24
030071058004	TA3	4H	32
030071058005	TA3	5H	40
030071058006	TA3	6H	48
030071058007	TA3	7H	56
030071058008	TA3	8H	64
Using 2.7mm Cortical Screw			

P/N	Material	Spec.	Length
030071059003	TA3	3H, L	32.4
030071059004	TA3	4H, L	40.4
030071059005	TA3	5H, L	48.4
030071059103	TA3	3H, R	32.4
030071059104	TA3	4H, R	40.4
030071059105	TA3	5H, R	48.4
Using 2.7mm Cortical Screw			