



# PRODUCT CATALOG

---

## 2023

**STRUmed Solutions Private Limited**

**Registered Office:** Plot No. 44, Ichhapore GIDC Industrial Estate, Bhatpore, Off Hazira Road, behind Tata Motors, Surat, Gujarat - 394510

**Division:** 1st Floor, BTCIF-06, No 5, TICEL Bio park Limited Ph-2, CSIR Road, Taramani, Chennai – 600113

**Contact: Tel:-** +91 6379454154 **Email:-** Customercare@STRUmed.com **Web:-**www.strumed.com



**CLOSED  
SYSTEM**

# UriQuick<sup>®</sup>

Urine Analyzer

ADVANCED, POINT-OF-CARE URINE ANALYZER



Introducing the

# Next-Gen Urine Analyzers

Indigenously Developed High-Performance Analyzer To Meet Your Urine Sample Testing With First In Class Features

Testing Range

Up to 14 Parameters



2.8" Colour Touch Screen

Throughput

Up to 120 Tests/Hr



WiFi & Bluetooth Enabled



Built In Thermal Printer



Wireless LIS



RFID Based Closed System

UriQuick Urine Analyzer



As Standard

## URINE ANALYZER TECHNICAL CHARACTERISTICS COMPARISON

FEATURES	STRUMED URIQUICK	MACHERY NAGEL URYXXON RELAX	SIEMENS CLINITEK STATUS +	URIT MEDICAL URIT-50
Reflectance Photometry Principle				
Complete Visible Spectrum Detection (400nm - 700nm)				
Wireless LIS				
Colour Display				
Unlimited Memory Storage via File Dump				
Internal Barcode Scanner				
Local Spare Parts Availability				
Replacement Warranty				
Low Cost Maintenance				
Indigenous Manufacturer				

Step into the Future with Wireless Connectivity

### STRUmed Solutions Private Limited

Registered Office: Plot No. 44, Ichhapore GIDC Industrial Estate, Bhatpore, Off Hazira Road, behind Tata Motors, Surat, Gujarat - 394510

Division: 1st Floor, BTCIF-06, No 5, TICEL Bio park Limited Ph-2, CSIR Road, Taramani, Chennai – 600113

Contact: Tel:- +91 6379454154 Email:- Customercare@STRUmed.com Web:-www.strumed.com