


# URINE ANALYZER



## SPECIFICATION AND FEATURES

<b>Features</b>	<b>Strumed Urine Analyzer</b>
<b>Principle</b>	Reflectance Colorimetry
<b>Colour Detection</b>	Complete Visible Spectrum (400nm - 700nm)
<b>Throughput</b>	Normal Mode: 65/Hour Continuous Mode: 130/Hour
<b>LIS</b>	Wired
<b>Display</b>	2.8" TFT Touch Screen Display with stylus
<b>Internal Thermal Printer</b>	Yes, 2" Thermal Printer
<b>Memory</b>	5000
<b>External Barcode Scanner</b>	Low Cost, Optional, USB Port
<b>Weight</b>	1.8 Kg
<b>Dimensions</b>	218x168x95 mm
<b>Parameters</b>	up to 14 Parameters
<b>Spare Parts &amp; Accessories</b>	Indigenous