


SEMI AUTOMATIC COAGULATION ANALYZER



FEATURES	SPECIFICATIONS
PRINCIPLE	TURBO-DENSITOMETRIC WITH OPTO-MECHANICAL SYSTEM
DISPLAY	2.8" TFT TOUCH SCREEN COLOR DISPLAY
TEST PARAMETERS	PROTHROMBIN TIME (PT), ACTIVATED THROMBOPLASTIN TIME (APTT), THROMBIN TIME (TT), FIBRINOGEN (FBG)
TEST SELECTION	RFID CARD BASED SELECTION WITH CLOSED SYSTEM
TESTING CHANNELS	1 CHANNEL
REAGENT PRE-WARMING CHANNELS	2 CHANNELS AT 37° C ± 1° C
CUVETTE HOLDING CHANNELS	20 CHANNELS AT 37° C ± 1° C
MEMORY	500 PATIENT RECORDS
BARCODE SCANNER	EXTERNAL BARCODE SCANNER
PRINTER	INTERNAL PRINTER – AUTOMATIC RESULT PRINT AVAILABLE
LIS CONNECTIVITY	WI-FI BASED WIRELESS LIS SUPPORT
POWER SUPPLY	INPUT: 150-285V AC, OUTPUT: 12V-3A
DIMENSION	300 X 300 X 80 MM (L X B X H)
WEIGHT	3 KG
MAKE IN INDIA	YES