

# uRIQUICK

Wireless | Compact | Advanced

Introducing  
the next-generation  
**Urine Analyzer**



**2.8" Colour  
Touch Screen**




**Wireless  
LIS**



**Auto  
Calibration**



**Closed  
System**



**Built In  
Thermal Printer**

# Step into the Future with Wireless Connectivity

## uRIQUICK Urine Analyzer

Testing Range  
Up to  
10 Parameters



RFID Based  
Closed  
System

Throughput  
Up to  
120 Tests/Hr



2.8" Colour  
Touch Screen



WiFi  
Enabled

**HL7**

Wireless LIS  
As per Standard

Urine Analyzer Technical Characteristics Comparison

FEATURES	STRUMED URIQUICK	EUROPEAN ANALYZERS	AMERICAN ANALYZERS	CHINESE ANALYZERS
Reflectance Photometry Principle				
Complete Visible Spectrum Detection (400nm - 700nm)				
Wireless LIS				
Colour Display				
Unlimited Memory Storage via File Dump				
Internal Barcode Scanner				
Local Spare Parts Availability				
Replacement Warranty				
Low Cost Maintenance				
Indigenous Manufacturer				

Uriquick Urine Strips - 100 Strips / Container

PARAMETERS	Leukocytes
	Nitrites
	Urobilinogen
	Protein
	pH
	Blood
	Specific Gravity
	Ketone
	Bilirubin
	Glucose