


iNSTA XPLOR P - PA001

NEPHELOMETER



Simplify your Screenings, Amplify your Accuracy



**Opto-
Mechanical
Technology**



Auto Calibration



**Internal
Incubation
System**



**PID Controlled
Heating System**

STRUMed Solutions Private Limited
Plot No. 44, Ichhapore GIDC Industrial Estate, Bhatpore,
Surat, Gujarat - 394510 | Contact:- +91 6379454154
Email:- customercare@strumed.com | Web:- www.strumed.com



iNSTA XPLOR P - PA001

NEPHELOMETER



FEATURES	SPECIFICATIONS
PRINCIPLE	TURBO-DENSITOMETRY AND NEPHELOMETRY WITH OPTO-MECHANICAL SYSTEM
TEST PARAMETERS	HbA1c, VITAMIN D, CYSTATIN C, ASO, RF, D-DIMER, hsCRP, CRP, MICROALBUMIN, IgG, IgA, IgE *UPCOMING - ANTI CCP, PCT, IgM, FERRITIN, C3, C4, ApoA-I, ApoB, ApoE
TEST SELECTION	RFID CARD BASED SELECTION WITH CLOSED SYSTEM
TESTING CHANNELS	1 CHANNEL
MEMORY	500 PATIENT RECORDS
LIS CONNECTIVITY	WI-FI BASED WIRELESS LIS SUPPORT
BARCODE SCANNER	EXTERNAL BARCODE SCANNER
PRINTER	INTERNAL PRINTER – AUTOMATIC RESULT PRINT AVAILABLE
POWER SUPPLY	INPUT: 150-240V, OUTPUT: 12V-3A
DISPLAY	2.8" TFT TOUCH SCREEN COLOR DISPLAY
DIMENSION	300 X 300 X 165 MM (L X B X H)
WEIGHT	3 KG

