



**iNSTA  
XPORT**  
(for COVID-19)

**iNSTA XPORT, STRUMed's Viral Transport Medium is suitable for collection, transport, maintenance and long-term freeze storage of clinical specimens containing viruses, including COVID-19.**

#### **Description:**

The InstaXport Viral Transport Medium comes in a plastic, screw cap tube and maintains organism viability for 48 hours at refrigerated temperature. The medium formulation includes antibiotics to inhibit Bacterial and Fungal growth. A phenol red pH indicator provides a visual gauge of the medium integrity. The product is suitable for samples collected using Sterile Nasal/ Throat swabs made of Polyester/ Nylon material.

#### **Product Advantages:**

- Phenol red indicator lends comfort to the users on the medium integrity throughout the storage and life of the product.
- The medium contains cryoprotectants to ensure viability of organisms through freezing & thawing.
- Can be stored at room temperature, and does not require refrigeration, prior to collection.
- Suitable for both Throat as well as Nasal Swab samples.

#### **Product Specification:**

- Medium Volume: 3ml
- Compatible swabs: Type - Throat or Nasal;  
Material - Polyester or Nylon
- Standards complied: CDC and WHO



#### **STRUMED SOLUTIONS PRIVATE LIMITED**

Unit - 1: No. 315, Luxuria Business Hub,  
Dumas Road, Surat - 395007, Gujarat, India.

Unit 2: 1st Floor, Module No. 06 - BTCIF,  
Ticel Bio Park - Phase 2, CSIR Road,  
Taramani, Chennai - 600113, Tamilnadu, India.