



**STRUMED**  
STRUCTURED INNOVATION

# iNSTAXPLOR<sup>®</sup>

COVID-19 ANTIGEN TEST KIT

Product code - STRUCOVAG



## PRODUCT HIGHLIGHTS

- High Sensitivity : 98%.
- High Specificity : 100%.
- Robust and highly Stable.
- Rapid results within 15 - 30 mins.
- Store at 2° C – 28° C.
- Longer Shelf life - 24 months.

## TECHNOLOGY

Lateral Flow

## SPECIMEN

Nasopharyngeal

## APPROVAL AND ACCREDITATION

