

# iNFLA QUICK

Erythrocyte Sedimentation Rate (ESR) Analyzer

**20 ESR RESULTS IN 30 MIN WITH MODIFIED  
WESTERGREN METHOD**



## TECHNICAL CHARACTERISTICS

Features	iNFLAQUICK ESR Analyzer
Principle	Infrared Photometry
Throughput	40 samples/hour. 20 samples random simultaneously measurement.
Accuracy	±1mm reading accuracy
Optical Reading Resolution	± 0.04mm Optical measurement precision
Temperature Compensation	Automatic temperature conversion to 18°C.(Manley table)
IOT Connectivity*	Yes, Cloud Based Closed / Subscription /Recharge & Open Model
LIS*	Wireless, Follows HL7 Protocol
Consumable Management System**	✓
Wifi-Connectivity	✓
Memory	1024 samples
Display	3.5” Colour TFT Touch Screen with stylus
Internal Thermal Printer	Optional, USB Port
External Barcode Scanner	Low Cost, Optional, USB Port
Weight	3.6 Kg
Dimensions	215 x 210 x 173 mm (L x W x H)
Power Supply	12Volt 2 Amp
Spare Parts & Accessories	Indigenous

\*Patent under progress. \*\* Minimum order quantity required.