


# GM Pro Peek Abutment



Temporary  
prosthesis

Biocompatible Peek of easy customization.



Consider in addition 1.5 - 2.0 mm for the restorative material

## ► Installation Sequence

GM Pro Peek Abutment <sup>1</sup>



	0.8 mm	1.5 mm	2.5 mm	3.5 mm	4.5 mm	5.5 mm
Ø 4.5	114.738	114.739	114.740	114.741	114.742	114.743
Ø 6.0	114.744	114.745	114.746	114.747	114.748	114.749



In mouth customization

1



Neo Screwdriver  
Torque Connection

+



Torque Wrench