

GM Universal Abutment



Single-unit
cement-retained
prosthesis

028



Consider in addition 1.5 - 2.0 mm for the restorative material
Minimum interocclusal space of 4.9 mm from the mucosa level

Installation Sequence

1
GM Exact Click
Universal Abutment



1
GM Exact Click
Universal Abutment 17°



1
GM Exact Click
Universal Abutment 30°



or

or

	0.8 mm	1.5 mm	2.5 mm	3.5 mm	4.5 mm	5.5 mm
4 mm	Ø 3.3	114.566	114.567	114.568	114.569	114.570
	Ø 4.5	114.578	114.579	114.580	114.581	114.582
6 mm	Ø 3.3	114.572	114.573	114.574	114.575	114.576
	Ø 4.5	114.584	114.585	114.586	114.587	114.588

	1.5 mm	2.5 mm	3.5 mm
4 mm	Ø 3.3	114.542	114.543
	Ø 4.5	114.548	114.549
6 mm	Ø 3.3	114.545	114.546
	Ø 4.5	114.551	114.552

	1.5 mm	2.5 mm	3.5 mm
4 mm	Ø 3.3	114.554	114.555
	Ø 4.5	114.560	114.561
6 mm	Ø 3.3	114.557	114.558
	Ø 4.5	114.563	114.564

Click Universal Abutment
Impression Coping



4 mm	Ø 3.3	108.172
	Ø 4.5	108.174
6 mm	Ø 3.3	108.173
	Ø 4.5	108.175

Click Universal Abutment
Provisional Coping



4 mm	Ø 3.3	118.304
	Ø 4.5	118.306
6 mm	Ø 3.3	118.305
	Ø 4.5	118.307

Universal Abutment Analog



4 mm	Ø 3.3	101.070	6 mm	Ø 3.3	101.071	Click (conventional)
	Ø 4.5	101.072		Ø 4.5	101.073	
4 mm	Ø 3.3	101.097	6 mm	Ø 3.3	101.098	Hybrid Repositionable (conventional/digital)
	Ø 4.5	101.099		Ø 4.5	101.100	

Universal Abutment
Burn-out Coping



4 mm	Ø 3.3	118.181
	Ø 4.5	118.183
6 mm	Ø 3.3	118.182
	Ø 4.5	118.184

1

Neo Screwdriver
Torque Connection

+

Torque Wrench