

Grand Morse™ Attachment TiN for Removable Prosthesis



Multiple-unit
overdenture
prosthesis



- Novaloc® compatible
- GM connection
- Straight and angled (15°)
- TiN coating
- Torque – 20 Ncm



Compatible with new
Neo Screwdriver

Articles

Description

102.148	GM Attachment TiN for Removable Prosthesis
102.149	GM Attachment TiN for Removable Prosthesis
102.150	GM Attachment TiN for Removable Prosthesis
102.151	GM Attachment TiN for Removable Prosthesis
102.152	GM Attachment TiN for Removable Prosthesis
102.153	GM Attachment TiN for Removable Prosthesis
102.154	GM Attachment TiN for Removable Prosthesis 15 Degrees
102.155	GM Attachment TiN for Removable Prosthesis 15 Degrees
102.156	GM Attachment TiN for Removable Prosthesis 15 Degrees
102.157	GM Attachment TiN for Removable Prosthesis 15 Degrees
102.158	GM Attachment TiN for Removable Prosthesis 15 Degrees
102.159	GM Attachment TiN for Removable Prosthesis 15 Degrees

Dimension

Material

0.8 mm	Titanium
1.5 mm	Titanium
2.5 mm	Titanium
3.5 mm	Titanium
4.5 mm	Titanium
5.5 mm	Titanium
0.8 mm	Titanium
1.5 mm	Titanium
2.5 mm	Titanium
3.5 mm	Titanium
4.5 mm	Titanium
5.5 mm	Titanium