



GM Titanium Base Angled Solution (AS)


With removable screw.



Single-unit screw-retained prosthesis

OR

Single-unit cement-retained prosthesis



034

Accessories

Replacement Sterile Screw



116.288

Screw for GM Titanium Base AS



Workflow Options

► Intraoral

► Model Scanning

GM Implant Intraoral Scanbody ²



108.183



GM Implant Analog



Ø 3.5/3.75	Ø 4.0/4.3	Ø 5.0/6.0	Hybrid Repositionable (conventional/digital)
101.103	101.089	101.090	



GM Implant Exact Impression Coping ² Closed and Open Tray



Regular	108.160	108.162	Regular
Long	108.161	108.163	Long



GM Implant Analog



Ø 3.5/3.75	Ø 4.0/4.3	Ø 5.0/6.0	Hybrid Repositionable (conventional/digital)
101.103	101.089	101.090	
	101.074	101.075	Conventional



GM Exact Implant Scanbody ²



108.181



GM Titanium Base Angled Solution (AS) ¹

			4 mm	20 N.cm	6 mm					
0.8 mm	1.5 mm	2.5 mm				0.8 mm	1.5 mm	2.5 mm	Ø 4.0	
Ø 4.0	135.327	135.328				135.329	135.330	135.331		135.332
Ø 4.5	135.333	135.334				135.335	135.336	135.337		135.338
Ø 5.5	135.339	135.340	135.341			135.342	135.343	135.344	Ø 5.5	

1

Short 105.150
Regular 105.151
Long 105.152

Angled Solution Screwdriver for Torque Wrench + Torque Wrench

OR

Short 105.147
Regular 105.148
Long 105.149

Angled Solution Screwdriver for Contra-angle + Contra-angle

2

Neo Screwdriver Torque Connection + Manual Screwdriver Torque