



# GM Titanium Base for Bridge

With removable screw.

Multiple-unit screw-retained prosthesis

OR

Multiple-unit cement-retained prosthesis

032



Cementable area: 4.0 mm for Ø 3.5  
4.5 mm for Ø 4.5 and Ø 5.5

## Accessories

Replacement Sterile  
Screws



Titanium 116.286



# Workflow Options

► Intraoral

► Model Scanning

GM Implant Intraoral Scanbody <sup>2</sup>



108.183



GM Implant Analog



Ø 3.5/3.75	Ø 4.0/4.3	Ø 5.0/6.0	Hybrid Repositionable (conventional/digital)
101.103	101.089	101.090	



GM Implant Impression Coping Open Tray (non-indexed) <sup>2</sup>



Regular 108.158 108.159 Long



GM Implant Analog



Ø 3.5/3.75	Ø 4.0/4.3	Ø 5.0/6.0	Hybrid Repositionable (conventional/digital)
101.103	101.089	101.090	
	101.074	101.075	Conventional



GM Exact Implant Scanbody <sup>2</sup>



108.181



GM Titanium Base for Bridge <sup>1</sup>

	0.8 mm	1.5 mm	2.5 mm	3.5 mm	4.5 mm
Ø 3.5	135.304	135.305	135.306	135.307	135.308
Ø 4.5	135.309	135.310	135.311	135.312	135.313
Ø 5.5	135.314	135.315	135.316	135.317	135.318



<sup>1</sup>

Neo Screwdriver Torque Connection + Torque Wrench

<sup>2</sup>

Neo Screwdriver Torque Connection + Manual Screwdriver Torque