

GM Mini Conical Abutment



071

Consider in addition 1.5 - 2.0 mm for the restorative material


Minimum interocclusal space of 4.5 mm from the mucosa level for straight abutments.

► Accessories

Mini Conical Abutment
Polishing Protector



Replacement
Coping Screw





Workflow Options

GM Mini Conical Abutment ¹

0.8 mm	1.5 mm	2.5 mm
115.243	115.244	115.245
3.5 mm	4.5 mm	5.5 mm
115.246	115.247	115.248



or



GM Exact Mini Conical Abutment 17°/30°/45°* ²

	1.5 mm	2.5 mm	3.5 mm
17°	115.249	115.250	115.251
30°	115.252	115.253	115.254
45°	115.267	115.268	

*The 45° Mini Conical Abutment is indicated for use only with Helix GM® Long and Zygoma GM™.

Intraoral ³

GM Mini Conical Abutment Scanbody



108.196

Model Scanning ³

Slim Mini Conical Abutment Open Tray Impression Coping



108.176

Conventional ³

Slim Mini Conical Abutment Open Tray Impression Coping ³

Slim Mini Conical Abutment Closed Tray Impression Coping ³



108.176



108.021

072

Mini Conical Abutment Hybrid Repositionable Analog



101.092

Mini Conical Abutment Hybrid Repositionable Analog



101.092

GM Mini Conical Abutment Scanbody ³



108.196

Neo Mini Conical Abutment One Step Hybrid Coping ²



118.330

Neo Mini Conical Abutment One Step Hybrid Coping ²



118.330

Neo Mini Conical Abutment Titanium Coping ²



118.302

Neo Mini Conical Abutment Protection Cylinder ³



106.220

Mini Conical Abutment Analog



101.092

101.020

Hybrid Repositionable (conventional/digital)
Conventional

Neo Mini Conical Abutment CoCr Coping ²



118.303

Neo Mini Conical Abutment Burn-out Coping ²



118.301

1

Hexagonal Prosthetic Driver + Torque Wrench

2

Neo Screwdriver Torque Connection + Torque Wrench


3

Neo Screwdriver Torque Connection + Manual Screwdriver Torque




Measurements GM Mini Conical Abutment

➤ 17°



➤ 30°



➤ 45°*



*The 45° Mini Conical Abutment is indicated for use only with Helix GM® Long and Zygoma GM™.