



Scan to download Signia App

Scan to Contact



Scan to Email



Find us on:



www.signia.net/en-in



www.youtube.com/signiahearing/



www.facebook.com/signiaindia/



www.instagram.com/signiahearingindia/

The information in this document contains general descriptions of the technical options available, which do not always have to be present in individual cases and are subject to change without prior notice. The required features should therefore be specified in each individual case at the time of conclusion of the respective contract.

Android and Google Play are trademarks of Google Inc. Apple App Store is a trademark of Apple Inc.



Be Brilliant™

Retail Price List

Effective from 20th November 2024





Be
Brilliant™

Signia is one of the leading hearing aid brands in the world. Our goal is to enhance human performance through iconic innovations that shape the prescription audio category. We refuse to accept hearing loss as a limitation. It's not about correcting a loss, it's about gaining an edge. Not just performing, but performing brilliantly. Not just hearing again, but hearing better than ever. Hearing every word, beat, and breath.

With Signia, you can Be Brilliant.

CONTENTS

1) Integrated Xperience (IX)	6
2) Orion C&G	16
3) Augmented Xperience (AX)	18
4) Xperience (X)	24
5) primax (px)	28
6) Other Hearing Aids	31
7) Accessories	33






























Terms & Conditions



1. Signia Hearing Aids carry warranty as indicated in the price list.
2. Kit Includes 2 units of hearing aids and a standard charger.
3. Accessories (Chargers, miniReceiver, Remote Control, Streamers, CROS / BiCROS) carry 1 year warranty.
4. Power adaptor for some accessories will have to be procured separately.
5. Custom hearing aids prices include shelling charges.
6. HP miniReceiver customization charges are included in the RIC pricing.
7. RIC / BTE instruments will be supplied in standard beige colour. Other colour housing options for RIC / BTE models are subject to a delivery lead time of 3 weeks.
8. Prices of color conversion kit may vary depending on hearing aid model.
9. Hearing Aids are exempted from GST.
10. GST applicable on accessories. Retail prices of accessories are inclusive of GST.
11. The pricing and product specifications are subject to change without prior notice.
12. Chargers for selective models will have to be purchased separately.

Sivantos India recommends all its valued consumers to seek a hearing care professional's opinion on the best suited hearing instrument. Online purchase may not necessarily give you the desired results.

Iconic innovations that empower people to hear what matters

Signia Hearing Aid Portfolio

	RIC				BTE			Instant Fit			Custom			
profound					BTE Li-ion	BTE		BTE						
severe	RIC Li-ion	RIC		RIC	Motion C&G Motion C&G P Motion C&G SP	Intuis 4 M Intuis 4 P Intuis 4 SP		Orion C&G Orion C&G P Orion C&G SP				Custom Li-ion	Custom	Custom
moderate to severe	Pure C&G Pure C&G T	Pure 312	Slim RIC Li-ion	Orion C&G								Insio ITC/ITE	Insio IIC/CIC ¹	Insio ITC/ITE ²
moderate			Styletto						ITE IF Li-ion Active ITC	CIC IF Micro Li-ion Silk C&G CIC ¹	CIC IF Micro Silk CIC ¹			
mild														
	 		 		 				 			 		
7	Pure Charge&Go IX	Pure 312 AX	Styletto IX /AX		Motion Charge&Go X	Intuis 4.7			Active Pro IX	Silk Charge&Go IX	Silk X	Insio Charge&Go AX	Insio IX	Insio Nx
5	Pure Charge&Go IX	Pure 312 AX	Styletto IX /AX		Motion Charge&Go X	Intuis 4.5				Silk Charge&Go IX	Silk X	Insio Charge&Go AX	Insio IX	Insio Nx
3	Pure Charge&Go IX	Pure 312 AX	Styletto IX /AX		Motion Charge&Go X	Intuis 4.3				Silk Charge&Go IX	Silk X	Insio Charge&Go AX	Insio IX	Insio Nx
2	Pure Charge&Go IX	Pure 312 AX	Styletto IX /AX	Orion Charge&Go	Motion Charge&Go X	Intuis 4.2		Orion Charge&Go	Active IX	Silk Charge&Go IX	Silk X	Insio Charge&Go AX	Insio IX	Insio Nx
1	Pure Charge&Go IX	Pure 312 AX	Styletto IX /AX	Orion Charge&Go	Motion Charge&Go X	Intuis 4.1		Orion Charge&Go		Silk Charge&Go IX	Silk X	Insio Charge&Go AX	Insio IX	Insio Nx

 Rechargeable
 Bluetooth

¹ Non-Bluetooth
² Bluetooth available in "top-up" model

Styletto IX, Pure Charge&Go IX, and Active Pro IX features

Features	Description	Performance levels						
		Styletto IX and Pure Charge&Go IX					Active Pro IX	
Conversation management		7 IX	5 IX	3 IX	2 IX	1 IX	Active Pro	Active
RealTime Conversation Enhancement	Provides uninterrupted engagement in noisy group conversations by seamlessly tracking the locations of all conversation partners and enhancing their speech.*	●●●●●	●●●●○	●●○○○			●●●●●	
Conversation Dynamics Analyzer	Analyzes conversation dynamics including turn-taking and timing, locations of pure noise without speech, and levels of background noise.	✓					✓	
Conversation Source Analyzer	Constantly scouts the acoustic environment to identify number of talkers and their location.	✓	✓	✓			✓	
Dynamic Focus Streams	Isolate and track speech independently.	3	1				3	
Conversation Booster	Dynamically enhances individual speakers.	●●●●●	●●●●○	●●○○○			●●●●●	
Augmented Focus	Separates the focus sound (speech) from surrounding sounds, creating a strong contrast.*	✓	✓	✓	✓	✓	✓	✓
Focus Stream	Enhances speech signals from the front hemisphere.	1	1	1	1	1	1	1
Surrounding Stream	Processes environmental sounds to preserve environmental awareness without masking speech.	1	1	1	1	1	1	1
Dynamic Soundscape Processing 3.0	Integrates inputs for precise hearing aid optimization, including (when available): acoustic sensors, motion sensors, OVP 2.0, noise detection, SNR analysis, and RealTime Conversation Enhancement.	●●●●●	●●●●○	●●●○○	●●○○○	●●○○○	●●●●●	●●○○○
Spatial SpeechFocus	Steers directional microphone beam to the dominant speech source (left, right, or behind the wearer automatically in a car situation or in a dedicated "Strolling" program).*	✓					✓	
SpeechFocus	Depending on the direction of the dominant speech source, Speech Focus steers the directional beam behind the wearer, or engages a surround-sound mode for speech from either side. Available for fittings with directional microphones*	✓	✓				✓	
Own Voice Processing 2.0 (OVP™ 2.0)	Uses real-time recognition of the wearer's voice to deliver a natural own-voice impression.*	✓	✓	✓				
Motion sensor	Detects whether the wearer is walking, running, or sitting still to guide how sound is optimized.	✓	✓	✓			✓	
Binaural Directionality	Provides scientifically proven enhanced speech understanding in noise.*	✓	✓	✓			✓	
Adaptive Directionality	Improves the signal to noise ratio in noisy environments by reducing dominant noise sources from the back and sides. The polar plot is adjusted continuously in multiple frequency channels.	✓	✓	✓	✓		✓	
Automatic Directionality	Improves the signal to noise ratio in noisy environments by automatically and seamlessly applying directionality in individual processing channels.	✓	✓	✓	✓	✓	✓	✓
TruEar 360	Reproduces benefits of the natural pinna effect by highlighting the difference between front and back sources of sound.	✓	✓	✓	✓	✓	✓	✓

*) Available for bilateral fittings with e2e wireless link

Signia Integrated Xperience, featuring RealTime Conversation Enhancement, is the heart of our innovative hearing aids. This technology empowers wearers to engage dynamically in conversations, ensuring every word is understood, every laugh appreciated, and every story fully experienced.

Features	Description	Performance levels						
		Styletto IX and Pure Charge&Go IX					Active Pro IX	
		7 IX	5 IX	3 IX	2 IX	1 IX	Active Pro	Active
SpeechStabilizer	Preserves the natural envelope of speech in any situation, including very noisy situations, which preserves acoustic cues that are important for speech understanding.	✓	✓	✓	✓	✓	✓	✓
Speech and Noise Management	Reduces listening effort using sophisticated sound processing.	✓	✓	✓	✓	✓	✓	✓
e2e wireless 4.0	Controls coupling, synchronization, CROS/BiCROS, and a range of sound features.*	✓	✓	✓	✓	✓	✓	✓
Listening Comfort								
Extended bandwidth	12 kHz available bandwidth for enhanced processing of high-frequency speech and environmental sounds.	✓					✓	
Auto EchoShield	Automatically reduces the effect of reverberations and sudden loud noises in reverberant environments.	✓					✓	
Reverberant Room Program	Dedicated program for reverberant environments.	✓	✓				✓	
SoundSmoothing™	Reduces the annoyance of sudden loud impact sound such as when handling glassware.	●●●●●	●●●●●	●●●●●	●●●○○	●○○○○	●●●●●	●○○○○
eWindScreen™	Reduces the interference of wind noise when outdoors.	●●●●●	●●●●●	●●●○○	●○○○○	●○○○○	●●●●●	●○○○○
Extended dynamic range	117 dB 24 bit processing for enhanced digital encoding of soft and loud sounds.	✓	✓	✓	✓	✓	✓	✓
Feedback Cancellation	High speed monitoring and control of feedback in individual processing channels.	✓	✓	✓	✓	✓	✓	✓
Wearer Empowerment								
Spatial Configurator	Allows the wearer to adjust the direction of the microphone beam via the Signia App.	✓	✓				✓	
My WellBeing	Accessed via the Signia App, wearers can get insights about how much they are using their hearing aids, their level of social interaction, their daily step count, and their overall activity level.	✓	✓	✓			✓	
Signia Assistant	Offers wearers instant 24/7 troubleshooting and customized sound adjustments via an intuitive AI-enabled chat. Provides detailed information in Connexx to address each client's fine-tuning needs more precisely.	✓	✓	✓	✓	✓	✓	✓
Signia App (iOS and Android)	Easy-to-use smartphone app enables remote control adjustments, audio streaming, access to My WellBeing, and support via Signia Assistant and TeleCare.	✓	✓	✓	✓	✓	✓	✓
Adaptive Streaming Volume	Speech intelligibility weighted adaptation of the streaming volume depending on the level of ambient noise.	✓	✓	✓	✓	✓	✓	✓
TeleCare	Provides remote services including the TeleCare Dashboard, Remote Fine Tuning, and CareChat (convenient text messages, phone calls, and video calls).	✓	✓	✓	✓	✓	✓	✓

*) Available for bilateral fittings with e2e wireless link

Styletto IX, Pure Charge&Go IX, and Active Pro IX features

Features	Description	Performance levels						
		Styletto IX and Pure Charge&Go IX					Active Pro IX	
		7 IX	5 IX	3 IX	2 IX	1 IX	Active Pro	Active
Connectivity								
Direct streaming for iOS (MFi) & Android (LE Audio and ASHA)	Hearing aids can be paired to Apple iOS devices (iPhones, iPads, iPods) and Android devices using ASHA or LE Audio for direct streaming of audio such as phone calls and music.	✓	✓	✓	✓	✓	✓	✓
HandsFree	Enables truly hands-free calling. The wearer's voice is picked up by the hearing aid microphone and delivered in high quality to the conversation partner, while the call audio is streamed directly to the wearer's hearing aids. Available for connections to iOS devices or to Android devices using LE Audio.	✓	✓	✓	✓	✓	✓	✓
CallControl	Calls can be answered, ended, and ignored with a simple press of the hearing aid button. Available for connections to iOS devices or to Android devices using LE Audio.	✓	✓	✓	✓	✓	✓	✓
Telecoil (Pure Charge&Go T IX only)	Enables direct connection to audio loop systems and compatible telephone signals.	✓	✓	✓	✓	✓		
StreamLine Mic compatibility	Stream high quality stereo sound directly from any Bluetooth device (PCs, older Apple devices, Android devices without ASHA or LE Audio, etc.). Can also be used as a remote mic.	✓	✓	✓	✓	✓	✓	✓
StreamLine TV compatibility	Stream high quality stereo sound directly from TV to hearing aids via Bluetooth.	✓	✓	✓	✓	✓	✓	✓
miniPocket compatibility	Optional keychain remote control.	✓	✓	✓	✓	✓	✓	✓
Tinnitus								
Ocean Waves Therapy Signal (presets)	Therapy signals that mimic the sound of the sea to provide a soothing listening experience and reduce the impact of tinnitus.	4	4	4	4		4	
Static Therapy Signal (channels / presets)	Soft, gentle therapy signals that reduce the impact of tinnitus.	20/5	18/5	16/5	12/1		20/5	
Notch Therapy	Signia's unique, patented tinnitus therapy which is proven especially effective against tonal tinnitus.	✓	✓	✓	✓		✓	
Tinnitus Noise Therapy Manager	A choice of tinnitus treatments based on traditional noise therapy signals, amplified notch therapy, or a combination of the two approaches.	✓	✓	✓	✓		✓	

*) Available for bilateral fittings with e2e wireless link

Features	Description	Performance levels						
		Styletto IX and Pure Charge&Go IX					Active Pro IX	
		7 IX	5 IX	3 IX	2 IX	1 IX	Active Pro	Active
Fitting Optimization								
Processing Channels / Gain & MPO handles	Number of processing channels/handles in Connexx	48/20	36/18	32/16	24/12	16/8	48/20	16/8
HD Music	Preset optimized music programs for enhanced enjoyment of non-streamed music.	3	3	1	1		3	
Hearing programs	Number of programs configurable in Connexx.	6	6	6	4	4	6	4
CROS/BICROS compatibility	A CROS device worn on an ear with unaidable hearing loss transmits sound to a hearing aid worn on the better ear.	✓	✓	✓	✓	✓	✓	✓
Frequency Compression	Improves the detection of speech sounds for wearers with reduced high frequency hearing thresholds. Frequency compression is selectively activated based on audiological criteria.	✓	✓	✓	✓	✓	✓	✓
Smart Optimizer and Data Logging	Tailored fine-tuning recommendations based on logging and analysis of the wearer's hearing aid use.	✓	✓	✓	✓	✓	✓	✓
Acclimatization Manager	Allows the wearer to benefit from enhanced audibility in the longer term, while enjoying optimal device acceptance when initially fitted.	✓	✓	✓	✓	✓	✓	✓
InSituGram	Allows quick and convenient measurement of hearing thresholds with the hearing aids in place.	✓	✓	✓	✓	✓	✓	✓
AutoFit™	Automatically provides an optimized real-ear insertion gain (REIG) setting without the need to manually match targets, enabling significantly faster real-ear verification.	✓	✓	✓	✓	✓	✓	✓

*) Available for bilateral fittings with e2e wireless link

Silk Charge&Go IX and Insio IX IIC/CIC features

Features	Description	Performance levels					Silk Charge&Go IX not available in 2IX and 1IX
Conversation management		7 IX	5 IX	3 IX	2 IX	1 IX	
Dynamic Soundscape Processing 3.0	Integrates inputs for precise hearing aid optimization, including (when available): acoustic sensors, noise detection and SNR analysis.	●●●●●	●●●●○	●●●○○	●●○○○	●●○○○	
Binaural OneMic Directionality 2.0	Microphones from each hearing aid work together to create a directional listening pattern preferentially amplifying speech from the front to support hearing in noisy environments.*	✓	✓	✓			
SpeechStabilizer	Preserves the natural envelope of speech in any situation, including very noisy situations, which preserves acoustic cues that are important for speech understanding.	✓	✓	✓	✓	✓	
Speech and Noise Management	Reduces listening effort using sophisticated sound processing.	✓	✓	✓	✓	✓	
e2e wireless 4.0	Controls coupling, synchronization, CROS/BiCROS, and a range of sound features.*	✓	✓	✓	✓	✓	
Listening Comfort							
Extended bandwidth	12 kHz available bandwidth for enhanced processing of high-frequency speech and environmental sounds.	✓					
Auto EchoShield	Automatically reduces the effect of reverberations and sudden loud noises in reverberant environments.	✓					
Reverberant Room Program	Dedicated program for reverberant environments.	✓	✓				
SoundSmoothing™	Reduces the annoyance of sudden loud impact sound such as when handling glassware.	●●●●●	●●●●●	●●●●●	●●●○○	●○○○○	
eWindScreen™	Reduces the interference of wind noise when outdoors.	●●●●●	●●●●●	●●●○○	●○○○○	●○○○○	
Extended dynamic range	117 dB 24 bit processing for enhanced digital encoding of soft and loud sounds.	✓	✓	✓	✓	✓	
Feedback Cancellation	High speed monitoring and control of feedback in individual processing channels.	✓	✓	✓	✓	✓	
Wearer Empowerment & Connectivity							
Signia App (iOS and Android)	Easy-to-use smartphone app enables remote control adjustments and TeleCare.	✓	✓	✓	✓	✓	
TeleCare	Provides remote services including the TeleCare Dashboard and CareChat (convenient text messages, phone calls, and video calls).	✓	✓	✓	✓	✓	
miniPocket compatibility	Optional keychain remote control.	✓	✓	✓	✓	✓	

*) Available for bilateral fittings with e2e wireless link

Features	Description	Performance levels					Silk Charge&Go IX not available in 2IX and 1IX
		7 IX	5 IX	3 IX	2 IX	1 IX	
Tinnitus							
Ocean Waves Therapy Signal (presets)	Therapy signals that mimic the sound of the sea to provide a soothing listening experience and reduce the impact of tinnitus.	4	4	4	4		
Static Therapy Signal (channels / presets)	Soft, gentle therapy signals that reduce the impact of tinnitus.	20/5	18/5	16/5	12/1		
Notch Therapy	Signia's unique, patented tinnitus therapy which is proven especially effective against tonal tinnitus.	✓	✓	✓	✓		
Tinnitus Noise Therapy Manager	A choice of tinnitus treatments based on traditional noise therapy signals, amplified notch therapy, or a combination of the two approaches.	✓	✓	✓	✓		
Fitting Optimization							
Processing Channels / Gain & MPO handles	Number of processing channels/handles in Connexx	48/20	36/18	32/16	24/12	16/8	
HD Music	Preset optimized music programs for enhanced enjoyment of non-streamed music.	3	3	1	1		
Hearing programs	Number of programs configurable in Connexx.	6	6	6	4	4	
CROS/BiCROS compatibility	A CROS device worn on an ear with unaidable hearing loss transmits sound to a hearing aid worn on the better ear.**	✓	✓	✓	✓	✓	
Frequency Compression	Improves the detection of speech sounds for wearers with reduced high frequency hearing thresholds. Frequency compression is selectively activated based on audiological criteria.	✓	✓	✓	✓	✓	
Smart Optimizer and Data Logging	Tailored fine-tuning recommendations based on logging and analysis of the wearer's hearing aid use.	✓	✓	✓	✓	✓	
Acclimatization Manager	Allows the wearer to benefit from enhanced audibility in the longer term, while enjoying optimal device acceptance when initially fitted.	✓	✓	✓	✓	✓	
InSituGram	Allows quick and convenient measurement of hearing thresholds with the hearing aids in place.	✓	✓	✓	✓	✓	
AutoFit™	Automatically provides an optimized real-ear insertion gain (REIG) setting without the need to manually match targets, enabling significantly faster real-ear verification.	✓	✓	✓	✓	✓	

*) Available for bilateral fittings with e2e wireless link
*) **) Requires a hearing aid with e2e wireless link

IX Power the conversation
with Signia Integrated Xperience



Product Name	Form Factor	Channels	Fitting Range (dB)	Warranty (Years)	Unit Price (₹)
INTEGRATED XPERIENCE					

KIT ACTIVE IX

Kit Active Pro IX ^{NEW}	Instant Fit ITE	48	10-80	4	2,99,990
Kit Active IX ^{NEW}	Instant Fit ITE	16	10-80	2	1,49,990

KIT STYLETTO IX

Kit Styletto 7IX	Slim-RIC	48	10 - 110	4	6,99,990
Kit Styletto 5IX	Slim-RIC	36	10 - 110	4	4,69,990
Kit Styletto 3IX	Slim-RIC	32	10 - 110	2	3,40,990
Kit Styletto 2IX ^{NEW}	Slim-RIC	24	10 - 110	2	1,49,990
Kit Styletto 1IX ^{NEW}	Slim-RIC	16	10 - 110	2	1,19,990

KIT PURE C&G IX / KIT PURE C&G T IX

Kit Pure C&G 7IX / T 7IX	RIC	48	10 - 115	4	7,45,990
Kit Pure C&G 5IX / T 5IX	RIC	36	10 - 115	4	4,63,990
Kit Pure C&G 3IX / T 3IX	RIC	32	10 - 115	2	3,35,990
Kit Pure C&G 2IX / T 2IX	RIC	24	10 - 115	2	1,59,990
Kit Pure C&G 1IX / T 1IX	RIC	16	10 - 115	2	1,25,990

PURE C&G IX / PURE C&G T IX

Pure C&G 7IX / T 7IX	RIC	48	10 - 115	4	3,79,990
Pure C&G 5IX / T 5IX	RIC	36	10 - 115	4	2,39,990
Pure C&G 3IX / T 3IX	RIC	32	10 - 115	2	1,75,990
Pure C&G 2IX / T 2IX	RIC	24	10 - 115	2	85,990
Pure C&G 1IX / T 1IX	RIC	16	10 - 115	2	67,990

Note:
The fitting range is completely dependent on the receiver matrix of the hearing aid.

IX Power the conversation
with Signia Integrated Xperience



Product Name	Form Factor	Channels	Fitting Range (dB)	Warranty (Years)	Unit Price (₹)
INTEGRATED XPERIENCE					

KIT SILK C&G IX

Kit Silk C&G 7IX	Instant Fit CIC	48	10 - 115	4	6,89,990
Kit Silk C&G 5IX	Instant Fit CIC	36	10 - 115	4	4,12,990
Kit Silk C&G 3IX	Instant Fit CIC	32	10 - 115	2	2,99,990
Kit Silk C&G 2IX ^{NEW}	Instant Fit CIC	24	10 - 115	2	1,45,990
Kit Silk C&G 1IX ^{NEW}	Instant Fit CIC	16	10 - 115	2	1,14,990

INSIO IX

Insio 7IX	CIC/IIC	48	10 - 110	4	3,09,990
Insio 5IX	CIC/IIC	36	10 - 110	4	1,84,990
Insio 3IX	CIC/IIC	32	10 - 110	2	1,09,990
Insio 2IX	CIC/IIC	24	10 - 110	2	62,990
Insio 1IX	CIC/IIC	16	10 - 110	2	45,990



Active Pro IX



Styletto IX









Pure Charge&Go IX











Silk Charge&Go IX



Insio IX

Product Name	Form Factor	Battery Size/Type	Warranty (Years)	Unit Price (₹)
CROS/BiCROS				
Kit Active Pro CROS IX 	Instant Fit ITE	Li-ion	4 yrs for Hearing Aid 1 yr for CROS	1,99,990
Kit Active CROS IX 	Instant Fit ITE	Li-ion	2 yrs for Hearing Aid 1 yr for CROS	1,24,990
Kit Styletto CROS 7IX	Slim-RIC	Li-ion	4 yrs for Hearing Aid 1 yr for CROS	3,89,990
Kit Styletto CROS 5IX	Slim-RIC	Li-ion	4 yrs for Hearing Aid 1 yr for CROS	2,74,990
Kit Styletto CROS 3IX	Slim-RIC	Li-ion	2 yrs for Hearing Aid 1 yr for CROS	2,10,990
Kit Styletto CROS 2IX 	Slim-RIC	Li-ion	2 yrs for Hearing Aid 1 yr for CROS	1,14,990
Kit Styletto CROS 1IX 	Slim-RIC	Li-ion	2 yrs for Hearing Aid 1 yr for CROS	99,990
Kit Silk C&G CROS 7IX	Instant Fit CIC	Li-ion	4 yrs for Hearing Aid 1 yr for CROS	3,94,990
Kit Silk C&G CROS 5IX	Instant Fit CIC	Li-ion	4 yrs for Hearing Aid 1 yr for CROS	2,55,990
Kit Silk C&G CROS 3IX	Instant Fit CIC	Li-ion	2 yrs for Hearing Aid 1 yr for CROS	1,99,990
Kit Silk C&G CROS 2IX 	Instant Fit CIC	Li-ion	2 yrs for Hearing Aid 1 yr for CROS	1,22,990
Kit Silk C&G CROS 1IX 	Instant Fit CIC	Li-ion	2 yrs for Hearing Aid 1 yr for CROS	1,07,990
CROS / BiCROS (Pure C&G IX Transmitter)	RIC	Li-ion	1	35,990



Product Name		Warranty (Year)	Unit Price (₹)
CHARGERS			
Active Pro IX Charger		1	11,990
Styletto IX Charger Case		1	11,990
Styletto Dry&Clean Charger		1	11,990
Pure Dry&Clean Charger		1	11,990
Silk Charger		1	11,990
Pure Portable Charger		1	7,990
Pure Charger		1	5,990
Multi-Charger		1	3,990

Orion Charge&Go feature overview

Features	Performance levels			
	200	100	75	50
Frequency channels / compression controls	16/8	16/8	12/6	8/4
Number of programs	4	4	4	4
48-channel first fit	●	●	●	●
Adaptive compression	●	●	●	●
Extended dynamic range	●	●	●	●
Signia app	●	●	●	●
miniPocket compatability	●	●	●	●
InSituGram	●	●	●	●
AutoFit	●	●	●	●
Smart optimizer and Data logging	●	●	●	●
Feedback cancellation	●	●	●	●
Speech and noise management	●	●	●	●
Microphone directionality	●	●	●	●
Automatic directionality	●	●	●	-
Adaptive directionality	●	-	-	-
TruEar	●	●	-	-
SoundBalance	●	●	-	-
Frequency compression	●	●	-	-
Automatic acclimitization manager	●	●	-	-
e2e wireless	●	●	-	-
User control coupling via e2e	●	●	-	-
SoundSmoothing™	●	-	-	-
eWindScreen™	●	-	-	-
Tinnitus static therapy signal	●	-	-	-

Effortless to use

- Rechargeability for all
Your clients can say goodbye to the hassle of small batteries and hello to rechargeability.
- Wireless fitting
Cumbersome cables are a thing of the past with wireless technology that streamlines the fitting process.
- First fit precision
Our 48 channel first fit calculation. ensures a precise and personalized fit from the start, reducing the need for fine-tuning and optimizing the client experience right from the start.

Elevated design

- Contemporary design
With an attractive premium finish, Orion Charge&Go's design is as sophisticated as its technology
- IP68 rating
Certified to resist dust and sweat, Orion Charge&Go provides dependable performance in all conditions.

Enhanced sound processing

- Clarity in quiet
Our TruEar technology emulates the effect of the natural ear shape by prioritizing speech frequencies from the front for optimal hearing support in quiet environments.
- Focus in noise
Adaptive directionality zeros in on the speech your client wants to hear in noisy environments.
- Versatility across soundscapes
Adaptive compression and an extended input dynamic range work in tandem to adjust amplification to the natural flow of sounds, ensuring impeccable sound quality everywhere.

Charge your expectations



Product Name	Form Factor	Channels	Fitting Range (dB)	Warranty (Years)	Unit Price (₹)
ORION C&G P/SP BTE					
Orion C&G P/SP 200	BTE	16	10 - 120	2	38,990
Orion C&G P/SP 100	BTE	16	10 - 120	2	36,990
Orion C&G P/SP 75	BTE	12	10 - 120	2	32,990
Orion C&G P/SP 50	BTE	8	10 - 120	2	28,990
ORION C&G P/SP T BTE					
Orion C&G P/SP T 200	BTE	16	10 - 120	2	38,990
Orion C&G P/SP T 100	BTE	16	10 - 120	2	36,990
ORION C&G RIC					
Orion C&G 200	RIC	16	20 - 115	2	38,990
Orion C&G 100	RIC	16	20 - 115	2	36,990
Orion C&G 75	RIC	12	20 - 115	2	32,990
Orion C&G 50	RIC	8	20 - 115	2	28,990

Product Name		Warranty (Year)	Unit Price (₹)
CHARGERS			
Multi-Charger		1	3,990

Note:
Chargers for Orion C&G will have to be purchased separately.
The fitting range is completely dependent on the receiver matrix of the hearing aid.

AX

Platform Features

Features	Information	Performance levels				
		7AX	5AX	3AX	2AX	1AX
Dynamic Soundscape Processing 2.0	Dynamic Soundscape Processing takes sound processing to a more precise level by including motion data for a more realistic sound adapted to everyday life.
Augmented Focus™	Separates the focus sound (speech) from surroundings, creating a clear contrast between the two, which makes it easier to follow conversations.	✓	✓	✓	✓	✓
Acoustic Sensors	The redesigned acoustic sensors can detect more variables of sound than ever before, providing more accurate information about the environment.	✓	✓	✓	✓	✓
e2e wireless 4.0	Control coupling, synchronization, binaural directionality, eWindscreen binaural, Spatial Speech Focus, CROS/ BiCROS.	✓	✓	✓	✓	✓
Motion Sensor	An integrated motion sensor detects whether the wearer is walking, running or sitting still. This motion data complements the acoustic data to optimize Dynamic Soundscape Processing 2.0.	✓	✓	✓		
OVP™ (Own Voice Processing)	Own Voice Processing (OVP™) utilizes real-time recognition of the wearer's voice to deliver a natural own-voice impression. Available for RIC/BTE devices with e2e wireless in binaural fittings.	✓	✓	✓		
Sound Clarity	Sound Clarity provides lifelike sound which covers the full spectrum of wearing situations.
Signal processing / Gain & MPO	Channels/handles in Connexx. For Signia AX, channels are per processing chain (Focus/Surroundings).	48/20	32/16	24/12	16/8	16/8
Hearing Programs	Number of programs configurable in Connexx.	6	6	6	4	4
Extended dynamic range	117 dB 24 bit processing for enhanced digital encoding of soft and loud sounds.	✓	✓	✓	✓	✓
Speech and noise management	Reduces listening effort using sophisticated sound processing	✓	✓	✓	✓	✓
SoundSmoothing™	Reduces the annoyance of sudden loud impact sound such as when handling glassware.	✓	✓	✓	✓	
Feedback cancellation	High speed monitoring and control of feedback in individual processing channels.	✓	✓	✓	✓	✓
HD Music	Preset programs for enhanced enjoyment of non-streamed music.	3	3	1	1	
eWindScreen	Reduces the annoyance of wind noise when outdoors.	✓	✓	✓	✓	
Extended bandwidth	12 kHz bandwidth for enhanced processing of high-frequency speech and environmental sounds.	✓				
EchoShield	Dedicated program for reverberant environments.	✓	✓			
Speech Quality	Speech Quality delivers improved speech understanding and communication ease in noise.
Binaural Directionality	Directionality provides scientifically proven enhanced speech understanding in noise. Available for bilateral fittings with directional microphones and e2e wireless link.	✓	✓	✓		

1) In Universal program: Automatic mic beam to the right / left / rear directions.
2) In Universal program: Automatic mic beam to the rear, surround mode for side speech.
In manual programs, microphone beam to the right / left / rear directions.

Features	Information	Performance levels				
		7AX	5AX	3AX	2AX	1AX
Wireless CROS/BiCROS	Requires CROS transmitters.	✓	✓	✓	✓	✓
Frequency compression	Improves the detection of speech sounds for wearers with reduced high frequency hearing thresholds. Frequency compression is selectively activated based on audiological criteria.	✓	✓	✓	✓	✓
Spatial SpeechFocus	Depending on the direction of the dominant speech source, Spatial SpeechFocus steers the directional beam to the left, right, or behind the wearer automatically in a car situation, or in a dedicated "Strolling" program. Available for bilateral fittings with directional microphones and e2e wireless link.	✓ ¹⁾	✓ ²⁾			
Wearer Interaction	Technically inclined wearers may appreciate having extra control over their hearing aids. Signia Augmented Xperience hearing aids offer these wearers a comprehensive range of adjustment options via the Signia App.
Signia Assistant¹	Signia Assistant offers wearers instant 24/7 support and troubleshooting around the clock. It provides detailed information in Connexx to address each client's fine-tuning needs more precisely	✓	✓	✓	✓	✓
Signia app (iOS and Android)	The new Signia app combines all functionalities of previous apps (touchControl, myControl, myHearing app). The functionality offered in the app is automatically customized depending on the connected hearing aid.	✓	✓	✓	✓	✓
Adaptive Streaming Volume	Speech intelligibility weighted adaptation of the streaming volume depending on ambient noise.	✓	✓	✓	✓	✓
Spatial Configurator	Spatial Configurator is part of the Signia app and allows the wearer to adjust the direction of the microphone beam.	✓	✓			
Direct streaming for iOS	Signia Augmented Xperience hearing aids can be connected directly to Apple devices for phone calls and audio streaming – without additional accessories.	✓	✓	✓	✓	✓
Direct streaming for Android	Direct audio streaming from ASHA compatible Android devices to Bluetooth enabled Augmented Xperience hearing aids.	✓	✓	✓	✓	✓
StreamLine TV accessory	Stream high quality stereo sound directly from a TV to Bluetooth equipped Augmented Xperience hearing aids.	✓	✓	✓	✓	✓
StreamLine Mic accessory	Stream high quality stereo sound directly from any Bluetooth device (Android devices, older Apple devices, other smartphones, PCs, laptops,etc.) to Bluetooth equipped Augmented Xperience hearing aids.	✓	✓	✓	✓	✓
Tinnitus	A choice of tinnitus treatments based on amplified notch therapy or traditional noise therapy signals. It is even possible to combine the two approaches.	✓	✓	✓	✓	
Notched AmplificationTherapy	Innovative tinnitus therapy delivered via everyday amplification.	✓	✓	✓	✓	
Smart Optimizer and Data Logging	Smart Optimizer (found under Data Logging) provides tailored fine-tuning recommendations based on an analysis of the wearer's usage of thehearing aid.	✓	✓	✓	✓	✓
Acclimatization manager	Allows the wearer to benefit from enhanced audibility in the longer term, while enjoying optimal device acceptance when initially fitted.	✓	✓	✓	✓	✓
InSituGram	Allows quick and convenient screening audiograms and checking of thresholds while the hearing aids are worn.	✓	✓	✓	✓	✓
AutoFit	Automatically provides an optimized real-ear insertion gain (REIG) setting without the need to manually match targets. The new IMC2 standard expands the range of REM equipment which supports AutoFit with Connexx and Noah 4.5.1 or higher.	✓	✓	✓	✓	✓
TeleCare	TeleCare provides remote services including additional tools to follow-up with patients. In addition, user engagement and autonomy can be increased via the Signia app.	✓	✓	✓	✓	✓
Remote Services	Includes TeleCare Dashboard, Remote Fine Tuning, CareChat (incl. text messages, phone calls and video calls), virtual follow-up.		✓	✓	✓	✓



Product Name	Form Factor	Channels	Fitting Range (dB)	Warranty (Years)	Unit Price (₹)
AUGMENTED XPERIENCE					

STYLETTO AX

Kit Styletetto 7AX	Slim-RIC	48	10 - 110	4	6,79,990
Kit Styletetto 5AX	Slim-RIC	32	10 - 110	4	4,45,990
Kit Styletetto 3AX	Slim-RIC	24	10 - 110	2	3,25,990
Kit Styletetto 2AX	Slim-RIC	16	10 - 110	2	1,39,990
Kit Styletetto 1AX	Slim-RIC	16	10 - 110	2	1,09,990

KIT PURE C&G AX / KIT PURE C&G T AX

Kit Pure C&G AX / Kit Pure C&G T 7AX	RIC	48	10 - 115	4	7,09,990
Kit Pure C&G AX / Kit Pure C&G T 5AX	RIC	32	10 - 115	4	4,41,990
Kit Pure C&G AX / Kit Pure C&G T 3AX	RIC	24	10 - 115	2	3,19,990
Kit Pure C&G AX / Kit Pure C&G T 2AX	RIC	16	10 - 115	2	1,49,990
Kit Pure C&G AX / Kit Pure C&G T 1AX	RIC	16	10 - 115	2	1,13,990

PURE C&G AX / PURE C&G T AX

Pure C&G AX / Pure C&G T 7AX	RIC	48	10 - 115	4	3,61,990
Pure C&G AX / Pure C&G T 5AX	RIC	32	10 - 115	4	2,27,990
Pure C&G AX / Pure C&G T 3AX	RIC	24	10 - 115	2	1,66,990
Pure C&G AX / Pure C&G T 2AX	RIC	16	10 - 115	2	81,990
Pure C&G AX / Pure C&G T 1AX	RIC	16	10 - 115	2	63,990

PURE 312 AX

Pure 312 7AX	RIC	48	10 - 115	4	3,14,990
Pure 312 5AX	RIC	32	10 - 115	4	1,99,990
Pure 312 3AX	RIC	24	10 - 115	2	1,41,990
Pure 312 2AX	RIC	16	10 - 115	2	64,990
Pure 312 1AX	RIC	16	10 - 115	2	46,990



Product Name	Form Factor	Channels	Fitting Range (dB)	Warranty (Years)	Unit Price (₹)
AUGMENTED XPERIENCE					

KIT INSIO C&G AX

Kit Insio C&G 7AX	ITC / ITE	48	10 - 110	4	7,09,990
Kit Insio C&G 5AX	ITC / ITE	32	10 - 110	4	4,41,990
Kit Insio C&G 3AX	ITC / ITE	24	10 - 110	2	3,19,990
Kit Insio C&G 2AX	ITC / ITE	16	10 - 110	2	1,49,990
Kit Insio C&G 1AX	ITC / ITE	16	10 - 110	2	1,13,990

INSIO C&G AX

Insio C&G 7AX	ITC / ITE	48	10 - 110	4	3,61,990
Insio C&G 5AX	ITC / ITE	32	10 - 110	4	2,27,990
Insio C&G 3AX	ITC / ITE	24	10 - 110	2	1,66,990
Insio C&G 2AX	ITC / ITE	16	10 - 110	2	81,990
Insio C&G 1AX	ITC / ITE	16	10 - 110	2	63,990

INTUIS 4

Intuis 4.7 (SP/P/M) ^{NEW}	BTE	48	10 - 120	4	2,04,990
Intuis 4.5 (SP/P/M) ^{NEW}	BTE	32	10 - 120	4	1,29,990
Intuis 4.3 (SP/P/M) ^{NEW}	BTE	24	10 - 120	2	74,990
Intuis 4.2 (SP/P/M)	BTE	16	10 - 120	2	41,990
Intuis 4.1 (SP/P/M)	BTE	16	10 - 120	2	35,990
Intuis 4.0 (SP/P/M) **	BTE	12	10 - 120	2	28,990



Styletetto AX



Pure Charge&Go AX



Insio Charge&Go AX



Pure 312 AX



Intuis 4

Note:
The fitting range is completely dependent on the receiver matrix of the hearing aid.
** Direct streaming option not available in Intuis 4.0.

AX

A revolution in augmented hearing


Product Name	Form Factor	Battery Size/Type	Warranty (Years)	Unit Price (₹)
CROS / BiCROS				
Kit Styletto CROS 7AX	RIC	Li-ion	4 years for Hearing Aid 1 year for CROS	3,69,990
Kit Styletto CROS 5AX	RIC	Li-ion	4 years for Hearing Aid 1 year for CROS	2,59,990
Kit Styletto CROS 3AX	RIC	Li-ion	2 years for Hearing Aid 1 year for CROS	1,99,990
Kit Styletto CROS 2AX	RIC	Li-ion	2 years for Hearing Aid 1 year for CROS	1,09,990
Kit Styletto CROS 1AX	RIC	Li-ion	2 years for Hearing Aid 1 year for CROS	94,990
CROS / BiCROS (Pure C&G AX Transmitter)	RIC	Li-ion	1	35,990
CROS / BiCROS (Pure 312 AX Transmitter)	RIC	312	1	29,490

Product Name		Warranty (Year)	Unit Price (₹)
CHARGERS			
Styletto AX Charger Case		1	11,990
Styletto Dry&Clean Charger		1	11,990
Insio AX Charger Case		1	11,990
Dry&Clean Charger RIC		1	11,990
Portable Charger RIC		1	7,990
Standard Charger RIC		1	5,990
Multi-Charger		1	3,990

EarWear 3.0 (Only for IX & AX Platform)

Better physical fit and on-ear attractiveness

Our new EarWear 3.0 offers a higher physical fit rate,improved cosmetics and better retention.



miniReceiver 3.0 Cable Bending
10/10 fit correctly-out of the box

- Your EarWear 3.0 advantages at a glance:
- Higher physical fit rate**
The combination of smaller receiver cross sections, softer eartips providing a smaller receiver front, and a new cable geometry results in major physical fit rate improvements for the miniReceiver 3.0 portfolio.
 - Cosmetically appealing**
Improved cable/tube bending and in-ear position ensures a closer fit on the ear, making it less visible. Every variant is available in the lengths 0-4.
 - Better retention**
A special bending of the cable/tube and its outlet locks the acoustic coupling reliably into the ear canal, with no lateral migration and no need for an anchor.
 - New eartips**
Compatible with miniReceiver 3.0 and ThinTube 3.0 and easier for your clients to exchange (compared to click domes).
 - Smaller hearing aids**
The introduction of a substantially smaller plug allows for smaller RIC designs overall.
 - Improved usability**
The proven design makes click domes much easier for consumers to exchange.

Acoustical and physical fitting performance of miniReceiver 3.0 portfolio

Physical fit rate	↗		↗		↑		→	
Acoustical performance 3.0	↗		↗		↘		↗	
	S		M		P		HP	
	2ccm coupler	Ear simulator	2ccm coupler	Ear simulator	2ccm coupler	Ear simulator	2ccm coupler	Ear simulator
OSPL 90 (Peak)	110 dBSPL	120 dBSPL	119 dBSPL	129dBSPL	122 dBSPL	131 dBSPL	131 dBSPL	138 dBSPL
138 dBSPL	46 dB	57 dB	60 dB	70 dB	65 dB	75 dB	75 dB	83 dB

Different values for Pure 312 AX

SigniaXperience Platform Features

Features	Information	Performance levels				
		7X	5X	3X	2X	1X
YourSound Technology	Redesigned acoustic sensors plus an integrated motion sensor provide detailed information about where the wearers are and what they are doing. Dynamic Soundscape Processing and Own Voice Processing (OVP™), process this information delivering the most natural and personalized sound.
Acoustic Sensors	The redesigned acoustic sensors can detect more variables of sound than ever before, providing more accurate information about the environment.	✓	✓	✓	✓	✓
Motion Sensor	An integrated motion sensor detects whether the wearer is walking, running or not. This information is used for Dynamic Soundscape Processing.	✓	✓	✓	-	-
Dynamic Soundscape Processing¹	Steers Sound Clarity and Sound Quality features for natural sound and speech in every situation from any direction, even when moving.
OVP™ (Own Voice Processing)¹	Own Voice Processing (OVP™) utilizes real-time recognition of the wearer's voice to deliver a natural own-voice impression. Available for RIC/BTE devices with e2e wireless in binaural fittings.	-	-
Sound Clarity	Sound Clarity provides lifelike sound which covers the full spectrum of wearingsituations.
Signal processing / Gain & MPO	Channels/handles in Connexx.	48/20	32/16	24/12	16/8	16/8
Hearing Programs	Number of programs configurable in Connexx	6	6	6	4	4
Extended dynamic range	24 bit processing for enhanced digital encoding of soft and loud sounds.	✓	✓	✓	✓	✓
Speech and noise management	Reduces listening effort using sophisticated sound processing.	✓	✓	✓	✓	✓
SoundSmoothing™	Reduces the annoyance of sudden loud impact sound such as when handling glassware.	✓	✓	✓	✓	-
Feedback cancellation	High speed monitoring and control of feedback in individual processing channels.	✓	✓	✓	✓	✓
HD Music	Preset programs for enhanced enjoyment of non-streamed music.	3	3	1	1	-
eWindScreen	Reduces the annoyance of wind noise when outdoors. Binaural e Windscreen requires e2e wirelessl link. 7X, - Binaural Auto in Universal program 5X, - Binaural 3X, 2X, - Monaural	✓	✓	✓	✓	-
Extended bandwidth	12 kHz bandwidth for enhanced processing of high-frequency speech and environmental sounds.	✓	-	-	-	-
EchoShield	Dedicated program for reverberant environments.	✓	-	-	-	-
Speech Quality	Speech Quality delivers improved speech understanding and communication ease in noise.
Directionality	Directionality provides scientifically proven enhanced speech understanding in noise. Available for bilateral fittings with directional microphones and e2e wireless link. 7X, 5X, 3X - Binaural processing Automatic / Adaptive 2X, 1X - Automatic / Adaptive	✓	✓	✓	✓	✓
Binaural OneMic Directionality	Advanced binaural beamforming technology for small in-the-ear instruments with onlyone microphone. Available for Silk X.	✓	✓	✓	-	-
Wireless CROS / BiCROS²	Requires CROS trasmitters: CROS Pure 312 X or CROS Pure Charge&Go X for Bluetooth enabled fittings CROS Silk X for fittings without Bluetooth	✓	✓	✓	-	-
Frequency compression	Improves the detection of speech sounds for wearers with reduced high frequency hearing thresholds. Frequency compression is selectively activated based on audiological criteria.	✓	✓	✓	✓	✓

Features	Information	Performance levels				
		7X	5X	3X	2X	1X
Spatial SpeechFocus	Depending on the direction of the dominant speech source, Spatial SpeechFocus 7X - In Universal program:automatic mic beam to the right / left / rear directions. 5X - In Universal program: automatic mic beam to the rear, surround mode for side speech. In manual programs, microphone beam to the right/left/ rear directions.	✓	✓	-	-	-
Wearer Interaction	Technically inclined wearers may appreciate having extra control over their hearing aids. Signia Xperience hearing aids offer these wearers a comprehensive range of adjustment options via the Signia App.
Signia Assistant¹	Signia Assistant offers the wearers instant 24/7 support and troubleshooting around the clock. It provides detailed information in Connexx to address each client's fine-tuning needs more precisely.	✓	✓	✓	✓	✓
Signia App (iOS and Android)	The new Signia App combines all functionalities of previous apps (touchControl, myControl, myHearing App). The functionality offered in the app is automatically customized depending on the connected hearing aid.	✓	✓	✓	✓	✓
Adaptive Streaming Volume	Speech intelligibility weighted adaptation of the streaming volume depending on ambient noise.	✓	✓	✓	✓	✓
Spatial Configurator	Spatial Configurator is part of the Signia App and allows the wearer to adjust the direction of the microphone beam.	✓	✓	-	-	-
Direct Streaming	Signia Xperience hearing aids can be connected directly to Apple devices for phonecalls and audio streaming – without additional accessories. For Android devices the StreamLine Mic accessory makes streaming easy. The Stream Line TV accessory sends high quality stereo sound from the wearer's TV directly to the Xperience hearing aids. Available for Bluetooth hearing aids.	✓	✓	✓	✓	✓
Made for iPhone	Direct audio streaming from compatible iOS devices to Bluetooth enabled Xperience hearing aids.	✓	✓	✓	✓	✓
StreamLine TV accessory	Stream high quality stereo sound directly from a TV to Bluetooth equipped Xperience hearing aids.	✓	✓	✓	✓	✓
StreamLine Mic accessory	Stream high quality stereo sound directly from any Bluetooth device (Android devices, older Apple devices, other smartphones, PCs, laptops, etc.) to Bluetooth equipped Xperience hearing aids.	✓	✓	✓	✓	✓
Tinnitus	A choice of tinnitus treatments based on amplified notch therapy or traditional noise therapy signals. It is even possible to combine the two approaches.	✓	✓	✓	✓	-
Notched Amplification Therapy	Innovative tinnitus therapy delivered via everyday amplification.	✓	✓	✓	✓	-
Tinnitus noise therapy signal	A set of fully configurable stationary noise signals and the modulated ocean wave sounds	✓	✓	✓	✓	-
Fitting	Connexx fitting software offers a number of advanced tools to suit every Hearing Care Professional's workflow.	✓	✓	✓	✓	✓
Smart Optimizer and Data Logging	Smart Optimizer (found under Data Logging) provides tailored fine-tuning recommendations based on an analysis of the wearer's usage of the hearing aid.	✓	✓	✓	✓	✓
Acclimatization manager	Allows the wearer to benefit from enhanced audibility in the longer term, while enjoying optimal device acceptance when initially fitted.	✓	✓	✓	✓	✓
InSituGram	Allows quick and convenient screening audiograms and checking of thresholds while the hearing aids are worn.	✓	✓	✓	✓	✓
AutoFit	Automatically provides an optimized real-ear insertion gain (REIG) setting without the need to manually match targets. The new IMC2 standard expands the range of REM equipment which supports AutoFit with Connexx and Noah 4.5.1 or higher.	✓	✓	✓	✓	✓
TeleCare	TeleCare provides remote services including additional tools to follow-up the patients. In addition, the user engagement and autonomy can be increased via the Signia app.	✓	✓	✓	✓	✓
Remote Services	Includes TeleCare Dashboard, Remote Fine Tuning, CareChat (incl. text messages, phone calls and video calls), virtual follow-up.	✓	✓	✓	✓	✓
Signia App	The Signia App provides TeleCare specific functionalities like Hearing Lessons, Daily Satisfaction Ratings, Tutorials and FAQs.	✓	✓	✓	✓	✓

¹ Not available with Silk X ² Not available with Styletto X



Product Name	Form Factor	Channels	Fitting Range (dB)	Warranty (Years)	Unit Price (₹)
XPERIENCE					

KIT MOTION C&G / KIT MOTION C&G P X / KIT MOTION C&G SP X

Kit Motion C&G 7X / P 7X / SP 7X	BTE	48/48/34	10 - 120	4	6,19,990
Kit Motion C&G 5X / P 5X / SP 5X	BTE	32/32/32	10 - 120	4	3,89,990
Kit Motion C&G 3X / P 3X / SP 3X	BTE	24/24/24	10 - 120	2	2,19,990
Kit Motion C&G 2X / P 2X / SP 2X	BTE	16/16/16	10 - 120	2	1,29,990
Kit Motion C&G 1X / P 1X / SP 1X	BTE	16/16/16	10 - 120	2	99,990

MOTION C&G / MOTION C&G P X / MOTION C&G SP X

Motion C&G 7X / P 7X / SP 7X	BTE	48/48/34	10 - 120	4	3,16,990
Motion C&G 5X / P 5X / SP 5X	BTE	32/32/32	10 - 120	4	2,01,990
Motion C&G 3X / P 3X / SP 3X	BTE	24/24/24	10 - 120	2	1,12,990
Motion C&G 2X / P 2X / SP 2X	BTE	16/16/16	10 - 120	2	71,990
Motion C&G 1X / P 1X / SP 1X	BTE	16/16/16	10 - 120	2	56,990

KIT PURE C&G X

Kit Pure C&G 2X	RIC	16	10 - 110	2	1,03,990
Kit Pure C&G 1X	RIC	16	10 - 110	2	89,990

PURE C&G X

Pure C&G 2X	RIC	16	10 - 110	2	54,990
Pure C&G 1X	RIC	16	10 - 110	2	47,990



Motion Charge&Go SP X



Pure Charge&Go X



Silk X

Note:
The fitting range is completely dependent on the receiver matrix of the hearing aid.



Product Name	Form Factor	Channels	Fitting Range (dB)	Warranty (Years)	Unit Price (₹)
XPERIENCE					

SILK X

Silk 7X	Instant Fit CIC	48	10-80	4	2,99,990
Silk 5X	Instant Fit CIC	32	10-80	4	1,79,990
Silk 3X	Instant Fit CIC	24	10-80	2	1,29,990
Silk 2X	Instant Fit CIC	16	10-80	2	59,990
Silk 1X	Instant Fit CIC	16	10-80	2	43,990

Product Name	Form Factor	Warranty (Year)	Unit Price (₹)
--------------	-------------	-----------------	----------------

CROS/BiCROS

CROS / BiCROS Pure Charge&Go X Transmitter	RIC	1	35,990
CROS / BiCROS Pure 312 X Transmitter	RIC	1	29,490
CROS / BiCROS Silk X Transmitter	Instant Fit CIC	1	29,490

Product Name	Warranty (Year)	Unit Price (₹)
--------------	-----------------	----------------

CHARGERS

Inductive Charger 2 (Pure C&G X and Motion C&G X) - 7,5,3		1	11,990
Dry&Clean Charger (Motion C&G X) (P/SP)		1	11,990
Inductive Charger E (Pure C&G X and Motion C&G X) - 2,1		1	9,990
Standard Charger (Motion C&G X) (P/SP)		1	5,990
Multi-Charger		1	3,990



Product Name	Form Factor	Channels	Fitting Range (dB)	Warranty (Years)	Unit Price (₹)
INSIO PX					
Insio 3PX	CIC & ITC/ITE	24	10 - 110	2	44,990
Insio 2PX	CIC & ITC/ITE	16	10 - 110	2	39,990
Insio 1PX	CIC & ITC/ITE	16	10 - 110	2	32,990



primax
All levels, all styles, all wireless.

Product Name	Form Factor	Channels	Fitting Range (dB)	Warranty (Years)	Unit Price (₹)
ESSENTIAL					
Intuis 3 Custom	ITE/ITC/CIC	12	10 - 110	2	26,990
Intuis RIC 312	RIC	12	10 - 115	2	25,990
Intuis 3 SP	BTE	12	20-120	2	25,990
Intuis 3 P	BTE	12	10 - 110	2	25,990
Intuis 3 M	BTE	12	10 - 90	2	25,990
Intuis 3 S	BTE	12	0 - 80	2	25,990
Intuis 3 Click	ITC/CIC	12	10 - 90	2	25,990




Intuis 3 (BTE)



Intuis 3 (ITE / ITC / CIC)



Intuis 3 (Click ITC)








Product Name	Form Factor	Channels	Fitting Range (dB)	Warranty (Years)	Unit Price (₹)
BASIC					
Prompt Custom	ITC/CIC	8	10 - 110	2	20,990
Prompt SP	BTE	8	20 - 120	2	19,990
Prompt P	BTE	8	10 - 90	2	19,990
Prompt S	BTE	8	0 - 80	2	19,990
Prompt Click	ITC/CIC	8	10 - 90	2	19,990
Run Click	ITC/CIC	8	10 - 90	2	17,990
Run SP	BTE	8	20 - 120	2	16,990
Run P	BTE	8	20 - 120	2	16,990
Fun SP	BTE	6	20 - 120	2	12,990
Fun P	BTE	6	10 - 90	2	12,990
Fast P	BTE	4	10 - 80	2	10,990



Product Name	Form Factor	Channels	Fitting Range (dB)	Warranty (Years)	Unit Price (₹)
OTHER HEARING AIDS					
Active Pro	Instant Fit ITE	48	10 - 80	4	2,99,990
Active	Instant Fit ITE	16	10 - 80	2	1,49,990
Kit Pure C&G 7X	RIC	48	10 - 110	4	6,49,990
Kit Pure C&G 5X	RIC	32	10 - 110	4	3,99,990
Kit Pure C&G 3X	RIC	24	10 - 110	2	2,91,990
Pure C&G 7X	RIC	48	10 - 110	4	3,31,990
Pure C&G 5X	RIC	32	10 - 110	4	2,06,990
Pure C&G 3X	RIC	24	10 - 110	2	1,52,990
Pure 312 7X	RIC	48	10 - 115	4	2,99,990
Pure 312 5X	RIC	32	10 - 115	4	1,79,990
Pure 312 3X	RIC	24	10 - 115	2	1,29,990
Pure 312 2X	RIC	16	10 - 115	2	45,990
Pure 312 1X	RIC	16	10 - 115	2	38,990
Insio 7NX	ITE/ITC*/CIC/IIC	48	10 - 110	4	2,99,990
Insio 5NX	ITE/ITC*/CIC/IIC	32	10 - 110	4	1,79,990
Insio 3NX	ITE/ITC*/CIC/IIC	24	10 - 110	2	1,04,990
Insio 2NX	ITE/ITC*/CIC/IIC	16	10 - 110	2	57,990
Insio 1NX	ITE/ITC*/CIC/IIC	16	10 - 110	2	41,990
Motion P/SP 7PX	BTE	24/48	10 - 120	2	2,49,990
Motion P/SP 5PX	BTE	24/32	10 - 120	2	1,54,990
Motion P/SP 3PX	BTE	24/24	10 - 120	2	44,990
Motion P/SP 2PX	BTE	16/16	10 - 120	2	39,990
Motion P/SP 1PX	BTE	16/16	10 - 120	2	32,990
Pure 3PX	RIC	24	10 - 115	2	44,990
Pure 2PX	RIC	16	10 - 115	2	39,990
Pure 1PX	RIC	16	10 - 115	2	32,990

Note:
 ITE & ITC – Available with and without Bluetooth.
 The fitting range is completely dependent on the receiver matrix of the hearing aids.

Product Name	Form Factor	Warranty (Years)	Unit Price (₹)
CROS/BiCROS			
CROS/BiCROS Pure PX Transmitter	RIC	1	24,490
CROS/BiCROS(Pure 312 NX Transmitter)	RIC	1	29,490

Product Name		Warranty (Year)	Unit Price (₹)
CHARGERS			
Active Charger Case		1	11,990
Dry&Clean Charger RIC		1	11,990
Portable Charger RIC		1	7,990
Standard Charger RIC		1	5,990
Multi-Charger		1	3,990



Accessories for Hearing Aids			
Product	Unit Price (₹)	Product	Unit Price (₹)

Remote Control and Accessories	
StreamLine Mic	22,000
StreamLine TV	19,470
miniPocket	7,670

Hearing Aid Zinc Air Batteries	
Signia # 675 Battery (Mercury Free)	51
Signia # 13 Battery (Mercury Free)	51
Signia # 312 Battery (Mercury Free)	51
Signia # 10 Battery (Mercury Free)	51

In-The-Ear (ITE) Fitting Accessories	
Quick Guard Kit	530

Receiver-In-Canal (RIC) Fitting Accessories	
Receiver Mold	1,168
miniReceiver Waxguard Set (10 pcs)	650

Behind-The-Ear (BTE) Fitting Accessories	
Comfort Coupling Set (for Prompt, Run, Fun)	236
Simple Coupling Set (for Fast)	177
BTE Eartip Assembly (Eartip + Elbow + Tube)	118
BTE Eartip Joint Elbow (L-Bend)	41
BTE Eartip (Any Size)	59
BTE Hucki Holder (One Pair)	148
Colour Exchange Kit	1,522

Miscellaneous Items	
UV Hard Pocket / BTE Earmold	900
Silicon Soft Pocket / BTE Earmold	1,100
Custom Swim Plug	1,100
Drying Kit	150
Cleaning Kit	590

CROS / BiCROS Transmitter	
Styletto IX Transmitter	39,990
Silk C&G IX Transmitter	39,990
Active IX Transmitter	39,990
Styletto AX Transmitter	39,990

Earwear 2.0	
Domes	
Click Domes Open (8mm, 10mm) (6pcs/blister)	932
Click Domes Closed (8mm, 10mm) (6pcs/blister)	932
Click Domes Double (8mm/10mm,,10/12 mm) (6pcs/blister)	932
Click Domes Semi Open (8mm/10mm) (6pcs/blister)	932

Sleeves	
Sleeves Vented (XS,S,M,L) (6pcs/blister)	1,168
Sleeves Closed (XS,S,M,L) (6pcs/blister)	1,168

miniReceiver 2.0	
S miniReceiver 2.0 (Sizes 0 - 3) Left / Right	5,890
M miniReceiver 2.0 (Sizes 0 - 4) Left / Right	5,890
P miniReceiver 2.0 (Sizes 0 - 3) Left / Right	5,890
HP miniReceiver 2.0 (Sizes 0 - 3) Left / Right	5,890

Earwear 3.0	
Sleeves 3.0	
Sleeves 3.0 Vented (XS, S, M, L) (6pcs/blister)	1,190
Sleeves 3.0 Closed (XS, S, M, L) (6pcs/blister)	1,190
Sleeves 3.0 Power (XS, S, M, L) (6pcs/blister)	1,190










Eartips 3.0	
Eartip 3.0 Open (5mm, 7mm, 10mm) (6pcs/blister)	790
Eartip 3.0 Tulip (8mm, 12mm) (6pcs/blister)	790

Thin Tube 3.0	
Set Thintube 3.0 (S0, S1, S2, S3, S4) (Set of Left & Right)	490
Set Thintube 3.0 P (S0, S1, S2, S3, S4) (Set of Left & Right)	490

miniReceiver 3.0	
S miniReceiver 3.0 (Sizes 0 - 4) Left / Right	5,890
M miniReceiver 3.0 (Sizes 0 - 4) Left / Right	5,890
P miniReceiver 3.0 (Sizes 0 - 4) Left / Right	5,890
HP miniReceiver 3.0 (Sizes 0 - 4) Left / Right	5,890

App, accessory, and charger compatibility

This overview shows you which apps and accessories work with our hearing aids. Depending on the solution your customers choose, they can even enhance their hearing experience and convenience with this well-matching portfolio.

Hearing aid	Signia app	CROS/ BiCROS	Streamers and transmitters		Remote Controls	Portable Charger	Inductive Charger II	Standard Charger	Dry&Clean Charger	Multi Charger
		CROS	StreamLine TV	StreamLine Mic	miniPocket					
										
Styletto IX	●	● ^P	●	●	●	●	-	-	●	-
Pure Charge&Go IX	●	● ^N	●	●	●	●	-	●	●	●
Pure Charge&Go T IX	●	● ^N	●	●	●	●	-	●	●	●
Active Pro IX	●	● ^R	●	●	●	●	-	-	-	-
Silk Charge&Go IX	●	● ^O	-	-	●	●	-	-	-	-
Insio IIC / CIC IX	●	● ^Q	-	-	●	-	-	-	-	-
Styletto AX	●	● ^I	●	●	●	●	-	-	●	-
Insio Charge&Go AX	●	● ^J	●	●	●	-	-	●	-	-
Pure Charge&Go AX	●	● ^J	●	●	●	●	-	●	●	●
Pure Charge&Go T AX	●	● ^J	●	●	●	●	-	●	●	●
Pure 312 AX	●	● ^J	●	●	●	-	-	-	-	-
Active Pro	●	-	●	●	●	●	-	-	-	-
Motion Charge&Go X	●	● ^K	●	●	●	-	●	-	-	-
Motion Charge&Go P X	●	● ^K	●	●	●	-	-	●	●	●
Motion Charge&Go SP X	●	● ^K	●	●	●	-	-	●	●	●
Styletto X	●	-	●	●	●	●	-	-	-	-
Silk X	●	● ^L	-	-	●	-	-	-	-	-
Pure Charge&Go X	●	● ^K	●	●	●	-	●	-	-	-
Pure 312 X	●	● ^K	●	●	●	-	-	-	-	-
Motion 13 Nx	●	● ^C	●	●	●	-	-	-	-	-
Motion 13P Nx	●	● ^C	●	●	●	-	-	-	-	-
Insio Nx	●	● ^{G,H}	● ^F	● ^F	●	-	-	-	-	-
Orion Charge&Go	●	-	-	-	●	-	-	-	-	●
Intuis 4	●	● ^M	● ^M	● ^M	●	-	-	-	-	-
Intuis 3 RIC312	●	-	-	-	●	-	-	-	-	-
Intuis 3	●	-	-	-	●	-	-	-	-	-
Prompt	●	-	-	-	●	-	-	-	-	-
Run	●	-	-	-	-	-	-	-	-	-

C) Compatible with CROS Pure 312 Nx or CROS Pure Charge&Go Nx

D) Compatible with the Pure Charge&Go/Motion Charge&Go/Cellion charger

E) Compatible with CROS Silk Nx

F) Compatible with Bluetooth-enabled Insio Nx ITC/ITE

G) Insio Nx IIC and CIC are compatible with CROS Silk Nx

H) Bluetooth-enabled Insio Nx ITC and ITE are compatible with CROS Pure 312 Nx

I) Compatible with CROS Styletto AX

J) Compatible with CROS Pure Charge&Go AX and CROS Pure 312 AX

K) Compatible with CROS Pure Charge&Go X and CROS Pure 312 X

L) Compatible with CROS Silk X

M) Excluding PL 4.0

N) Compatible with CROS Pure Charge&Go IX

O) Compatible with CROS Silk Charge&Go IX

P) Compatible with CROS Styletto IX

Q) Requires wireless faceplate option, Compatible with CROS Insio IIC /CIC IX

R) Compatible with CROS Active IX

