




**URSUS™**  
special stroller



Ursus™ seat is a modern, safe and functional solution for patients, which prefers a high comfort and stabilization. **TILT-IQ®** function, **Tilt-in-space ClickClack®** are some of many practical solutions of **URSUS SYSTEM™**


- **Outdoor frame Ursus™** Means tours and travels in very comfortable conditions either for a child and caregivers. It provides a convenient ride in varied terrain. Swivel front wheels ensure smooth changing of direction.
- **Hi-Low frame Ursus Home™** Is a comfortable and functional solution in home conditions. Adjustable seat height greatly facilitates contact with peers, and push handles make maneuvering, even in tight spaces, much easier.



one seat  
many possibilities

- **Active frame Ursus Active™** is a mode that gives your child a freedom of movement. Thoughtful solution enables independent movement of the user or with assistance of a caregiver. Now, therapy, learning, and fun will be even more enjoyable and effective.
- **Pneumatic Hi-Low frame Quark™** Will convert **URSUS™** seat into a comfortable rehabilitation chair. Innovative system of height adjustment ensures smooth and easy operation. Included push handle allows easy relocation of the device.

Choose  
your  
upholstery



### BREATHABLE UPHOLSTERY

- ☑ Lightweight
- ☑ Hypoallergenic
- ☑ Natural protection against mold
- ☑ The antibacterial properties
- ☑ Easy to clean
- ☑ Easy decontamination
- ☑ Wide range of hardness and elasticity
- ☑ Ecological
- ☑ Easy to recycle
- ☑ Resistant to crushing



### WATERPROOF UPHOLSTERY

- ☑ Lightweight
- ☑ Easy to clean
- ☑ Soft
- ☑ Waterproof
- ☑ Easy maintenance
- ☑ Ecological
- ☑ Durable
- ☑ Easy decontamination



Ultralight stroller **URSUS™** is designed for daily rehabilitation of children with locomotor disabilities. As far as its construction is concerned, the frame and seat are made of high-strength aluminum, and the chassis is additionally cushioned. **Ursus™** is equipped with **ClickClack®** feature, which means that you can easily swap the seats and adjust the stroller to the most suitable mode. The device has a wide range of adjustments such as width, seat depth, backrest height, backrest tilt angle, infinitely adjustable tilt-in-space function – **TILT-IQ®**, reversible seat. **Ursus™** has also been equipped with spatial, breathable upholstery that prevents excessive sweating as well as bedsores, and more to the point, positively influences the comfort of use. Moreover, the stroller can easily be converted into a lounger so that the child can sleep in it.

The product is available in **3** different sizes with a very wide range of accessories that will help you to meet the individual needs of user.

## IT IS RECOMMENDED IN THE FOLLOWING CASES:

- Infantile cerebral palsy – different forms.
- Meningomyelocele.
- Brain and spinal cord injuries.
- Muscle diseases connected with paresis and paralysis.
- Genetic disorders connected with paresis.
- Other diseases connected with paresis, paralysis and injuries regarding locomotor system.
- After traumatic brain injury, spine with spinal cord and lower limbs injury.



Stroller is equipped with Light Grey rims. However there is possibility to change the rim's colour to black.



Reversible seat!



## STANDARD EQUIPMENT:

Light aluminium frame | Fully foldable | Wheels with glass fiber, front PU swivel with direction | blockage, Rear inflatable with 2 separate wheel brace | Front shock absorbers | Ergonomic folding system | Seat unit with ClickClack® system for easy plug in/out | 1 piece style height adjustable push handle | Reflective safety marks | Seat unit that can be used forward or rearward facing | Fully reclining backrest to make the stroller as a bed | Height adjustable backrest | Angle adjustable backrest | 5-points safety belts with height adjustment | Adjustable headrest | Anatomical lateral supports | Depth and width adjustable seat | Pneumatic tilt-in-space system TILT-IQ® | Height and angle adjustable footrest | 5 points feet stabilizing belts with heel support and foot | abduction-adduction adjustments | Angle adjustable footplate | Soft upholstery with spine load reducing and ventilating | Soft front bar with reflective safety marks |

*Our passion – your independence!*

Colors may differ from the presented images. The manufacturer reserves the right to make changes in products without updating the directory. Current offer on our website.

Akces-MED sp. z o.o.  
Jasionka 955B  
36-002 JASIONKA

export@akces-med.com  
+48 (17) 864 04 77  
www.akces-med.com





ADVANTAGES:



**Stable structure**

The device ensures adequate stability and safety of Patient.



**Reflective elements**

The product has reflective elements to ensure your safety.



**Lightweight construction**

The product is made of lightweight and durable materials.

24

**24 months warranty**

We offer 24 months warranty on the product. We also provide after sales service.



**Easy to use**

Intuitive way of adjustment and easy assembly.



**Movable seat**

You can turn the seat easily so that the child may be facing towards caregiver or movement direction.



**Tilt-in-space**

It is possible to change the tilt angle of the seat.



Five-point straps for keeping safety.



Front swivel wheels for easy maneuverability.



The stylish and modern upholstery offers the highest comfort.

STANDARD COLORS:

BREATHABLE UPHOLSTERY:



Grass  
Size 1



Sea  
Size 2



Latte  
Size 3

WATERPROOF UPHOLSTERY:



City Grass  
Size 1



City Rose  
Size 2



City Coffee  
Size 3

You can change upholstery colours between sizes without additional charges. Production time up to 14 days.

OPTIONAL COLORS:



Black\*



Black&White\*



Red\*



City Carbon\*

\*On special request. Production time ca. 14 days.


*Our passion - your independence!*

Colors may differ from the presented images. The manufacturer reserves the right to make changes in products without updating the directory. Current offer on our website.

Akces-MED sp. z o.o.  
Jasionka 955B  
36-002 JASIONKA

export@akces-med.com  
+48 (17) 864 04 77  
www.akces-med.com





Symbol	Measurement	Unit	Size 1	Size 2	Size 3
A	Seat width	[cm]	19-25	22-30	24-34
B	Seat depth	[cm]	21-26	27-35	34-42
C	Backrest height	[cm]	44-55	49-59	55-67
E	Backrest width	[cm]	23	27	32
F	Footrests length	[cm]	13-29	16-40	16-43
	Footer length	[cm]	20	20	20
	Push handle height	[cm]	100	100	100
X	Seat tilt angle	[°]	90-180	90-180	90-180
Y	Backrest tilt angle (forward-facing)	[°]	3-20	3-25	3-30
Y1	Backrest tilt angle (rearward facing)	[°]	3-22	3-17	3-20
Z	Footrests tilt angle	[°]	90-180	90-180	90-180
	Max user weight	[kg]	35	35	50

Product dimensions and weight		Size 1	Size 2	Size 3
Width	[cm]	59	66	67
Width after folding	[cm]	59	66	67
Length	[cm]	90	90	90
Length after folding	[cm]	90	90	92
Height	[cm]	100	100	100
Height after folding	[cm]	45	48	53
Weight	[kg]	19,5	21	22,5

Note: Measurement tolerance range is +/- 1 cm / 1/2 in.



Safety



**Soft cover for tray**  
Soft cover protects the user from hitting the tray. Indicated for patients with impaired limbs coordination.



**Led lights (2 pieces)**  
Front and rear LED lights for increased safety.




**Lesser side covers for head**  
They provide additional head protection from unwanted contact with stroller's frame. Mounted around the push handle.




**Patient lift belt**  
Robust and ergonomic design as well as high durability facilitates daily care of the patient.



Stabilization



**H vest**  
It acts as a trunk stabilizer. HOLD&PULL© system makes it easier to adjust the vest to the patient body. You can quickly detach the vest if necessary.



**4 point vest**  
It acts as a trunk stabilizer. Wide range of adjustment helps to fit it to the user's needs. Equipped with HOLD&PULL© system.



**6 point vest**  
Wide range adjustable vest that acts as a trunk stabilizer. HOLD&PULL© system makes it easier to adjust the vest to the patient body.



**Thigh abduction belts**  
They can be easily attached to the chair by use of Velcro straps. Their use allows to keep legs in abduction and also prevents the patient from slipping off of the device.



**Thigh abduction belts**  
They can be easily attached to the chair by use of Velcro straps. Their use allows to keep legs in abduction and also prevents the patient from slipping off of the device.



**Pelvic belt**  
Wide belt equipped with "HOLD AND PULL" system. It provides better stabilization in the device.



**Firm headrest**  
Made of poliurethane foam. Provides with basic support of the head




**Head-neck support**  
Perfectly stabilizes the head. Indicated for people with limited head control. You can mount it at different height.



**Lumbar support**

It allows better stabilization of the lumbar back.




**Trunk stabilizing belt**

It acts as an additional trunk protection which helps to keep a proper sitting position. The belt is fastened by use of Velcro straps and its size can be adjusted to the user.



**Trunk belt**

Universal solution applicable to many of our products. It secures the user from falling out of the device. Anatomical shape and Velcro straps allow for a comfortable fit.



**Abduction block**

Shaped abduction block provides with correct seating position. Can be easily removed.



**Lateral supports**

Width adjustable supports that provide additional side stabilization of the trunk. You can mount it at different height.



**Seat narrowing cushions 3cm**

Made of elastic foam trimmed with soft material. They are attached to the seat upholstery from both sides. Thanks to them you can narrow the seat by 6cm.



**Seat narrowing cushions 5cm**

Made of elastic foam trimmed with soft material. They are attached to the seat upholstery from both sides. Thanks to them you can narrow the seat by 10cm.



**BodyMap Stabilization**



**BodyMap A seat cushion**

Vacuum seat cushion allows for stabilisation and proper positioning of hips.



**BodyMap B Backrest cushion**

The vacuum back cushion provides comfort and correction for patients with scoliosis and poor posture. Cushion maintains and controls required posture in the stroller.




**Vismemo cover for BodyMap cushion**

Made of thermo-active foam for extra bed sore prevention. Allows for longer and more comfortable use of the cushion.



**Cotton cover for BodyMap cushion**

Made of soft material it offers protection and absorbs moisture.



**BodyMap spacer**

Aerating-massaging pad to strengthen the prophylactic anti-bed sore function in patients using BodyMap. It can only be used together with the terry cover. The effect of being 'suspended in the air' creates a skin-friendly microclimate.






Comfort



**Umbrella**

Attachable to the stroller's frame, perfectly protects from the sunlight or rain.



**Tray**

Made of plastic. It supports the therapy and allows you to play, eat, and work.



**Folding canopy without side panels**

Protects against both sun and rain.



**Folding canopy with side panels**

Side panels protect not only against sun and rain but also the wind.



**Cup holder**

An universal holder mounted to the frame. It enables easy and quick access to a bottle of drink.



**Headrest cover**

It helps to keep the headrest clean. Easy to assemble/disassemble, it can be washed.



**Summer sleeping bag**

Protects legs and trunk. This esthetic body sheet fits the height of the user covering also the footrest of the pushchair.



**Winter sleeping bag**

Protects legs and the torso area. A very tasteful cover of the body, adjusting to the height of the user, covering the footrest of the chair. Insulated lining provides extra protection against the cold.



**Winter hood**

Can be attached to the sleeping bag. Protects head against frost, wind and moisture.



**Plexiglas® Tray**

Made of Plexiglas®. It supports the therapy and allows you to play, eat, and work.



**Rain Cover**

An excellent protection from the rain, snow, and wind either for the user and stroller.



**Mosquito Net**

Excellent protection against insects



**Wheel covers (4 pcs.)**

Use it on the dirt tires when you want to take the products indoor or transport in the trunk.



**Decorated spokes guard**

Spokes guards 20", made OF polycarbonate 1,5mm, outsider diameter 415mm, insi de diameter 100mm



Storage



**Storage basket**

Enables you to carry small objects while you stroll. It is attached to the base of the stroller.



**Women' bag**

This capacious bag designed for mothers will be perfect for carrying all the necessities.



**Men's bag**

Capacious men-style bag for carrying all the necessary things.



**DeLux bag**

Large, durable and stylish bag useful for carrying necessary things. Made from high-quality materials.



**Oxygen tank case**

Oxygen tank case can be attached to the stroller by use of Velcro straps with adjustable length. Made to measure to perfectly fit the patient's oxygen tank.



**Stroller cover 'n carry bag**

Perfect solution to protect folded stroller against damage (e.g. during transport). Easily fits the shape and the size of a stroller.



Build-up



**Kid-sit buggypod IO**

For convenient and safe transportation of a second child without having to buy a double stroller. The buggypod IO has been designed so that within a few seconds, it can be folded, unfolded, or completely disassembled. 5 point belts seat and folding footrest provide the child with comfort and safety. Using the conversion kit, you are able to convert the kid-sit into an independent three-wheeled stroller. Load capacity up to 20kg.



**Conversion kit (for The buggypod kit-sits)**

Set of elements thanks to which the Buggypod kit-sit can be converted into a fully functional stroller.



**Rain cover for the kid-sits**

An excellent protection from the rain, snow, and wind either for the user and stroller



Buggy pod PERLE



Buggy board MAXI



**Extension bar**

An adapter enabling to mount the Buggy pod kid-sits to the strollers of other manufacturers.



**Mounting set**

It enables easy assembly of the Buggy pod kit-sits to other strollers whose frame diameter is larger than 40mm.



**Pneumatic Hi-Low frame Quark**

Will convert URSUS™ seat into a comfortable rehabilitation chair. Innovative system of height adjustment en-sures smooth and easy operation. Included push handle allows easy relocation of the device.



**Hi-Low frame Ursus Home**

Is a comfortable and functional solution in home conditions. Adjustable seat height greatly facilitates contact with peers, and push handles make maneuvering, even in tight spaces, much easier.



**Outdoor frame Ursus**

Means tours and travels in very comfortable conditions either for a child and caregivers. It provides a convenient ride in varied terrain. Swivel front wheels ensure smooth changing of direction.



**Active frame Ursus Active**

is a mode that gives your child a freedom of movement. Thoughtful solution enables independent movement of the user or with assistance of a caregiver. Now, therapy, learning, and fun will be even more enjoyable and effective.





Spare parts



**Spacer fabric**  
Breathable and easy to wash. Usefull especially during hot days and high temperature climate.



**Wheel small, full silver metallized (Instead of / in addition)**  
Wheel 12" x 2" 1/2



**Wheel duże, full, silver metallized (Instead of / in addition)**  
Wheel 12" x 2" 1/2 srebrna felga



**Wheel duże, pompowane, silver metallized (Instead of / in addition)**  
Wheel 12" x 2" 1/2 srebrna felga. Możliwość zamówienia również samej dętki lub opony.



**Wheel small full**  
Wheel 12" x 2" 1/2



**Small wheel**  
12" x 2" 1/2 wheel.



**Large wheel (instead of / in addition)**  
12" x 2" 1/2 wheel.

