








Ursus[™] seat is a modern, safe and functional solution for patients, which prefers a high comfort and stabilization. TILT-IQ© function, Tilt-in-space ClickClack© are some of many practical solutions of URSUS SYSTEM[™]

- Outdoor frame Ursus[™] Means tours and travels in very comfortable conditions either for a child and caregivers. It provides a convenient ride in varied terrain. Swivel front wheels ensure smooth changing of direction.
- Hi-Low frame Ursus Home[™] Is a comfortable and functional solution in home conditions. Adjustable seat height greatly facilitates contact with peers, and push handles make maneuvering, even in tight spaces, much easier.



- Active frame Ursus Active[™] is a mode that gives your child a freedom of movement. Thoughtful solution enables independent movement of the user or with assistance of a caregiver. Now, therapy, learning, and fun will be even more enjoyable and effective.
- Pneaumatic Hi-Low frame Quark[™] Will convert URSUS[™] seat into a comfortable rehabilitation chair. Innovative system of height adjustment en¬sures smooth and easy operation. Included push handle allows easy relocation of the device.



Our passion - your independence!

Akces-MED sp. z o.o. Jasionka 955B 36-002 JASIONKA export@akces-med.com +48 (17) 864 04 77 www.akces-med.com





Ultralight stroller **URSUS™** is designed for daily rehabilitation of children with locomotor disabilities. As far as its construction is concerned, the frame and seat are made of high-strength aluminum, and the chassis is additionally cushioned. **Ursus™** is equipped with **ClickClack**© feature, which means that you can easily swap the seats and adjust the stroller to the most suitable mode. The device has a wide range of adjustments such as width, seat depth, backrest height, backrest tilt angle, infinitely adjustable tilt-in-space function – **TILT-IQ**©, reversible seat. **Ursus™** has also been equipped with spatial, breathable upholstery that prevents excessive sweating as well as bedsores, and more to the point, positively influences the comfort of use. Moreover, the stroller can easily be converted into a lounger so that the child can sleep in it.

The product is available in **3** different sizes with a very wide range of accessories that will help you to meet the individual needs of user.

IT IS RECOMMENDED IN THE FOLLOWING CASES:

- Infantile cerebral palsy different forms.
- Meningomyelocele.
- Brain and spinal cord injuries.
- Muscle diseases connected with paresis and paralysis.
- Genetic disorders connected with paresis.
- Other diseases connected with paresis, paralysis and injuries regarding locomotor system.
- After traumatic brain injury, spine with spinal cord and lower limbs injury.



Stroller is equipped with Light Grey rims. However there is possibility to change the rim's colour to black.



STANDARD EQUIPMENT:

Light aluminium frame | Fully foldable | Wheels with glass fiber, front PU swivel with direction | blockade, Rear inflatable with 2 separate wheel brace | Front shock absorbers | Ergonomic folding system | Seat unit with ClickClack© system for easy plug in/out | 1 piece style height adjustable push handle | Reflective safety marks | Seat unit that can be used forward or rearward facing | Fully reclining backrest to make the stroller as a bed | Height adjustable backrest | Angle adjustable backrest | 5-points safety belts with height adjustment | Adjustable headrest | Anatomical lateral supports | Depth and width adjustable seat | Pneumatic tilt-in-space system TILT-IQ© | Height and angle adjustable footrest | 5 points feet stabilizing belts with heel support and foot | abduction-adduction adjustments | Angle adjustable footplate | Soft upholstery with spine load reducing and ventilating | Soft front bar with reflective safety marks |

Our passion - your independence!

Akces-MED sp. z o.o. Jasionka 955B 36-002 JASIONKA export@akces-med.com +48 (17) 864 04 77 www.akces-med.com





ADVANTAGES:



Stable structure The device ensures adequate stability and safety of Patient.

Reflective elements The product has reflective elements to ensure your safety.



Lightweight construction

The product is made of lightweight and durable materials.



24 months warranty

We offer 24 months warranty on the product. We also provide after sales service.



Five-point straps for keeping safety.



Front swivel wheels for easy maneuverability.



Easy to use

Intuitive way of adjustment and easy assembly.



Movable seat

You can turn the seat easily so that the child may be facing towards caregiver or movement direction.



Tilt-in-space

It is possible to change the tilt angle of the seat.



The stylish and modern upholstery offers the highest comfort.

STANDARD COLORS:

BREATHABLE UPHOLSTERY:





WATERPROOF UPHOLSTERY:



Akces-MED sp. z o.o.

Jasionka 955B 36-002 JASIONKA

Size 1



City Coffee Size 3

OPTIONAL COLORS:



You can change upholstery colours between sizes without additional charges. Production time up to 14 days.

*On special request. Production time ca. 14 days.




Colors may be differ from the presented images. The manufacturer reserves the right to make changes in products without updating the directory. Current offer on our website.

export@akces-med.com








| Symbol | Measurement | Unit | Size 1 | Size 2 | Size 3 |
|--------|---------------------------------------|------|--------|--------|--------|
| А | Seat width | [cm] | 19-25 | 22-30 | 24-34 |
| В | Seat depth | [cm] | 21-26 | 27-35 | 34-42 |
| С | Backrest height | [cm] | 44-55 | 49-59 | 55-67 |
| E | Backrest width | [cm] | 23 | 27 | 32 |
| F | Footrests length | [cm] | 13-29 | 16-40 | 16-43 |
| | Footer lenght | [cm] | 20 | 20 | 20 |
| | Push handle height | [cm] | 100 | 100 | 100 |
| Х | Seat tilt angle | [°] | 90-180 | 90-180 | 90-180 |
| Y | Backrest tilt angle (forward-facing) | [°] | 3-20 | 3-25 | 3-30 |
| Y1 | Backrest tilt angle (rearward facing) | [°] | 3-22 | 3-17 | 3-20 |
| Z | Footrests tilt angle | [°] | 90-180 | 90-180 | 90-180 |
| | Max user weight | [kg] | 35 | 35 | 50 |
| | | | | | |
| | Product dimensions and weight | | Size 1 | Size 2 | Size 3 |
| | Width | [cm] | 59 | 66 | 67 |
| | Widthafter folding | [cm] | 59 | 66 | 67 |
| | Lenght | [cm] | 90 | 90 | 90 |
| | Lenght after folding | [cm] | 90 | 90 | 92 |
| | Height | [cm] | 100 | 100 | 100 |
| | Height after folding | [cm] | 45 | 48 | 53 |
| | Weight | [kg] | 19,5 | 21 | 22,5 |
| | | | | | |

Note: Measurement tolerance range is +/- 1 cm / ½ in.





URSUS HOME[™]

special stroller



URSUS HOME™ is the perfect stroller designed for daily rehabilitation of children with locomotor dysfunction.

Stable construction and an extensive system of supports and belts makes children feel very confidently and safely in **URSUS HOME™**.

Consequently, the small patient is able to take a meal, play, learn and perform the therapy exercises in seating or half-lying position. Caregivers can easily change the seat height. Angle adjustment of the backrest, seat, and footrest allows your child to rest in supine position.

A particular advantage that will be appreciated by therapist is infinitely adjustable tilt-in-space function.

Rich equipment and wide range of adjustment makes **URSUS HOME™** fit easily to the needs and dysfunctions of the child.

The compact construction of the stroller is perfect for small and tight spaces.

Thanks to the ClickClack© system we can quickly detach the seat and attach it to the other products that are in our offer **(Ursus™, Ursus Home™, Ursus Active™)**. The URSUS™ seat is available in **3** sizes.






IT IS RECOMMENDED IN THE FOLLOWING CASES:

- Cerebral palsy various types and forms.
- Brain and spinal cord injuries with lower limbs paralysis and weakness, balance and motor coordination disorders with different etiology.
- Demyelinating disease (e.g. multiple sclerosis SM).
- Meningocele.
- Muscles diseases with weakness and paralysis (dystrophy, myasthenia gravis, myositis).
- Genetic disorders.
- Degeneration and neurological diseases.
- Other condition connected with muscles weakness, paralysis and musculoskeletal disorders.
- After head and spinal cord injuries and spine trauma with paralysis and paresis.

STANDARD EQUIPMENT:

Adjustable headrest | Anatomical lateral supports | 5-points safety belts with height adjustment | Height adjustable backrest | Angle adjustable backrest | Depth and width adjustable seat | Pneumatic tilt-in-space system tilt-iq© | Height and angle adjustable footrest | 5 points feet stabilizing belts with heel support and foot | abduction-adduction adjustments | Soft upholstery with spine load reducing and ventilating | function | Soft front bar with reflective safety marks | Hi-low frame equipped with castors with brakes and push | handle



Akces-MED sp. z o.o. Jasionka 955B 36-002 JASIONKA export@akces-med.com +48 (17) 864 04 77 www.akces-med.com





ADVANTAGES:



Stable structure The device ensures adequate stability and safety of Patient.

Reflective elements The product has refle

The product has reflective elements to ensure your safety.



Lightweight construction

The product is made of lightweight and durable materials.



24 months warranty

We offer 24 months warranty on the product. We also provide after sales service.



Easy to use

Intuitive way of adjustment and easy assembly.



Movable seat

You can turn the seat easily so that the child may be facing towards caregiver or movement direction.



Tilt-in-space It is possible to change the tilt angle of the seat.



Five-point straps for keeping safety.



Easy height adjustment with a lever.



The stylish and modern upholstery offers the highest comfort.

STANDARD COLORS:

BREATHABLE UPHOLSTERY:



You can change upholstery colours between sizes without additional charges. Production time up to 14 days.

WATERPROOF UPHOLSTERY:



City Grass Size 1



Size 2



City Coffee Size 3

OPTIONAL COLORS:









City Carbon*

Akces-MED sp. z o.o.

Jasionka 955B 36-002 JASIONKA

*On special request. Production time ca. 14 days.


Our passion - your independence!

Colors may be differ from the presented images. The manufacturer reserves the right to make changes in products without updating the directory. Current offer on our website.

export@akces-med.com








| Symbol | Measurement | Unit | Size 1 | Size 2 | Size 3 |
|--------|---|------|--------|--------|--------|
| A | Seat width | [cm] | 19-25 | 22-30 | 24-34 |
| В | Seat depth | [cm] | 21-26 | 27-35 | 34-42 |
| С | Backrest height (without he- adrest) | [cm] | 44-55 | 49-59 | 55-67 |
| E | Backrest width | [cm] | 23 | 27 | 32 |
| F | Footrests length | [cm] | 13-29 | 16-40 | 16-43 |
| | Footer lenght | [cm] | 20 | 20 | 20 |
| | Push handle height | [cm] | 86-96 | 86-96 | 86-96 |
| Х | Seat tilt angle | [°] | 90-180 | 90-180 | 90-180 |
| Y | Backrest tilt angle | [°] | 0-25 | 0-25 | 0-25 |
| Z | Footrests tilt angle | [°] | 90-180 | 90-180 | 90-180 |
| | Max user weight | [kg] | 35 | 35 | 50 |
| | Product dimensions and weight | | Size 1 | Size 2 | Size 3 |
| | Width | [cm] | 55 | 55 | 55 |
| | Lenght | [cm] | 78 | 78 | 78 |
| | Height (maximum) | [cm] | 105 | 110 | 115 |
| | Weight | [kg] | 17 | 18,5 | 19,5 |

Note: Measurement tolerance range is +/- 1 cm / $\ensuremath{{1\!\!\!/}}$ in.



URSUS ACTIVE[™]

active special stroller



URSUS ACTIVE™ is a special wheelchair designed for daily rehabilitation of children with locomotor dysfunctions. Furthermore, this wheelchair may be very useful concerning orthopaedic disorders and lower limbs relief needed after surgery.

The wheelchair is equipped with anti-tippers wheels that prevent the patients to tip over backwards. Optionally, **Ursus Active^M** can also be equipped with push handles for a caregiver

Tilt-in-space function allows you to smoothly adjust the tilt angle of the seat which helps you to relief hips by moving the weight of body from pelvis area to the back.

Thanks to the ClickClack© system we can quickly detach the seat and attach it to other products that are in our offer (Ursus™, Ursus Home™, Ursus Active™, Quark™).





IT IS RECOMMENDED IN THE FOLLOWING CASES:

- Cerebral palsy various types and forms.
- Brain and spinal cord injuries with lower limbs paralysis and weakness, balance and motor coordination disorders with different etiology.
- Demyelinating disease (e.g. multiple sclerosis SM).
- Meningocele.
- Muscles diseases with weakness and paralysis (dystrophy, myasthenia gravis, myositis).
- Genetic disorders.
- Degeneration and neurological diseases.
- Other condition connected with muscles weakness, paralysis and musculoskeletal disorders.
- After head and spinal cord injuries and spine trauma with paralysis and paresis.

STANDARD EQUIPMENT:

Light aluminium frame | Large, inflatable rear wheels in sizes from 22" to 26" | equipped with quick-release couplings and strings | 2 independent rear brakes | Front shock absorbers | Seat unit with system ClickClackTM system for easy stick in/out | Fully reclining backrest to make the stroller as a bed (couchb function) | Adjustable backrest | Anatomical lateral supports | 5-points safety belts with height adjustment | Height adjustable backrest | Angle adjustable backrest | Reflective safety marks | Pneumatic tilt-in-space system TILT-IQ© | Height and angle adjustable footrest | 5 points feet stabilizing belts with heel support and foot | abduction-adduction adjustments | Angle adjustable footplate | Breathable upholstery | Soft front bar with reflective safety marks

| our passion - your | ·independence! |
|--------------------|----------------|
|--------------------|----------------|

Akces-MED sp. z o.o. Jasionka 955B 36-002 JASIONKA export@akces-med.com +48 (17) 864 04 77 www.akces-med.com





ADVANTAGES:



Stable structure The device ensures adequate stability and safety of Patient.



Reflective elements

The product has reflective elements to ensure your safety.



Lightweight construction

The product is made of lightweight and durable materials.



24 months warranty

We offer 24 months warranty on the product. We also provide after sales service.



Ergonomic push handle with height adjustment



Podpory przeciwywrotne zapobiegają przed przewróceniem się wózka.



Easy to use

Intuitive way of adjustment and easy assembly.



Movable seat You can turn the seat easily so that the child may be facing towards caregiver or movement direction.



Tilt-in-space It is possible to change the tilt angle of the seat.



Small front casters helps to manouver even in small rooms.

COLORS

BREATHABLE UPHOLSTERY:



Black

Standard



Optional

Black & White





Optional

Grass Optional



Sea Optional

Akces-MED sp. z o.o.

Jasionka 955B 36-002 JASIONKA




Latte Optional

You can change upholstery colours between sizes without additional charges. Production time up to 14 days.

CITY UPHOLSTERY:





Colors may be differ from the presented images. The manufacturer reserves the right to make changes in products without updating the directory. Current offer on our website.


export@akces-med.com

+48 (17) 864 04 77 www.akces-med.com

URSUS ACTIVE


size table





| Symbol | Measurement | Unit | Size 1 | Size 2 | Size 3 |
|------------|--|------|--------|--------|--------|
| А | Seat width | [cm] | 19-25 | 22-30 | 24-34 |
| В | Seat depth | [cm] | 21-26 | 27-35 | 34-42 |
| С | Backrest height (without he- adrest) | [cm] | 44-55 | 49-59 | 55-67 |
| E | Backrest width | [cm] | 23 | 27 | 32 |
| F | Footrests length | [cm] | 13-29 | 16-40 | 16-43 |
| | Footer lenght | [cm] | 20 | 20 | 20 |
| | Push handle height | [cm] | 85 | 85 | 85 |
| Х | Backrest tilt angle | [°] | 90-180 | 90-180 | 90-180 |
| Y | Seat tilt angle | [°] | 0-25 | 0-25 | 0-25 |
| Z | Footrests tilt angle | [°] | 90-180 | 90-180 | 90-180 |
| | Max user weight | [kg] | 35 | 35 | 50 |
| | Product dimensions and weight | | Size 1 | Size 2 | Size 3 |
| | Width | [cm] | 62 | 66 | 70 |
| | Lenght | [cm] | 100 | 100 | 100 |
| | Height | [cm] | 87 | 92 | 87 |
| | Weight | [kg] | 18,5 | 20 | 22,5 |
| Noto: Moor | uromont toloranco rango is 1/ 1 cm / 1/ in | | | | |

Note: Measurement tolerance range is +/- 1 cm / $\ensuremath{\mathscr{V}}_2$ in.





accessories



► ← Safety



Soft cover for tray Soft cover protects the user from hitting the tray. Indicated for patients with impaired limbs coordination.



Led lights (2 pieces) Front and rear LED lights for increased safety.



Lesser side covers for head They provide additional head protection from unwanted contact with stroller's frame. Mounted around the push handle.



Patient lift belt

Robust and ergonomic design as well as high durability facilitates daily care of the patient.



Stabilization



H vest

It acts as a trunk stabilizer. HOLD&PULL© system makes it easier to adjust the vest to the patient body. You can quickly detach the vest if necessary.



4 point vest

It acts as a trunk stabilizer. Wide range of adjustment helps to fit it to the user's needs. Equipped with HOLD&PULL© system.



6 point vest

Wide range adjustable vest that acts as a trunk stabilizer. HOLD&PULL© system makes it easier to adjust the vest to the patient body.



Thigh abduction belts

They can be easily attached to the chair by use of Velcro straps. Their use allows to keep legs in abduction and also prevents the patient from slipping off of the device.



Thigh abduction belts They can be easily attached to the chair by use of Velcro straps. Their use allows to keep legs in abduction and also prevents the patient from slipping off of the device.



Pelvic belt

Wide belt equipped with "HOLD AND PULL" system. It provides better stabilization in the device.



Firm headrest Made of poliurethane foam. Provides with basic support of the head



Head-neck support Perfectly stabilizes the head. Indicated for people with limited head control. You can mount it at different height.

URSUS™

Our passion - your independence!

I URSUS HOME™

A URSUS ACTIVE™

Akces-MED sp. z o.o. Jasionka 955B 36-002 JASIONKA export@akces-med.com +48 (17) 864 04 77 www.akces-med.com



accessories

AKCESmed



Lumbar support It allows better stabilization of the lumbar back.



Trunk stabilizing belt It acts as an additional trunk protection which helps to keep a proper sitting position. The belt is fastened by use of Velcro straps and its size can be adjusted to the user.



Trunk belt

Universal solution applicable to many of our products. It secures the user from falling out of the device. Anatomical shape and Velcro straps allow for a comfortable fit.



Abduction block Shaped abduction block provides with correct seating position. Can be easily removed.



Latteral supports

Width adjustable supports that provide additional side stabilization of the trunk. You can mount it at different height.



Seat narrowing cushions 3cm Made of elastic foam trimmed with soft material. They are attached to the seat upholstery from both sides. Thanks to them you can narrow the seat by 6cm.



Seat narrowing cushions 5cm Made of elastic foam trimmed with soft material. They are attached to the seat upholstery from both sides. Thanks to them you can narrow the seat by 10cm.



BodyMap Stabilization



BodyMap A seat cushion Vacuum seat cushion allows for stabilisation and proper positioning of hips.



BodyMap B Backrest cushion The vacuum back cushion provides comfort and correction for patients with scoliosis and poor posture. Cushion maintains and controls required posture in the stroller.



Vismemo cover for BodyMap cushion Made of thermo-active foam for extra bed sore prevention. Allows for longer and more comfortable use of the cushion.




Cotton cover for BodyMap cushion Made of soft material it offers protection and absorbs moisture.



BodyMap spacer

Aerating-massaging pad to strengthen the prophylactic anti-bed sore function in patients using BodyMap. It can only be used together with the terry cover. The effect of being 'suspended in the air' creates a skin-friendly microclimate.



H URSUS HOME™



Akces-MED sp. z o.o. Jasionka 955B 36-002 JASIONKA export@akces-med.com +48 (17) 864 04 77 www.akces-med.com



accessories





Comfort



Umbrella

Attachable to the stroller's frame, perfectly protects from the sunlight or rain.



Made of plastic. It supports the therapy and allows you to play, eat, and work.



Folding canopy without side panels Protects against both sun and rain.



Folding canopy with side panels Side panels protect not only against sun and rain but also the wind.



Cup holder

An universal holder mounted to the frame. It enables easy and quick access to a bottle of drink.



Headrest cover

It helps to keep the headrest clean. Easy to assemble/disassemble, it can be washed.



Summer sleeping bag Protects legs and trunk. This esthetic body sheet fits the height of the user

body sheet fits the height of the user covering also the footrest of the pushchair.



Winter sleeping bag

Protects legs and the torso area. A very tasteful cover of the body, adjusting to the height of the user, covering the footrest of the chair. Insulated lining provides extra protection against the cold.



Winter hood Can be attached to the sleeping bag. Protects head against frost, wind and moisture.



Plexiglas® Tray

Made of Plexiglas®. It supports the therapy and allows you to play, eat, and work.



Rain Cover An excellent protection from the rain, snow, and wind either for the user and stroller.



Mosquito Net Excellent protection against insects



Wheel covers (4 pcs.) Use it on the dirt tires when you want to take the products indoor or transport in the trunk.

URSUS™

Our passion - your independence!



Decorated spokes guard

Spokes guards 20", made OF polycarbonate 1,5mm, outsider diameter 415mm, insi de diameter 100mm

A URSUS ACTIVE™

Akces-MED sp. z o.o. Jasionka 955B 36-002 JASIONKA export@akces-med.com +48 (17) 864 04 77 www.akces-med.com

I URSUS HOME™



Storage





Storage basket

Enables you to carry small objects while you stroll. It is attached to the base of the stroller.



Women' bag This capacious bag designed for mothers will be perfect for carrying all the necessities.

Stroller cover 'n carry bag

of a stroller.

Perfect solution to protect folded stroller against damage (e.g. during transport). Easily fits the shape and the size



Men's bag Capacious men-style bag for carrying all the necessary things.



DeLux bag

Large, durable and stylish bag useful for carrying necessary things. Made from high-quality materials.



Oxygen tank case

Oxygen tank case can be attached to the stroller by use of Velcro straps with adjustable length. Made to measure to perfectly fit the patient's oxygen tank.





Kid-sit buggypod IO

For convenient and safe transportation of a second child without having to buy a double stroller. The buggypod IO has been designed so that within a few seconds, it can be folded, unfolded, or completely disassembled. 5 point belts seat and folding footrest provide the child with comfort and safety. Using the conversion kit, you are able to convert the kid-sit into an independent three--wheeled stroller. Load capacity up to 20kg.



Conversion kit (for The buggypod kit-sits)

Set of elements thanks to which the Buggypod kit-sit can be converted into a fully functional stroller.



Rain cover for the kid-sits An excellent protection from the rain, snow, and wind either for the user and stroller

Akces-MED sp. z o.o.

Jasionka 955B 36-002 JASIONKA

URSUS™

H URSUS HOME™

A URSUS ACTIVE™

Our passion - your independence!

Colors may be differ from the presented images. The manufacturer reserves the right to make changes in products without updating the directory. Current offer on our website.

export@akces-med.com +48 (17) 864 04 77 www.akces-med.com



accessories





Buggypod PERLE

Buggyboard MAXI



Extension bar

An adapter enabling to mount the Buggypod kid-sits to the strollers of other manufacturers.



Pneaumatic Hi-Low frame Quark Will convert URSUS" seat into a comfortable rehabilitation chair. Innovative system of height adjustment en-sures smooth and easy operation. Included push handle allows easy relocation of the device.



Mounting set

It enables easy assembly of the Buggypod kit-sits to other strollers whose frame diameter is larger than 40mm.



Hi-Low frame Ursus Home

Is a comfortable and functional solution in home conditions. Adjustable seat height greatly facilitates contact with peers, and push handles make maneuvering, even in tight spaces, much easier.



Outdoor frame Ursus

Means tours and travels in very comfortable conditions either for a child and caregivers. It provides a convenient ride in varied terrain. Swivel front wheels ensure smooth changing of direction.



Active frame Ursus Active

is a mode that gives your child a freedom of movement. Thoughtful solution enables independent movement of the user or with assistance of a caregiver. Now, therapy, learning, and fun will be even more enjoyable and effective.

URSUS™

Our passion - your independence!

I URSUS HOME™



Akces-MED sp. z o.o. Jasionka 955B 36-002 JASIONKA export@akces-med.com +48 (17) 864 04 77 www.akces-med.com



accessories





Spare parts



Spacer fabric Breathable and easy to wash. Usefull especially during hot days and high temperature climate.



Wheel small, full silver metallized (Instead of / in addition) Wheel 12" x 2" 1/2



Wheel duże, full, silver metallized (Instead of / in addition) Wheel 12" x 2" 1/2 srebrna felga



Wheel duże, pompowane, silver metallized (Instead of / in addition) Wheel 12" x 2" 1/2 srebrna felga. Możliwość zamówienia również samej dętki lub opony.



Wheel small full Wheel 12" x 2" 1/2



Small wheel 12" x 2" ½ wheel.



Large wheel (instead of / in addition) 12" \times 2" % wheel.



URSUS™

Our passion - your independence!

URSUS HOME™

▲ URSUS ACTIVE[™]

Akces-MED sp. z o.o. Jasionka 955B 36-002 JASIONKA

export@akces-med.com +48 (17) 864 04 77 www.akces-med.com