

## **SMART™** Upright Standing Frame

Represented in Canada by: East West Medical Products Ltd. www.eastwestmedical.ca

Email: info@eastwestmedical.ca

To order / Distributed in Canada by: Assistability
Unit G 7803 – 35 street SE

Calgary AB T2C 1V3

http://www.assistability.ca/

Email: info@assistability.ca


		Tel: 1-800-403	9019/ 403-629-7697 / Fax: 403.800.9900			
Order #:		Date:				
Name: Phone #:		E	mail:			
	Bil	l to:				
Account #:		Company name:				
Street:		Street 2:				
City:	Province:	Postal:				
	Shi	p to:				
Account #:		Company name:				
Street:		Street 2:				
City:	Province:	P	ostal:			
Name:	Phone #:	E	mail:			
STA	NDARD EQUIPMENT (C	ircle size and options bel	ow)			
	Frame with casto	rs and brakes	SIZE 0A: \$2,089			
	<ul><li>Foot plate with fo</li><li>Adjustable knees</li></ul>		SIZE 1: \$2,102			
	<ul><li>Chest brace with chest belt (adjusts</li><li>Pelvic brace with pelvic belt (adjusts</li><li>Tray</li></ul>		'   3			
			SIZE 2: \$2,169			
O O			SIZE 2A: \$2,271			
			SIZE 3: \$2,338			
0			SIZE 3A: \$2,405			
TU VI			SIZE 1: \$2 436			

		ACCESSORIES			
CATEGORY	SYMBOL	EQUIPMENT	PIECES	QTY	
SAFETY	SMD_002	Belt for moving patients AM-P  - Robust and ergonomic design as well as high durability facilitates daily care of the patient.	1	\$312	
	SMD_003	Soft cover for tray - Soft cover protects the user from being hit by the tray. Indicated for patients with impaired limbs coordination.	1	\$276	
	rango wheels are noise	eless, easy to maintain, and do not leave marks or corrode.			
	SMD_015	TANGO Wheel (75 mm) (I)	1	\$126	
	31115_013	TANGO Wheel (75 mm) (A)		\$186	
	SMD_016	TANGO Wheel (125 mm) (I)		\$178	
		TANGO Wheel (125 mm) (A)	1	\$219	
	SMD_020	TANGO Wheel (100 mm) (I)	1	\$152	
	31010_020	TANGO Wheel (100 mm) (A)	1	\$201	
POSITIONING	SMD_111	Independently adjustable knee supports  - Adjustable and ergonomically shaped knee supports perfectly stabilize the lower limbs.	2	\$465	
	Elisures all adequate is	evel of stability and safety for the patient for whom the standard be Longer pelvic belt (I)	eit is too short.	\$276	
	SMD_116	Longer pelvic belt (I)  Longer pelvic belt (A)	1 1	\$474	
	Ensures an adequate le	Ensures an adequate level of stability and safety for the patient for whom the standard be			
			\$276		
	SMD_159	Longer chest belt (I) Longer chest belt (A)	1 1	\$474	
COMFORT	Tray with edge (I)  Made of wood. It allows you to play and perform therapy exercises. Perfect for working or taking meals. Tray edges prevent items from sliding off.		1	\$156	
		(A) ADDITIONAL/(I) INSTEAD			

## **AVAILABLE SIZES**







Cumbal	Measurement	Unit	Size							
Symbol		Unit	0a	1	1a	2	2a	3	3a	4
w	Indicative user height**	[cm]	50 - 70	70 - 100	70 - 100	90 - 130	90 - 130	110 - 170	110 - 170	140 - 200
Α	Pelvic support width	[cm]	22	22	28	28	37	37	44	44
С	Pelvic support height	[cm]	18 - 44	24 - 68	24 - 68	28 - 86	28 - 86	30 - 107	30 - 107	33 - 129
A1	Chest support width	[cm]	20	20	26	26	35	35	42	42
C1	Chest support height	[cm]	41 - 56	56 - 76	56 - 76	71 - 98	71 - 98	86 - 121	86-121	106 - 145
	Chest support perimeter	[cm]	50 - 64	50 - 64	64 - 76	64 - 76	80 - 105	80 - 105	95 - 125	95 - 125
G	Knee supports width	[cm]	8 - 10,5	8 - 10,5	11 - 14	11 - 14	14 – 20	14 - 20	16,5 - 25	16,5 - 25
Н	Knee supports height*	[cm]	11 - 35	11 - 55	12 – 57	13 - 72	14 - 78	14 - 88	14 - 96	14 - 98
К	Footer length	[cm]	17	17	20	20	24	24	30	30
	Max. user weight	[kg]	35	35	55	55	75	75	95	95

Note: Measurement tolerance range is +/- 1 cm /  $\frac{1}{2}$  in.

<sup>\*</sup>Knee supports height depends on tray height.

<sup>\*\*</sup>The value cannot be main factor for selection of product size.