



JORDI™
positioning chair



JORDI™ is a therapeutic chair that is perfect for therapy and fun as well as other daily activities such as learning or having meals. The chair can be used at home, school, kindergarten or in a group of peers. **JORDI™** turns out to be an indispensable tool in the rehabilitation process, both in terms of therapy and social integration.

This therapeutic chair provides the child with proper support and safety, moreover, it can be easily adjusted to meet your child's needs and dysfunctions. Pelvic belt, headrest, height adjustable seat and armrests further enhance comfort and safety of the user.

It is made of user-friendly materials. Available in **4** sizes.



IT IS RECOMMENDED IN THE FOLLOWING CASES:

- Cerebral palsy – various types and forms.
- Brain and spinal cord injuries with lower limbs paralysis and weakness, balance and motor coordination disorders with different etiology.
- Demyelinating disease (e.g. multiple sclerosis – SM).
- Meningocele.
- Muscles diseases with weakness and paralysis (dystrophy, myasthenia gravis, myositis).
- Genetic disorders.
- Degeneration and neurological diseases.
- Other condition connected with muscles weakness, paralysis and musculoskeletal disorders.
- After head and spinal cord injuries and spine trauma with paralysis and paresis.

STANDARD EQUIPMENT:

Adjustable depth of the seat | Adjustable backrest tilt angle | Pelvic belt | Medical upholstery | Adjustable armrests | Lightweight and simple construction with adjustable height



ADVANTAGES:



Stable structure

The product ensures adequate stability and safety of a Patient.



Eco-friendly materials

We use only human and environment friendly materials.



24 months warranty

We offer 24 months warranty on the product. We also provide after sales service.



Easy to use

Intuitive way of adjustment and easy assembly.



Hold and pull©

Our harnesses and pelvic belts are equipped with Hold and pull system that provides fast and comfortable support for your child in the chair.



HexagonPro© Headrest has been designed to create the most comfort for user by supporting and protecting their head.



Medical upholstery used in this product is vapor permeable, waterproof and easy to disinfect.

COLOURS

MEDICUS™ UPHOLSTERY



WATERPROOF



EASY TO CLEAN



STEAM PERMEABLE



ECOLOGICAL




FIRE RETARDANT



MED Sea Standard



MED Blue *Optional



MED Skin *Optional

EXTREME-MED™ UPHOLSTERY



PHTHALATE FREE



ABRASION



DURABLE



ECOLOGICAL



CRIB7



ANTIMICROBIAL & ANTIFUNGAL



STAIN RESISTANT



Citrus Green *Optional



Candy *Optional



Pillarbox *Optional


*Optional colors available for special request. Lead time ca. 14 days

Our passion - your independence!

Colors may differ from the presented images. The manufacturer reserves the right to make changes in products without updating the directory. Current offer on our website.

Akces-MED sp. z o.o.
Jasionka 955B
36-002 JASIONKA

export@akces-med.com
+48 (17) 864 04 77
www.akces-med.com



Symbol	Measurement	Unit	Size 1	Size 2	Size 3	Size 4
A	Seat width	[cm]	25	30	37	45
L	Seat depth	[cm]	23-29	26-34	30-40	35-45
C	Seat height to floor	[cm]	26-30	31-39	36-48	41-55
F	Seat height to footplate	[cm]	10-27	10-36	10-45	10-52
	Back height	[°]	33	39	46	55
	Length of armrest	[cm]	23	25	29	34
	Headrest height	[cm]	37-53	43-59	50-66	59-75
	Max user weight	[kg]	50	65	80	100
	Product weight	[kg]	11	13	16	19

Note: Measurement tolerance range is +/- 1 cm / ½ in.



Comfort



Tray
It helps to keep the headrest clean. Easy to assemble/disassemble, it can be washed. Tray is available in gray/black colour.



Push handles
Helps to transfer the Jordi from place to place. Recommended with JRI_423



Castors with brakes (4 pcs)
You can easily transfer Jordi from one place to another. The castors increase the height of the chair by 6-10cm (it depends on the product size)



JOY™ Table
School table that can be used with Jordi or any other school chair. Height adjustable, easy to clean and very durable. JOY™ table is available in various options.
JOY™ table is available in following options:


- 436/1 – Cut-Out Board
- 436/2 – Plain Board
- 436/3 – Tilt Cut-Out Board
- 436/4 – Tilt Plain-Board



Feet Cover
Comfortable pad, protecting the feet from direct contact with the footplate



Shaped seat cushion
Seat cushion with anatomic cuts for thighs which provides a comfortable position while using the chair.



Shaped backrest cushion
Backrest cushion with an ergonomic shape which has a significant impact on the patient's comfort while using the chair.



Positioning



Hexagon Pro™ Headrest
Individual "ladyfingers" on the Hexagon Pro™ are adjusted to cushion and support any head position. It has been designed to create the most comfort for user by supporting and protecting their head.



Knee abductor
Adjustable and easy to remove. Prevent in legs crossing. Made of soft foam, upholstered with easy to wash medical fabric. Collision with: abduction block



Abduction block
Shaped abduction block provides with correct seating position. Can be easily removed. Collision with: Knee abductor. Restriction: tray height



Lateral supports
Width adjustable supports that provide additional side stabilization of the trunk. You can mount them at different height.



Hip guide
Gives great stabilization in hip area. Can narrow the seat width. Made of MEDICUS™ upholstery.



4 points safety vest
It acts as a trunk stabilizer. Wide range of adjustment helps to fit it to the user's needs. Equipped with HOLD&PULL© system.



H Harness
It acts as a trunk stabilizer. HOLD&PULL© system makes it easier to adjust to the patient's body. You can quickly detach the vest if necessary.



Positioning



2 points foot stabilizing belts
Adjustable straps with Velcro, perfectly stabilize the feet. Made of durable material.



5 points foot stabilizing belts
Adjustable straps with Velcro, perfectly stabilize the feet and ankle. Made of durable material.



Build-up



Folded Footrest
Easy to fold when not needed. Made of durable PU board.