

JUMBO BODYMAP Positioning Chair

Represented in Canada by: East West Medical Products Ltd. www.eastwestmedical.ca Email: info@eastwestmedical.ca

> To order / Distributed in Canada by: Assistability Unit G 7803 - 35 street SE

Calgary AB T2C 1V3

http://www.assistability.ca/

Email: info@assistability.ca

		Tel: 1-800-403-9	9019/ 403-629-769	7 / Fax: 403.800.9900
Order #:		Date:		
Name:	Phone #:		Email:	
	Bill to:			
Account #:		Company name:		
Street:		Street 2:		
City:	Province:		Postal:	
	Ship to:			
Account #:		Company name:		
Street:		Street 2:		
City:	Province:		Postal:	
Name:	Phone #:		Email:	
	PMENT (Circle size a		elow)	
ETHANAS	Frame with castors and brakesFootrest with feet stabilizing belts		C: 1.	
	Seat			Size 1:
	Abduction block Backrest Backrest extension			\$5,341
ATMANNE .			45,511	
·E				
	• 4 points safety vest ('hold & pull' system)		,	Size 2:
	 BODYMAP® A cushion – abduction block shaped cushion BODYMAP® B+ cushion – backrest cushion 		n block	ĆF F20
			\$5,539	
			st cusnion	
The same of the sa	with side stabilization • BODYMAP® DZ cushion – triangular headrest		C: 2.	
	• Reversible electric pump – mobile pump			Size 3:
Aug.	lating and deflatir			\$7,122
	ower limb stabiliz	•		7,,

		ACCESSORIES			
CATEGORY	SYMBOL	EQUIPMENT	PIECES	QTY	
SAFETY	SLB_002	Belt for moving patients AM-P - A solid and ergonomic design as well as high durability make it easier for the caregiver to take care of the patient every day.		\$312	
	SLB_003	Soft cover for tray - It cushions the user from the hard surface of the tray.		\$327	
		TANGO Wheels are noiseless, maintenance-free, non-marking and corrosion free.			
	CLD 045	TANGO Wheel with brake (75 mm) (I) TANGO Wheel with brake (75 mm) (A)	1 -	\$126	
	SLB_015			\$186	
	SLB_016	TANGO Wheel with brake (125 mm) (I) TANGO Wheel with brake (125 mm) (A)	1 -	\$178	
	3LB_010			\$219	
	SLB_020	TANGO Wheel with brake (100 mm) (I)	1	\$152	
	315_020	TANGO Wheel with brake (100 mm) (A)		\$201	
		s as a trunk stabilizer. HOLD&PULL© system makes it easier to adjust the vest to the patient body.			
		vest if necessary.			
POSITIONING SLB_114	SLB_114	H Harness ('hold & pull' system) (I)		\$186	
	H	H Harness ('hold & pull' system) (A)	1	\$238	
BODYMAP POSITIONING	andylfian byetem co	*Cotton Covers: In order to minimize the risk of bedsores, we recommend using a cotton cover along with each BodyMap cushion if the permanent duration of use is longer than 2 hours.			
		**VisMemo Covers: For people using BodyMap cushions for more than 8 hours it is recommended to use a thermo-elastic VISMEMO cover.			
	Control of the Contro	***Aerating-Massaging Pad: It is recommended use of the BodyMap cushion with aerating-massaging pad for people who use BodyMap cushions non-stop longer than 8 hours.			

		Cotton cover for Bodymap® A+ - Sz 1*		\$160
	SLB_306	Cotton cover for Bodymap® A+ - Sz 2*	1	\$178
	_	Cotton cover for Bodymap® A+ - Sz 3*		\$204
		VisMemo cover for Bodymap® A+ - Sz 1**		\$516
	SLB_307	VisMemo cover for Bodymap® A+ - Sz 2**	1	\$595
	_	VisMemo cover for Bodymap® A+ - Sz 3**	_	\$665
		Aerating-massaging pad for Bodymap® A+ - Sz 1***		\$182
	SLB_308	Aerating-massaging pad for Bodymap® A+ - Sz 2***	1	\$204
	_	Aerating-massaging pad for Bodymap® A+ - Sz 3***		\$223
	SLB_318	Cotton cover for Bodymap® B+ - Sz 1*		\$175
		Cotton cover for Bodymap® B+ - Sz 2*	1	\$189
		Cotton cover for Bodymap® B+ - Sz 3*		\$245
		VisMemo cover for Bodymap® B+ - Sz 1**		\$609
	SLB_319	VisMemo cover for Bodymap® B+ - Sz 2**	1	\$717
		VisMemo cover for Bodymap® B+ - Sz 3**		\$873
		Aerating-massaging pad for Bodymap® B+ - Sz 1***		\$245
	SLB_320	Aerating-massaging pad for Bodymap® B+ - Sz 2***	1	\$282
		Aerating-massaging pad for Bodymap® B+ - Sz 3***		\$323
		Cotton cover for Bodymap® - Sz 1*		\$119
	SLB_330	Cotton cover for Bodymap® - Sz 2*	1	\$137
		Cotton cover for Bodymap® - Sz 3*		\$167
		VisMemo cover for Bodymap® DZ - Sz 1**		\$279
	SLB_331	VisMemo cover for Bodymap® DZ - Sz 2**	1	\$301
		VisMemo cover for Bodymap® DZ - Sz 3**		\$346
COMFORT	SLB_401	Tray for manual therapy (1) - Allows for active exercises assisting upper limbs with a particular focus on palm function.	1	\$826
	SLB_403	Tray - Tray made of ecological wood. Great for therapy, daily activities and meals.		\$396
	SLB_407	Tray for nonverbal communication (1) - Set consists of 38 cards that may be put in random order into 24 places of the tray. It is covered with transparent and aesthetic plate made of Plexiglas®. Parents and child can create their own symbols they believe are significant for their communication.	1	\$1,247
STORAGE	SLB_507	Oxygen tank holder - Oxygen tank case can be attached to the stroller by use of Velcro straps with adjustable length. Made to measure to perfectly fit the patient's oxygen tank.	1	Price depends on tank size

SLB_002 BELT FOR MOVING THE PATIENTS AM-P

Size	Waist circumference	How to measure	Thigh circumference	How to measure
	[cm]		[cm]	
S	65 - 75		40 - 44	
М	75 - 85		44 - 48	
L	85 - 97	75	48 - 52	6
XL.	97 - 110		52 - 56	\/
XXL	110 - 125		56 - 60	13
X000L	125 - 145		60 - 64	

Size is determined when ordering.

More Information:

JUMBO positioning chair is designed for maintaining initial upright position of children when sitting. Despite its simple construction, it is a very functional and easy to operate device. It is safe for spastic as well as flaccid children with decreased muscle tone. Use of BodyMap system helps in correcting postural problems and, moreover, it may be very useful regarding orthopedic disorders. The device is equipped with a safety vest with innovative "hold & pull" © system thanks to which the child is both quickly and precisely secured in the chair. Size 3 is equipped with a folding footrest which is very helpful in placing a child in the chair.

Concerning the people and environment, the product has been made of environmentally friendly materials.

The BodyMap system includes:

- BodyMap A+ seat cushion with pelvis stabilization and abduction block – thanks to its appropriate shape, sitting in correct position is even easier and profiled abduction block prevents crossing of the legs.



BodyMap A+ cushion.

- BodyMap B+ backrest cushion with side stabilization of the trunk – it perfectly adapts to the shape of the spine ensuring additional stabilization.



BodyMap B+ backrest.

- BodyMap DZ headrest – provides correct and comfortable support for the head. It easily adapts to the head shape and thanks to the special granules, filling the cushion, you can increase the stabilization from every side.



BodyMap DZ headrest.

- Leg separator with belt – it helps to keep the lower limbs in proper position and stabilize the feet. Furthermore, it prevents crossing of the legs.



- Dual-purpose electric pump — quick and simple way to inflate the cushions, more importantly, you can form the shape of the cushion without removing the patient from the chair or stroller.



Eclectic pump.