

<b>Order #:</b>	<b>Date:</b>
-----------------	--------------

<b>Name:</b>	<b>Phone #:</b>	<b>Email:</b>
--------------	-----------------	---------------

**Bill to:**

<b>Account #:</b>	<b>Company name:</b>
-------------------	----------------------

<b>Street:</b>	<b>Street 2:</b>
----------------	------------------

<b>City:</b>	<b>Province:</b>	<b>Postal:</b>
--------------	------------------	----------------

**Ship to:**

<b>Account #:</b>	<b>Company name:</b>
-------------------	----------------------


<b>Street:</b>	<b>Street 2:</b>
----------------	------------------

<b>City:</b>	<b>Province:</b>	<b>Postal:</b>
--------------	------------------	----------------




<b>Name:</b>	<b>Phone #:</b>	<b>Email:</b>
--------------	-----------------	---------------

<b>Notes:</b>
---------------

**STANDARD EQUIPMENT (Circle size and choose options below)**

	<ul style="list-style-type: none"> <li>• Adjustable seat depth</li> <li>• Waterproof upholstery made of nylon mesh</li> <li>• Adjustable head support</li> <li>• Adjustable trunk belt</li> <li>• Independently adjustable leg belts</li> <li>• Water-resistant adjustable frame equipped with foam covers to avoid scratches on the bathtub</li> <li>• DRY FLOOR TECH™. In all of our bath products the frame is filled with polyurethane foam which prevents water from getting inside.</li> </ul>	<b>Size 0: \$1,072</b>
		<b>Size 1: \$1,072</b>
		<b>Size 2: \$1,105</b>
		<b>Size 3: \$1,141</b>


## ACCESSORIES

CATEGORY	SYMBOL	EQUIPMENT	PIECES	QTY
	The BodyMap cushion can be used generally or during baths as all BodyMap cushions are waterproof! The use of BodyMap cushion makes stabilization of the patient more efficient during bath/shower.			
BODYMAP POSITIONING	AKS_325	Fixing kit for BodyMap cushion	1	\$206
		BodyMap cushion K/F0 (40x120cm)	1	\$2,336
		BodyMap cushion K/F1 (40x130 cm)		\$2,336
		BodyMap cushion K/F2 (50x145 cm)		\$2,491
		BodyMap cushion K/F3 (55x180 cm)		\$2,647
BUILDUP	AKS_804 	AKVO Active waterproof platform for bath chairs - Active waterproof platform for bath devices for disabled children	1	\$951
POSITIONING	AKS_805 	Pelvic suport (Z) - Pelvic support can be equipped instead Pelvic-femoral belt.	1	\$144



An adjustable seat depth helps to fully adapt the product to the patient's height, which provides a comfortable and stable sitting position

## SIZE TABLES



Symbol	Measurement	Unit	Size 0	Size 1	Size 2	Size 3
A	Product width	[cm]	37	42	45	50
B	Seat depth	[cm]	23 - 28	28 - 33	32 - 40	37 - 50
C	Headrest height	[cm]	16	16	19	23
C1	Backrest height	[cm]	31	40	44	50
C2	Lower leg length	[cm]	26	31	36	41
C3	Seat height (from the floor)	[cm]	11 - 44	11 - 44	13 - 49	15 - 54
		[°]	90 - 180	90 - 180	90 - 180	90 - 180
		[°]	90 - 180	90 - 180	90 - 180	90 - 180
		[°]	205	210	210	205
	Max. user weight	[kg]	25	30	50	75
	Indicative user height*	[cm]	< 100	100 - 120	120 - 140	130 - 170

**Note:** Measurement tolerance range is +/- 1 cm / ½ in.

\*The value cannot be main factor for selection of product size.