



# AKVO waterproof bath chair platform

Represented in Canada by: East West Medical Products Ltd.

www.eastwestmedical.ca

Email: info@eastwestmedical.ca

To order / Distributed in Canada by: Assistability

Unit G 7803 – 35 street SE

Calgary AB T2C 1V3

http://www.assistability.ca/

Email: info@assistability.ca

Tel: 1-800-403-9019/ 403-629-7697 / Fax: 403.800.9900

<b>Order #:</b>		<b>Date:</b>	
<b>Name:</b>	<b>Phone #:</b>	<b>Email:</b>	
<b>Bill to:</b>			
<b>Account #:</b>		<b>Company name:</b>	
<b>Street:</b>		<b>Street 2:</b>	
<b>City:</b>	<b>Province:</b>	<b>Postal:</b>	
<b>Ship to:</b>			
<b>Account #:</b>		<b>Company name:</b>	
<b>Street:</b>		<b>Street 2:</b>	
<b>City:</b>	<b>Province:</b>	<b>Postal:</b>	
<b>Name:</b>	<b>Phone #:</b>	<b>Email:</b>	
<b>Notes:</b>			

## STANDARD EQUIPMENT (Choose quantity to order below)



With the lightweight, aluminum and water-resistant base AKVO™ you can install the AKVOSEGO or AKVOLITO bath chairs in just a few seconds. Bath chair installed on the base is a great aid which you can use during the shower or in bath tub. No long is it required to lift and carry child into the bath. Just put child on our base AKVO™ in room, go to bath, take a shower, dry the child with a towel and back to room. The AKVO™ is useful as an mobile mover. The bath chair or deckchair, assembled on the base, is placed higher lessening excessive bending.

**1 Size:  
\$951**

**Qty:**

SPECIAL FEATURES




Product can be adjusted to fit all sizes of AKVOSEGO™ and AKVOLITO™ bath chairs



The belt system allows easy installation of the product on the basis ensuring safety during use



Easy to adjust length



Symbol	Measurement	Unit	Size
A	Base width	[cm]	64
A1	Frame width	[cm]	43
B	Base length	[cm]	93 - 124
B1	Frame length	[cm]	55 - 85
C	Frame height	[cm]	38
	Maximum user weight	[kg]	83
	Weight net	[kg]	4,4
	Weight gross	[kg]	6,5

