

HIPPO+ Special needs stroller

Represented in Canada by: East West Medical Products Ltd. www.eastwestmedical.ca Email: info@eastwestmedical.ca

To order / Distributed in Canada by: Assistability

Unit G 7803 – 35 street SE

Calgary AB T2C 1V3 http://www.assistability.ca/

Email: info@assistability.ca

Tel: 1-800-403-9019/ 403-629-7697 / Fax: 403.800.9900					
Order #:		Date:			
Name:	Phone #:		Email:		
	Bill to:				
Account #:		Company name	:		
Street:		Street 2:			
City:	Province:		Postal:		
	Ship to:				
Account #:		Company name	:		
Street:		Street 2:			
City:	Province:		Postal:		
Name:	Phone #:		Email:		
Notes:					
STANDARD EQUIPME	NT (Circle size and ch	noose options bel	ow)		
	Foldable frame with wheels and rear brake Swivel front wheels 360° Drive system with shock absorbers Under seat storage basket Adjustable footrest with supporting belt Removable abduction block Reclining backrest Adjustable lateral supports 5 points safety belts Adjustable headrest Adjustable folding canopy Detachable and reversible seat unit Adjustable push handle made of ECO leather Leg cover Pumped wheels (front/rear) Front rail Rain cover Bag Pump Reduction set consisting of additional seat size 1) Lumbar pad (size 1)		brake	SIZE 1: \$2,642	
				Size 2: \$2,796	

ACCESSORIES					
CATEGORY	SYMBOL	EQUIPMENT	PIECES	QTY	
SAFETY	HPO_001	LED lights - Front and rear LED lights for increased safety.	1	\$119	
	HPO_007	Infant car seat 0-10kg - Size 1 only	1	\$1,078	
	HPO_101	Thigh abduction belts - Can be easily attached to the chair by use of Velcro straps. Assists to keep legs in abduction and prevent sliding forward.	1	\$256	
		tabilizer. HOLD&PULL© system makes it easier to adjust the ver	st to the patie	nt body. You	
POSITIONING	can quickly detach the vest if necessary. HPO_114 H Harness ('hold & pull' system) (I)		1	\$186	
		H Harness ('hold & pull' system) (A)		\$238	
	HPO_124	Lumbar support (size 2) - It allows better stabilization of the lumbar region of the back.	1	\$152	
	The 4 points safety vest acts as a trunk stabilizer. Wide range of adjustment helps to fit it to the user's needs. Equipped with HOLD&PULL© system. SIZE Chest measurements (CM): S-42 – 56 / M-53 – 64 / L-60 – 70 XL -66 – 82 XXL-80 – 97				
	HPO_130	, ,		\$197	
		4 points safety vest (A)	1	\$256	
	HPO_137	Reduction set (with additional seat - size2)	1	\$275	
	HPO_138	Headrest with belt - Mounted in place of the standard headrest. It can be precisely adjusted to the user's needs. Soft head belt additionally stabilizes the forehead.	1	\$516	

		Death and decreased a con-		
	HPO_151	Depth reducing pad 3 cm - Reduce seat depth by 3 cm. Made of foam covered with upholstery.	1	\$320
	HPO_152	Depth reducing pad 5 cm - Reduce seat depth by 5 cm. Made of foam covered with upholstery.	1	\$346
	HPO_402	Umbrella - Attachable to the stroller's frame adding protection from sunlight and rain.	1	\$375
	HPO_409	Mosquito net - Excellent protection against insects	1	\$130
	HPO_410	Headrest cotton cover - Helps keep the headrest clean/ can be washed. Easy to		\$30
		add/remove.	2	\$52
COMFORT	HPO_411	Cup holder - A universal holder mounted to the frame. It enables easy and quick access to drinking bottle.	1	\$163
	HPO_414	Tray Plexiglas - Supports therapy and allows you to play, eat, and perform daily activities.	1	\$702
		Covers dirty tires when taking stroller indoors or transporting in trunk or in vehicle.		
	HPO_415	Rear wheel covers	2	\$26
	HPO_442	Front wheel covers	2	\$19
	HPO_422	Carrycot - Designed to carry the child in a lying position. It can be mounted to the frame without additional adapters.	1	\$936
	HPO_429	Soft footrest cover	1	\$230

	HPO_601	Kid-sit Buggypod IO - For convenient and safe transportation of a second child without having to buy a double stroller. The buggypod IO has been designed so that within a few seconds, it can be folded, unfolded, or completely disassembled. 5 point belts seat and folding footrest provide the child with comfort and safety. Using the conversion kit, you are able to convert the kid-sit into an independent three-wheeled stroller. Load capacity up to 20kg.	1	\$1,717
BUILD-UP	HPO_603	Extension bar for Buggypod - An adapter enabling to mount the Buggypod kid-sits to the strollers of other manufacturers.	1	\$275
	HPO_604	Mounting set 70 for Buggypod - An adapter enabling to mount the Buggypod kid-sits to the strollers of other manufacturers.	1	\$275
	HPO_605	Rain cover for Buggypod IO - An excellent protection from the rain, snow, and wind either for the user and stroller	1	\$312
	HPO_606	Adaptation set for Buggypod IO - Set which allows the Buggypod kit-sit to be converted into a fully functional stroller.	1	\$587
STORAGE	HPO_506	Carrier bag for the stroller - An excellent solution to protect the folded stroller from damages (e.g. during the transportation). It perfectly fits the stroller shape and dimensions.	1	\$516
	HPO_509	Organizer bag - A handy bag attached to the handle of the stroller. Specially designed compartments make it easy to organize items, keeping the most essential items at hand. The organizer has two useful spaces for liquid containers.	1	\$234


	HPO_714	Rear wheel	1	\$141
	HPO_715	Rear tire	1	\$59
	HPO_716	Rear innertube	1	\$48
SPARE PARTS	HPO_717	Front wheel	1	\$141
	HPO_718	Front tire	1	\$59
	HPO_719	Front innertube	1	\$48
	HPO_731	Front rim	1	\$56
	HPO_732	Rear rim	1	\$56

COLOURS



^{*}Due to a previous information about discontinue Mint and Chilli colours for Hippo+ buggy we would like to inform that for size 2 these colours are not available anymore. For size 1 are still available while stocks last

SIZE TABLES



Symbol	Measurement	Unit	Size 1	Size 2
w	Indicative user height**	[cm]	70 - 100	90 - 120
Α	Seat width	[cm]	20*/ 31	34
L	Seat depth	[cm]	28	32
С	Headrest height	[cm]	40 - 51*/ 44 - 54	55 - 64
E	Backrest width	[cm]	31	34
F	Footrests length	[cm]	20*, 24*/16, 20	22, 26, 30
	Backrest height	[cm]	40*/ 45	55
L1	Seat depth to the abduction block	[cm]	18	22
	Footplate length	[cm]	14	18,5
Х	Backrest tilt angle	[°]	90 – 140	90 – 140
Υ	Seat tilt angle	[°]	10	10
Z	Footrests tilt angle	[°]	90 – 180	90 – 180
	Maximum user weight	[kg]	20	30

Note: Measurement tolerance range is +/- 1 cm / ½ in.

^{*}Measurement starts from narrowing pad.

^{**}The value cannot be main factor for selection of product size.