

ULISES EVO Special Needs Stroller

Represented in Canada by: East West Medical Products Ltd.
www.eastwestmedical.ca
Email: info@eastwestmedical.ca

To order / Distributed in Canada by: Assistability

Unit G 7803 – 35 street SE

Calgary AB T2C 1V3

http://www.assistability.ca/

		Tel: 1-800-403-	9019/ 403-629-7697 / Fax: 403.800.9900
Order #:		Date:	30137 313 7
Name:	Phone #:	1	Email:
	Bill to:		
Account #:		Company name:	
Street:		Street 2:	
City:	Province:		Postal:
	Ship to:		
Account #:		Company name:	
Street:		Street 2:	
City:	Province:		Postal:
Name:	Phone #:		Email:
Notes:			
STANDARD E	QUIPMENT (Circle size	and choose options b	below)
	5 points safety belts	th 2 points feet stabili	Size 1: \$2,374
	Bill to: Company name: Street 2: Post		
	BodyPrint seat© belt sension of the seat uph	system adjusting the olstery and backrest	Size 2: \$2,572
	Backrest extension Headrest	ne seat	Size 2A: \$2,769
	SIZE 1 GREEN GRA SIZE 1A HOT RED, SIZE 2 BLUE SEA, SIZE 2A BLUE WAV SIZE 3 GREY VULC	Έ,	Size 3: \$2,835

	ACCESSORIES						
CATEGORY	SYMBOL	EQUIPMENT	PIECES	QTY			
	ULE_001	LED lights - Front and rear LED lights for increased safety.	1	\$119			
	ULE_002	Belt for moving patients AM-P - Robust and ergonomic design as well as high durability facilitates daily care of the patient.	1	\$312			
SAFETY	ULE_003	Soft cover for tray - Soft cover protects the user from being hit by the tray. Indicated for patients with impaired limbs coordination.	1	\$305			
	ULE_013	Safety rail - Makes the user safe in the stroller. Protects against falling out.	1	\$312			
	ULE_021	Knee protectors - Protects the knee and crus from hit and/or chafes. Made of soft foam covered with spacer, breathable fabric.	2	\$149			
DOSITIONING	ULE_101	Thigh abduction belts - They can be easily attached to the chair by use of Velcro straps. Their use maintains legs in abduction and also prevents the patient from slipping off of the device.	1	\$256			
POSITIONING	ULE_107	Pelvic brace with pelvic belt (adjustable) - Wide belt equipped with "HOLD AND PULL" system. It provides better stabilization in the device.	1	\$234			

ULE_113	Trunk stabilizing belt -It acts as an additional trunk protection which helps to keep a proper sitting position. The belt is fastened by use of Velcro straps and its size can be adjusted to the user.	1	\$561
ULE_114	The H Harnes acts as a trunk stabilizer. HOLD&PULL© system make the patient body. You can quickly detach the vest if necessary.	s it easier to a	djust the vest to
ŭ	H Harness ('hold & pull' system) (I)	1	\$186
	H Harness ('hold & pull' system) (A)		\$238
ULE_115	Calf belt - It provides additional protection for lower limbs. Mounted on the front side limits the legs extension. Mounted on the rear side limits the legs bending.	1	\$212
ULE_121	Head-neck support - Stabilizes head. Can be mounted at different heights on the base.	1	\$327
ULE_124	Lumbar support - It allows better stabilization of the lumbar back region.	1	\$305
ULE_125	The 6 points safety vest is a wide range adjustable vest that acts as HOLD&PULL© system makes it easier to adjust the vest to the patients.		zer.
	SIZE Chest measurements (CM) S 42 – 56 / M 53 – 64 / L 60 – 70) / XL 66 – 82 ,	
	6 points safety vest (I)	1	\$334
ULE_126	6 points safety vest (A)		\$446
OLL_120	Trunk belt - Universal solution applicable to many of our products. It secures the user from falling out of the device. Anatomical shape and Velcro straps allow for a comfortable fit.	1	\$327
ULE_128	Adjustable abduction block (3) - ELEMENTS ULE_136, ULE_128 AND ULE_129 CANNOT BE INSTALLED TOGETHER	1	\$256

ULE_129	Abduction-adduction thigh belts - Allow positioning of legs with adjustable length straps for range of desired fit.	1	\$516
ULE_130	It acts as a trunk stabilizer. Wide range of adjustment helps to fit it with HOLD&PULL© system.	to the user's n	eeds. Equipped
	4 points safety vest (I)		\$197
ni res	4 points safety vest (A)	1	\$256
ULE_132	Elastico headrest - Covered with thermo-active foam easily fits the shape of a head. Offers more comfortable support and better stabilisation. Attached in the place of a standard headrest.	1	\$212
ULE_133	Firm foam headrest - Covering the entire width of the backrest. Increases the safety of user.	1	\$223
ULE_134	Narrowing pads 6 cm (2) - Made of elastic foam trimmed with soft material. They are attached to the seat upholstery from both sides. Narrow the seat by 6 cm total.	2	\$353
ULE_136	Thigh belts - Helps maintain a proper sitting position. Velcro-adjustable straps as well as ergonomic design allow for a better fit for the patient. (elements ule_136, ule_128 and ule_129 cannot be installed together. elements ule_136 & ule_129 are mutually exclusive with ule_412 element.)	1	\$353
ULE_137	Narrowing pads 10 cm (2) - Made of elastic foam trimmed with soft material. They are attached to the seat upholstery from both sides. Narrow the seat by 10 cm total.	2	\$398

	ULE_147	5 points feet stabilizing belts (I)	2	\$171
	ULE_148	High lateral supports	2	\$624
	Protects the hip area 3 cm on each side.	a on both sides of the trolley. To the surface of the material is mount	ed padding nar	rowing the seat,
	ULE_149	Rigid pelvic supports with narrowing pads 6 cm	2	\$624
	ULE_150	Rigid pelvic supports with narrowing pads 10 cm	2	\$702
	ULE_153	Low lateral supports	2	\$468
	ULE_402	Umbrella - Attachable to the stroller's frame, protects from sunlight or rain.	1	\$375
COMFORT	ULE_403	Tray - Easy to attach/detach and can be washed.	1	\$430
	ULE_405	Folding canopy - Protects against sun and rain.	1	\$516

ULE_408	Rain cover - Protection from the rain, snow, and wind for the user and product. Elements ule_408, ule_409 could be used only with folded canopy ule_405.	1	\$398
ULE_409	Mosquito net - Excellent protection against insects Elements ule_408, ule_409 could be used only with folded canopy ule_405.	1	\$93
ULE_410	Headrest cotton cover	1	\$30
	- It helps to keep the headrest clean. Easy to assemble/disassemble, it can be washed.	2	\$52
ULE_411	Cup holder - Universal holder mounted to the frame enables quick access to a bottle or cup.	1	\$163
ULE_412	Elastico seat cushion - It precisely adapts to the patient's body. Made of memory foam covered with breathable and spatial fabric. It is indicated in antibedsore prophylaxis. (elements ule_136 & ule_129 are mutually exclusive with ule_412 element)	1	\$282
ULE_413	Elastico backrest cushion - It perfectly adapts to the patient's body. Made of dense shape memory foam covered with breathable and spatial fabric. It is indicated in anti-bedsore prophylaxis.	1	\$375
ULE_415	Rear wheel covers - For use on dirt tires to take stroller indoors or transport in trunk.	2	\$26

ULE_416	Summer footmuff - Protects legs and trunk. This esthetic body sheet fits the height of the user and covering footrest of the pushchair.	1	\$587
ULE_417	Winter footmuff - Protects legs and the torso area. A very tasteful body cover, adjusts to the height of the user, covering the footrest of the chair. Insulated lining provides protection against the cold.	1	\$843
ULE_418	Hood for winter footmuff - Can be attached to the winter footmuff. Protects head against frost, wind and moisture.	1	\$305
ULE_421	Shaped seat cushion (pommel shape) - Made of dense shape memory foam covered with breathable and spatial fabric. It is indicated in anti-bedsore prophylaxis (Elements ule_136 & ule_129 are mutually exclusive with ule_412 element)	1	\$375
ULE_419	Shaped seat cushion (thigh shape) - Anatomical notch on the thighs improves sitting comfort. (elements ule_136 & ule_129 are mutually exclusive with ule_412 element)	1	\$234
ULE_429	Soft footrest cover - Comfortable pad that protects the feet from direct contact with the footplate	1	\$230
ULE_432	Footrest angle adjustment - The adaptive set allows you to incline the footplate to the desired angle.	2	\$186
ULE_433	Tip assist (lifting) bar - Easier to tilt the wheelchair, entrances to the curbs, steps and other elements of the vertical surface.	1	\$156
_	f the footrest length. The installation of the front wheels blockade is a	required to av	oid the collision

between the footrest and the front wheels while turning the stroller.

		Extended footrest (size 1 - 6cm; size 1a,2 - 8 cm; size 2a, 3 - 10 cm) with front swivel wheels brake (I)		\$334
	ULE_437	Extended footrest (size 1 - 6cm; size 1a,2 - 8 cm; size 2a, 3 - 10 cm) with front swivel wheels brake (A)	1	\$464
	ULE_440	Directional lock for front wheels - Used instead of a standard fork. It allows the front wheels to be locked straight ahead.	2	\$186
	ULE_439	Hand brake(4) - The hand brake makes it possible to stop the stroller quickly and securely. It also eases handling the stroller. The hand brake is equipped with a four-level blocking which functions as a parking brake as well.	1	\$669
	ULE_442	Front wheel covers - Easy to assemble and practical to use. They protect against transferring dirt from tires to indoor spaces.	1	\$19
	ULE_502	Bag DeLux - A large, durable, stylish bag will fit a lot of useful things. Made of high quality materials.	1	\$702
STORAGE	ULE_503	Lady's style bag - A spacious bag made for mums, it will fit a lot of useful things.	1	\$327
STORAGE	ULE_504	Man's style bag - The spacious bag made in men's style will fit a lot of useful trinkets.	1	\$327
	ULE_505	Under seat storage basket	1	\$420

	_	,		
	ULE_506	Carrier bag for the stroller - The perfect solution to protect the folded trolley against damage (eg during transport). Adapted to the size and shape of the trolley.	1	\$516
	ULE_507	Oxygen tank holder - Fastened with adjustable length straps. It is made to measure, so that it perfectly matches the oxygen tank dimensions.	1	Price depends on tank size
	ULE_509	Organizer bag - A handy bag attached to the handle of the trolley. Specially designed compartments make it easier to stay organized, keeping the most essential items at hand. The organizer has two useful openings for a cup or bottle.	1	\$234
	ULE_723	Upholstery comfort (A) - Breathable, air-permeable upholstery with spatial material.	1	\$1,011
	ULE_701	Rear inflatable wheel	1	\$186
	SPARE PARTS ULE_704	Rear PU wheel (I)		\$186
SPARE PARTS		Rear PU wheel (A)		\$282
	ULE_707	Front inflatable wheel (I)	1	\$186
		Front inflatable wheel (A)	1	\$282
	ULE_711	Front PU wheel	1	\$282

Attention, some accessories mutually exclude each other. Please refer to the table below.

	ULISES EVO														
ELEMENT NUMBER	013	128	129	134	136	137	148	149	150	153	403	405	408	409	412
013															
128															
129															
134							1			1					
136															
137							2			2					
148				1		2									
149															
150															
153				1		2									
403															
405													3	3	
408												3			
409												3			
412															

ACCESSORY CAN BE MOUNTED UNDER SPECIFIC CONDITIONS

1. ELEMENT ULE_134 AVAILABLE ONLY WITH ULE_148 OR ULE_153

2. ULE_134 AND ULE_137 ARE AVAILABLE ONLY WITH ULE_148 OR ULE_153

3. ULE_408 AND ULE_409 AVAILABLE ONLY WITH ULE_405

SIZE TABLES



Symbol	Measurement	Unit	Size 1	Size 1a	Size 2	Size 2a	Size 3
w	Indicative user height*	[cm]	80 - 120	100 - 150	100 - 150	130 - 180	130 - 180
Α	Seat width	[cm]	31	31	37	37	43
L	Seat depth	[cm]	22-29	29-37	29-37	35-45	35-45
С	Backrest height	[cm]	40-70	43-80	43-80	51-90	51-90
E	Backrest width	[cm]	28	28	34	34	40
F	Footrests length	[cm]	19 - 25	24 - 32	24 - 32	30 - 40	30 - 40
C1	Backrest height (without backrest extension)	[cm]	46	54	54	62	62
C2	Backrest height (with backrest extension)	[cm]	61 – 70	69 – 80	69 – 80	77 – 90	77 – 90
	Footplate length	[cm]	19	19	21	21	25
	Push handle height	[cm]	102	102	102	102	102
Υ	Seat tilt angle	[°]	13	13	13	13	13
X	Backrest tilt angle	[°]	90 – 150	90 – 150	90 – 150	90 – 150	90 – 150
Z	Footrests tilt angle	[°]	90 – 180	90 – 180	90 – 180	90 – 180	90 – 180
	Maximum user weight	[kg]	55	55	75	75	75

Note: Measurement tolerance range is +/- 1 cm / $\frac{1}{2}$ in. *The value cannot be main factor for selection of product size.