

URSUS ACTIVE transferrable seat wheelchair

Represented in Canada by: East West Medical Products Ltd. www.eastwestmedical.ca Email: info@eastwestmedical.ca

To order / Distributed in Canada by: Assistability Unit G 7803 - 35 street SE

Calgary AB T2C 1V3

http://www.assistability.ca/

Email: info@assistability.ca

		Tel: 1-800-403	-9019/ 403-629-76	9//Fax: 403.800.9900	
Order #:		Date:			
Name:	Phone #:		Email:		
	Bil	l to:			
Account #:		Company name:			
Street:		Street 2:			
City:	Province:		Postal:		
	Shi	p to:			
Account #:		Company name:			
Street:		Street 2:			
City:	Province: Postal:		Postal:		
Name:	Phone #:		Email:		
Notes:					
-		size and choose options	below)		
	 Light aluminium frame Large, inflatable rear wheels in sizes from 22" to 26" with quick-release axles 2 independent rear brakes Front shock absorbers Seat unit with ClickClack™ system for easy in/out Fully reclining backrest 		Size 1: \$6,462		
	ckrest teral supports y belts with height adjustment able backrest ble backrest ety marks		Size 2: \$6,858		
	 Pneumatic tilt Height and ar 5 points feet s foot abduction- Breathable up BLACK&WHIT 	-in space system TILT-l ngle adjustable footrest stabilizing belts with hee adduction adjustments pholstery ATMO™BLAC	el support and	Size 3: \$7,517	


		ACCESSORIES		
CATEGORY	SYMBOL	EQUIPMENT	PIECES	QTY
USA_001		LED lights - Front and rear LED lights for increased safety.	1	\$119
SAFETY	USA_003 USA_101	Soft cover for tray - Soft cover protects the user from being hit by the tray. Indicated for patients with impaired limbs coordination.	1	\$327
	USA_101	Thigh abduction belts - Can be easily attached to the chair by use of Velcro straps. Keeps legs in abduction and prevents sliding.	1	\$256
	USA_107	Pelvic brace with pelvic belt (adjustable) - Wide belt equipped with "HOLD AND PULL" system. It provides better stabilization in the device.	1	\$256
POSITIONING	USA_113	Trunk stabilizing belt - Acts as an additional trunk stabilizer helping assist the sitting position. The belt is fastened by use of Velcro straps and its size can be adjusted to the user.	1	\$516
	USA_114	The H Harness acts as a trunk stabilizer. HOLD&PULL© system makes the patient body. You can quickly detach the vest if necessary.	it easier to a	djust the vest to
		H Harness ('hold & pull' system) (I)	1	\$186
		H Harness ('hold & pull' system) (A)	1	\$238
USA_121 USA_124	USA_121	Head-neck support - Stabilizes the head. Indicated for people with limited head control. You can mount it at different heights.	1	\$327
	USA_124	Lumbar support - Allows better stabilization of the lumbar back area.	1	\$420

USA_125	Dints safety vest (A) nk belt hiversal solution applicable to many of our products. It secures his puser from falling out of the device. Anatomical shape and his straps allow for a comfortable fit. In the straps allow for a comfortable fit.		
	6 points safety vest (I)	1	\$334
Co y les	6 points safety vest (A)	_	\$446
USA_126	Trunk belt - Universal solution applicable to many of our products. It secures the user from falling out of the device. Anatomical shape and Velcro straps allow for a comfortable fit.	1	\$327
USA_127	Chest laterals - Stabilizes the user's upper body. The height can be easily adjusted.	2	\$420
USA_128	Abduction block (3) - Shaped abduction block provides with correct seating position. Can be easily removed. USS_128 and USS_129 cannot be installed together	1	\$256
USA_130	The 4 points safety vest cts as a trunk stabilizer. Wide range of adjustiuser's needs. Equipped with HOLD&PULL© system.	ment helps to	o fit it to the
	4 points safety vest (I)		\$197
	4 points safety vest (A)	g position. be installed 1 \$256 range of adjustment helps to fit it to the 1 \$197 1 \$256	
USA_131	Firm foam headrest - Made of polyurethane foam. Provides support of the head across a wide area.	1	\$632
USA_134	Narrowing pads 6 cm - Made of elastic foam trimmed with soft material. They are attached to the seat upholstery from both sides. Narrow the seat by 6 cm total.	2	\$294
USA_137	Narrowing pads 10 cm - Made of elastic foam trimmed with soft material. They are attached to the seat upholstery from both sides. Narrow the seat by 10 cm total.	2	\$294

			T	
	USA_403	Tray - Easy to assemble/disassemble; can be washed.	2 \$52 1 \$163 1 \$653 is adapted to the size of the wheel 2 \$1,092 2 \$1,170 \$1,326 ith 1 \$230	
	USA_410	Headrest cotton cover (4)	1	\$30
		- Helps to keep headrest clean. Easy to assemble/disassemble and can be washed.	2	\$52
	USA_411	Cup holder - Universal holder mounted to the frame. It enables easy and quick access to bottle or cup.	1	\$163
COMFORT	USA_414	Tray Plexiglass® (2) - Easy to add/remove and can be washed.	1	\$653
	USA_426	Wheel spoke cover, 1.5 mm polycarbonate material. The cover is ad for the size of the unit.	apted to the s	ize of the wheel
		Decorated spokes guards – Sz 1 (2)		\$1,092
		Decorated spokes guards – Sz 2 (2)	2	\$1,170
		Decorated spokes guards – Sz 3 (2)		\$1,326
	USA_429	Soft footrest cover - Comfortable pad that protects the feet from direct contact with the footplate.	1	\$230
	USA_502	Bag DeLux - Large, durable and stylish bag useful for carrying necessary things. Made from high-quality materials.	1	\$702
STORAGE	USA_503	Lady's style bag - This spacious bag designed for mothers to carry all necessities.	1	\$327

USA_504 Man's style bag - Spacious men-style bag for carrying all necessary items. 1 \$327	
Oxygen tank holder - Oxygen tank case can be attached to the stroller by use of Velcro straps with adjustable length. Made to measure to perfectly fit the patient's oxygen tank. 1 depends o tank size	
Organizer bag - A handy bag attached to the stroller bar. Specially designed compartments make it easy to keep items in order and the most essential items at hand. The organizer has two openings for a cup or bottle. \$234	
URSUS Hi-Low Frame (10) - Comfortable and functional solution for in-home use. Adjustable seat height facilitates contact with peers, and push handles make maneuvering, even in tight spaces, much easier.	
URSUS Outdoor Frame - Switch seat easily from indoor to outdoor units. Travel outside in comfort and ease for a child and caregivers. It provides a convenient ride in varied terrain. Swivel front wheels ensure smooth directional changes.	
URSUS seat - The URSUS ™ seat is a modern safe and functional solution ensuring stability and comfort during use. The TILT-IQ © bucket function , the ClickClack © interleaved seat is one of many practical solutions of the URSUS ™ system .	
USA_623 QUARK Hi-Low frame (10) Turns the URSUS ™ seat into a comfortable rehabilitation system. Innovative height adjustment system guarantees smooth and easy operation. The handle attached to the base enables easy movement of unit. Waterproof upholstery Hydro™ is made of soft and comfortable material resistant to water and moisture. This upholstery is very	

	durable and easy to clean and disinfect. This makes it perfect for indoor use when eating and drinking fluids, and outdoors during unfavorable weather conditions such as rain or snow. Hydro™ is made of environmental friendly ecological materials.				
unia					
	USA 700	Waterproof upholstery HYDRO™ (I)	1	\$0	
	U3A_700	Waterproof upholstery HYDRO™ (A)	1	\$991	
The	The breathable upholstery Atmo™ is made of elastic, soft material. It is characterized by high durability and is resistant to stretching.				
Usin	Using three-dimensional material allows air to flow, which is crucial in case of pressure sores prevention. Also, it protects against				
bact	bacterias and mould. Atmo™ is easy to clean with mild cleaning agents. Upholstery is made of environmental friendly materials.				
Atm	Atmo™ is a high level of comfort and safe for the skin.				
	USA 720	Breathable upholstery ATMO ™ (I) (1)	1	\$0	
	U3A_720	Breathable upholstery ATMO ™ (A) (1)	1	\$991	



STANDARD COLORS

You can read more about our upholstery <u>here</u>

ATMO ™ BREATHABLE UPHOLSTERY











Black

OPTIONAL COLORS

ATMO ™ BREATHABLE UPHOLSTERY



Red **

HYDRO ™ WATERPROOF UPHOLSTERY











/ Simple_tooltip]








City Grass *


City Rose *

City Coffee *

Color upholstery shades may differ from those shown in the pictures.

- * There is a possibility of changing colors between sizes without additional payment. Waiting time 14 days.
- ** Optional colors are available on special order, waiting time up to 30 days.





Symbol	Description	Unit	Size 1	Size 2	Size 3
IN	approximate user height *	[Cm]	80 - 100	100 - 130	130 - 150
AND	seat width	[Cm]	19-25	22 - 30	24 - 34
L	seat depth	[Cm]	21 - 26	27 - 35	34 - 42
c	backrest height (without headrest)	[Cm]	44 - 55	49 - <mark>5</mark> 9	55 - 67
E	backrest width	[Cm]	23	27	32
F	the length of the footrests	[Cm]	13 - 29	16 - 40	16 - 43
	the length of the foot	[Cm]	20	20	20
	height of the guide handle	[Cm]	85	85	85
	Diameter of wheels	[Cm]	22 "	24 "	26 "
×	reclining back angle	[°]	90 - 180	90 - 180	90 - 180
Υ	seat inclination angle	[°]	0-25	0-25	0-25
FROM	the angle of tilt of the footrests	[°]	0 - 90	0 - 90	0 - 90
	max. user's weight	[Kg]	35	35	50

Note: The tolerance of dimensions contained in the table is +/- 1 cm.

^{*} Value should not be the basis for selecting the device.