


ENDOCH



MEDICAL

SURGERY



YOU WILL SEE THINGS DIFFERENTLY

REUSABLE CONSUMABLES

Confidentiality:

This document contains information that is confidential and proprietary property of Endoch Medical.

Neither this document nor the information therein may be reproduced, used or distributed to or for the benefit of any third party without prior written consent of Endoch Medical.

HANDLES (LAPAROSCOPY)

Tip



Shaft



STANDARD HANDLE



WORKING LENGTH	HANDLE
220 mm	EAI-0395
330 mm	EAI-0395
450 mm	EAI-0395

With Ratchet

WORKING LENGTH	HANDLE
220 mm	EAI-0395R
330 mm	EAI-0395R
450 mm	EAI-0395R

PREMIUM HANDLE




WORKING LENGTH	HANDLE
220 mm	EAI-0396
330 mm	EAI-0396
450 mm	EAI-0396

With Ratchet


WORKING LENGTH	HANDLE
220 mm	EAI-0396R
330 mm	EAI-0396R
450 mm	EAI-0396R

5mm GRASPING AND DISSECTING FORCEPS (LAPAROSCOPY)




**GRASPING FORCEPS,
POINTED**

WORKING LENGTH	COMPLETE
220 mm	EAI-0456R
330 mm	EAI-0457R
450 mm	EAI-0458R




KELLY GRASPING FORCEPS

WORKING LENGTH	COMPLETE
220 mm	EAI-0492R
330 mm	EAI-0493R
450 mm	EAI-0494R




**LONG DELPHIN GRASPING
FORCEPS**

WORKING LENGTH	COMPLETE
220 mm	EAI-0459R
330 mm	EAI-0460R
450 mm	EAI-0461R




**LONG ATRAUMATIC
GRASPING FORCEPS**

WORKING LENGTH	COMPLETE
220 mm	EAI-0495R
330 mm	EAI-0496R
450 mm	EAI-0497R




**SHORT DELPHIN GRASPING
FORCEPS**

WORKING LENGTH	COMPLETE
220 mm	EAI-0462R
330 mm	EAI-0463R
450 mm	EAI-0464R




**LONG ATRAUMATIC CONCAVE
GRASPING FORCEPS**

WORKING LENGTH	COMPLETE
220 mm	EAI-0498R
330 mm	EAI-0499R
450 mm	EAI-0500R




**LONG GRASPING FORCEPS
WITH SPOON**

WORKING LENGTH	COMPLETE
220 mm	EAI-0465R
330 mm	EAI-0466R
450 mm	EAI-0467R




**FENESTRATED GRASPING
FORCEPS**

WORKING LENGTH	COMPLETE
220 mm	EAI-0501R
330 mm	EAI-0502R
450 mm	EAI-0503R



**DISS.- AND GRASPING
FORCEPS, HEAVY PATTERN**


WORKING LENGTH	COMPLETE
220 mm	EAI-0468R
330 mm	EAI-0469R
450 mm	EAI-0470R



**MICRO-FENESTRATED
GRASPING FORCEPS**


WORKING LENGTH	COMPLETE
220 mm	EAI-0504R
330 mm	EAI-0505R
450 mm	EAI-0506R

5mm GRASPING AND DISSECTING FORCEPS (LAPAROSCOPY)




MICRO GRASPING FORCEPS

WORKING LENGTH	COMPLETE
220 mm	EAI-0507R
330 mm	EAI-0508R
450 mm	EAI-0509R




ALLIGATOR GRASPING FORCEPS

WORKING LENGTH	COMPLETE
220 mm	EAI-0540R
330 mm	EAI-0541R
450 mm	EAI-0542R




MICRO-MARYLAND GRASPING FORCEPS

WORKING LENGTH	COMPLETE
220 mm	EAI-0528R
330 mm	EAI-0529R
450 mm	EAI-0530R




MIXER GRASPING FORCEPS 60°

WORKING LENGTH	COMPLETE
220 mm	EAI-0543R
330 mm	EAI-0544R
450 mm	EAI-0545R




LONG ATRAUMATIC GRASPING FORCEPS

WORKING LENGTH	COMPLETE
220 mm	EAI-0531R
330 mm	EAI-0532R
450 mm	EAI-0533R




MIXER GRASPING FORCEPS 90°

WORKING LENGTH	COMPLETE
220 mm	EAI-0564R
330 mm	EAI-0565R
450 mm	EAI-0566R




GRASPING FORCEPS, STRAIGHT

WORKING LENGTH	COMPLETE
220 mm	EAI-0534
330 mm	EAI-0535
450 mm	EAI-0536




ALLIGATOR GRASPING FORCEPS, WITH SPOON

WORKING LENGTH	COMPLETE
220 mm	EAI-0567R
330 mm	EAI-0568R
450 mm	EAI-0569R



GRASPING FORCEPS, CURVED


WORKING LENGTH	COMPLETE
220 mm	EAI-0537R
330 mm	EAI-0538R
450 mm	EAI-0539R



DA VINCI AGGRESSIVE GRASPING FORCEPS


WORKING LENGTH	COMPLETE
220 mm	EAI-0570R
330 mm	EAI-0571R
450 mm	EAI-0572R

5mm GRASPING AND DISSECTING FORCEPS (LAPAROSCOPY)




MAXI ATRAUMATIC GRASPING FORCEPS

WORKING LENGTH	COMPLETE
220 mm	EAI-0573R
330 mm	EAI-0574R
450 mm	EAI-0575R




BIOPSY SPOON FORCEPS, SINGLE ACTION

WORKING LENGTH	COMPLETE
220 mm	EAI-0606R
330 mm	EAI-0607R
450 mm	EAI-0608R




OVIDUCT GRASPING FORCEPS

WORKING LENGTH	COMPLETE
220 mm	EAI-0576R
330 mm	EAI-0577R
450 mm	EAI-0578R




BIOPSY SPOON FORCEPS WITH 2 TEETH, SINGLE ACTION

WORKING LENGTH	COMPLETE
220 mm	EAI-0609R
330 mm	EAI-0610R
450 mm	EAI-0611R




DISSECTING AND GRASPING FORCEPS, FINE PATTERN

WORKING LENGTH	COMPLETE
220 mm	EAI-0579R
330 mm	EAI-0580R
450 mm	EAI-0581R




BLAKESLEY BIOPSY FORCEPS

WORKING LENGTH	COMPLETE
220 mm	EAI-0612R
330 mm	EAI-0613R
450 mm	EAI-0614R




BIOPSY SPOON FORCEPS, DOUBLE ACTION

WORKING LENGTH	COMPLETE
220 mm	EAI-0600R
330 mm	EAI-0601R
450 mm	EAI-0602R




BIOPSY PUNCH

WORKING LENGTH	COMPLETE
220 mm	EAI-0615R
330 mm	EAI-0616R
450 mm	EAI-0617R



BIOPSY SPOON FORCEPS WITH 2 TEETH, DOUBLE ACTION


WORKING LENGTH	COMPLETE
220 mm	EAI-0603R
330 mm	EAI-0604R
450 mm	EAI-0605R



MULTIPLE TEETH GRASPING FORCEPS


WORKING LENGTH	COMPLETE
220 mm	EAI-0627R
330 mm	EAI-0628R
450 mm	EAI-0629R

5mm GRASPING AND DISSECTING FORCEPS (LAPAROSCOPY)




RAT-TOOTH GRASPING FORCEPS

WORKING LENGTH	COMPLETE
220 mm	EAI-0630R
330 mm	EAI-0631R
450 mm	EAI-0632R




BOWEL GRASPING FORCEPS, FENESTRATED

WORKING LENGTH	COMPLETE
220 mm	EAI-0657R
330 mm	EAI-0658R
450 mm	EAI-0659R




CLAW GRASPING FORCEPS

WORKING LENGTH	COMPLETE
220 mm	EAI-0633R
330 mm	EAI-0634R
450 mm	EAI-0635R



ALLIS GRASPING FORCEPS

WORKING LENGTH	COMPLETE
220 mm	EAI-0672R
330 mm	EAI-0673R
450 mm	EAI-0674R




MARYLAND GRASPING FORCEPS

WORKING LENGTH	COMPLETE
220 mm	EAI-0471R
330 mm	EAI-0472R
450 mm	EAI-0473R




BAKEY ATRAUMATIC GRASPING FORCEPS

WORKING LENGTH	COMPLETE
220 mm	EAI-0678R
330 mm	EAI-0679R
450 mm	EAI-0680R




CROCODILE GRASPING FORCEPS

WORKING LENGTH	COMPLETE
220 mm	EAI-0645R
330 mm	EAI-0646R
450 mm	EAI-0647R



BABCOCK GRASPING FORCEPS


WORKING LENGTH	COMPLETE
220 mm	EAI-0684R
330 mm	EAI-0685R
450 mm	EAI-0686R



RETRACTION GRASPING FORCEPS


WORKING LENGTH	COMPLETE
220 mm	EAI-0651R
330 mm	EAI-0652R
450 mm	EAI-0653R

10mm GRASPING AND DISSECTING FORCEPS (LAPAROSCOPY)




MARYLAND GRASPING FORCEPS

WORKING LENGTH	COMPLETE
220 mm	EAI-0702R
330 mm	EAI-0703R
450 mm	EAI-0704R




ALLIS ATRAUMATIC GRASPING FORCEPS

WORKING LENGTH	COMPLETE
220 mm	EAI-0729R
330 mm	EAI-0730R
450 mm	EAI-0731R




KELLY GRASPING FORCEPS

WORKING LENGTH	COMPLETE
220 mm	EAI-0705R
330 mm	EAI-0706R
450 mm	EAI-0707R




BABCOCK GRASPING FORCEPS

WORKING LENGTH	COMPLETE
220 mm	EAI-0735R
330 mm	EAI-0736R
450 mm	EAI-0737R




ALLIGATOR GRASPING FORCEPS

WORKING LENGTH	COMPLETE
220 mm	EAI-0708R
330 mm	EAI-0709R
450 mm	EAI-0710R




SPOON GRASPING FORCEPS

WORKING LENGTH	COMPLETE
220 mm	EAI-0750R
330 mm	EAI-0751R
450 mm	EAI-0752R




MIXER GRASPING FORCEPS 90°

WORKING LENGTH	COMPLETE
220 mm	EAI-0711R
330 mm	EAI-0712R
450 mm	EAI-0713R




CLAW GRASPING FORCEPS, SINGLE ACTION

WORKING LENGTH	COMPLETE
220 mm	EAI-0753R
330 mm	EAI-0754R
450 mm	EAI-0755R



ALLIS GRASPING FORCEPS


WORKING LENGTH	COMPLETE
220 mm	EAI-0723R
330 mm	EAI-0724R
450 mm	EAI-0725R



CLAW GRASPING FORCEPS, DOUBLE ACTION


WORKING LENGTH	COMPLETE
220 mm	EAI-0756R
330 mm	EAI-0757R
450 mm	EAI-0758R

5mm SCISSORS (LAPAROSCOPY)




**PERITONEAL SCISSORS,
SERRATED**

WORKING LENGTH	COMPLETE
220 mm	EAI-0777R
330 mm	EAI-0778R
450 mm	EAI-0779R




**METZENBAUM SCISSORS,
STRAIGHT, HEAVY BLADES
PATTER**

WORKING LENGTH	COMPLETE
220 mm	EAI-0807R
330 mm	EAI-0808R
450 mm	EAI-0809R




HOOK SCISSORS

WORKING LENGTH	COMPLETE
220 mm	EAI-0780R
330 mm	EAI-0781R
450 mm	EAI-0782R




**METZENBAUM SCISSORS,
CURVED, HEAVY BLADES
PATTERN**

WORKING LENGTH	COMPLETE
220 mm	EAI-0810R
330 mm	EAI-0811R
450 mm	EAI-0812R




**BABY-METZENBAUM SCIS-
SORS, CURVED**

WORKING LENGTH	COMPLETE
220 mm	EAI-0783R
330 mm	EAI-0784R
450 mm	EAI-0785R




Micro Scissors, curved

WORKING LENGTH	COMPLETE
220 mm	EAI-0816R
330 mm	EAI-0817R
450 mm	EAI-0818R



**METZENBAUM SCISSORS,
CURVED**

WORKING LENGTH	COMPLETE
220 mm	EAI-0786R
330 mm	EAI-0787R
450 mm	EAI-0788R




ENDO SCISSORS, CURVED

WORKING LENGTH	COMPLETE
220 mm	EAI-0792R
330 mm	EAI-0793R
450 mm	EAI-0794R


NEEDLE HANDLE (LAPAROSCOPY)



**METAL HANDLE
WITH SPRING
AND
CLEANING PORT,
NOT ROTATING**




**AXIAL HANDLE WITH RATCHET
AND CLEANING PORT, NOT ROTATING**




TUNGSTEN CARBIDE, 3 MM

WORKING LENGTH	METAL HANDLE	AXIAL HANDLE
220 mm	EAI-0819R	EAI-0840R
330 mm	EAI-0820R	EAI-0841R
450 mm	EAI-0821R	EAI-0842R




"FLAMINGO", 5 MM

WORKING LENGTH	METAL HANDLE	AXIAL HANDLE
220 mm	EAI-0831R	EAI-0852R
330 mm	EAI-0832R	EAI-0853R
450 mm	EAI-0833R	EAI-0854R




SERRATED JAW, 3 MM

WORKING LENGTH	METAL HANDLE	AXIAL HANDLE
220 mm	EAI-0822R	EAI-0843R
330 mm	EAI-0823R	EAI-0844R
450 mm	EAI-0824R	EAI-0845R




"PARROT", 5 MM

WORKING LENGTH	METAL HANDLE	AXIAL HANDLE
220 mm	EAI-0834R	EAI-0855R
330 mm	EAI-0835R	EAI-0856R
450 mm	EAI-0836R	EAI-0857R




TUNGSTEN CARBIDE, 5 MM

WORKING LENGTH	METAL HANDLE	AXIAL HANDLE
220 mm	EAI-0825R	EAI-0846R
330 mm	EAI-0826R	EAI-0847R
450 mm	EAI-0827R	EAI-0848R




**TUNGSTEN CARBIDE,
CURVED LEFT, 5 MM**

WORKING LENGTH	METAL HANDLE	AXIAL HANDLE
220 mm	EAI-0837R	EAI-0858R
330 mm	EAI-0838R	EAI-0859R
450 mm	EAI-0839R	EAI-0860R



SERRATED JAW, 5 MM


WORKING LENGTH	METAL HANDLE	AXIAL HANDLE
220 mm	EAI-0828R	EAI-0849R
330 mm	EAI-0829R	EAI-0850R
450 mm	EAI-0830R	EAI-0851R



**TUNGSTEN CARBIDE,
CURVED RIGHT, 5MM**

WORKING LENGTH	METAL HANDLE	AXIAL HANDLE
220 mm	EAI-1775R	EAI-1778R
330 mm	EAI-1776R	EAI-1779R
450 mm	EAI-1777R	EAI-1780R


EASYKNOT SUTURING FORCEPS (LAPAROSCOPY)



EASY KNOT

EAI-0861R

For easy intracorporeal
suturing and knot making




1. Take the stitch with standard needle holder leaving a 1" - 2" short tail on the needle entry side. Grasp the needle end of the suture with an assisting forceps about 3" - 5" from the operative site to fix the suture. Trap the suture under the upper jaw of the Easyknot device about 1.5" from the point of suture fixation.

2. Use the assisting forceps to wrap the suture around the Easyknot device twice, distal to the upper jaw (preferably going first under the device than over it).

3. Bring both instruments closer to the short tail and hold the tail with the lower jaw of the Easyknot device. Open the upper jaw to release trapped suture and pull the device (holding the suture tail) and assisting forceps (holding the needle suture end) in opposite directions (relative to initial positions) to establish a flat half hitch. Cinch and tighten as needed.

4. Repeat the procedure to form additional half hitches making sure to reverse the direction of pulling on the suture ends to establish a square knot. For instance, if the short end starts on the right side of the incision, it should end at the left side after the first half hitch and again on the right side after the second half hitch.



ENDOCH



— MEDICAL —

1- 12 Bradwick Drive,
Concord, Ontario,
Canada L4K 3P6

Tel: (+1) 905-669-6999
Fax: (+1) 905-669-1448
Toll Free: (+1) 888-4-ENDOCH
Email: info@endoch.com
Web: www.endoch.com