M-FINGERS by Partial Hand SOLUTIONS, LLC



M-Fingers are simple to operate and controlled

by the user's own motion, allowing the patient to perform tasks more effectively, more responsively, and with optimized precision. They conform around objects to provide a secure and confident grip for the user. M-Fingers were designed with high strength, injected molded composite materials to deliver the durability and







X-Large M-Finger

M-FINGER SIZING

Total 2.4 in Medium M-Finger


1.375" (34.9mm)


.55












(41.2mm)





Partial Hand M-FINGER PARTS SOLUTIONS, LLC



M-FINGER PART NUMBER



* (TAN) Tan, (BRN) Brown, or (BLK) Black. Example part number: ULTRA-SM-MF-BLK



PARTIAL M-FINGER PART NUMBER

COLOR*

* (TAN) Tan, (BRN) Brown, or (BLK) Black. Example part number: SM-PART-BLK







THUMB PART NUMBER

COLOR*

TAN, BRN, BLK

*(TAN) Tan, (BRN) Brown, or (BLK) Black. Example part number: THUMB-BLK

2



Partial Hand M-FINGER PARTS SOLUTIONS, LLC

MOUNTING PLATES





 $4\,\mathsf{Finger}\,\mathsf{M}\text{-}\mathsf{Finger}\,\mathsf{System}\,\mathsf{with}\,\mathsf{4}\text{-}\mathsf{Finger}\,\mathsf{Mounting}\,\mathsf{Plate},\mathsf{DCM},\mathsf{and}\,\mathsf{BOA}\,\mathsf{mechanism}$

4 FINGER SYSTEM PART NUMBER

4FS - KIT .

COLOR*

TAN, BRN, BLK

SM, MED, LG

* (TAN) Tan, (BRN) Brown, or (BLK) Black. Example part number: 4FS-KIT-BLK-SM



4 Finger M-Finger System with 4-Finger Mounting Plate, Thumb, DCM, and BOA mechanism

5 FINGER SYSTEM PART NUMBER

5FS - KIT -

TAN, BRN, BLK

SM, MED, LG

* (TAN) Tan, (BRN) Brown, or (BLK) Black. Example part number: 5FS-KIT-BLK-SM

HINGES/CONNECTORS

| | PRODUCT NUMBER | DESCRIPTION |
|--------|----------------|--|
| 11 | GAFFNEY | Gaffney Joint/Hinge (Colt) |
| | ВОА | BOA Fastener for M-Fingers |
| 0. 0 0 | DCM1 | 1-Finger Direct Connect Mechanism (DCM) Black |
| | DCM 2 | 2-Finger Direct Connect Mechanism (DCM) Black |
| | DCM 3 | 3-Finger Direct Connect Mechanism (DCM) Black |
| | DCM 4 | 4-Finger Direct Connect Mechanism (DCM) Black |

GRIP TIPS

| PRODUCT NUMBER | DESCRIPTION |
|----------------|---------------------------------------|
| GRIPTIP | Silicone protective covering , 5 pack |

REPLACEMENT PARTS


| REPLACEMEN | II PARTS |
|----------------|-------------------------------|
| PRODUCT NUMBER | DESCRIPTION |
| PMF-BASE-BLK | Partial M-Finger base (Black) |
| PMF-BASE-BRN | Partial M-Finger base (Brown) |
| PMF-BASE-TAN | Partial M-Finger base (Tan) |
| TB-TIP-BLK | Thumb Tip (Black) |
| TB-TIP-BRN | Thumb Tip (Brown) |
| TB-TIP-TAN | Thumb Tip (Tan) |
| TB-BASE-BLK | Base of thumb (Black) |
| TB-BASE-BRN | Base of thumb (Brown) |
| TB-BASE-TAN | Base of thumb (Tan) |
| SM-MJ-BLK | Small Middle joint (Black) |
| SM-MJ-BRN | Small Middle Joint (Brown) |
| SM-MJ-TAN | Small Middle Joint (Tan) |
| MED-MJ-BLK | Medium Middle joint (Black) |
| MED-MJ-BRN | Medium Middle Joint (Brown) |
| MED-MJ-TAN | Medium Middle Joint (Tan) |
| LG-MJ-BLK | Large Middle joint (Black) |
| LG-MJ-BRN | Large Middle Joint (Brown) |
| LG-MJ-TAN | Large Middle Joint (Tan) |
| FB-BLK | Finger Base Joint (Black) |
| FB-BRN | Finger Base Joint (Brown) |
| FB-TAN | Finger Base Joint (Tan) |
| SM-FT-BLK | Small Finger Tip (Black) |
| SM-FT-BRN | Small Finger Tip (Brown) |
| SM-FT-TAN | Small Finger Tip (Tan) |
| MED-FT-BLK | Medium Finger Tip (Black) |
| MED-FT-BRN | Medium Finger Tip (Brown) |
| MED-FT-TAN | Medium Finger Tip (Tan) |
| LG-FT-BLK | Large Finger Tip (Black) |
| LG-FT-BRN | Large Finger Tip (Brown) |
| LG-FT-TAN | Large Finger Tip (Tan) |
| | |

4


TITAN SERIES by Partial Hand SOLUTIONS, LLC



The newest innovation from Partial Hand Solutions, the titanium Titan series fingers are designed for the most active heavy duty user, ideal for those with occupational requirements. They are suitable for someone with a full or partial finger amputation and are available in a variety of sizes. The ratchet design allows the patient to manually position the joints in the most functional position, while providing a secure grasp. Various bases allow for the device to be fit as a thumb with full rotation and adjustable friction.



TITAN SIZING



Partial Hand TITAN FINGER PARTS



TITAN FINGER PART NUMBER

†(1) Small: length size 3.35", (2) Medium: length size 3.52", or (3) Large: length size 3.62". Example part number: FULL-TITAN3



TITAN PARTIAL FINGER PART NUMBER

PART - TITAN

†(1) Small: length size 2.25", (2) Medium: length size 2.5", or (3) Large: length size 2.62".



TITAN THUMB PART NUMBER

THUMB - TITAN

Example part number: THUMB-TITAN3



M-FINGER ORDER FORM

| CUSTOMER # | P.O.# | | DATE |
|---------------------------|---------------------------------------|----------------------|-----------------------|
| PROSTHETIST INFORM | ATION | | |
| BILLING | | SHIPPING (LEAVE BLAN | K IF SAME AS BILLING) |
| FACILITY/ATTN: | | FACILITY/ATTN: | |
| ADDRESS | | ADDRESS | |
| CITY | STATE/ZIP | CITY | STATE/ZIP |
| COUNTRY | | COUNTRY | |
| PHONE | FAX | PHONE | FAX |
| CARRIER* UPS OTH | HER | DATE REQUIRED | TIME |
| *ORDERS SHIP UPS GROUND C | N THE FOLLOWING DAY UNLESS SPECIFIED. | EMAIL | |
| PATIENT INFORMATION | ١ | | |
| PROSTHETIST NAME | | REQUISITIONER | |
| PATIENT ID / NAME | | | |
| NOTES | | | |

For sizing guides, visit college-park.com/fingers





| Ultra-Small | ULTRA-SM |
|-------------|----------|
| Small | SM |
| Medium | MED |
| Large | LG |
| Extra Large | XL |

Required Accessories:

- ☐ Gaffney Joint ☐ BOA Fastener
- ☐ Mounting Plate (see below)

MOUNTING PLATES



PARTIAL M-FINGER ALL PARTS INCLUDED





EMAIL COMPLETED ORDER FORM TO: orders@college-park.com **US TOLL FREE** 800.728.7950 **PHONE** 586.294.7950 **FAX** 586.294.0067 **ON THE WEB** www.college-park.com (college park

M-FINGER ORDER FORM

| CUSTOMER # | P.O. # | | DATE |
|-------------------------|--|----------------------|-----------------------|
| PROSTHETIST INFORM | MATION | | |
| BILLING | | SHIPPING (LEAVE BLAN | K IF SAME AS BILLING) |
| FACILITY/ATTN: | | FACILITY/ATTN: | |
| ADDRESS | | ADDRESS | |
| CITY | STATE/ZIP | CITY | STATE/ZIP |
| COUNTRY | | COUNTRY | |
| PHONE | FAX | PHONE | FAX |
| CARRIER* □ UPS □ OT | HER | DATE REQUIRED | TIME |
| *ORDERS SHIP UPS GROUND | ON THE FOLLOWING DAY UNLESS SPECIFIED. | EMAIL | |
| PATIENT INFORMATIO | N | | |
| PROSTHETIST NAME | | REQUISITIONER | |
| PATIENT ID / NAME | | | |
| NOTES | | | |

For sizing guides, visit college-park.com/fingers





BRN

BLK

ADDITIONAL PARTS

THUMB

| DESCRIPTION | PART ID |
|--|---------|
| GRIP TIPS SILICONE PROTECTIVE COVERING, 5 PACK | GRIPTIP |
| GAFFNEY JOINT/HINGE (COLT) | GAFFNEY |
| 1-FINGER DIRECT CONNECT MECHANISM | DCM1 |
| 2-FINGER DIRECT CONNECT MECHANISM | DCM 2 |
| 3-FINGER DIRECT CONNECT MECHANISM | DCM 3 |
| 4-FINGER DIRECT CONNECT MECHANISM | DCM 4 |
| BOA FASTENER FOR M-FINGER | ВОА |

Brown

Black




EMAIL COMPLETED ORDER FORM TO: orders@college-park.com **US TOLL FREE** 800.728.7950 **PHONE** 586.294.7950 FAX 586.294.0067 ON THE WEB www.college-park.com



M-FINGER ORDER FORM

PAGE 3

| CUSTOMER # | P.O.# | | DATE |
|---------------------------|---------------------------------------|----------------------|--|
| PROSTHETIST INFORM. | ATION | | |
| BILLING | | SHIPPING (LEAVE BLAN | IK IF SAME AS BILLING) |
| FACILITY/ATTN: | | FACILITY/ATTN: | |
| ADDRESS | | ADDRESS | |
| CITY | STATE/ZIP | CITY | STATE/ZIP |
| COUNTRY | | COUNTRY | |
| PHONE | FAX | PHONE | FAX |
| CARRIER* □ UPS □ OTH | HER | DATE REQUIRED | TIME |
| *ORDERS SHIP UPS GROUND C | N THE FOLLOWING DAY UNLESS SPECIFIED. | EMAIL | |
| PATIENT INFORMATION | N | | |
| PROSTHETIST NAME | | REQUISITIONER | |
| PATIENT ID / NAME | | | |
| NOTES | | | |
| | | For sizing | guides, visit college-park.com/fingers |



REPLACEMENT PARTS

| PART NUMBER | QUANTITY |
|-------------|----------|
| | |
| | |
| | |
| | |
| | |
| | |

| PART NUMBER | QUANTITY |
|-------------|----------|
| | |
| | |
| | |
| | |
| | |
| | |
| | |

Small

Medium

Large

TITAN1

TITAN2

TITAN3



EMAIL COMPLETED ORDER FORM TO: orders@college-park.com **US TOLL FREE** 800.728.7950 **PHONE** 586.294.7950 FAX 586.294.0067 ON THE WEB www.college-park.com