



*7106 Garden Park Rd.  
Cañon City, CO 81212*

*(719) 371-7707*

*eliteequine.net*



The NeedleView arthroscope helps give our clients the information they need to make good decisions about how to proceed with their horse. At 1.2mm, the NeedleView scope is only the size of an 18-gauge needle. We can scope with just a local block and sedation in the standing horse. It helps us realize at the time of the procedure what paths to consider. For

example, if the animal has been treated with biologic agents or regenerative therapies and doesn't seem to be responding, NeedleView lets us easily take a look in the effected joint to see what's going on. NeedleView can also help you give the owner the information and options – for example, to proceed with surgery, debride to speed healing, or – if the injury is too severe, particularly with performance horses – to choose the right next step.