

SPINE SYSTEM





POINTER

OF

CONTENT

Sequence No.	Implants / Instrumentations	Page No.
1	MONOAXIAL SCREW	01
2	POLYAXIAL SCREW	02
3	MONOAXIAL REDUCTION SCREW	03
4	POLYAXIAL REDUCTION SCREW	04
5	SPINAL RODS	05
6	AXIAL CONNECTOR	05
7	PARALLEL CONNECTOR	05
8	TRANSVERSE CONNECTOR	05
9	TRANSVERSE CONNECTOR PLATE	06
10	MULTIAXIAL CONNECTOR	06
11	STAPLE	06
12	PEDICLE HOOK	06
13	WIDE BLADE HOOK	07
14	NARROW BLADE HOOK	07
15	SUPRA LAMINAR HOOK	07
16	ANGLED BLADE HOOK - LEFT & RIGHT	07
17	ANTERIOR CERVICAL PLATE	08
18	Ø 4.3MM REVISION BONE SCREW SELF TAPPING FOR CERVICAL PLATE	08
19	Ø 4.0 MM BONE SCREW SELF TAPPING FOR CERVICAL PLATE	08
20	CERVICAL DEO CAGE	09
21	CERVICAL FUSION CAGE	09
22	PLIF CAGE (BULLET CAGE)	10
23	TLIF CAGE (BANANA CAGE) TYPE - 2	11
24	TLIF CAGE (KIDNEY CAGE) TYPE - 1	12
25	MESH CAGE	13
26	EXPANDABLE CAGE	14
27	SPINE INSTRUMENTATIONS	16
28	CERVICAL INSTRUMENTATIONS	18
29	PLIF / TLIF INSTRUMENTATIONS	20



Strengthening Your Grip



MONOAXIAL SCREW

REFERENCE NO. (Titanium)	SIZE: Ø 4.5
2020-TT-4525	25 MM
2020-TT-4530	30 MM
2020-TT-4535	35 MM
2020-TT-4540	40 MM
2020-TT-4545	45 MM
2020-TT-4550	50 MM



REFERENCE NO. (Titanium)	SIZE: Ø 5.5
2020-TT-5525	25 MM
2020-TT-5530	30 MM
2020-TT-5535	35 MM
2020-TT-5540	40 MM
2020-TT-5545	45 MM
2020-TT-5550	50 MM
2020-TT-5555	55 MM

REFERENCE NO. (Titanium)	SIZE: Ø 6.5
2020-TT-6525	25 MM
2020-TT-6530	30 MM
2020-TT-6535	35 MM
2020-TT-6540	40 MM
2020-TT-6545	45 MM
2020-TT-6550	50 MM
2020-TT-6555	55 MM



REFERENCE NO. (Titanium)	SIZE: Ø 7.0
2020-TT-7025	25 MM
2020-TT-7030	30 MM
2020-TT-7035	35 MM
2020-TT-7040	40 MM
2020-TT-7045	45 MM
2020-TT-7050	50 MM
2020-TT-7055	55 MM



POLYAXIAL SCREW



REFERENCE NO. (Titanium)	SIZE: Ø 4.5
2021-TT-4525	25 MM
2021-TT-4530	30 MM
2021-TT-4535	35 MM
2021-TT-4540	40 MM
2021-TT-4545	45 MM
2021-TT-4550	50 MM

REFERENCE NO. (Titanium)	SIZE: Ø 5.5
2021-TT-5525	25 MM
2021-TT-5530	30 MM
2021-TT-5535	35 MM
2021-TT-5540	40 MM
2021-TT-5545	45 MM
2021-TT-5550	50 MM
2021-TT-5555	55 MM

REFERENCE NO. (Titanium)	SIZE: Ø 6.5
2021-TT-6525	25 MM
2021-TT-6530	30 MM
2021-TT-6535	35 MM
2021-TT-6540	40 MM
2021-TT-6545	45 MM
2021-TT-6550	50 MM
2021-TT-6555	55 MM



REFERENCE NO. (Titanium)	SIZE: Ø 7.0
2021-TT-7025	25 MM
2021-TT-7030	30 MM
2021-TT-7035	35 MM
2021-TT-7040	40 MM
2021-TT-7045	45 MM
2021-TT-7050	50 MM
2021-TT-7055	55 MM



MONOAXIAL REDUCTION SCREW

REFERENCE NO. (Titanium)	SIZE: Ø 4.5
2039-TT-4525	25 MM
2039-TT-4530	30 MM
2039-TT-4535	35 MM
2039-TT-4540	40 MM
2039-TT-4545	45 MM
2039-TT-4550	50 MM



REFERENCE NO. (Titanium)	SIZE: Ø 5.5
2039-TT-5525	25 MM
2039-TT-5530	30 MM
2039-TT-5535	35 MM
2039-TT-5540	40 MM
2039-TT-5545	45 MM
2039-TT-5550	50 MM
2039-TT-5555	55 MM

REFERENCE NO. (Titanium)	SIZE: Ø 6.5
2039-TT-6525	25 MM
2039-TT-6530	30 MM
2039-TT-6535	35 MM
2039-TT-6540	40 MM
2039-TT-6545	45 MM
2039-TT-6550	50 MM
2039-TT-6555	55 MM



REFERENCE NO. (Titanium)	SIZE: Ø 7.0
2039-TT-7025	25 MM
2039-TT-7030	30 MM
2039-TT-7035	35 MM
2039-TT-7040	40 MM
2039-TT-7045	45 MM
2039-TT-7050	50 MM
2039-TT-7055	55 MM



POLYAXIAL REDUCTION SCREW



REFERENCE NO. (Titanium)	SIZE: Ø 4.5
2022-TT-4525	25 MM
2022-TT-4530	30 MM
2022-TT-4535	35 MM
2022-TT-4540	40 MM
2022-TT-4545	45 MM
2022-TT-4550	50 MM

REFERENCE NO. (Titanium)	SIZE: Ø 5.5
2022-TT-5525	25 MM
2022-TT-5530	30 MM
2022-TT-5535	35 MM
2022-TT-5540	40 MM
2022-TT-5545	45 MM
2022-TT-5550	50 MM
2022-TT-5555	55 MM

REFERENCE NO. (Titanium)	SIZE: Ø 6.5
2022-TT-6525	25 MM
2022-TT-6530	30 MM
2022-TT-6535	35 MM
2022-TT-6540	40 MM
2022-TT-6545	45 MM
2022-TT-6550	50 MM
2022-TT-6555	55 MM



REFERENCE NO. (Titanium)	SIZE: Ø 7.0
2022-TT-7025	25 MM
2022-TT-7030	30 MM
2022-TT-7035	35 MM
2022-TT-7040	40 MM
2022-TT-7045	45 MM
2022-TT-7050	50 MM
2022-TT-7055	55 MM



SPINAL RODS

REFERENCE NO. (Titanium)	SIZE: Ø 5.5
2028-TT-55050	50 MM
2028-TT-55060	60 MM
2028-TT-55070	70 MM
2028-TT-55080	80 MM
2028-TT-55090	90 MM
2028-TT-55120	120 MM
2028-TT-55140	140 MM
2028-TT-55160	160 MM
2028-TT-55180	180 MM
2028-TT-55200	200 MM
2028-TT-55250	250 MM
2028-TT-55300	300 MM
2028-TT-55350	350 MM
2028-TT-55500	500MM



AXIAL CONNECTOR

REFERENCE NO. (Titanium)
2033-TT-0001



PARALLEL CONNECTOR

REFERENCE NO. (Titanium)
2033-TT-0002



TRANSVERSE CONNECTOR

REFERENCE NO. (Titanium)
2035-TT-5500






TRANSVERSE CONNECTOR PLATE




REFERENCE NO. (Titanium)	SIZE: Ø 5.5
2035-TT-55030	30 MM
2035-TT-55040	40 MM
2035-TT-55050	50 MM
2035-TT-55060	60 MM
2035-TT-55070	70 MM
2035-TT-55080	80 MM
2035-TT-55090	90 MM
2035-TT-55100	100 MM

MULTIAXIAL CONNECTOR



REFERENCE NO. (Titanium)	SIZE
2036-TT-3035	30 MM X 35 MM
2036-TT-3540	35 MM X 40 MM
2036-TT-4050	40 MM X 50 MM
2036-TT-5070	50 MM X 70 MM

STAPLE



REFERENCE NO. (Titanium)
2040-TT-0001

PEDICLE HOOK



REFERENCE NO. (Titanium)
2040-TT-0002



WIDE BLADE HOOK

REFERENCE NO. (Titanium)

2040-TT-0003



NARROW BLADE HOOK

REFERENCE NO. (Titanium)

2040-TT-0004



SUPRA LAMINAR HOOK

REFERENCE NO. (Titanium)

2040-TT-0005



ANGLED BLADE HOOK

REFERENCE NO. (Titanium) - LEFT

2040-TT-0006

REFERENCE NO. (Titanium) - RIGHT

2040-TT-0007





ANTERIOR CERVICAL PLATE



REFERENCE NO. (Titanium)	SIZE
2024-TT-0021	21 MM
2024-TT-0023	23 MM
2024-TT-0025	25 MM
2024-TT-0275	27.5 MM
2024-TT-0030	30 MM
2024-TT-0325	32.5 MM
2024-TT-0035	35 MM
2024-TT-0375	37.5 MM
2024-TT-0040	40 MM
2024-TT-0425	42.5 MM
2024-TT-0045	45 MM
2024-TT-0475	47.5 MM
2024-TT-0050	50 MM
2024-TT-0525	52.5 MM
2024-TT-0055	55 MM
2024-TT-0575	57.5 MM
2024-TT-0060	60 MM
2024-TT-0625	62.5 MM
2024-TT-0065	65 MM
2024-TT-0675	67.5 MM
2024-TT-0070	70 MM
2024-TT-0725	72.5 MM
2024-TT-0075	75 MM
2024-TT-0775	77.5 MM
2024-TT-0080	80

Ø 4.3MM REVISION BONE SCREW SELF TAPPING FOR CERVICAL PLATE



REFERENCE NO. (Titanium)	SIZE
2025-TT-4310	10 MM
2025-TT-4312	12 MM
2025-TT-4314	14 MM
2025-TT-4316	16 MM
2025-TT-4318	18 MM
2025-TT-4320	20 MM

Ø 4.0 MM BONE SCREW SELF TAPPING FOR CERVICAL PLATE



REFERENCE NO. (Titanium)	SIZE
2025-TT-4010	10 MM
2025-TT-4012	12 MM
2025-TT-4014	14 MM
2025-TT-4016	16 MM
2025-TT-4018	18 MM
2025-TT-4020	20 MM



CERVICAL DEO CAGE

REFERENCE NO. (Titanium)	SIZE: 12 MM
2042-TT-1205	5 MM
2042-TT-1206	6 MM
2042-TT-1207	7 MM
2042-TT-1208	8 MM

REFERENCE NO. (Titanium)	SIZE: 14 MM
2042-TT-1405	5 MM
2042-TT-1406	6 MM
2042-TT-1407	7 MM
2042-TT-1408	8 MM



REFERENCE NO. (Titanium)	SIZE: 16 MM
2042-TT-1605	5 MM
2042-TT-1606	6 MM
2042-TT-1607	7 MM
2042-TT-1608	8 MM

CERVICAL FUSION CAGE

REFERENCE NO. (Titanium)	SIZE: 14 MM
2043-TT-1406	6 MM
2043-TT-1407	7 MM
2043-TT-1408	8 MM
2043-TT-1409	9 MM

REFERENCE NO. (Titanium)	SIZE: 16 MM
2043-TT-1606	6 MM
2043-TT-1607	7 MM
2043-TT-1608	8 MM
2043-TT-1609	9 MM





PLIF CAGE (BULLET CAGE)



REFERENCE NO. (Titanium)	SIZE: 0°
2032-TT-00820	8 x 20 MM
2032-TT-00823	8 x 23 MM
2032-TT-00825	8 x 25 MM
2032-TT-00920	9 x 20 MM
2032-TT-00923	9 x 23 MM
2032-TT-00925	9 x 25 MM
2032-TT-01020	10 x 20 MM
2032-TT-01023	10 x 23 MM
2032-TT-01025	10 x 25 MM
2032-TT-01120	11 x 20 MM
2032-TT-01123	11 x 23 MM
2032-TT-01125	11 x 25 MM
2032-TT-01220	12 x 20 MM
2032-TT-01223	12 x 23 MM
2032-TT-01225	12 x 25 MM
2032-TT-01320	13 x 20 MM
2032-TT-01323	13 x 23 MM
2032-TT-01325	13 x 25 MM
2032-TT-01420	14 x 20 MM
2032-TT-01423	14 x 23 MM
2032-TT-01425	14 x 25 MM

REFERENCE NO. (Titanium)	SIZE: 4°
2032-TT-40820	8 x 20 MM
2032-TT-40823	8 x 23 MM
2032-TT-40825	8 x 25 MM
2032-TT-40920	9 x 20 MM
2032-TT-40923	9 x 23 MM
2032-TT-40925	9 x 25 MM
2032-TT-41020	10 x 20 MM
2032-TT-41023	10 x 23 MM
2032-TT-41025	10 x 25 MM
2032-TT-41120	11 x 20 MM
2032-TT-41123	11 x 23 MM
2032-TT-41125	11 x 25 MM
2032-TT-41220	12 x 20 MM
2032-TT-41223	12 x 23 MM
2032-TT-41225	12 x 25 MM
2032-TT-41320	13 x 20 MM
2032-TT-41323	13 x 23 MM
2032-TT-41325	13 x 25 MM
2032-TT-41420	14 x 20 MM
2032-TT-41423	14 x 23 MM
2032-TT-41425	14 x 25 MM



TLIF CAGE (BANANA CAGE) TYPE - 2

REFERENCE NO. (Titanium)	SIZE: 0°
2031-TT-00825	8 x 25 MM
2031-TT-00830	8 x 30 MM
2031-TT-00925	9 x 25 MM
2031-TT-00930	9 x 30 MM
2031-TT-01025	10 x 25 MM
2031-TT-01030	10 x 30 MM
2031-TT-01125	11 x 25 MM
2031-TT-01130	11 x 30 MM
2031-TT-01225	12 x 25 MM
2031-TT-01230	12 x 30 MM
2031-TT-01325	13 x 25 MM
2031-TT-01330	13 x 30 MM
2031-TT-01425	14 x 25 MM
2031-TT-01430	14 x 30 MM



REFERENCE NO. (Titanium)	SIZE: 4°
2031-TT-40825	8 x 25 MM
2031-TT-40830	8 x 30 MM
2031-TT-40925	9 x 25 MM
2031-TT-40930	9 x 30 MM
2031-TT-41025	10 x 25 MM
2031-TT-41030	10 x 30 MM
2031-TT-41125	11 x 25 MM
2031-TT-41130	11 x 30 MM
2031-TT-41225	12 x 25 MM
2031-TT-41230	12 x 30 MM
2031-TT-41325	13 x 25 MM
2031-TT-41330	13 x 30 MM
2031-TT-41425	14 x 25 MM
2031-TT-41430	14 x 30 MM



TLIF CAGE (KIDNEY CAGE) TYPE - 1



REFERENCE NO. (Titanium)	SIZE: 0°
2041-TT-00825	8 x 25 MM
2041-TT-00830	8 x 30 MM
2041-TT-00925	9 x 25 MM
2041-TT-00930	9 x 30 MM
2041-TT-01025	10 x 25 MM
2041-TT-01030	10 x 30 MM
2041-TT-01125	11 x 25 MM
2041-TT-01130	11 x 30 MM
2041-TT-01225	12 x 25 MM
2041-TT-01230	12 x 30 MM
2041-TT-01325	13 x 25 MM
2041-TT-01330	13 x 30 MM
2041-TT-01425	14 x 25 MM
2041-TT-01430	14 x 30 MM

REFERENCE NO. (Titanium)	SIZE: 4°
2041-TT-40825	8 x 25 MM
2041-TT-40830	8 x 30 MM
2041-TT-40925	9 x 25 MM
2041-TT-40930	9 x 30 MM
2041-TT-41025	10 x 25 MM
2041-TT-41030	10 x 30 MM
2041-TT-41125	11 x 25 MM
2041-TT-41130	11 x 30 MM
2041-TT-41225	12 x 25 MM
2041-TT-41230	12 x 30 MM
2041-TT-41325	13 x 25 MM
2041-TT-41330	13 x 30 MM
2041-TT-41425	14 x 25 MM
2041-TT-41430	14 x 30 MM



MESH CAGE

REFERENCE NO. (Titanium)	SIZE: Ø 12
2029-TT-12015	15 MM
2029-TT-12020	20 MM
2029-TT-12025	25 MM
2029-TT-12030	30 MM
2029-TT-12035	35 MM
2029-TT-12100	100 MM

REFERENCE NO. (Titanium)	SIZE: Ø 14
2029-TT-14020	20 MM
2029-TT-14025	25 MM
2029-TT-14030	30 MM
2029-TT-14035	35 MM
2029-TT-14040	40 MM
2029-TT-14045	45 MM
2029-TT-14050	50 MM
2029-TT-14100	100 MM

REFERENCE NO. (Titanium)	SIZE: Ø 16
2029-TT-16020	20 MM
2029-TT-16025	25 MM
2029-TT-16030	30 MM
2029-TT-16035	35 MM
2029-TT-16040	40 MM
2029-TT-16045	45 MM
2029-TT-16050	50 MM
2029-TT-16060	60 MM
2029-TT-16070	70 MM
2029-TT-16100	100 MM

REFERENCE NO. (Titanium)	SIZE: Ø 18
2029-TT-18020	20 MM
2029-TT-18025	25 MM
2029-TT-18030	30 MM
2029-TT-18035	35 MM
2029-TT-18040	40 MM
2029-TT-18045	45 MM
2029-TT-18050	50 MM
2029-TT-18060	60 MM
2029-TT-18070	70 MM
2029-TT-18100	100 MM





EXPANDABLE CAGE



REFERENCE NO. (Titanium)	SIZE: Ø 10
2030-TT-1025	(20 MM TO 25 MM)
2030-TT-1035	(25 MM TO 35 MM)
2030-TT-1055	(35 MM TO 55 MM)
2030-TT-1075	(55 MM TO 75 MM)

REFERENCE NO. (Titanium)	SIZE: Ø 12
2030-TT-1225	(20 MM TO 25 MM)
2030-TT-1235	(25 MM TO 35 MM)
2030-TT-1255	(35 MM TO 55 MM)
2030-TT-1275	(55 MM TO 75 MM)

REFERENCE NO. (Titanium)	SIZE: Ø 14
2030-TT-1425	(20 MM TO 25 MM)
2030-TT-1435	(25 MM TO 35 MM)
2030-TT-1455	(35 MM TO 55 MM)
2030-TT-1475	(55 MM TO 75 MM)



REFERENCE NO. (Titanium)	SIZE: Ø 16
2030-TT-1625	(20 MM TO 25 MM)
2030-TT-1635	(25 MM TO 35 MM)
2030-TT-1655	(35 MM TO 55 MM)
2030-TT-1675	(55 MM TO 75 MM)

REFERENCE NO. (Titanium)	SIZE: Ø 18
2030-TT-1825	(20 MM TO 25 MM)
2030-TT-1835	(25 MM TO 35 MM)
2030-TT-1855	(35 MM TO 55 MM)
2030-TT-1875	(55 MM TO 75 MM)



Spine

Instrumentations



SPINE INSTRUMENTATIONS




TRAY - 1



TRAY - 2



TRAY - 3




BASE




REFERENCE NO.	INSTRUMENTATIONS
7600-INS-0001	PEDICLE PROBE STRAIGHT
7600-INS-0002	PEDICLE PROBE CURVED
7600-INS-0003	PEDICLE SOUND STRAIGHT
7600-INS-0004	PEDICLE SOUND CURVED
7600-INS-0005	PEDICLE FEELER
7600-INS-0006	4.5 MM BONE TAP WITH QUICK COUPLING
7600-INS-0007	5.5 MM BONE TAP WITH QUICK COUPLING
7600-INS-0008	6.5 MM BONE TAP WITH QUICK COUPLING
7600-INS-0009	MONOAXIAL SCREW DRIVER WITH QUICK COUPLING
7600-INS-0011	POLYAXIAL SCREW DRIVER QUICK COUPLING
7600-INS-0010	ROD TEMPLATE 100 mm
7600-INS-0020	ROD TEMPLATE 200 mm
7600-INS-0035	ROD TEMPLATE 350 mm
7600-INS-0012	UNIVERSAL ROD BENDER (ROLLER TYPE)
7600-INS-0013	COMPRESSOR
7600-INS-0014	DISTRACTOR
7600-INS-0015	PEDICLE AWL WITH STOPPER
7600-INS-0016	ROD PUSHER STRAIGHT
7600-INS-0017	ROD PUSHER CURVED
7600-INS-0018	QUICK COUPLING T - HANDLE
7600-INS-0019	ROD PERSUADER
7600-INS-0021	ROD HOLDING FORCEPS
7600-INS-0022	ROD ROCKING FORCEPS (FORK TYPE)
7600-INS-0023	WISE GRIP
7600-INS-0024	BENDING IRON(PAIR)
7600-INS-0025	INNER HOLDER
7600-INS-0026	FINAL TIGHTNER
7600-INS-0027	ANTI TORQUE KEY
7600-INS-0028	ROD HOLDING FORCEP MULTI PURPOSE
7600-INS-0029	PEDICLE CENTRALISER (4 PCS)
7600-INS-0030	CONNECTOR FINAL TIGHTNER
7600-INS-0031	MEASURING GAUGE FOR CONNECTOR
7600-INS-0032	QUICK COUPLING HANDLE STRAIGHT
7600-INS-0033	CONNECTOR HOLDING FORCEPS
7600-INS-0034	SPINE INSTRUMENT SET WITH BOX
7600-INS-0036	CONNECTOR INNER HOLDER
7600-INS-0037	REDUCTION BREAK OFF (PAIR)
7600-INS-0038	4.0 MM BONE TAP WITH COUPLING
7600-INS-0039	5.0 MM BONE TAP WITH COUPLING
7600-INS-0040	7.5 MM BONE TAP WITH COUPLING
7600-INS-0041	PEDICLE ELEVATOR
7600-INS-0042	WIDE BLADE ELEVATOR
7600-INS-0043	NARROW BLADE ELEVATOR
7600-INS-0044	SUPRA LAMINAR ELEVATOR
7600-INS-0045	ANGLE HOOK ELEVATOR LEFT
7600-INS-0046	ANGLE HOOK ELEVATOR RIGHT
7600-INS-0047	HOOK PUSHER
7600-INS-0048	STAPLE INTRODUCER
7600-INS-0049	SPINE IMPLANT SET WITH BOX (MONO-POLY)
7600-INS-0050	SPINE IMPLANT SET WITH BOX (REDUCTION MONO-POLY)



CERVICAL INSTRUMENTATIONS



TRAY - 1



BASE




REFERENCE NO.	INSTRUMENTATIONS
7800-INS-0001	DOUBLE DRILL & TAP SLEEVE
7800-INS-0002	DRILL BIT 2.5 WITH QUICK COUPLING
7800-INS-0003	QUICK COUPLING HANDLE
7800-INS-0004	SCREW HOLDER
7800-INS-0005	SCREW DRIVER
7800-INS-0006	PLATE HOLDER
7800-INS-0007	CERVICAL BODY DISTRACTOR
7800-INS-0008	DISTRACTOR PIN
7800-INS-0009	CERVICAL RETRACTOR(LARGE)
7800-INS-0010	CERVICAL RETRACTOR(SMALL)
7800-INS-0011	LOCK DRIVER
7800-INS-0012	PIN DRIVER
7800-INS-0013	PIN
7800-INS-0014	TAP WITH T - HANDLE
7800-INS-0015	UNIVERSAL PLATE BENDER
7000-INS-0028	ANTERIOR CERVICAL PLATE INSTRUMENT BOX WITH IMPLANT & SCREW TRAY
7800-INS-0016	BONE CURETTES
7800-INS-0017	AWL
7800-INS-0018	RONGER FOR CERVICAL



TRAY - 1



TRAY - 2




BASE



REFERENCE NO.	INSTRUMENTATIONS
8200-INS-0001	BONE GRAFTING BLOCK
8200-INS-0002	BONE GRAFTING PUNCH
8200-INS-0003	CAGE INTRODUCER (PLIF)
8200-INS-0004	CAGE INTRODUCER (TLIF)
8200-INS-0005	NERVE ROOT RETRACTOR (LARGE)
8200-INS-0006	NERVE ROOT RETRACTOR (SMALL)
8200-INS-0007	REAMER/SHAVER WITH QUICK COUPLING 8.0 MM
8200-INS-0008	REAMER/SHAVER WITH QUICK COUPLING 9.0 MM
8200-INS-0009	REAMER/SHAVER WITH QUICK COUPLING 10.0 MM
8200-INS-0010	REAMER/SHAVER WITH QUICK COUPLING 11.0 MM
8200-INS-0011	REAMER/SHAVER WITH QUICK COUPLING 12.0 MM
8200-INS-0012	REAMER/SHAVER WITH QUICK COUPLING 13.0 MM
8200-INS-0013	QUICK COUPLING HANDLE STRAIGHT
8200-INS-0014	TRIAL FOR PLIF CAGE 8.0 MM
8200-INS-0015	TRIAL FOR PLIF CAGE 9.0 MM
8200-INS-0016	TRIAL FOR PLIF CAGE 10.0 MM
8200-INS-0017	TRIAL FOR PLIF CAGE 11.0 MM
8200-INS-0018	TRIAL FOR PLIF CAGE 12.0 MM
8200-INS-0019	TRIAL FOR PLIF CAGE 13.0 MM
8200-INS-0020	TRIAL FOR TLIF CAGE 8.0 MM
8200-INS-0021	TRIAL FOR TLIF CAGE 9.0 MM
8200-INS-0022	TRIAL FOR TLIF CAGE 10.0 MM
8200-INS-0023	TRIAL FOR TLIF CAGE 11.0 MM
8200-INS-0024	TRIAL FOR TLIF CAGE 12.0 MM
8200-INS-0025	TRIAL FOR TLIF CAGE 13.0 MM
8200-INS-0026	CAGE PUSHER CURVED
8200-INS-0027	CUP CURETTE CURVED LEFT
8200-INS-0028	CUP CURETTE CURVED RIGHT
8200-INS-0029	CUP CURETTE STRAIGHT
8200-INS-0030	BONE RASP STRAIGHT
8200-INS-0031	BONE RASP CURVED
8200-INS-0032	BOX CURETTE STRAIGHT
8200-INS-0033	BOX CURETTE CURVED
8200-INS-0034	RONGER FOR SPINE
8200-INS-0035	CAGE HOLDING FORCEPS
7000-INS-0029	TLIF/PLIF CAGE INSTRUMENT AND IMPLANT BOX

CERTIFICATES



Notified Body 1023
INSTITUTE FOR TESTING AND CERTIFICATION, Inc.,
trida Tomáše Bati 280, Louky, 753 02 Zlín, Czech Republic

EC Certificate - Full Quality Assurance System
No. 18 0309 QS/NB


The quality system of manufacturer


Griportho Surgicals Pvt. Ltd
R.S. NO. 99/2, P-2, VILLAGE & POST KARMASIYA KHEDA,
RASULABAD VITTOJ ROAD, TALUKA WAGHODIA, DIST.
VADODARA – 391510 GUJARAT, INDIA

has been certified as meeting the requirements of
Directive 93/42/EEC
on medical devices, Annex II excluding (4)
for the following product category(ies):
Orthopaedic Implants

The Notified Body No. 1023 declares that the aforementioned manufacturer has implemented a quality assurance system for design, manufacture and final inspection of the respective devices / device categories in accordance with MDD Annex II. This quality assurance system conforms to the requirements of this Directive and is subjected to periodical surveillance. For placing on the market of Class III devices covered by this certificate, an EC Design-Examination Certificate according to Annex II (Section 4) is required.

Valid from: 2018-09-11
Valid until: 2023-09-10
First issued: 2018-09-11
Revision: -


Date: 2018-09-11


Mr. Jit Heš
Representative of the Notified Body No. 1023

Page 1/1

| evaluate | assess | certify |  assessors of quality

Certificate No. MD-QMS/91/R/2395

This is to certify that

Griportho Surgicals Pvt. Ltd.

Corporate Office / Head Office : 310, Sanket Avenue, Opp. Ambe Vidyalaya,
Sama - Savli Road, Vadodara - 390 024 Gujarat, India
(a list of certified workplaces is contained in Annex which is inseparable part of this Certificate)

has been found to conform to the requirements of
Medical Devices - Quality Management System Standard

ISO 13485:2016

This certificate is valid for the following scope :

Manufacturer, Exporter and Supplier of Orthopedic Implants, Instruments and related Accessories.

Initial Certification : 20/05/2017
Valid until : 19/05/2020

  
ACCREDITED
Management Systems
Certification Body
Authorised Signatory

The authority of this document may be affirmed by writing to certification@zqa.com. List of Effluent of conditions as set out in the Certification Agreement may be used for this purpose. Origin or indication of the product or appearance of this document is intended to be permanent in the latter part of the document.

Accreditation Body: International Accreditation Service, Inc.
2800 River Street, Suite 100, Brea, California 92617-1722 U.S.A.
www.iasinc.org

Issuing Authority: Zentib Quality Assessors Pvt. Ltd. - MSC0991
100, 99/2, P-2, Karmasiya Kheda, Rasulabad - Vitoj Road, Taluka Waghodia, Dist. Vadodara, Gujarat - 391510, India
www.zentib.com

FORM 28 FRESH LICENCE
(See Rule 76)

Licence to manufacture for sale (or for distribution) of drugs specified in
Schedule C, C(I) excluding those specified in schedule X

Number of licence and date of issue No **G/28/1500 Dated -04/02/2015**

1. **DIRECTORS** of **M/s. GRIPORTHO SURGICALS PVT LTD** is hereby granted licence to manufacture at the premises situated at the **R.S. NO : 99/2, P-2 - OPP: SHRI KABIR SAHEB VIDYA MANDIR, VILLAGE: & POST: KARMASIYA KHEDA, RASULABAD - VITTOJ ROAD, TA: WAGHODIA, DIST. - VADODARA-391 510.**

the following drugs being drugs specified in schedules C & C(I) Excluding those specified in sch. X to the Drugs and Cosmetics Rules 1945.


Names of Drugs : **As Per List Approved & Annexed**

2. Names of approved competent technical staff: **As Per List Approved & Annexed**


3. The licence authorizes the sale by way of wholesale dealing by the licensee and storage for sale by the licensee of the drugs manufactured under the licence, subject to the conditions applicable to licences for sale.

4. The licence will be in force from **04/02/2015 To 03/02/2020**

5. The licence is subject to the conditions stated below and to such other conditions as may be specified in the rules for the time being in force under the Drugs and Cosmetics Act, 1940.


Commissioner
Food & Drugs Control Administration
Gujarat State.

Date: **04/02/2015**


Central Licensing Approving Authority
Dr. S. S. P. P. P.
Drug Control Administration (India)
Dist. General of
Ministry of Health
P.O. Bhavnagar, Kolia Road, T.V.T.

Conditions of Licence

- This licence and any certificate issued in force shall be kept on the approved premises and shall be produced at the request of an Inspector appointed under the Drugs and Cosmetics Act, 1940.
- If the licensee wishes to undertake during the currency of the licence the manufacture of any drug specified in Schedule C and C(I) excluding those specified in Schedule X not included above, he should apply to the Licensing Authority for the necessary endorsement as provided in Rule 75(3). This licence will be deemed to extend to the items so endorsed.
- Any change in the (competent technical staff) shall be forthwith reported to the Licensing Authority.
- The licensee shall inform the Licensing Authority in writing in the event of any change in the constitution of the firm operating under the licence. If there is any change in the constitution of the firm takes place, the current licence shall be deemed to be valid for a maximum period of three months from the date on which the change takes place unless, in the meantime, a fresh licence has been taken from the Licensing Authority in the name of the firm with the changed Constitution.

FDCA/GUJARAT/28219

IMPORTANT NOTES



IMPORTANT NOTES







Corporate Office

310, Sanket Avenue, Opp. Ambe Vidyalaya
Sama - Savli Road, Vadodara - 390024
Gujarat - India.
Tel: +91 - 265 2711162 / + 91 99740 88425
Email: marketing@griportho.com

Registered Office & Works

R.S No. 99/2, P-2, Village & Post
Karmasia Kheda Rasulabad - Vitoj Road,
Taluka Waghodia, Dist. Vadodara - 391510
Gujarat - India.