

TS,TC,PD,SK,DB

Sample ID: BIA251201S0054
Strain: HL-14
Harvest Lot:
Matrix: Plant
Type: Flower - Cured
Sample Size:
Lot#:

Produced:
Collected:
Received: 12/03/2025
Completed: 12/12/2025
Batch#:

Client
Green Castle
Lic. #
853 RT 15W
Johnson, VT 05656



Summary

| Test | Date Tested | Result |
|------------|-------------|------------|
| Sample | | Complete |
| Moisture | 12/03/2025 | Not Tested |
| Pesticides | 12/08/2025 | Complete |




Luke Emerson-Mason
Laboratory Director
12/12/2025

Confident LIMS
All Rights Reserved
coa.support@confidentlims.com
(866) 506-5866
www.confidentlims.com



TS,TC,PD,SK,DB

 Sample ID: BIA251201S0054
 Strain: HL-14
 Harvest Lot:
 Matrix: Plant
 Type: Flower - Cured
 Sample Size:
 Lot#:

 Produced:
 Collected:
 Received: 12/03/2025
 Completed: 12/12/2025
 Batch#:

 Client
Green Castle
 Lic. #
 853 RT 15W
 Johnson, VT 05656

Pesticides

Completed

| Category 1 Pesticides | LOD | LOQ | Results |
|-----------------------|---------|--------|---------|
| | PPM | PPM | PPM |
| Chlorpyrifos | 0.0003 | 0.0010 | ND |
| Imazalil | 0.0003 | 0.0010 | ND |
| Category 2 Pesticides | LOD | LOQ | Results |
| | PPM | PPM | PPM |
| Abamectin | 0.0003 | 0.0010 | ND |
| Acephate | 0.001 | 0.0050 | ND |
| Acequinocyl | 0.0003 | 0.0010 | ND |
| Azoxystrobin | 0.00005 | 0.0010 | ND |
| Bifenazate | 0.0001 | 0.0010 | ND |
| Bifenthrin | 0.0001 | 0.0010 | ND |
| Carbaryl | 0.0001 | 0.0010 | ND |
| Cypermethrin | 0.001 | 0.0050 | ND |
| Etoxazole | 0.0001 | 0.0010 | ND |
| Imidacloprid | 0.00005 | 0.0010 | ND |
| Myclobutanil | 0.0001 | 0.0010 | ND |
| Pyrethrins | 0.001 | 0.0050 | ND |
| Spinosyn A | 0.0001 | 0.0010 | ND |
| Spinosyn D | 0.0003 | 0.0010 | ND |

Analyst: 056

Pesticides Methodology: Liquid Chromatography with Tandem Mass Spectrometry using PerkinElme QSight® LX50 UHPLC and QSight 220 Mass Spectrometer

LOQ = The lowest quantity this method can reliably quantify. Any pesticides or mycotoxins that were not quantifiable are less than the stated LOQ (<LOQ).

ppm = parts per million

All moisture and water activity analysis is determined by dewpoint measurement using an AQUALAB water activity meter.

ND = Not Detected (<LOD)




 Luke Emerson-Mason
 Laboratory Director
 12/12/2025

 Confident LIMS
 All Rights Reserved
coa.support@confidentlims.com
 (866) 506-5866
www.confidentlims.com
