



# Measuring for Carts

(No need to fill out this page, this is just an example)

## Your Pet's Measurements

Fold of the Flank to Ground

  
inches ▼

Approximate Weight

  
pounds ▼


Main Breed

Pembroke Welsh Corgi ▼

## How to Measure

To measure the height of a dog for our Walkin' Wheels® Dog Wheelchair, with the dog lying down, measure from the ground to the **Fold of the Flank** as shown here. The Fold of the Flank (FoF) is a thin piece of skin that connects the dog's upper leg to his body. Look at the webbing between your thumb and index finger. This skin is very similar to the fold of the flank (FoF).

\* When measuring the fold of the flank of your dog's leg, do not pull the leg tight. Leave a little natural bend.



*Helpful Tip: You can also use a piece of cardboard, mark the cardboard and then measure your marks on the cardboard*

1. Measure the fold of flank of your pet laying down to the toe pad, as shown here. **Do not pull the leg tight; leave some natural bend.**

The Fold of the Flank (FOF) is a thin piece of skin that connects the dog's upper rear leg to his body. Look at the webbing between your thumb and index finger. This skin is very similar to the fold of flank (FOF).

2. Measure your pet's armpit to rump (ATR).

