




Repositioning

with Centrobed

MADE IN UK

Preventing injuries to both users and carers.

Treating pressure ulcers costs the NHS more than £1.4 million everyday.



The right equipment at the right time saves

Improves quality of life.

Centrobed: Life changing turning beds



Our unique turning beds are designed to change the position of a client in a relaxed and secure way. Using our bed can help with posture related pain, spinal surgery recovery as well as using the bed for torso mobility to reduce the effects of chest infections. Used in conjunction with our Auto Box you can reduce the cost of turning client's day and night.

Considerations when choosing a bed



Prescribing a bed for a disabled person is a personalized process that should consider the client individual needs and circumstances. It's essential to involve healthcare experts who specialise in care and assistive devices to make informed recommendations.

Medical condition

Bed Frame

Mobility & Transfers

Long-term Needs

Mattress Selection

Width & Length

Reduce Care Support & Financial Outcomes



- 24hr Automatic re-positioning using the auto box.
- Personalised programmable turning.
- Improved comfort using multi positional Arctic turning bed.
- Postural Support
- Hands free turning
- Improved pressure care

Justifications Evidence Based Practice



- Determine the need and individuals capability.
- Consider and trial bed owned by the organisation
- Research alternatives bed and Journal reviews.
- Work with Colleagues, individual and family to evaluate the bed.
- Does the bed improve care and improve financial outcomes?
- Does the bed reduce the carer support?
- Complete your Clinical Reasoning form.
- Trial the bed

Clinical Reasoning (Balanced decision making)



- Identify the underlying risk of non-provision of the bed.
- Identify the control measures you have put in place.
- Can you reduce the risk to reasonably practicable level?
- Have you identified how the equipment will meet individuals need?
- Have you included your cost analysis and justification in the document?




Book your
Presentation / Training Day today
Call us on 01233 635353
or email sales@centrobed.com

Prescribing a bed for an Individual



Re-positioning with Centrobed.




Cost Saving Calculation

Hours Saved	Weekly Savings	Annual Savings
1	£161	£8,372
2	£322	£16,744
3	£483	£25,116
4	£644	£33,488

Based on £23 per hour

Amazon – Adult Size Re-positioning bed



Amazon will turn someone from side to side throughout the night, the Amazon turning bed also has the added advantage of growing in length, 78", 84", 90". Irrespective of the height of the person the Amazon turning bed will fit and profile in the correct places.

Aegean Amputee Bed



The Aegean Amputee turning bed has two section profiling with the addition of turning. The bespoke design prevents the user from sliding down the bed whilst being able to help with pressure relief using the turning feature.

Caribbean Children Turning Bed



Fully profiling bed that turns, has independent leg movement and will grow with the child. Wide range of movements available that can help children positioning and find comfort with minimal intervention.

Arctic Turning Bed



Multi-functional turning & full profiling bed that enables hands-free turning with its automatic turning feature, will turn someone from side to side throughout the night, the programme has the ability to accommodate individual needs.

Tungsten Elite '80 Stone Bed'



A Bariatric bed with unparalleled strength and support, the Tungsten Elite is a 80 stone turning bed that keeps the client stable and safe while providing pressure relief and reduce the risk of injury, with the Tungsten Elite there will be no more broken casters or bent framework.

Hire - Weighing Service - Trial - Bed Assessments

www.centrobed.com

t: 01233 635353

e: sales@centrobed.com

How can Centrobed play a role in bed provision?

- Work with prescriber to assess for the appropriate bed
- Seamless provision
- Working with the current equipment
- Rental provision
- Individual or group training
- Plan for ongoing needs
- Centrobed are flexible and work to overcome barriers
- Trials and home assessments
- Help to provide dignity enabled care
- Work with prescriber to provide the right equipment at the right time

Equipment Risk Assessment

Underpins all other decisions and actions that are taken to provide a reduced risk environment. It will involve either a single person, or multi-agency/collaborative agency approach. It will determine the equipment needs.




- Determining patient weight
- Communication between all parties including bed supplier
- Are adaptations required before the bed is installed?
- Number of personnel involved
- Enables patient led interventions
- What are the manual handling risks and how can they be minimised?
- Does the task of delivering the bed involve a 90-degree manoeuvre, if so, is there an alternative?
- If the road is a main road and/or near a school, you need to minimise the exploitation of the patient (dignity)
- What other equipment is required and is it compatible to the bed being delivered.

Manual Handling Regulation state that you need to:

- Avoid** - Hazardous Manual Handling
- Access** - Make suitable and sufficient assessments of hazards using the ergonomic approach.
- Consider** – Nature of the task, THE LOAD, working environment, the individual’s capability and re-positioning needs.
- Reduce the risk of injury** - Risk Assessments, provision and maintenance of equipment, Training, supervision and changes to the environment.
- Regular reviews** - Risk Assessments should be regularly reviewed when changes are made or happen.


Specification

Product	Width	Length	Start Height	Maximum Patient Weight
Arctic Turning Bed	36", 42", 48"	78" or bespoke	16"	158kg
Amazon Re-positioning Bed	36", 42", 48", 54", 60"	78", 84", 92"	16"	158kg
Aegean Amputee Bed	36"	48" or bespoke	16"	158kg
Caribbean Children Bed	36"	60", 66", 72"	16"	114kg
Tungsten Elite Turning Bed	54"	78"	16"	510kg




Pressure Areas

Regular re-positioning can reduce pressure sores.
Areas at risk of pressure sores.



Centrobed Turning Beds



Cost Saving

Hours Saved	Daily Cost	Weekly Savings	Annual Savings
2	£46	£322	£16,744
3	£69	£483	£25,116
4	£92	£644	£33,488
6	£138	£966	£50,232
8	£184	£1288	£66,976
10	£230	£1610	£83,720
12	£276	£1932	£100,464



Case Study

The Caribbean Turning Bed has literally transformed our lives. Previously we were getting up multiple times during the night to re-position him, and that was severely impacting both his and our well being. From the first time he slept in the Caribbean, he slept all night with no requests to be moved, and has continued to do so.

It's given him some independence back in terms of being able to manually re-position himself by using the handheld controls, as well as peace for us at night as we're able to pre-set the automatic turning function at bedtime.

Your kindness and understanding really was appreciated - as well as your patience as we went though the funding process for this essential piece of equipment for our son. I have recommended you (and the bed) to all the families I know who have a child with a similar condition.



Re-positioning helps to prevent Pressure Ulcers

Re-positioning with ease



Help with re-positioning at ease in one of our CentroBed.

Re-positioning with ease with patients, involves the careful and smooth adjustment of their position to ensure their comfort and safety. With our range of our beds that allows turning / re-positioning as well as profiling. An automatic box can be fitted to turn the bed at regular intervals.


These beds can be used for hospitals, clinics, nursing homes, and out in the home community where patients may have limited mobility or be bedridden.

- Reduces Friction

- Reduces Pressure Sores

- Comfort


- Improve Patient Experience



Fact

- Treating pressure ulcers costs the NHS more than £1.4 million every day.
- 180,000 new patients with pressure sores every year
- Over 700,000 patients are affected
- High risk patients are people over 70 or people confined to a bed for long periods of time.

1



Probable Causes

- Lack of mobility
- Lack of fluid (dehydration)
- Moist skin
- Sliding down the bed (Shearing)
- Weight loss
- Poor diet
- Friction or rubbing
- Other medical conditions such as diabetes

2

How can Centrobed help?

- Regular re-positioning
- Correct products and equipment at the correct time
- Keep patient moving
- Availability of products and equipment
- Cost Saving equipment
- Option to hire

4

Prevention

- Regular inspection and monitoring
- Keep patient moving
- Educate client and carer
- Awareness to all
- Correct products and equipment at the correct time
- Availability of staff, product and equipment

3

Long Term

- Cost Saving
- Independence
- Improve quality of life for patient & the family
- Reduced pressure sores
- Reduced sleepless nights
- Reduced hospital stays
- Reduction in medication & carers*

5

Hire - Weighing Service - Trial - Bed Assesments