

CaviWipes™

Easy as

1 2 3



Improved lid design dispenses wipes with precision. Easy to open and close securely.

CaviWipes HP

Features (S.A.F.E.)

Simple

1-minute universal contact time.
1-step cleaning/disinfection.¹

Alcohol-free

Hydrogen peroxide suitable for sensitive patient populations.

Future-proof

Fully qualified for the EPA Emerging Viral Pathogen Claim for all virus types (enveloped, large and small non-enveloped viruses) to meet current and potentially future infection prevention needs.²

Effective

Disinfection against 67 pathogens including multidrug-resistant organisms (MDRO).

Part No.	Description
16-1100	CaviWipes HP (6" x 6.75") 160 towelettes per canister
16-1150	CaviWipes HP XL (9" x 12") 65 towelettes per canister

CaviWipes 2.0

Features

Simple

2-minute universal contact time.
1-step cleaning/disinfection.¹

Agreeable

Neutral pH, low alcohol/quat formulation agreeable with healthcare surfaces and materials.

Future-proof

Fully qualified for the EPA Emerging Viral Pathogen Claim for all virus types (enveloped, large and small non-enveloped viruses) to meet current and potentially future infection prevention needs.²

Effective

Disinfection against 60 pathogens including MDRO.

Part No.	Description
14-1100	CaviWipes 2.0 (6" x 6.75") 160 towelettes per canister
14-1150	CaviWipes 2.0 XL (9" x 12") 65 towelettes per canister

CaviWipes Bleach

Features

Simple

3-minute universal contact time.
1-step cleaning/disinfection.¹

Alcohol-free

1:10 bleach dilution ratio recommended by the CDC.

Future-proof

Fully qualified for the EPA Emerging Viral Pathogen Claim for all virus types (enveloped, large and small non-enveloped viruses) to meet current and potentially future infection prevention needs.²

Effective

Disinfection against 43 pathogens including MDRO and *Clostridioides difficile* spores.

Part No.	Description
13-9100	CaviWipes Bleach (6"x10.5") 90 towelettes per canister
13-1175	CaviWipes Bracket Fits all CaviWipes canisters

To learn more scan the QR code or visit [metrex.com](https://www.metrex.com)



1. Always follow product label for special precleaning and PPE instructions.
2. <https://www.epa.gov/pesticide-registration/emerging-viral-pathogen-guidance-and-status-antimicrobial-pesticides>. Accessed 02.23.23.