

ABOUT SWIPESENSE

SwipeSense was started by two people whose families were impacted by hospital acquired infections. Their mission is to help reduce these infections in all hospitals.

Questions? Feedback?

We want to hear it! Email us at support@swipesense.com.

HOW DOES THE SYSTEM KNOW I CLEANED MY HANDS?



The badge is linked with your employee information. When you clean your hands the drip-guard hygiene sensor recognizes your badge and a green light will flash.

WHERE IS THE IDEAL PLACE TO WEAR MY BADGE? HOW DO I KNOW IF IT'S WORKING?



Wear the badge above waist level and visible at all times. To avoid damaging your badge, avoid high frequency hospital equipment (MRIs) and excess water. If you are concerned that your badge is no longer working, talk to your manager.

WHAT IF I'M WEARING GLOVES? OR IF THERE IS NO PATIENT IN THE ROOM?



You are still expected to clean in / clean out to ensure that the gloves are touched by clean hands and that the patient room remains clean. If you leave and re-enter the same patient room, be sure to wash your hands again.

CAN I USE ANY DISPENSER?




Yes, any dispenser can be used. Soap, sanitizer, or gel are all acceptable.




Hand Hygiene. Redefined.

Quick Start



HOW DOES THE SWIPESENSE SYSTEM WORK?


Badges communicate to the SwipeSense® System when a user enters and exits a patient room and records whether or not the user washed his or her hands.



WHAT IS EXPECTED?

Clean in / Clean out. Always wash or sanitize your hands upon entering and exiting a patient room, even if a patient is not present.

The system provides a time window of **30 seconds** upon entering and 30 seconds upon exiting a patient room for you to clean your hands.



WHERE TO CLEAN?

You can clean before entering a patient room or once inside. The same rule applies for exiting. Any dispenser can be used.

WHAT IS THE SYSTEM TRACKING?

Only hand hygiene events and opportunities associated with a patient room visit are reported. No other data is reported.

+100,000

lives could be saved every year by eliminating hospital acquired infections.

THE END GOAL.



You are the primary defense between hospital acquired infections and the patient. The most effective method to eliminate these infections is frequent hand washing. The SwipeSense® System helps hospitals quantify how hand hygiene is impacting infections.