



PRADA TEST SIEVES

Sieve Construction:

The frame is made of spun Stainless Steel frame without any joint in it. A clear rigid design of minimum weight and free from corners where the material could be lodged. Stainless Steel 18/8 quality sieve cloth is used for highest performance. An essential feature of Test Sieves is that the mesh should be flat, and uniform and with permanent tension round the whole circumferences. The same is fully observed in PRADA Test Sieves

Micro Test Sieves:

Micro Test Sieves are available in one piece brass frame of 3, 4 and 6 inch diameter with suitable Lid and Receiver.

Coarse and Fine Test Sieves:

Coarse and Fine Test Sieves 8, 12 and 18 Inch in diameter with Lid and Receiver assembly are also available.

Wet Washing Sieves:

For samples that need to be separated with the help of wet washing. Available in 8 Inch diameter by 4 or 8 Inch deep with brass frames. A complete range of aperture sizes are available. There is no national or international specifications for wet washing sieves PRADA sieves will nest with sieves of similar diameter and can be fitted with mesh to any desired specification.

Lid and Receivers:

Lids and receiving pans are available in stainless steel or brass in 8 and 12 Inch diameter.

B.S.S (410/1986)	
mm	microns
125	850
106	710
90	600
75.0	500
63.0	425
53.0	355
50.0	300
45.0	250
37.5	212
31.5	180
28.0	150
26.5	125
25.0	106
22.4	90
20.0	75
19.0	63
16.0	53
14.0	45
13.2	38
12.5	32
11.2	25
10.0	20
9.5	
8.00	
6.70	
6.30	
5.60	
5.00	
4.75	
4.00	
3.35	
3.15	
2.80	
2.36	
2.00	
1.70	
1.40	
1.18	
1.00	

SYARIKAT SAINTIFIK JAYA

NO 69 Jalan Penerbit UI/43,
 Temasya Industrial Park Glenmarie
 40150 Shah Alam
 Tel : 603- 5569 0758/0759
 HP : 018 274 0955
 Fax : 603- 5569 0763



SYARIKAT SAINTIFIK JAYA

<https://sites.google.com/view/syarikatsaintifikjaya>



PARADA Test Sieves Available Sizes:

mm

125mm	31.5mm	13.2mm
106mm	28mm	12.5mm
90mm	26.5mm	11.2mm
75mm	25mm	10mm
63mm	22.4mm	9.5mm
53mm	20mm	8mm
50mm	19mm	6.7mm
45mm	16mm	6.3mm
37.5mm	14mm	5.6 mm (No.3 ½)

Micro Test Sieves:

5 mm	2.36mm (No.8)
4.75mm (No.4)	2mm (No.10)
4mm (No.5)	1.7mm (No.12)
3.35 mm (No.6)	1.4mm (No.14)
3.15 mm	1.18mm (No.16)
2.8 mm (No.7)	1mm (No.18)

Coarse and Fine Test Sieves:

µm

850µm (No.20)	212µm (No.70)	53µm (No.270)
710µm (No.25)	180µm (No.80)	45µm (No.325)
600µm (No.30)	150µm (No.100)	38µm (No.400)
500µm (No.35)	125µm (No.120)	32µm (No.450)
425µm (No.40)	106µm (No.140)	25µm (No.500)
355µm (No.45)	90 µm (No.170)	20µm (No.635)
300µm (No.50)	75 µm (No.200)	15µm
250µm (No.60)	63 µm (No.230)	10µm

B.S.S (410/1986)	
mm	microns
125	850
106	710
90	600
75.0	500
63.0	425
53.0	355
50.0	300
45.0	250
37.5	212
31.5	180
28.0	150
26.5	125
25.0	106
22.4	90
20.0	75
19.0	63
16.0	53
14.0	45
13.2	38
12.5	32
11.2	25
10.0	20
9.5	
8.00	
6.70	
6.30	
5.60	
5.00	
4.75	
4.00	
3.35	
3.15	
2.80	
2.36	
2.00	
1.70	
1.40	
1.18	
1.00	

SYARIKAT SAINTIFIK JAYA

NO 69 Jalan Penerbit UI/43,
 Temasya Industrial Park Glenmarie
 40150 Shah Alam
 Tel : 603- 5569 0758/0759
 HP : 018 274 0955
 Fax : 603- 5569 0763



SYARIKAT SAINTIFIK JAYA

<https://sites.google.com/view/syarikatsaintifikjaya>