

A new and advanced amalgam separator just in time for the EPA Deadline

2 MONTHS LEFT
for EPA Amalgam compliance



The ECO II+ Amalgam Separator

- Higher flow rate supports 1-10 chairs (10-20 chairs in tandem configuration)
- Compatible with dry and wet vacs and without disruptive suction loss
- Easier, quicker install
- No maintenance required & lifetime warranty
- Larger inner/outer outlets with multi-adaptors 1.5" and 2"
- Free compliance management software
- Dimensions (HxWxD): 12.2 x 7.7 x 8.4 in
- Collection container volume: 4.5 lb. (700 ml) amalgam waste
- Visit pureway.com/ecoui-amalgam-separator for more info



51001 - Single Unit
52001 - Tandem Unit
53001 - Recycling Unit



Amalgam separator is installed on the main vacuum line.



Solids are collected. When full, replacement canister is ordered.



Full canister is placed in pre-paid return shipping box and picked up via UPS.



Mercury is extracted and recycled at an approved facility.




Recycling certificate is automatically sent to the dental office.



(727) 735.1895 or info@coynemed.com

Manage all your amalgam waste disposal and recycling with just one click



**50%
SAVINGS
OVER COMPETING
SYSTEMS**

The EPA amalgam rule doesn't just require the installation and maintenance of an amalgam separator. It also requires the disposal and proper recycling of solid amalgam, monthly checks on your amalgam separator and proper record keeping. Thankfully, Pureway's proprietary management software

**drastically simplifies compliance
while reducing overhead costs.**

This electronic compliance platform documents, tracks and communicates key compliance alerts and sends automated recycling certificates as required by state and federal regulations. With PureWay, office managers can now stay compliant with the click of a button.

